

TECHNICAL INFORMATION
SD-37

Characteristics and use of

Charge amplifier




Table of contents

1. General description	3
2. Principle of operation	3
3. Gain	3
3-1 Amplifier	3
3-2 Amplifier with detector	4
4. Characteristics	4
4-1 Open-loop gain	4
4-2 Noise	5
5. Applications	7
5-1 Gamma ray measurement (Direct detection using a PIN photodiode)	7
5-2 Power and stability measurements from lasers	7
6. Specifications	8
7. Precautions for handling charge amplifier	10

Characteristics and use of Charge amplifier

1. General description


When a semiconductor detector such as Si is used for the measurement of soft X-rays and low to high-energy gamma rays, the output signal is a weak charge pulse having a pulse width of several tens of nanoseconds. As the detector element itself is a capacitive device, its impedance is very high. Therefore, the performance of the preamplifier to be connected, must be taken into consideration when amplifying this output signal.

In such applications, operational amplifier mode integrators using feedback capacitance are commonly used. As these amplifiers have high input impedance, they integrate weak charge pulses and convert them into voltage pulses for amplification then provide a low-impedance output.

Because of this operation, this type of amplifier is called a "charge amplifier". The first stage of a charge amplifier is usually a low-noise FET and its open-loop gain is set sufficiently high so that the amplification is not influenced by the detector capacitance. The output stage is a low-impedance buffer so as to drive an external stage which is connected using a long cable.

2. Principle of operation

Figure 2-1 Principle of operation



When soft X-rays or gamma rays strike for example a Si semiconductor detector, signal charge pulses Q_s are generated, with an amplitude according to the particle energy. Due to this charge generation, the input-end potential of the charge amplifier rises and at the same time, a potential with reverse polarity appears at the output end. However because the amplifier's open-loop gain is sufficiently large, the output-end potential works through the feedback loop so as to make the input-end potential zero instantaneously.

As a result, the signal charge pulses Q_s are all integrated to the feedback capacitance C_f and then output as voltage pulses $e_{out}(t)$. At this point, since the feedback resistance R_f for direct current is connected in parallel to the feedback capacitance C_f , the output becomes voltage pulses that slowly discharge with the time constant determined by $\tau=C_f \cdot R_f$.

If a detector provides a constant charge generation over a time interval $t=0$ to t_0 , the output signal charge Q_s is given by the following equation using the Laplace transform.

$$Q_s(S) = Q_s \left(\frac{1}{S} - \frac{e^{-St_0}}{S} \right) \quad \dots \dots (2-1)$$

Similarly, the transmission coefficient $T(S)$ is given by

$$T(S) = -\frac{1}{C_f} \cdot \frac{1}{S + \frac{1}{\tau}} \quad \dots \dots (2-2)$$

Thus the output voltage $V(S)$ is expressed using Eqs (2-1) and (2-2) as follows:

$$\begin{aligned} V(S) &= Q_s(S) \cdot T(S) = Q_s \left(\frac{1}{S} - \frac{e^{-St_0}}{S} \right) \cdot \frac{1}{C_f} \cdot \frac{1}{S + \frac{1}{\tau}} \\ &= -\frac{Q_s}{C_f} \left(\frac{1}{S} \cdot \frac{1}{S + \frac{1}{\tau}} - \frac{e^{-St_0}}{S} \cdot \frac{1}{S + \frac{1}{\tau}} \right) \end{aligned} \quad \dots \dots (2-3)$$

As a result, the output voltage pulse $e_{out}(t)$ is given by

$$\begin{aligned} e_{out}(t) &= -\frac{Q_s}{C_f} \cdot \frac{1 - e^{-t/\tau}}{t_0/\tau} \quad 0 \leq t < t_0 \\ &= -\frac{Q_s}{C_f} \cdot \frac{(e^{t_0/\tau} - 1)}{t_0/\tau} e^{-t/\tau} \quad t_0 \leq t \end{aligned} \quad \dots \dots (2-4)$$

Because generally $t_0 \ll \tau$, Eq (2-4) can be simplified as follows:

$$e_{out}(t) = -\frac{Q}{C_f} e^{-t/\tau} \quad \dots \dots (2-5)$$

As can be seen from Eq (2-5), the signal charge pulses Q_s are converted into voltage pulses with amplitude $V_{out} = -\frac{Q_s}{C_f}$, which is damped with time constant $\tau=C_f \cdot R_f$.

3. Gain

The gain of charge amplifier is given in one of two ways: the gain for amplifier or the gain for a detector/amplifier combination.

3-1 Amplifier

The gain of amplifier G_c , referred to also as "charge gain", is given by the following equation:

$$G_c = \frac{V_{out}}{Q_s} \left(= \frac{1}{C_f} \right) \quad [V/\text{coulomb}] \text{ or } [V/\text{pico coulomb}] \quad \dots \dots (3-1)$$

3-2 Amplifier with detector

In this case we usually use the term called “sensitivity” rather than “gain”. Sensitivity is expressed in the output voltage (mV) per one MeV of particle energy irradiated onto a detector. The amplitude of the signal charge obtained with a semiconductor detector is determined by the input particle energy such as soft X-rays and gamma rays and also by the material of the semiconductor.

$$Q_s = \frac{E \cdot e^-}{\epsilon} \text{ (coulomb) or (pico coulomb)} \quad \dots \dots (3-2)$$

E: Particle energy (MeV)

e^- : Elementary charge 1.6×10^{-19} (coulomb)

ϵ : Energy required to create one electron/hole pair.

For example with silicon, Q_s ranges from 3.62 eV (at 300 K) to 3.71 eV (at 77 K).

Thus, from Eqs (3-1) and (3-2), sensitivity is given by

$$\begin{aligned} R_s &= \frac{V_{out}}{E} = \frac{\frac{Q_s}{C_f}}{Q_s \cdot \frac{\epsilon}{e^-}} \\ &= \frac{e^-}{C_f} \cdot \frac{1}{\epsilon} \text{ (mV/MeV)} \quad \dots \dots (3-3) \end{aligned}$$

For example when using a Si detector and an H4083 charge amplifier ($C_f=2$ pF), the sensitivity at room temperatures R_s becomes

$$\begin{aligned} R_s &= \frac{e^-}{C_f} \cdot \frac{1}{\epsilon} \\ &= \frac{1.6 \times 10^{-19}}{2 \times 10^{-12}} \cdot \frac{1}{3.62} \\ &= 2.2 \times 10^{-8} \text{ (V/eV)} \\ &= 22 \text{ (mV/MeV)} \quad \dots \dots (3-4) \end{aligned}$$

4. Characteristics

In general, the following characteristics are required of charge amplifier used for the detection of soft X-rays and low to high-energy gamma rays.


- High gain
- Low noise
- Excellent integration linearity
- High-speed rise time
- High temperature stability, etc.

The following sections discuss major characteristics of charge amplifier.

4-1 Open-loop gain

There are various semiconductor detectors used for the detection of soft X-rays and gamma rays. Even among Si detectors for example, a variety of types are used to match the application, which have different active areas and depletion layer thicknesses. Furthermore, detectors also differ in regards to capacitance. However, when the same Si detector is used for the detection of soft X-rays and gamma rays, the amount of generated charge must be the same if the particle energy of the soft X-rays and gamma rays is equivalent. Therefore, charge amplifier must provide a constant gain regardless of the capacitance value. In fact as shown in Eq (2-5) in “Principle of operation”, the output from the detector is independent of the junction capacitance C_j . This is because the open-loop gain of the charge amplifier is very high.

Figure 4-1 Equivalent circuit



C_j : Semiconductor detector capacitance
 C_f : Feedback capacitance
 AOL : Open-loop gain of amplifier

KACCC0016EA

When a charge amplifier is connected with a Si detector, its equivalent circuit is like that shown in Figure 4-1.

In this equivalent circuit, when seen from the amplifier's input, the input impedance Z_{in} is given by

$$Z_{in} = \frac{1}{j\omega C_f} \cdot \frac{1}{1 + AOL(j\omega)} \quad \dots \dots (4-1)$$

If signal charge Q_s is generated in the Si detector, the voltage e_{in} at the amplifier's input becomes

$$e_{in} = \frac{Q_s}{j\omega C_j + \{1 + AOL(j\omega)\} \cdot j\omega C_f} \quad \dots \dots (4-2)$$

Thus the output voltage e_{out} is expressed using Eq (4-2), as follows:


$$\begin{aligned} e_{out} &= AOL(j\omega) \cdot e_{in} \\ &= AOL(j\omega) \cdot \frac{Q_s}{j\omega C_j + \{1 + AOL(j\omega)\} \cdot j\omega C_f} \\ &= \frac{Q_s}{j\omega C_f + \frac{j\omega}{AOL(j\omega)} (C_f + C_j)} \quad \dots \dots (4-3) \end{aligned}$$

Here, assuming that $AOL \gg 0$, in other words, the open-loop gain of the amplifier is very large, then e_{out} can be simplified as follows:

$$e_{out} = \frac{Q_s}{j\omega C_f} \dots\dots (4-4)$$

As discussed above, the output voltage e_{out} of charge amplifier is not dependent on the capacitance of Si detectors.

Figure 4-2 Open-loop gain (H4083)




KACCB0022EA

4-2 Noise

Figure 4-3 shows the noise equivalent circuit of a charge amplifier.

Figure 4-3 Noise equivalent circuit of the charge amplifier



C_j : Capacitance of semiconductor detector
 C_s : Input capacitance of charge-sensitive amplifier
 C_f : Feedback capacitance
 R_f : Feedback resistance

KACCC0017EA

Noise in charge amplifier comes from the following three major sources:

◆Thermal noise of first-stage FET

Thermal noise of the first-stage FET, en_1 , is given by

$$en_1 = \sqrt{\frac{8}{3} \frac{KT}{gm}} \left(V/\sqrt{Hz} \right) \dots\dots (4-5)$$

- K : Boltzmann constant
- T : Absolute temperature
- gm: Mutual conductance of first-stage FET

◆Shot noise caused by gate current of first-stage FET and dark current of detector

The shot noise in is given by

$$in = \sqrt{2q (I_G + I_D)} \left(A/\sqrt{Hz} \right) \dots\dots (4-6)$$

- q : Elementary charge
- I_G : Gate leakage current of first-stage FET
- I_D : Dark current of detector

◆Thermal noise caused by feedback resistance

The thermal noise en_2 caused by the feedback resistance R_f is given by

$$en_2 = \sqrt{4KTRf} \left(V/\sqrt{Hz} \right) \dots\dots (4-7)$$


where R_f is the feedback resistance.

From Eqs (4-5), (4-6) and (4-7), the total noise $ent(j\omega)$ becomes as follows:

$$ent^2(j\omega) = en_1^2 \cdot \left(1 + \frac{C_{in}}{C_f} \right)^2 + \left\{ in^2 + \left(\frac{en_2}{R_f} \right)^2 \right\} \frac{1}{(j\omega C_f)^2} \dots\dots (4-8)$$


In Eq (4-8) above, the first term component is constant over the entire frequency range and amplified by the noise gain $(1+C_{in}/C_f)$ determined by the input capacitance C_f . The second term component is constant regardless of the input capacitance C_{in} , but decreases with increasing frequency. Figure 4-4 and 4-5 show noise characteristics of the H4083 charge amplifier.

Figure 4-4 Noise spectrum



KACCB0024EA


Figure 4-5 Input capacitance vs. number of noise electrons



KACCB0023EA

Generally, the noise measurement and evaluation for charge amplifier are carried out using a measurement system like that shown in Figure 4-6.

Figure 4-6 Noise measurement system




KACCC0018EA

In this system, a charge Q is supplied from the pulse generator via the capacitance connected to the input end of a charge amplifier.

The output from the charge amplifier is amplified once with the pulse-shaping amplifier and input to the pulse height analyzer. Then the pulse height distribution is measured to obtain the noise based on the half width.


The pulse-shaping amplifier improves the S/N of the charge amplifier and also serves as a filter. As various types of pulse-shaping amplifiers are available depending on the circuit configurations, choose the best amplifier for your application. One of the more popular pulse-shaping amplifiers is a Gaussian shaping circuit. Its typical frequency characteristic and output waveform are shown in Figure 4-7 and 4-8, respectively.

Figure 4-7 Frequency characteristic



KACCB0025EA

Figure 4-8 Output waveform




Since the output signal from a charge amplifier is as small as several tens of millivolts even after being integrated. The pulse-shaping amplifier is also used to amplify this signal to a level that matches the input range (0 to 10 V) of the pulse height analyzer.

5. Applications

5-1 Gamma ray measurement (Direct detection using Si PIN photodiode)

Figure 5-1 shows a typical pulse height distribution when the H4083 charge amplifier combined with S3590-05 Si PIN photodiode is irradiated with gamma rays from an ²⁴¹Am radiation source.


Figure 5-1



5-2 Power and stability measurements from lasers

Figure 5-2 shows a typical pulse height distribution when the H4083 charge amplifier combined with S3590-01 Si PIN photodiode is irradiated from a laser having a pulse width of about 100 ns and a wavelength at 830 nm. The peak channel of the pulse height distribution indicates the average laser power and the half width represents the fluctuation in the power level.


Figure 5-2



6. Specifications (H4083)

Charge amplifier H4083 is a hybrid low-noise amplifier that can be used in a wide range of application, including soft X-ray and low to high energy gamma-ray spectrometers. H4083 is used in conjunction with semiconductor detectors such as Si detectors. Even when combined with a detector having a relatively large terminal capacitance, H4083 can make measurements with good energy resolution.

For optimum operation it is recommended that H4083 be used with Hamamatsu Si PIN photodiodes (S3590/S3204 series, etc.). In particular, S3590 series eliminates the risk of trouble from increased capacitance because it solders directly onto the rear side of H4083 thus reducing wire length. The increased compactness and light weight offered by this configuration make it easy to incorporate in with other test equipment, etc.



■ Specifications

Parameter	Specification
Amplification method	Charge-sensitive type
Input/output polarity	Inverted
Charge gain	0.5 V/pC
	22 mV/MeV (Si)
Noise characteristic	550 electrons/FWHM
Negative feedback constant	50 M Ω /2 pF
Power supply	\pm 12 V
Power consumption	150 mW
Configuration	9-pin, single line type
Dimensional outline	24 (W) \times 19 (H) \times 4 (T) mm

■ Pin connections

Pin No.	Symbol	Content
①	IN	Input terminal
②	GAIN	Feedback constant adjustment terminal
③	GND *	Input ground terminal
④	CAL	Last pulse input terminal
⑤	GND *	Power and output ground terminal
⑥	-12 V	Power terminal
⑦	+12 V	Power terminal
⑧	GND *	Output ground terminal
⑨	OUT	Output terminal
Through-hole A	P	Anode connection terminal
Through-hole B	N	Cathode connection terminal

* GND is internally connected.

The left end pin is designated No. 1 when viewed from the component side with the pins facing downwards.

● Using pin No. 2

Use pin No. 2 when changing the feedback constant which is typically 50 M Ω /2 pF. Connect feedback resistance and capacitance between pin No. 2 and No. 9. Note that this connection is made parallel to the internal feedback resistance and capacitance.

● Changing the feedback resistance and capacitance

If you want to change the feedback constant without making


parallel connections, directly replace the resistor (Rf) and capacitor (Cf) shown at the right, which serve as feedback resistance and capacitance.

Replacement components should be 2012 size.

● Using through-holes A and B


When H4083 is not covered with protective coating, a photodiode (such as S3590 series) can be directly soldered on the board. In this case, insert the photodiode leads into the through-holes from the solder side.

■ Component side diagram



KACCC0138EA


■ Equivalent circuit




KACCC0019EA

■ Circuit example 1 (external connection to photodiode)

■ Circuit example 2 (direct connection to photodiode)





KACCC0020EA



KACCC0021EA

■ Dimensional outline (unit: mm)



KACCA0017EB

7. Precautions for handling charge amplifier

- Because charge amplifier H4083 is designed so that Si PIN photodiode of S3590 series can be directly soldered on it, no protective coating is provided. Avoid touching any part with bare hands when handling H4083.
- For proper storage of the charge amplifier, insert its leads into a non-conductive sponge material or wrap it with aluminum foil to protect it against static electricity because the impedance at the input end is very high.
- When soldering Si PIN photodiode S3590 series onto a charge amplifier or the amplifier onto a PC board, carry out the soldering as quickly as possible to prevent the amplifier and other components from being damaged by excessive heat.
- When configuring a measurement system using Si PIN photodiode and a charge amplifier, the following points should be observed.
 - When mounting Si PIN photodiode on a charge amplifier, take care to minimize stray capacitance.
 - The output end of the amplifier should be kept away from the input end as much as possible.
 - In the power supply line connect a ceramic capacitor of about 0.1 μF between the amplifier and the ground.
 - In the bias line for the photodiode, use a bias resistor and capacitor which can adequately withstand the bias voltage applied.
 - Si PIN photodiodes are capacitive elements with high impedance. Provide appropriate shielding to prevent external noise from interfering with the signal.
 - As a power supply for the system, we recommend using a series power supply with minimum noise characteristics.
- The following points must also be observed during operation:
 - To increase the bias voltage applied to Si PIN photodiode, slowly increase it while monitoring the output so that the output amplitude does not swing past the optimum level.
 - Because Si PIN photodiodes also have high sensitivity in the visible range, provide sufficient light-shielding against visible light when measuring soft X-rays or gamma rays.

HAMAMATSU

Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2001 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81) 53-434-3311, Fax: (81) 53-434-5184, <http://www.hamamatsu.com>

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P.O.Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-1218

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49) 08152-3750, Fax: (49) 08152-2658

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: 33-(1) 69 53 71 00, Fax: 33-(1) 69 53 71 10

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44) 1707-294888, Fax: (44) 1707-325777

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171 41 Solna, Sweden, Telephone: (46) 8-509-031-00, Fax: (46) 8-509-031-01

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39) 02-935-81-733, Fax: (39) 02-935-81-741