

# Amptek ASCII Configuration File Format

The Amptek ASCII Configuration File Format describes how Amptek applications save and read ASCII command configurations. ASCII commands are sent in the 'Text Configuration' Request Packet (PID1=0x20, PID2=2). The commands can be sent singly, but most often, they are grouped together into a single "Text Configuration" packet so that the entire configuration can be sent at once.

Amptek Digital Pulse Processor (DPP) and Multichannel Analyzer (MCA) devices using ASCII commands include DP5(FW6), X-123(FW6), PX5, DPG, MCA8000D, and GAMMA-RAD5.

The number of bytes sent/received in a "Text Configuration" packet must not exceed the maximum size of the corresponding buffer. If the number of bytes exceeds the buffer size the "Text Configuration" packet must be split into multiple packets. (See the DP5 Programmer's Guide.)

## Configuration File General Rules

- No whitespace before commands
- Commands must start at left column
- Commands must end with semicolon
- Anything after semicolon on a line is ignored (i.e. comments)
- All commands are in upper case

# Configuration File Format

## ASCII Commands and Values

Every ASCII command has a name and a value, delimited by an equals sign (=). The command name appears to the left of the equals sign. ASCII commands are always the 4-character Mnemonic in upper case that represents the ASCII command. If no value is present between the equals sign and the semicolon terminator the command is ignored.

*Example ASCII command and value:*

MCAC=1024; MCAC - Selects the Number of MCA Channels

## Sections

Commands are grouped into named sections. The section name appears on a line by itself, in square brackets ([ and ]). All keys after the section declaration are associated with that section. There is no explicit "end of section" delimiter. Sections end at the next section declaration, or the end of the file. Sections may not be nested. Section names may be case sensitive depending upon the implementation.

*Example Section Names:*

[DP5 Configuration File]

[DP5 Configuration Values]

[DP5 SCA Configuration]

## Comments

Semicolons (;) at the beginning of the line indicate a comment line. Comment lines are ignored. Semicolons (;) at the end of the command indicates where the command ends and optional comments may be placed. Comments are ignored from after the semicolon up to the end of the line.

*Example commented line:*

; comment line text

## Blank Lines

Blank lines are ignored with the exception of the "Advanced Programming" example application.

## Repeated Section Names or Commands

All section names are unique and declared only once. All commands that belong in a section will only be recognized if it is in the corresponding section. Commands are only allowed once per section, with the exception of the "Advanced Programming" example application.

# Configuration File Sections

There are three configuration file sections:

1. [DP5 Configuration File] – Stores ASCII commands.
2. [DP5 Configuration Values] – Stores ASCII command alternate values.
3. [DP5 SCA Configuration] – Stores SCA register setup values.

## DP5 Configuration File Section

This section is required. This is the main ASCII command storage section. All commands are stored here except for indexed SCA settings.

## DP5 Configuration Values Section

This section is optional. This section stores values used to setup configuration dialogs. Some values are predefined (discrete) while other values are selected from a range (continuous). Some hybrid commands have both discrete and continuous values. For hybrid commands it is useful to store the continuous values. This allows for seamless switching between discrete and continuous values.

## DP5 SCA Configuration Section

This section is required only if SCA counters are being used. This section is only used by DppMCA and DP5 SDK applications.

## Using/Setting DP5 SCA Configuration Section Values

There are four SCA ASCII commands that do not adhere to "INI File" acceptable rules. The same commands are repeated when sent to an Amptek device. Alternative (INI File compatible) storage must be used to take advantage of widely used (cross platform) "INI File" libraries. The four commands are:

- SCAH - Set SCA High Threshold
- SCAL - Set SCA Low Threshold
- SCAO - Select SCA Output Level
- SCAI - Set SCA Index

SCAI is used in conjunction with the SCAL, SCAH and SCAO commands to specify which SCA these commands will apply to. When the SCA index is set to 1 through 8, it also enables the specified SCA output and sets it to active high, so the SCAO command doesn't need to be specified unless a different output option is needed.

The SCA index is stored with the SCA configuration command in the DP5 SCA Configuration section.

**So:**

SCAI=4;SCAO=OFF;SCAL=1;SCAH=8192;

**Becomes:**

SCAO4=OFF;

SCAL4=1;

SCAH4=8192;

**And:**

SCAO6=OFF;

SCAL6=1;

SCAH6=8192;

**Becomes:**

SCAI=6;SCAO=OFF;SCAL=1;SCAH=8192;

## Using Configurations with the "Advanced Programming" Visual Basic Application

The "Advanced Programming" example application only uses the " DP5 Configuration File" section. When using the "Advanced Programming" example application, other sections must be removed.

The "Advanced Programming" example application uses repeated SCAI, SCAL, SCAH and SCAO commands in the same order that they are sent out. (This is not compatible with many "INI File" libraries).

**So:**

SCAI=4;SCAO=OFF;SCAL=1;SCAH=8192;

**Becomes:**

SCAI=4;

SCAO=OFF;

SCAL=1;

SCAH=8192;

**And:**

SCAI=6;

SCAO=OFF;

SCAL=1;

SCAH=8192;

**Becomes:**

SCAI=6;SCAO=OFF;SCAL=1;SCAH=8192;

# Example DppMCA ASCII Configuration File

[DP5 Configuration File]

RESC=YES; Reset Configuration  
CLCK=20; 20MHz/80MHz  
TPEA=12.800; Peaking Time  
GAIF=0.9562; Fine Gain  
GAIN=66.738; Total Gain (Analog \* Fine)  
RESL=3276; Detector Reset Lockout  
TFLA=0.200; Flat Top  
TPFA=1600; Fast Channel Peaking Time  
PURE=OFF; PUR Interval On/Off  
RTDE=OFF; RTD On/Off  
MCAS=NORM; MCA Source  
MCAC=1024; MCA/MCS Channels  
SOFF=OFF; Set Spectrum Offset  
AINP=NEG; Analog Input Pos/Neg  
INOF=DEF; Input Offset  
GAIA=20; Analog Gain Index  
CUSP=0; Non-Trapezoidal Shaping  
PDMD=NORM; Peak Detect Mode (Min/Max)  
THSL=3.613; Slow Threshold  
TLLD=OFF; LLD Threshold  
THFA=99.81; Fast Threshold  
DACO=SHAPED; DAC Output  
DACF=50; DAC Offset  
RTDS=603; RTD Sensitivity  
RTDT=8.98; RTD Threshold  
BLRM=1; BLR Mode  
BLRD=3; BLR Down Correction  
BLRU=0; BLR Up Correction  
GATE=OFF; Gate Control  
AUO1=ICR; AUX\_OUT Selection  
PRET=OFF; Preset Time  
PRER=OFF; Preset Real Time  
PREC=OFF; Preset Counts  
PRCL=1; Preset Counts Low Threshold  
PRCH=8191; Preset Counts High Threshold  
HVSE=160; HV Set  
TECS=220; TEC Set  
PAPZ=OFF; Pole-Zero  
PAPS=8.5; Preamp 8.5/5 (N/A)  
SCOE=FALLING; Scope Trigger Edge  
SCOT=12; Scope Trigger Position  
SCOG=1; Digital Scope Gain  
MCSL=1; MCS Low Threshold  
MCSH=8191; MCS High Threshold  
MCST=0.00; MCS Timebase  
AUO2=ICR; AUX\_OUT2 Selection  
TPMO=OFF; Test Pulser On/Off  
GPED=RISING; G.P. Counter Edge  
GPIN=AUX1; G.P. Counter Input  
GPME=OFF; G.P. Counter Uses MCA\_EN?  
GPGA=OFF; G.P. Counter Uses GATE?  
GPMC=OFF; G.P. Counter Cleared With MCA Counters?  
MCAE=ON; MCA/MCS Enable  
VOLU=OFF; Speaker On/Off  
CON1=DAC; Connector 1  
CON2=AUXIN2; Connector 2  
BOOT=ON; Turn Supplies On/Off At Power Up  
SCAW=100; SCA Pulse Width (Not Indexed - SCA1-8)

[DP5 Configuration Values]

TPEA=12.800;  
GAIF=0.9562;  
GAIN=66.738;  
RESL=3276;  
TFLA=0.200;  
SOFF=;  
INOF=;  
CUSP=0;  
THSL=3.613;  
TLLD=;  
THFA=99.81;  
DACO=;  
DACF=50;  
RTDS=603;  
RTDT=8.98;  
BLRD=3;  
BLRU=0;  
AUO1=;  
PRET=;  
PRER=;  
PREC=;  
PRCL=1;  
PRCH=8191;  
HVSE=160;  
TECS=220;  
PAPZ=;  
MCSL=1;  
MCSH=8191;  
MCST=0.00;  
AUO2=;  
GPIN=;

[DP5 SCA Configuration]

SCAO1=OFF;  
SCAL1=0;  
SCAH1=1023;  
SCAO2=OFF;  
SCAL2=0;  
SCAH2=1023;  
SCAO3=OFF;  
SCAL3=0;  
SCAH3=1023;  
SCAO4=OFF;  
SCAL4=0;  
SCAH4=1023;  
SCAO5=OFF;  
SCAL5=0;  
SCAH5=1023;  
SCAO6=OFF;  
SCAL6=0;  
SCAH6=1023;  
SCAO7=OFF;  
SCAL7=0;  
SCAH7=1023;  
SCAO8=OFF;  
SCAL8=0;  
SCAH8=1023;