

## Efficiency and Attenuation in CdTe Detectors

Bob Redus, May 25, 2010

Amptek's XR-100T-CdTe is a high performance x-ray and gamma ray detector system. Like Amptek's other XR100 products, a detector element and preamplifier components are mounted on a thermoelectric cooler. The CdTe products replace the Si diodes used in the other XR100 products with CdTe, a wide bandgap, compound semiconductors as the detector element. The primary advantage of CdTe is its much greater efficiency, due to its higher atomic number, Z. The photoelectric cross-section scales as  $Z^5$ . For Si,  $Z=14$ , while for CdTe,  $Z=50$ . The efficiency of Amptek's 500  $\mu\text{m}$  Si detectors begins to fall above 10 keV, while for 1 mm CdTe, efficiency is high to 100 keV.

The detection efficiency is a very important consideration, but due to charge transport effects, defining it is somewhat subtle. The purpose of this note is to provide general information on efficiency, for estimating system performance, and to recommend procedures for measuring the actual efficiency of a given detector. The charge transport is discussed in more detail in another Amptek Application Note, "Charge Trapping in XR-100T-CdTe and -CZT Detectors" AN-CZT-002<sup>1</sup>.

### 1 INTRODUCTION

As is well known<sup>1</sup>, when a beam of energetic photons, X-rays or  $\gamma$ -rays, passes through a material the result is a simple exponential attenuation of the primary beam. Each of the possible interaction processes can be characterized by a probability of occurrence per unit path length in the absorber. The sum of the probabilities for the individual processes is the total probability per unit length that the photon is removed from the beam. This is termed the linear attenuation coefficient, is denoted  $\mu$ , and has units of inverse length ( $\text{cm}^{-1}$ ). The number of primary photons transmitted through a thickness  $t$  is

$$I_{trans} = I_0 e^{-\mu t}$$

where  $I_0$  is the flux of incident photons,  $t$  is the thickness of the attenuator,  $\mu$  is the linear attenuation coefficient, and  $I_{trans}$  is the flux of transmitted primary photons. The number of primary photons interacting in a thickness  $t$  is obviously

$$I_{in} = I_0 [1 - e^{-\mu t}]$$

The linear attenuation coefficient obviously depends strongly on energy, since the interaction mechanisms are energy dependent. The attenuation is often described using the mass attenuation coefficient,  $\mu_m = \mu/\rho$  where  $\rho$  is the density of the medium. This can be written in units of  $\text{cm}^2/\text{g}$ , with the density in  $\text{g}/\text{cm}^3$ , or in units of barns ( $1 \text{ barn} = 10^{-24} \text{ cm}^2$ ) with the density in atoms per  $\text{cm}^3$ .

There are several different processes by which photons interact. In the energy range most often measured with the XR-100T-CdTe, the most important processes are the photoelectric interaction and Compton scattering. In a photoelectric interaction, the entire incident energy of the interacting photon is deposited in the detector, while in Compton scattering, only a portion of the incident energy will generally be deposited in the detector. Photoelectric interactions contribute to the full energy, which is usually of primary interest. The probability of a photoelectric interaction is usually of primary interest.

### 2 APPLICATION TO THE XR-100T-CDTE

Amptek's standard XR-100T-CdTe consists of a 1 mm thick CdTe detector located behind a 4 mil Be window. The probability of a photon interaction somewhere in the detector thickness is the product of (1) the probability of transmission through the Be window, dead layer, and contacts, and (2) the probability of interaction in the material,

$$P = \left( e^{-\mu_{Be} t_{Be}} \right) \left( e^{-\mu_{dead} t_{dead}} \right) \left( e^{-\mu_c t_c} \right) \left( 1 - e^{-\mu_{CdTe} t_{CdTe}} \right)$$

<sup>1</sup> This note focuses on CdTe, which is Amptek's current product available for higher energies. The  $\text{CdZn}_x\text{Te}_{1-x}$  (CZT) detectors used previously have very similar properties.




Figure 1 and Figure 2 show results which have been computed for Amptek's standard 1 mm CdTe detector. Also shown, for comparison, are results for a 0.5 mm Si detector and for a CdTe stack detector, 2.25 mm thick. The calculation includes the effects of transmission through the Be window and of stopping in the detector, but neglect the effects of trapping and hole tailing. Both total and photoelectric probabilities are given.




Figure 1. Linear plot of interaction probability.





Figure 2. Log-log plot of interaction probability between 1 keV and 1 MeV.

The data plotted in Figure 1 and Figure 2 are listed in ASCII text files on Amptek's web site, [www.amptek.com/anczt1.html](http://www.amptek.com/anczt1.html). The user may download these files for detailed calculations.



### 3 CONSEQUENCES OF TRAPPING AND HOLE TAILING

CdTe is a wide bandgap, high-Z, compound semiconductor material. It is used for X-ray and  $\gamma$ -ray spectroscopy because it has a very high linear attenuation coefficient, permitting high efficiency in a small volume, and low leakage currents, permitting low electronic noise without cryogenic cooling<sup>ii</sup>. However, like other compound semiconductors, it exhibits significant spectral distortions due to hole trapping. As is discussed elsewhere<sup>iii</sup>, the trapping length of holes in CdTe is smaller than the linear dimensions of the detector. For interaction occurring near the anode, virtually all of the signal is due to electrons and so the full charge is collected. For interactions occurring near the cathode, virtually all of the signal is due to holes and so a smaller charge is collected.

The result is that the measured signal, the measured “energy,” depends upon the depth of the interaction on the detector, decreasing with increasing depth. In the output spectrum, one observes a tail of counts towards lower amplitudes, an effect known as “hole tailing.” Figure 5 is a plot of the pulse height as a function of depth in the detector, computed from the Hecht<sup>iv</sup> relation for two different values for the hole lifetime. Figure 6 is a plot of a measured spectrum demonstrating hole tailing with a 2 mm CZT detector. This was a measurement of a <sup>57</sup>Co source, which emits primarily at 122 keV. The blue trace in Figure 6 is the raw spectrum. The tail of counts extending to low values is due to hole tailing.





Figure 5. Plot showing the induced signal size as a function of depth, computed for two different values of the hole lifetime.

Note that in Figure 5, for the blue curve, about 40% of the depth of the detector produces a full signal. The rest of the depth produces a smaller signal. The effective volume of the detector, which contributes to the primary, full energy peak, is only 40% of the physical volume of the detector. For the red curve, about 20% of the detector volume contributes to the full energy peak.

For CdTe and other detectors in which charge trapping is important, it is critical to distinguish between the total geometric efficiency and the efficiency of the full energy peak. The total geometric efficiency, which is due to the total physical volume of the detector, can be used to compute the total rate of counts in the detector. The efficiency of the full energy peak, which is due to the volume of the detector leading to “full charge collection”, can be used to compute the total rate of counts in the primary peak. The phrase “full charge collection” is in quotation marks because it is not well defined. The charge collection efficiency decreases smoothly with increasing depth. Different users may define the full energy peak differently, depending on the specific application.

In Amptek’s XR-100T-CdTe, risetime discrimination (RTD) may be used to minimize spectral distortions due to hole tailing. Figure 5 showed that the induced signal size is correlated with depth in the detector. The risetime of the pulse from the preamp is also well correlated with depth in the detector. Therefore, Amptek’s signal processors measure the risetime of the pulses and reject those with a long risetime. This leads to a significant improvement in the quality of the spectrum<sup>v</sup>.

When RTD is used, the detection efficiency is significantly lower than would be anticipated from the physical dimensions of the detector. The effective depth depends upon the charge transport properties of



the material. These are not well controlled by the manufacturer, so significant variations exist from one detector to the next. The effective depth also depends upon the RTD setting.




Figure 7. <sup>57</sup>Co spectrum with RTD on and off for the 1 mm thick CdTe.

#### 4 MEASURING THE EFFICIENCY

In many applications it is important to know the detection efficiency at a particular energy. Because there exists a wide variation in the effective depth of Amptek's XR-100T-CdTe detector, if a user requires precise knowledge, the best solution is to measure the actual efficiency at the energy of interest from some well-known standard. Since this can be difficult, another approach is to measure the effective depth of the particular detector. There are two approaches to this measurement.

If one has a calibrated source, with a very well known strength and an energy high enough to be only partially absorbed in the detector material, then the effective depth can be readily computed by inverting the relation above:

$$P = 1 - e^{-\mu t} \Rightarrow t = -\ln(1 - P) / \mu$$


For example, from Figure 1 the 122 keV line from <sup>57</sup>Co should be detected with 44% efficiency and the linear attenuation coefficient is 5.85 cm<sup>-1</sup>. Assume that a lab measurement shows that this line is detected with 35% efficiency. This gives t=0.74 mm. This effective depth can be used to compute the efficiency at any other energy (above the energy where attenuation in the Be window is significant).

In the absence of a calibrated source, then a single source which emits  $\gamma$ -rays at two distinct lines with a well known ratio can be used, if at least one of the lines is high enough in energy to be detected with <100% efficiency and both are above the energy where the Be window is significant. The source must not attenuate either line significantly. For example, <sup>57</sup>Co emits at 14.4 keV with 9.8% efficiency and 122 keV with 85.6% efficiency. We define P<sub>1</sub> and P<sub>2</sub> as the probabilities of emission of the two lines, N<sub>1</sub> and N<sub>2</sub> and the measured counts of the two lines, and  $\mu_1$  and  $\mu_2$  are the linear attenuation coefficients of the two lines. The ratio of the measured counts will obviously be

$$\frac{N_1}{N_2} = \frac{P_1(1 - e^{-\mu_1 t})}{P_2(1 - e^{-\mu_2 t})}$$

Applying some algebra, this implies

$$t = \frac{\ln(N_1 P_2 - N_2 P_1) - [\ln(N_1 P_2) - \ln(N_2 P_1)]}{(\mu_1 - \mu_2)}$$



- <sup>i</sup> a) Knoll, Glenn F., Radiation Detection and Measurement, John Wiley & Sons, New York. 1989.
- b) Tsoulfanidis, Nicholas, Measurement and Detection of Radiation, Hemphire Publishing Corporation, New York, 1983.
- <sup>ii</sup> Jordanov, V.T., J.A. Pantazis, and A.C. Huber, "Thermoelectrically-Cooled Cadmium Zinc Telluride Detectors (CZT) for X-Ray and Gamma-Ray Detection," Radiation, Vol. 43, No. 1, July 1996.
- <sup>iii</sup> a) Charge Trapping in XR-100T-CdTe/CZT Detectors, Amptek Application Note (ANCZT-2) by Bob Redus, 2002.
- b) Semiconductors and Semimetals vol. 43, Semiconductors for Room Temperature Nuclear Detector Applications, section on Characterization and Quantification of Detector Performance by Vernon M. Gerrish, Volume Editors T.E. Schlesinger and R.B James, Academic Press, San Diego, 1995.
- <sup>iv</sup> Hecht (1932)
- b) Semiconductors and Semimetals vol. 43, Semiconductors for Room Temperature Nuclear Detector Applications, section on Characterization and Quantification of Detector Performance by Vernon M. Gerrish discusses Hecht relation, Volume Editors T.E. Schlesinger and R.B James, Academic Press, San Diego, 1995.
- <sup>v</sup> Squillante, M.R., Presentation at 11th International Workshop on Room Temperature Semiconductor X-Ray and Gamma-Ray Detectors and Associated Electronics, 1999 Vienna Austria.