
High Performance Image Control System

HiPic users manual

By Hamamatsu Photonics Deutschland GmbH

This manual was produced using *ComponentOne Doc-To-Help*.™

Contents

Introduction	1
Legal terms	1
Security/Precautions	2
Precautions of individual devices	2
Overexposure	2
Image Intensifier	2
Help / Info	2
How do I get more information about this program?	2
In which circumstances should I read the information in this documentation?	3
How do I use this help system?	3
Help	3
About	3
Features	4
General	4
What is new in version 8.3.0?	4
Hints for users familiar with version 7 or before	5
Supported hardware and hardware setup	6
Required components	6
Internal driver and DCAM driver	6
Supported cameras/grabber by internal driver	7
Camera/grabber combination table by internal driver	7
Cameras/Grabbers exclusively supported by Internal driver	8
Cameras/Grabbers exclusively supported by DCAM	9
Installation	9
Setup launcher	9
Framegrabber / Camera / DCAM selection	11
Installing individual frame grabber driver etc.	11
Hardware installation	14
Changing/Selecting different hardware from the intro screen.	15
Setting a different file as *INI file.	18
Basic operations	19
Getting started	19
Start the program	19
Acquiring images from a camera	21
LUT	23
Change the display size of the image (zooming)	24
Setting camera parameters	24
Acquisition/Options dialog	25
Saving and loading images	26
Options	27
Quick profile	28
The elements of the user I/F	29
The menu	29
The Toolbar	30
The status bar	30
Child windows and System modal windows	30
Options, Acquisition and Tool dialog windows	30
Arranging windows/Activating dedicated windows	31
Image display windows	32
LUT	32
ROIs	32
Context sensitive menus	33

Using commands from the menu	34
The main menu	34
File Menu	35
Tools menu	35
The acquisition menu	35
The Corrections menu	36
The analysis menu	36
The display menu	36
The Processing menu	36
The Window menu	37
The Info menu	37
Using commands from the toolbar	37
CCD Camera control	38
Acquisition/Options Dialog	38
List of all camera dialogs	40
List of camera parameters	60
Acquisition modes	67
Live	69
Acquire Mode	69
Difference between Live and Acquire mode	69
Analog Integration	70
Photon Counting	70
Camera data types and data storage	70
Corrections	72
Background subtraction correction	72
Shading correction	74
Defect pixel correction	75
Rolling average	76
Saving and loading images	76
Saving images	76
Loading images	77
Corrections	79
Overview	79
Background Subtraction	79
Background Subtraction from Camera	80
Background Subtraction from File	80
Constant	82
Real time corrections	83
Real time background correction	83
Real time shading correction	83
Shading Correction	84
Get a Shading Reference Image	84
Declare an Image as Shading Reference Image	84
Algorithm of Shading Correction	85
Defect Pixel Correction	86
Defect pixel correction tool	86
Enabling defect pixel correction	89
Offline Defect Pixel Correction	89
Background/Shading correction assistant	89
Performing Background subtraction	90
Performing Background subtraction and Shading correction	93
Reserved Filenames	96
Other Options	97
Image Display & LUT	103
Image display windows	103
Caption	104
Zoom	104
Rulers	105
Quickprofile	105
FWHM	108
ROIs	109

Create ROI	112
Mouse Coordinates	113
ROI I/F	113
Image Status	114
Context sensitive menu	115
Print Images	115
LUT	116
Selecting measurement parameters	121
Options	121
General options	121
CCD camera	122
Acquisition	122
Calibration	124
Sequence	125
Background	126
Shading	128
Defect pixel	129
Images	130
Quick profile	132
Analysis profiles	133
Saving and loading parameters	134
Permanent parameters	134
Workfiles	134
Image positions	136
Calibration & Analysis	137
Calibration	137
What is Calibration?	137
Linear/Table	138
System / Image / Profile Calibration	139
Using Different Calibration Methods	140
Profiles	142
Quick profiles	143
Analysis profiles	143
Histogram	154
3D Data	155
Extract 3D Data from an Image	155
Using the Clipboard or a 3D Data DDE Link	157
Special acquisition modes & special setups	159
Line sensors	159
Introduction	159
Line sensor Options	160
Image Acquisition	165
PIV mode	169
Introduction	169
Setting the camera to PIV mode	169
Recording two PIV images in sequence mode	170
X-ray Flat Panels	171
Options	171
Live Mode	173
Acquire mode	174
TDI mode for ORCA	176
Introduction	176
Selecting TDI mode	176
Live and Acquire mode	177
Display a quick profile	177
TDI Mode for C10000	177
Area Mode	178
TDI mode	179
Normal and scrolling LIVE display.	180
Exposure time in TDI mode	182
Using Sequences	183

Forward/Reverse scanning	184
Generic camera	185
Selecting Generic Camera	185
Generic Camera Options Dialog	185
Generic camera acquisition dialog	186
Sequence acquisition	189
Selecting acquisition parameters	190
Acquisition mode and camera parameters	191
Cycle definition	192
Other information on the acquisition tab	192
Start sequence acquisition	193
Data storage options	196
Storing images	196
Differences between storage modes	197
Additional features and limitations of fast hard disc recording	198
Specifying file names	199
Storing profiles	200
Image Sequence display (normal images)	200
Image Sequence display (Linesensor and TDI data)	201
Sequence image file naming conventions	202
Profile sequences	203
Profile sequence acquisition	203
Profile sequence display	204
Processing sequences	205
Averaging	205
Photon Counting.....	207
Introduction	207
Required Hardware	207
Conditions for Photon Counting Mode	208
Sensitivity to single photons	209
Low light level	209
Basic principle of photon counting	209
Threshold	209
Reduce the spot information to a single count	209
Coordinate determination	209
Photon counting algorithm and implementation	209
Background subtraction and threshold calculation	209
Coordinate calculation	211
Advanced processing	211
Data Storage	212
Counting errors and maximum count rate	212
Correct threshold	212
Overlapping photon spots	212
Maximum count rate	213
Operating Photon Counting	213
Photon Counting Setup	213
Photon counting acquisition	214
Options related to photon counting	215
Dynamic Photon Counting (time resolved 2-D photon counting)	215
Acquisition	216
Analysis	216
Offline photon counting	217
Reserved Filename	218
Image processing.....	219
Arithmetic	219
Superimpose	221
Technical Explanation of Superimpose	221
The Superimpose Control	222
Average	223
Creating data for special purposes / Tools	225
Calibration file editor	225

Defect pixel correction tool	227
Programming techniques	229
User function	229
General	229
Function Parameters	230
Memory and Display parameters	230
User Function Dialog	232
Prototype	233
User function in Sequence mode	236
Remote control programming	237
RemoteEx	237
Script programming	238
Accessing the CCD camera directly	238
Appendix	241
Image file formats	241
ITEX image (*.img format).	241
TIFF format	242
Display to TIFF image	242
Data to TIFF image	242
ITEX Sequence	243
TIFF Sequence	243
ASCII image	244
ASCII Sequence	244
HIS image sequence file	245
Profile File Format	250
Calibration File Format	251
DPC File Format	251
Status String Format	252
Application section	253
Camera section	253
Acquisition section	265
Grabber Section	265
Display LUT Section	265
Calibration Section	266
Comment Section	266
Constants	266
Details of sequence mode	269
LUT	270
Linear LUT (default)	270
Gamma LUT	271
Sigmoid LUT	271
Error handler	272
Glossary of Terms	273
Index	277

Introduction

Legal terms

Please read this license carefully before installing or using the software. By using the software, you are agreeing to be bound by the terms of this license. If you do not agree to the terms of this license, do not install or use the software and erase it from your computer.

COPYRIGHT

The HiPic software (herein called the Software) and all accompanying documentation (herein called the Documentation) are copyrighted by Hamamatsu Photonics Deutschland GmbH (herein abbreviated as HPD) with all rights reserved. It is illegal to copy (except as granted under section "License Agreement" below), duplicate, sell, or otherwise distribute the Software or the Documentation or any parts of them without written permission from HPD.

Any attempts to "reverse-engineer", "disassemble", "re-source" the Software or to modify the Software in any way are violations of HPD's copyright and may be legally prosecuted.

LICENSE AGREEMENT

HPD gives you the permission to copy the software for your own use as long as it is absolutely necessary for using the Software for its intended purpose. This includes your right of loading the Software into the computer's memory and making copies of the installation disks for the purpose of backup. You may also run multiple copies of the Software on several computers for the purpose of data inspection and analysis. Under no circumstances, however, are you allowed to use the Software to drive data acquisition hardware on more than one computer simultaneously at any given time.

WARRANTY

All efforts have been made to make the Software and the Documentation as accurate and free of "bugs" as possible. However, **HPD makes no warranties, express or implied, and specifically disclaims any warranty of merchantability or fitness for a particular purpose. In no event is HPD liable for direct, indirect or consequential damages of product or property, personal damage, inconvenience, loss of time, data or profits, which may result from the use of the Software, even if advised the possibility thereof. The entire risk as to the results and performance of the Software is assumed by the purchaser.**

Specifically, where life or health is involved, results obtained with the Software should be examined carefully before any conclusions are made. The Software and the Documentation are not designed and tested intending to ensure a level of reliability suitable for use in diagnosis and treatment of humans.

TRADEMARKS

Product names and company names mentioned in the Software or the Documentation are trademarks or trade names of their respective holders.

UPDATES

Update conditions may vary depending on where you purchased this product. Please contact your local Hamamatsu subsidiary or retailer to obtain information about the update conditions valid for you.

HPD is committed to a policy of customer satisfaction and product reliability. If you find any errors, "bugs" or other inconveniences in the Software or the Documentation, please inform HPD or your local retailer so that countermeasures are possible with the next update.

Security/Precautions

This chapter deals with the precautions you should observe when handling and operation the devices used in your system. This is of course primarily to avoid risk of dangers related to your live and health (like electrical shock) but also to prevent any damage to the used device (like damage due to overexposure). Please read the following chapters carefully and don't operate the devices unless you are sure that you are operating them on a save basis.

Precautions of individual devices

At first please read the safety instructions of the individual devices. Every device has its own handling manual where its safety precautions are described properly.

Overexposure

One risk for an operated device is that it might be overexposed, which in certain cases may lead to a damage or even the destruction of the device. Devices which may be damaged during overexposure are Image intensifiers, other type of tube cameras and also some kind of CCD cameras at least if they are exposed under special conditions (UV, very strong light etc.).

Image Intensifier

An image intensifier is a device which is very sensitive to illumination and can be easily damaged. This is especially true when the high voltage is applied to the photocathode. If a system which includes an image intensifier is operated illumination should very low in the beginning and stepwise increased by continuously monitoring the output on the computers screen with the maximum available contrast enhancement appropriate for the given case.

Help / Info

HiPic designed to offer a user-friendly operation system for these cameras as well as all basic image processing functions needed for routine work. In many cases the user can just guess by the command names and the arrangement of the user I/F controls how a desired function works.

How do I get more information about this program?

The following possibilities are available to get more information about this program, what it can do and how this can be achieved.

- A windows help system explains divided in small easy to understand topics the program and all its parts
- The help system is context sensitive, which means that you get the information about a specific topic just by pressing F1 while the program is in this context
- For those how like to study the feature of this with the help of a printed manual also such a manual is available.

In which circumstances should I read the information in this documentation?


- To get started with the digital image processing system.
- To get detailed information about complex operation methods and procedures. Since HiPic is a "look and feel" program, you may not need this manual for your daily work. However there are several functions which may not be understood easily by just using the program. We have put priority to explain such items in this manual.

How do I use this help system?

To use this help system simply press F1 from the program or double-click to the Help file from the explorer. The Help system will open a dialog. If you call the help system from the program in many cases it directly leads you to the topic related to the current circumstances.


Help

The help system can also be invoked by executing the Help... menu entry from the Info main menu.



About

Executing the About menu entry from the Info main menu shows the about Dialog which gives you information about the current program version and other version related information.



Features

General

HiPic is an image processing software for high performance CCD cameras. It can be used in combination with most Hamamatsu digital cameras like the ORCA, ORCA II, C4880, C8800, C9100, C9300 cameras, standard analog video cameras (CCIR or EIA standard) CMOS X-ray flat panel sensors and X-Ray line sensors.

What is new in version 8.3.0?

The following features are new or improved in version 8.3 compared to version 8.2:

- All functions dealing with corrections are transferred to a DLL (CorrLib.Dll). There are several advantages for this:
 - Higher speed in most cases.
 - Better maintainability.
 - DLL can be sold to customer for his development, then having the same functionality as under HiPic.
- The Defect Pixel Correction Tool has additional functions:
 - Can detect overflow lines and columns.
 - Calculates the number of pixels which cannot be corrected with the current method (using neighbors).
 - Can show the location of defects on the original image as overflow.

- Calculating of the defect coordinates runs many times faster.
- Frame trigger for line sensor / TDI sensor in combination with PC-DIG, PCCamLink, Phoenix LVDS, Phoenix CamLink, X64 LVDS, X64 CamLink
- Phoenix Dig 36 supports all Flat Panel Sensors.
- Flat Panels with special 2 TAP readout improved: C9311, C9312, C9252DK, C10502D.
- 13 Bit mode supported
- Recursive Filter for all area cameras (n=2...255)
- Generic camera control extended for Mode control, trigger output, alternate trigger, max image area, data valid.
- Grabber options dialog for most of the grabbers: IC-PCI+AMDIG, IC-PCI+AMVS, PCDIG, PCCamLink, PCVision, Phoenix LVDS, Phoenix CL, X64 LVDS, X64 CL.

Hints for users familiar with version 7 or before

Even though version 8 has been designed to give maximum user friendliness users familiar with older versions (like version 7 or 6) may find several features confusing at the first glance. This is due to the fact that some features have been improved but needed to be modified. Reading the following small chapter may help to resolve this confusing situation.

First of all such users may desperately search for options at well defined locations and cannot find them any more. The solution is simple: All options including the former camera setup parameters have been collected in the options dialog and can be called with the main menu entry **file – options** or with a context sensitive menu. Please see also the chapter Basic operations - Options and Selecting Measurement parameters – Options in the document.

Another change was to introduce context sensitive menus with right mouse click. As a consequence the right mouse click which is used to zoom out images is no longer available. There are two solutions for this: Generally speaking all right mouse clicks are replaced by the shift key. In other words zoom out can be done by Shift Key + Left Mouse click. The second solution is to use the context sensitive menu to select the zoom factor directly. Please see also the chapter Context sensitive menus and zoom in this document.

The third topic which is very different from earlier versions is how to draw and handle ROIs. The biggest difference is that a new ROI can only be draw when a ROI type is selected in the toolbar. After selection the ROI this ROI type is no longer valid thus the ROI cannot be easily overwritten. To draw more then one ROI just select a ROI type again. To select one ROI just click to it and to delete press the delete key. The new ROI types Full horizontal and Full vertical may satisfy the request from many users in the past. Please see also the chapter The

Elements of the user I/F - ROIs and Image display Windows – ROIs in this document.

Supported hardware and hardware setup

Required components

The HiPic requires the following system components:

- Cameras with RS 422 or LVDS interface:
Frame grabber: **IC-PCI+AM-DIG** or **PCDig**, or **PCI 1422** or **PCI 1424** or **AS-PHX-D36-PCI32** or **AS-PHX-D36-PCI64**
- Cameras with CameraLink interface:
Frame grabber: **PcCamLink**, **PCI1426** or **PCI1428** or **PCI1429e** or **Active Silicon Camera Link board** or **X64 Full** or **X64 iProLight**
- Cameras with IEEE1394 interface:
IEEE1394 adapter
- Analog cameras
Frame grabber: **IC-PCI+AM-VS** or **PC-Vision** or **PC2-Vision**
- Computer environment:
An industry standard PC with at least 512 MB RAM, one free serial port (not required for cameras with IEEE1394 interface), one PCI slot, Windows '98, Windows ME, Windows NT, Windows 2000 or Windows XP is required.
- The HiPic
- Camera:
A Hamamatsu digital camera C4742-95, -98 (black and white version), C4880, C8800, C7190, C7300-10, C8000, C9100, C9300, C10000, X-ray Flat panel sensor, X-Ray line sensor or a camera with standard analog video signal (e.g. Hamamatsu C3077).
- Accessories (if C4880 or C4742-98 is used): Depending on the camera type, a water cooler may be necessary.

Internal driver and DCAM driver

To understand which cameras or more precisely speaking which combinations of cameras and frame grabbers are supported we have to know that there are two methods of interfacing the camera/frame grabber:

- Internal driver (using software modules written together with the main application)
- DCAM API (using software modules provided by the manufacturer, Hamamatsu KK Japan)

Many combinations of cameras/frame grabber are supported using both methods. However there are combinations of cameras/grabbers which are only supported by one of these methods. The following is an overview of hardware supported by internal drivers. To find out which cameras/grabbers are supported by DCAM please see the release notes on the DCAM section of your distribution CD or consult the following website.

<http://www.dcamapi.com/>

Please note that color cameras are not supported at all by this program.

Supported cameras/grabber by internal driver

This is a list of supported hardware by internal drivers:

Frame grabber: IC-PCI+AM-DIG, IC-PCI+AM-VS, PCVision, PCDig, PCCamLink, PCI1422, PCI1424, PCI1426, PCI1428, PCI1429e (single tap mode only), PC2Vision, X64 Full, X64 iProLight, All Phoenix grabbers

Cameras: C4880, C4880-8X, C4742-95(ORCA 100), C4742-98 (ORCA II), ORCA II ER, ORCA II BT, C5987, C6918-05, C7300, C3077, C5405, ORCA ER, C8000-10, C8000-20, C8484-01, FlatPanel (C7942, C7943, C7921, C7920, C9728/C9730 (USB, internal driver) C9732 (USB, with DCAM), C9250 (30Hz, no binning), C9311, C9312, C9321, C10013, C9252, C10322, C10500D, C10502D), ORCA-HR, C8800, C9100, C9300, C7390, C8133, C8750, C9133, ORCA ER TDI, C9750, C10000-201, C9750,

Camera/grabber combination table by internal driver

This is a detailed table which cameras can be operated with which frame grabbers by internal drivers (Note: Some cameras are available with different interfaces. Example: The C8484 is available with RS422, IEEE1394 and Cameralink. If the camera is equipped with Cameralink you have to operate it with a camera link board)

Analog video output

Camera \ Grabber	IC-PCI+ AM-VS (Coreco)	PC-VISION (Coreco)	PC2-Vision (Dalsa)
CCIR or EIA analog camera	X	X	X
C8000-20	X	X	X

RS422 / LVDS

Camera \ Grabber	IC-PCI+A M-DIG (Coreco)	PC-DIG (Coreco)	AS-PHX-D36-PCI32 or AS-PHX-D36-PCI64 (Active Silicon)	PCI1422 (National Instruments)	PCI1424 (National Instruments)

C4880 standard	X	X	X	X	X
C4742-95 (ORCA)	X	X	X	X	X
C4880-8X	X	X	X	X	X
C4742-98 (ORCA II)	X	X	X	X	X
C7300	X	X	X	X	X
C8000-10	X	X	X	X	X
C8484	X	X	X	X	X
ORCA HR	-	X	X	X	X
Single Tap Flat Panel (C7920, C7921, C7942, C7943, C9250DP, C9321, C10322D)	-	X	X	X	X
Dual Tap Flat Panel (C9311SK, C9312SK, C9252DK, C10502D)	-	-	-	-	X
LineSensor	-	X	-	X	X
ORCA TDI	X	X	X	X	X

Camera Link

Camera	PC-Cam-Link	PCI1426 / PCI 1428	PCI 1429e	Phoenix Camera Link D24	Phoenix Camera Link D48	X64 Full	X64 iProLight
C8484	X	X	X	X	X	X	X
C8800	X	X	X	X	X	X	X
C9100	X	X	X	X	X	X	X
C9300	X	X	-	X	X	X	X
C10000-201/501	-	-	-	-	X	X	-
C10000-401	-	-	-	-	X	X	X

USB:

The C9266 dental X-ray sensor and the USB Flat panels C9728/C9730 are supported by internal driver.

Cameras/Grabbers exclusively supported by Internal driver

The following cameras/interfaces/grabber are exclusively supported by internal drivers:

- The Coreco IC-PCI board
- The Coreco board PCVision
- The Dalsa board PC2Vision
- The Dalsa board X64 Full and X64 IProLight
- The National Instruments boards PCI 1422 / PCI 1424, PCI 1426, PCI 1428, PCI 1429e
- The USB Flat panels C9728/C9730
- The C9266 dental X-ray sensor
- All RS422/LVDS Flat panels
- The C7300 camera.

Cameras/Grabbers exclusively supported by DCAM

The following cameras/interfaces/grabber are exclusively supported by DCAM:

- All fire wire (IEEE 1394) cameras
- All USB Cameras (except the C9266 dental X-ray sensor)

Installation

This chapter describes the hardware configuration, set-up and software installation of the HiPic system. If your system was already pre-installed by Hamamatsu you do not need to read this chapter.

Note: Please install the software at first and install the hardware in a second step.

Setup launcher

HiPic software is delivered on a CD.

Note: You should install the frame grabber board or IEEE1394 adapter into your computer **after** the software installation is completed.



Start Windows 98SE, Windows ME, WindowsNT 4.0 SP6, Windows 2000 or Windows XP.

Note: Cameras with IEEE1394 interface are not supported under Windows NT.

Insert the CD into your CD drive. The installation program will normally be started automatically by a "autostart" routine.

If it does not run automatically, you have to select "**Run**" from the "**Start**" menu and type **D:\SetupHi** then click "**OK**" (assuming that your CD drive has the name D).

The set-up program will start displaying the setup launcher dialog.



The setup launcher is a collection of setup programs from Hamamatsu and third parties to install the main program, several drivers and other utilities.

Notes: Please make sure that older versions of all programs/drivers and utilities are uninstalled before the new version is installed.

Please collect information about the camera type and interface type (grabber, IEEE1394 etc.) before you start the installation. Press **Help on selecting cameras...** to get some assistance about camera and interface selections.


Start with the software installation by pressing the **Install the HiPic** button. Follow the instructions displayed on the screen.

Framegrabber / Camera / DCAM selection

Depending on the camera model and camera interface type you have to decide which **setup type** you want to use.

Select **DCamAPI modules** for configurations supported by a Hamamatsu DCAM driver (e.g. cameras connected by IEEE 1394 interface).


Select **Internal drivers** for cameras not supported by DCAM.



Please ask your camera supplier if you are not sure which setup type you should use. Many cameras are supported by both drivers and the customer can freely select between the two drivers. If in any case a problem arises with one of the drivers use the other one could be of help.

Installing individual frame grabber driver etc.


At first, please select the camera which you want to use.




Camera selection dialogue

Then, please select the frame grabber which you intend to use in the next dialogue.


In case of a RS422/LVDS camera please select as follows:



In case of a Camera Link camera please select as follows:



In case of an analog camera please select as follows:



At the end of the installation you will be advised to install additional program components such as a hardware lock driver or a frame grabber driver and you will be asked if you want to reboot your computer. If you decide to do it later, you have to make sure that you do not start the program before you have restarted the computer. If these drivers are already installed on the computer and you are sure that they are operating you do not need to install it again.

Please select the appropriate pushbuttons on the setup launcher to start the installation of the driver/utility programs.

Typically this is the installation of the hardware lock driver.

The next installation step is the installation of a hardware interface driver. You have to choose one out of the four alternative drivers (ITEX, IFC, SAPERA NI-IMAQ or DCAM) shown in the setup launcher dialog. Please choose the proper driver for your actual installation.

The final installation step is the installation of the Acrobat Reader which is required to read the help file.

After completion of the installation you have to re-boot your system.

Hardware installation

Proceed with following steps to install the hardware of your HiPic system:

Step 1: Install the frame grabber or IEEE1394 interface board in your computer.

Install the frame grabber board in a free PCI slot of your computer. Following frame grabber types are supported: IC-PCI with AM-DIG and AM-VS module, PCVision, PC2Vision, PCDig, NI PCI 1422, NI PCI 1424, NI PCI 1426, NI PCI 1428, NI PCI 1429e, Active Silicon Phoenix, X64 Full, X64 iPro Light

If you use a camera with IEEE1394 (FireWire) interface, an interface adapter is needed instead of a frame grabber board.

Step 2: Install the software protection.

In order to work with HiPic you have to install a hardware lock.

If this is not installed, HiPic works in demonstration mode only.

Connect the hardware lock (this is a 25 pin connector labeled with the name of the program or a USB stick type) to the parallel port (LPT1) respective a USB port of your computer. The hardware lock is delivered with the HiPic program disk.

Step 3: Make the cable connections

- Connect all cables of your computer (AC cable, keyboard, mouse, monitor etc.)
- Connect all AC line cables of your other devices.
- Connect the camera to the input connector of the frame grabber board.

If you use a Hamamatsu digital camera:

- Connect the camera head to the camera controller.
- Connect the serial cable to a **serial interface** port (default: COM 1) of your computer. You can use any COM port, but you have to indicate this port during start-up. If you use a camera with CameraLink interface all control

commands are sent via the CameraLink cable, so no additional serial interface cable is required.


- Connect the water cooler (in case of a water cooled camera type).

Changing/Selecting different hardware from the intro screen.

Normally the correct access mode (internal driver/D-CAM API), the frame grabber and the camera selection are done by the setup program and/or by a Hamamatsu engineer during installation. Only in seldom cases - like if you selected a wrong camera model or if you want to change it later - you have to change it after you started the program but before you start the program finally by clicking OK.

Changing CCD camera setting


The first selection you have to do is to select between Internal drivers and DCAM API modules.



The screenshot shows a dialog box titled "CCD camera access". At the top, there are two radio buttons: "Internal drivers" (which is selected) and "DCam API modules". Below this, there is a section for "Internal Drivers" containing two dropdown menus: "Frame Grabber" set to "PC-CamLink" and "Camera" set to "C8800". To the right of these dropdowns are "OK" and "Cancel" buttons. At the bottom, there is a "Configuration file:" label followed by a text box containing "C:\Programme\HiPic\HiPic800\ifc-hpk.txt" and a "Get" button. Below the text box is a "Board No." label followed by a text box containing "0".

Select **Internal drivers** when a camera which is supported by HPD-TAs internal camera drivers shall be used. (Most cameras using a frame grabber use this method).

Select **DCAM API modules** when a camera which is supported by a DCAM API module shall be used (e.g. cameras with IEEE1394 interface).



When **Internal drivers** is selected, you have to make following settings:

Frame Grabber

Please select the frame grabber which is installed in your system.

- If no frame grabber is installed you can select the type **None**.
- Select **IC-PCI+AM-DIG** if you have the IC-PCI board installed with a digital input module.
- Select **IC-PCI+AM-VS** if you have the IC-PCI board installed with an analog input module.
- Select **PCDig** if you have the PCDig board installed which has a digital input module.
- Select **PCI 1422/1424** if you have a National Instruments board installed which has a digital input module.
- Select **PCVision** if you have the PCVision board installed which has an analog input module.
- Select **Phoenix/Cameralink** if you have any Phoenix/Cameralink board installed.
- Select **Phoenix/LVDS** if you have any Phoenix/LVDS board installed.
- Select **X64/CameraLink** if you have a X64 Full or X64 iProLight board installed.
- Select **PC2Vision** if you have the PC2Vision board installed which has an analog input module.

If you have installed more than one frame grabber of the same type in your computer you have to select the number of the grabber which shall be used during the actual session in the **Board No.:** field. Please select 0 if only one grabber is installed.

Camera Type

Please select the camera type which is installed in your system.

Configuration files

The frame grabber will be initialised by use of special configuration files. Normally a pointer to the correct file is automatically established during software installation.

During installation of the program the user is already asked for the camera he wants to use and an entry in the initialization file is already set specifying the configuration file and camera.

If you want to select another configuration file you should enter the correct file name for the configuration file in the text box **Configuration file..** Press **Get** in order to open a dialogue for the file selection.

Note: The program will save your settings automatically. When you start it the next time you will find that the previous settings are already set.

If all settings are finished you have to click **OK** in order to start the application.

Camera configuration files

The following camera configuration files are delivered with the HiPic. These files are installed and specified at installation time. Only if you change the frame grabber (or in certain cases also the camera) you need to change this specification manually.

Configuration file..	...with frame grabber...	... for the cameras
DIGITAL.CNF	IC-PCI+AM DIG	All digital cameras like C4742 (ORCA) series and C4880 series
CCIR.CNF	IC-PCI+AM VS	For analog video cameras with CCIR standard
EIA.CNF	IC-PCI+AM VS	For analog video cameras with EIA standard
PCVSCCIR.CNF	PCVision	For analog video cameras with CCIR standard
PCVS_EIA640.CNF	PCVision	For analog video cameras with EIA standard
PCDig.txt	PCDig	All digital cameras
Ifc-hpk.txt	PCDig/PCCamLink	All digital cameras
Hamamatsu_DCAM.pcf	Phoenix grabber	All digital cameras
Hamamatsu_DCAM.ccf	X64 Full, X64 iProLight	All digital cameras
P2V_CCIR.ccf	PC2-Vision	CCIR cameras

The National Instruments frame grabber PCI 1422, PCI 1424, NI PCI 1426, PCI 1428 and PCI 1429e do not require such file.

Setting a different file as *INI file.

The HiPic uses the file hipic32r.ini in the windows directory as the default to store all parameters.

To specify a different file as the file which stores permanent parameters a command line argument can be used. The command line can be specified in the link to the program.



Instead of the default target:

“C:\Program Files\Hipic\HiPic830\HiPic32ab.exe”

you can use the following statement:

“C:\Program Files\Hipic\HiPic830\HiPic32ab.exe /ini=myhipic.ini”

If the /ini=file statement contains a full path statement this path is used otherwise the file is searched in the windows directory.


Basic operations

This topic describes how to perform basic operations of the program, how to get started and how to get familiar with it in short time.


Getting started

This topic describes how to get started with the program.

Start the program

To start the program simply click to the icon  displayed on the desktop or select the appropriate command within the program tree: Start – Programs – HiPic 8 – HiPic 8.3.0. . If you are familiar with the windows explorer and with starting programs by clicking on the executable files you can also locate the application directory and double-click to hipic32ab.exe hpdta32ab.


The introduction screen will appear. While this screen is open you can select access mode (internal driver/DCAM), frame grabber, camera and in certain cases also COM port and baud rate of the serial communication.



All these settings have been done during setup and/or installation by a Hamamatsu engineer and are kept permanently. We therefore assume that we do not need to make any modifications here and we can go on starting the program by hitting Enter or clicking to the OK pushbutton. If you like to Cancel the starting process click to Cancel now.

Serial communication problems during startup

Most CCD cameras are controlled by a serial interface (sometimes this serial interface is operated through the cable which also transfers the digital data. This is the case for Camera Link cameras). If a problem during serial I/F communication happens the RS232 Timeout dialogue will appear.



RS232 Timeout dialogue

Several reasons can lead to such a timeout. In most cases either the RS232 cable is not connected or the device simply not switched on. Another reason can be that the wrong baud rate (clock speed) is selected in the startup screen. Most cameras use 9600 baud, however there camera you are using may need a different setting. Correct the error here and go on by selecting “Retry”. Sometimes only a certain serial command cannot be understood by the camera. In such cases Ignore continues the program without sending the command again. Abort aborts the communication process closes the communication and continues the program without any communication.

Acquiring images from a camera

We assume that the CCD camera is already switched on.

If the camera is equipped with water cooling switch on the cooling water now. Some older cooled cameras have a vacuum pump to avoid condensation on chip. If so, please switch on this vacuum pump and open all necessary valves as has been explained during installation of the system. Most modern cameras however do not have such peripheral equipment as they are permanently sealed and there is nothing to do here.

Images can be acquired using the camera acquisition dialog. Even though the acquisition dialogs for different cameras look somehow different you will find the pushbuttons for Live mode and Acquire mode on every camera acquisition dialog. Most dialogs have also a Single image pushbutton which acquires an image with the parameters selected for live mode. Any image acquisition can only be started if the camera acquisition dialog is visible. This is to avoid acquisition with wrong or unexpected parameters. If the camera acquisition dialog is not visible, the first attempt to acquire an image will place the acquisition dialog on screen and only the next trial will start the acquisition.


Live mode/Single image

Live mode is a mode to acquire and display images continuously. This mode is available for all cameras. In case of C4880 and ORCA II type cameras, the camera is switched to the fast readout mode. If a user wants to get only one single image with the currently selected parameters most camera acquisition dialogs have a single image function.

To start Live mode proceed as follows:


1.) Show the acquisition dialog

Select Live from the Acquisition menu




- or -

Click to the Live Toolbar button



The acquisition dialog will appear on screen



2.) Start Live mode

(If the acquisition dialog was already visible Live mode will start with step 1.)

To Start the LIVE mode

Click to the Live pushbutton on the Acquisition dialog


- or -

Click to the Live Toolbar button

- or -

Select Live from the Acquisition menu

The live image will appear on screen and show a continuously updated image.



To stop LIVE mode

Click to the Freeze pushbutton on the Acquisition dialog

- or -

Click to the Freeze Toolbar button


- or -

Select Freeze from the Acquisition menu

To execute single exposure Click to Single Exposure on the acquisition dialog.

LUT

If you acquire images you will find that the images are frequently displayed with too low or too strong contrast. If you acquire images under low light level conditions the display may be too dark. Use the LUT tool to adapt the image display to the desired contrast.



Lut tool.

The LUT tool has two cursors (two colored lines, one is blue and one is red) which defines the intensity limits of the currently displayed intensity range. To find the correct range the easiest way is to click on the small asterisk * on the right bottom side of the LUT tool. A similar pushbutton can be found on the

toolbar:

Please see also the chapter LUT, LUT Parameters and the Appendix LUT in this document.

Change the display size of the image (zooming)

Important note for users of previous versions: While older versions (before version 8) of this programs use the right mouse button to zoom out, this version uses the left mouse button + Shift key to zoom out.

Once an image is acquired it will be displayed on screen within a window. Depending on the number of pixels in the image and the resolution of the computers screen the image may be either too small or too large to be seen correctly. Therefore the image can be zoomed with a factor larger than one (if the image is too small) or a factor smaller than one (if the image is too large).

To change the zooming factor of an image


- 1.) Select zooming on the toolbar



- 2.) Click on the image with the left mouse button to enlarge the image and click with the left mouse while holding down the Shift key to reduce the image.

- OR -

Select the zooming factor with the context sensitive menu by clicking on the image with the right mouse button.



See also the chapter Image Display Windows in the document.

Setting camera parameters

If the image you get by this procedure is not satisfactory there may be several reasons. In case of a cooled camera you may need to switch on the peltier cooler now. It may be necessary to focus or adjust the optics and it may also be necessary to adjust the exposure time here. The most important camera parameter is the exposure time. With very few exceptions all CCD cameras have an exposure time control on the camera acquisition dialog.



In general it is a sliderbar where you can increase the exposure time by clicking on the right arrow and decrease it by clicking on the left arrow. Normally also an edit box allows to enter the exposure time directly. In this case the exposure time should be activated by pressing the Tab key.


If the image you get is too faint increase the exposure time, if it is too bright decrease the intensity. If portions of the image became red (or white in the case the color setting is rainbow) the CCD camera is saturated, which means the exposure time is too high. If the intensity is still too high you may need to decrease the light level of your sample or close the iris of your optics.

See also the chapter CCD camera control and CCD camera in the document.


Acquisition/Options dialog

The exposure time however is only one parameter of the CCD camera which can be set by the software. Depending on the type of CCD camera many other parameters can be set now. It is a good idea to consult the CCD cameras hardware or operations manual and find out which features and parameters are provided by the model. All these features can be controlled from the software by using two main dialogs: the camera acquisition dialog and the options dialog. The parameters are separated in a way that parameters which are normally used and changed very often can be found in the camera acquisition dialog. To avoid that the camera acquisition dialog is burdened with many - often unchanged - parameters, these less needed parameters are collected in the camera options dialog.

The camera acquisition dialog is always placed on screen when you execute an acquisition command. The simplest way to show the camera options dialog is to right click on the camera acquisition dialog and select Camera options.



The camera options dialog will then appear on screen.




Once the camera options dialog is displayed you can make all modifications of camera parameters. When clicking OK these modifications will be applied, when clicking to Cancel the old settings are restored. See also the chapter List of all camera dialogs for details.

Saving and loading images

Once you have acquired useful image data you may want to save them to disk.

To save an image to disk proceed as follows:

Select the File – Save As menu command.




- or -

Click to the Save As toolbar button



A dialog appears to allow specifying file name directory and filing type of the file. Additionally you can type a comment which is then saved into the file header. There are several file types which you can use. At start we recommend to use

either ITEX (*.img) or Data2Tiff (*.tif) as these file types allows to store the full dynamic range in an image.




Click to Save if you want to save the image now or Cancel if you want to cancel the save operation.

See also the chapter Saving and loading images in this document.

Options

The exact behavior of the program is specified in different options dialogs. To access the options either


Click the File – options menu



- or -

Click with the right mouse button on the dialog in which circumstance you want to specify options.

The options dialog will appear



Quick profile


During image acquisition one or more profiles could be shown on the same window where the image data is displayed. Because it is immediately updated when a new image appears it is also called quick profile.

To display a quick profile:

- 1.) Select the type of quick profile you want to get on the tool bar. We assume that we want to display a vertical profile over the full size of the image.




- 2.) Click to the image with the left mouse button and move the mouse while holding down the left mouse button. A new full size vertical ROI will be created.



- 3.) Click to the PRF button on the toolbar.



The quickprofile appears.




Clicking to PRF button on the toolbar again hides the quickprofile.

To move the ROI click on the ROI with the left mouse button and move the mouse.

To delete the ROI make sure that the image dialog is the active window and press the Del key on your keyboard. See also the chapter Image Display and LUT - Quick profile Selecting Measurement Parameters – Quickprofile for details.

The elements of the user I/F

This chapter explains the basic elements of the user interface and its general functionality. Every function is explained in more detail in further chapters if necessary. The following screenshot shows the main window after program start.



The menu

The menu of the program allows easy access to the most important commands. It consists of the main menu and a submenu for every entry in the main menu.

The Toolbar

The toolbar is a selection of pushbuttons which provides easy access to many important commands.



The status bar

On the bottom part of the program window status information is displayed. The content depends on the actual conditions. E.g. during LIVE mode the actual frame rate will be displayed

While the mouse cursor is inside of an image window, the co-ordinates and the intensity value at the current mouse position are displayed.



Child windows and System modal windows

Some windows used in the software are placed in the client area of the main window and cannot be larger than the main window. These windows are called child windows. Most of the windows in this software are child windows like the acquisition dialog or the image display windows.

Other windows are used to make special setting and it is useful that the process of setting is finished when the user wants to continue his routine work. These windows are so called modal (or system-modal) windows. Other windows of the software e.g. the main window cannot be accessed unless the modal window has been closed. Normally these modal windows have to pushbutton on the bottom:

- A pushbutton OK to apply the settings and close the dialog
- A pushbutton Cancel to restore the previous settings and close the dialog.

Options, Acquisition and Tool dialog windows

There are three different types of dialogs within this software:

- Option dialogs
- Acquisition dialogs
- Tool dialogs

Option dialogs

Option dialogs allow setting values of certain parameters which specify the behavior of the software. These parameters are static settings, which are normally restored from one session to the next. Option dialogs have to be closed after the settings are done to continue with the work.

Acquisition dialogs

Acquisition dialogs are dedicated to create image and related data. They are placed in the client area of the main window and can coexist with other windows especially image display windows.

Tool dialogs


Tool dialogs are dedicated to create other type of data like calibration data or defect pixel data which are further used for certain system settings. They are also placed in the client area of the main window and can coexist with other windows especially image display windows. These dialogs are normally not used during routine work

Arranging windows/Activating dedicated windows

To perform several tasks normally several dialogs have to be on the screen. This can sooner or later lead to an untidy and somehow messy screen. Even though the user still may need several dialogs at once there are two measures to handle this problem

Restore window positions

One can arrange the windows freely at the screen and there is an option to restore a window at the same location where it has been previously. There is an option “Restore window positions” in the general options dialog.



If this option is checked the position of every window will be remembered and the window placed at the same location where it has been. So it is easy to arrange the dialogs in a fixed way where they can be find easily.

Find windows which are hidden under other windows

On the status bar there are some pushbuttons which makes it very easy to get the desired windows on top which are mostly used:




These pushbuttons can be used to show the dialogs if they are already on screen but hidden under another dialog (if they are not on screen nothing will happen).

- IMG Shows all image display windows
- PRF Shows all dialogs related to profile
- SEQ Shows the sequence dialog
- LUT Shows the LUT tool


Image display windows

Image data is displayed on separate image display windows. A maximum of 20 image display windows can be displayed at once. Once an image window is created an entry is made in the menu windows. If more than one window is displayed at a time one of these is always the current (or selected) one. Every function which refers to image data is always processed on the current window. As an example, if you want to save one specific image to file click on the image windows caption (or to another part of the image). This makes the image the current one. Then click to the Save As function to save the current image.



LUT

The LUT tool is used to control the image display by manipulating its brightness and contrast.



The LUT Tool

If you acquire images you will find that the images are frequently displayed with too low or too strong contrast. If you acquire images under low light level conditions the display may be too dark. Use the LUT tool to adapt the image display to the desired contrast.

ROIs

There are many image processing functions which refer to image data like saving the data to file or extracting profile data. Most of these functions however can not

only be applied to the full image but also to a subset of the image data. For this purpose Regions of Interest (shortly called ROI) can be used. There are different types of ROIs which can be used for different purposes. The types of ROIs are:

- Point ROI
- ↘ Arbitrary line ROI
- ▭ Rectangle ROI
- ▬ Rectangle ROI with full size in horizontal direction
- ▮ Rectangle ROI with full size in vertical direction

To create a ROI first of all select the ROI type on the toolbar:




Then draw the ROI with clicking to the starting point with the left mouse button, then draw the mouse to the end point of the ROI and release it. Afterwards you can move the ROI by clicking to it with the left mouse button and moving the mouse to the desired position. You can select up to 10 different ROIs on one image. A more detailed description of how to create, select, delete and use ROIs can be found in the chapter ROIs

Context sensitive menus

From any given dialog the user can reach the most related dialogs and commands by just right click to the dialog. A context sensitive menu will appear and show commands which may be important in this circumstance. We want to look at two different examples to show this feature:

If you click with the right mouse button on the sequence dialog you will get the following menu:




If you click with the right mouse button on the image display dialog you will get the following menu:



Using commands from the menu

The menu is a way to easily access the most important commands or dialogs. It consists of the main menu which can be seen on the top of the main window and several submenus. To select a submenu click to the main menu entry and select the submenu entry. For example if you want to select the arithmetic menu first click to processing main menu, then select Arithmetic and click to it.



The main menu

The menu of the program consists of the following main entries:

Command	Description
File	Contains a set of commands related to image files and other general properties
Tools	Contains commands to create special data
Acquisition	Contains commands to acquire image data
Corrections	Contains commands to perform different corrections on images
Analysis	Contains commands to analyze image data
Display	Contains commands to control and optimize image display
Processing	Contains image processing commands
Window	Contains commands to select and close image display windows
Info	Contains commands to get information about the program a description of commands and the connected camera

File Menu

The file menu contains the following entries

Open	Loads an image or image sequence in IMG, TIFF or other formats.
Save As	Saves the current image or image sequence in IMG, TIFF or other formats.
Save ROI As	Saves the specified ROI of the current image or image sequence in IMG, TIFF or other formats.
Print	Print the current image
Options	Shows a dialog which displays and allows to modify all available options
Workfiles	Allows to save and load all current options and other settings to a file
Image status	Shows the current image status
Exit	Shuts down the program

Tools menu

The tools menu contains the following entries

Calibration files	Shows a dialog which allows to create, modify and save calibration files
Defect pixels	Shows a dialog which allows to analyze images containing defect pixel

The acquisition menu

The acquisition menu contains the following entries

Live	If the Live mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Live mode dialogue is open, image acquisition in Live mode will be started. The acquisition parameters as shown in the Live mode dialogue are used.
Freeze	Live mode stops. The last image will be stored.
Acquire	If the Acquire mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Acquire mode dialogue is open, an image acquisition in Acquire mode will be executed. The acquisition parameters as shown in the Acquisition mode dialogue are used
Analog Integration	If the Analog integration mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Analog integration mode dialogue is open, an image acquisition in Analog integration mode will be executed. The acquisition parameters as shown in the Analog integration mode dialogue are used.

Photon Counting	If the Photon counting mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Photon counting mode dialogue is open, an image acquisition in Photon Counting mode will be executed. The acquisition parameters as shown in the Photon counting mode dialogue are used. See.
Clear Image	The Clear Image command clears the image data
Sequence	Acquisition of series of images, storage on computer memory (RAM) or hard disk. Sequence mode can be used if a sequence of images shall be recorded with a high frequency. This command includes also a sequence replay function.

The Corrections menu

The acquisition menu contains the following entries

Background subtraction	This function allows subtracting a background image. This can be used to subtract the camera dark current
Shading correction	This functions allows to correct the shading of images

The analysis menu

The analysis menu contains the following entries

Profile	This function allows to display and analyze profiles
Histogram	This function allows analyzing the intensity distribution (histogram) of an image.
3D-Data	This function allows displaying image data in a numerical (table) form.

The display menu

The display menu contains the following entry:

LUT...	This function shows the LUT parameters dialog on screen. The LUT parameters dialog allows the user to change parameters like to color of the display.
--------	---

The Processing menu

The Processing menu contains the following entries:

Arithmetic...	This function allows to perform arithmetic operations on image data
User function...	This function allows to operate user function commands
Superimpose...	This function allows to superimpose the image data of two images

The Window menu

The Window menu displays all the currently available image windows. It also contains the “Close All” command which allows closing all image windows with one command.

The Info menu

The Info menu contains the following entry:

- Help... This function calls the help dialog
- RS232... This function allows to send and receive commands from the connected camera (if it is connected by RS232)
- About... This function displays the About dialog which tells the user details about the version of the program














Using commands from the toolbar










The toolbar allows the user – similar to the menu – to access the most important commands in a very easy way. To execute a command by the toolbar simply click to the icon on the toolbar.



Most of these icons can be understood without an explanation.

The following is an explanation of all toolbar icons used on the toolbar:

	Open	Opens an image file
	Save As	Saves the currently selected image to file
	Live	Starts live mode (or shows the camera acquisition is not on screen)
	Freeze	Freezes live mode
	Acquire	Acquires an image in Acquire mode
	Analog Integration	Acquires an image in Analog Integration mode
	Photon Counting	Acquires an image in Photon Counting mode
	Background Subtraction	Applies Background Subtraction to the currently selected image
	Shading Correction	Applies Shading Correction to the currently selected image
	Profile	Shows or hides the profile dialog
	LUT dialog	Shows or hides the LUT dialog
	About	Shows the About dialog
	Pointer	Selects Pointer mode. This allows selecting and

		modifying ROIs.
	Zoom	Selects zooming mode. This allows to change the zooming factor of the image
	Create Point ROI	Selects Point ROI mode. This allows creating a Point ROI.
	Create Line ROI	Selects Line ROI mode. This allows creating a line ROI.
	Create general rectangle ROI	Selects Rectangle ROI mode. This allows creating a rectangle ROI of any size.
	Create horizontal rectangle ROI (full image width)	Selects horizontal rectangle ROI mode. This allows creating a rectangle ROI with full horizontal size.
	Create vertical rectangle ROI (full image height)	Selects vertical rectangle ROI mode. This allows creating a rectangle ROI with full vertical size.
	ROI Interface	Shows or hides the ROI Interface dialog
	Quick profile	Shows or hides the quick profile
	Unzoom current	Sets the current image to the default zooming factor

CCD Camera control



Hamamatsu CCD cameras and other imaging devices are sophisticated state of the art technical devices which in general have many different operation modes. To get familiar with these operating modes you should consult your operations manual of the CCD camera. Depending on your purchase contract there may also be an installation and introduction done by a Hamamatsu engineer on your site. Also this introduction is a good opportunity to get familiar by the features and operations modes of the camera. Most of the features of your CCD camera can be controlled by software (There may be some features which can only be controlled by hardware devices like potentiometers or the like. In this case you cannot control the features by software).

Acquisition/Options Dialog

All features can be controlled by two dialogs:

The camera acquisition dialog and the camera options dialog. In very few cases there is a third dialog which controls the subarray feature of the camera. The camera acquisition dialog is a dialog which is placed on the client area of the main window and can stay there without a limitation. The camera options dialog is showed as a modal dialog. This means that the options dialog has to be closed before the user can continue his routine work.

As the features of different cameras are different also the camera control dialogs look different. As an example we show the camera acquisition dialog and camera options dialog of the C4742-95 (ORCA).



Often/Seldom used parameters

The parameters to control a camera are distributed in a way the seldom used parameters are placed in the camera options dialog, whereas parameters which are likely that they are used more often are placed in the camera acquisition dialog.

Camera acquisition dialog: Smaller/larger

The camera acquisition dialog can be made smaller or larger with two small arrows on the upper right side. Depending on the number of acquisition parameters the user wants to see he can make the dialog smaller or larger



Small camera dialog



Large camera dialog

Acquisition modes

Different acquisition modes can be used to acquire image data from a camera. The parameters are stored together with the acquisition mode and reactivated when the user reselects the acquisition mode.

Camera/Acquisition parameters

There are parameters which directly modify settings inside the camera. We call these parameters camera parameters.

There are parameters which do not change any setting inside of the camera but still influence the result of the acquisition, like the exposure number of an Analog Integration. We call these parameters acquisition parameters.

Value range of parameters



Depending on the camera and the setting of other parameters the value range of certain parameters can vary. For example if the Binning mode is set to 1 x 1 for an Orca camera, subarray settings for X and Y can range from 8 to 1024. If the same Orca camera is operated in 2 x 2 binning the subarray settings can only range from 8 to 512. Many parameters are dependent on each other. If possible parameters are arranged in a way the parameters from the options dialog influence parameters in the acquisition dialog but not vice versa.

List of all camera dialogs


This topic includes screenshots of all supported cameras and a list of supported parameters by these cameras. Please note that for many of these cameras there are subversions and optional features which leads to that fact that certain parameters are only available for certain cameras. If you are not certain please verify the precise camera type number, consult your camera documentation, get information from our website (www.hamamatsu.de or www.hamamatsu.com) or contact your local Hamamatsu dealer. Both the Camera acquisition dialog and the Camera options dialog are shown.

Analog camera



This camera selection includes all cameras which output a standard video signal (either CCIR or EIA). Among others these are the Hamamatsu C2400, C3077 and C5405.



C4742-95 (ORCA) + C8484





C4742-98 (ORCA II)





C4880



The C4880 standard model comprises all sub types of the C4880 except the C4880-8X. It includes the special versions which have incorporated an image intensifier C4880-91, C4880-92, C4880-92, C5987 and C6918-05.





C4880-8X





C7300





C8000-10




C8000-20





C8800



C9100



C9266 X-Ray dental sensor



ImageEM (C9100-13)

C9100-13 acquisition control

Live
 Acquire
 Analog Integration
 Photon counting

Action

Exposure time

31 ms

Gain & Offset

Gain 0

EMG 0

DSP Commands

Recursive filter Number: 2

Background Sub Offset: 0

Storage: Memory 1

Shading correction

Storage: Memory 1

Real time background subtraction (RTBS)

RT backsub:

RT shading corr:

All Options: <C9100-13 options>

Worksettings

- General
- Camera
- Acquisition
- Calibration
- Sequence
- Corrections
- Images
- Profiles

Timing mode: Scan Speed:

Trigger mode: OnChipAmp:

Trigger source: Cooling Fan

Trigger polarity: Cooler:

Trigger number: Target Temperature:

Cycle number: Ext. Output Polarity:



Scan mode: Ext. Output Delay [s]:

Ext. Output Width [s]:



Photon imaging mode:

—Camera information —

C9300





C10000




Cameras operated by DCAM module


Most cameras can be controlled also by a programming interface called DCAM. DCAM modules are provided by Hamamatsu and comprise a combination of camera and interface board. Even though the cameras can have different parameters the user interface when using DCAM is similar and look like the following screens.



Cameras operated by DCAM module (Version 3.0)


Recently a new version of the DCAM interface has been developed. Until now only a small number of cameras are supported by this interface, but the number of these cameras will increase quickly. The User Interface looks as follows:






X-Ray Flat Panel

There is a variety of Hamamatsu X-Ray flat panel sensors. The following models are currently supported: C7942, C7943, C7921, C9720, C7930, C9728DK (USB, internal driver), C9250, C9312).










X-Ray Line Sensor

The HiPic also supports Hamamatsu X-Ray line sensors. The following models are currently supported:

C7390, C8133, C8750, C9750, C9133.







List of camera parameters

The following is a complete list of all camera parameters, its range of values and an explanation of the meaning of the parameter.

Parameter	Range	Description
Binning	1x1 / 2x2 / 4x4 / 8x8 / 16x16 / different factors in two directions	A special readout mode of Hamamatsu digital cameras. Already in the camera the charges from a predefined number of pixel is accumulated. The accumulated charge signal is then digitized and read out. Binning allows getting higher camera sensitivity and reducing the amount of

		data. Of course this function works in expense of spatial resolution.
Hor. Binning	Integer number	Binning factor in horizontal direction
Ver. Binning	Integer number	Binning factor in vertical direction
Trigger polarity	Positive / Negative	Defines the active level of an electrical signal starting a camera exposure. If Trigger polarity is negative a falling edge (transition from high to low) starts image exposure (provided that the camera is switched to edge trigger mode).
Scan Mode	Normal/ Superpixel / Subarray / Interlace / Outline / Extended	<p>Scan mode defines how the content of a CCD chip is read out.</p> <p>Normal: The whole CCD chip is read out pixel per pixel</p> <p>Superpixel: A special binning mode where a square number of pixel is binned (e.g. 2 x 2, 4 x 4). This function can be used to increase camera sensitivity, increase the readout rate (frame rate) and reduce the amount of data.</p> <p>Subarray: Only a user-predefined area of the CCD is read out. This function can be used to reduce the amount of data and speed up the readout, if only a part of the image is necessary for analysis.</p> <p>Interlace: A mode where every second line is read out.</p> <p>Outline: A mode where the full chip is read out with low resolution and high frame rate</p>
Light mode	Low/High Or 0-1	Two different settings of the camera can be selected. In Light mode “ Low ” (=“0”) the camera works with reduced sensitivity but maximum dynamic range, while in Light mode “ High ” (=“1”) the camera works with maximum sensitivity at a slightly lower dynamic range.
Scan speed	High / Low or High / Low / medium or 0 / 1 (0=fast, 1=slow) or (single tapped / dual	Speed of CCD readout. This is mainly determined by the pixel clock of the camera. Fast readout has the advantage of high frame rate, whereas slow readout results in high dynamic range. Sometimes slow readout is called high precision mode. Cameras with slow readout are also called slow-scan cameras.

	tapped mode)	
Trigger Mode	Edge / Level / PIV	<p>This parameter defines how a camera reacts on an external trigger.</p> <p>Edge: Starts exposure with a (rising or falling according to polarity) edge of an electrical signal (Normally a TTL signal)</p> <p>Level: Exposure starts with an edge (falling or rising according to polarity) and stops with the inverse edge of an electrical signal (Normally a TTL signal)</p> <p>PIV: PIV stands Particle image velocitometry. With an edge (falling or rising according to polarity) of an electrical signal (Normally a TTL signal) two images are acquired and readout within short time.</p>
Timing mode	<p>Internal / External or Internal / External Grb / External Lemo or Internal / External Edge 2 / External Level 1 / External Level 2 / External Edge 1</p>	<p>The timing or timing mode defines whether start and stop of an exposure is defined internally (in the camera) or by an external electrical signal</p> <p>Internal defines the timing inside the camera where External starts/stops exposure by an electrical signal.</p> <p>External-Grb defines the timing by the connected frame grabber</p> <p>External Lemo defines the timing by an electrical signal inputted by a Lemo connector</p> <p>External Edge 1 or 2 and External Level 1 or 2 are special timing modes dedicated to the C9266 X-ray dental sensor (see C9266_InstructionManual.pdf delivered with the C9266 for details)</p>
H-Offset	Integer number	Horizontal offset
H-Width	Integer number	Horizontal width (Setting used only in subarray sacn mode)
V-Offset	Integer number	Vertical offset (Setting used only in subarray scan mode)
V-Width	Integer number	vertical width (Setting used only in subarray scan mode)
X direction (pixel)	Integer number	Size of the acquired image in X-direction
Scan Area	Two integer numbers	Number of pixel in both directions (x and y)
Trigger source	BNC / D-Sub / Digital I/F / Multiport	BNC: A BNC connector normally at the rear side of the camera controller

		D-SUB: A 9 pin D-Sub connector at the rear side of the camera controller
		Digital I/F: An electrical signal normally connected from the digital data connector of the camera to the frame grabber
		Multiport: A 4 pin Hirose type connector at the camera head, which eventually has other features on other pins.
Exposure time	Time in h, m, s, ms, us, ns / External / Controlled by PC	The exposure time is the time the CCD camera is sensitive to light. This timing is dependent on timing and trigger mode and electrical signals applied to the camera controller or head. In internal mode and in special other externally triggered modes like the edge trigger mode this time is just specified as a time in physical units (e.g. seconds). In other cases like level trigger modes it is defined by the length of the electrical pulses. In again other cases it is defined by the PCs clock or by the timing of external events.
Gate delay	Time in h, m, s, ms, us, ns	The delay time after a gate pulse. Specified in physical units.
Gain/Amp gain	Low / high / super-high or 0-255 or 0 / 1 or 0 / 1 / 2	Gain or Amp gain specifies the setting of an analog amplifier inside the CCD camera controller or head. Sometimes the range of settings depends on other parameters like scan speed.
MCP (MCP gain)	Integer number	MCP gain specifies the setting of a high voltage applied to a photocathode inside an image intensifier placed on the front of a CCD camera.
EMG	Integer number	EMG specifies the gain of an electron multiplying device inside a CCD camera.
Offset	Integer number	Offset specifies an analog which is subtracted to the analog signal of the CCD chip. It is normally used in combination with gain of amp gain.
Panel switch	On / Off	Specifies whether a small panel on the front side of the CCD camera controller is active or not. If Off the parameters can be specified only by software.
Shading control	On / Off	Specifies whether an analog shading corrector is active on the front side of the camera.

Contrast enhancement	Off / Potentiometer / Computer control	Specifies whether and how an analog contrast enhancement is active and whether the parameters can be specified by potentiometers or by computer control.
Target temperature	Real number, normally in Steps of 5 degrees.	Sets the temperature the peltier cooler tries to reach at the location of the chip.
Actual temperature	Real number normally with a precision of 0.1 degree.	Displays the actual CCD temperature
ROM version	String	Displays the ROM version
Chip type	String	Displays the Type of CCD chip
Vol. Offs.	Integer number	Displays the potentiometer setting for offset
Vol. gain.	Integer number	Displays the potentiometer setting for gain
Shutter control	"Auto" / "Open" / "Close" / "Close on readout"	Control how a mechanical shutter which is mounted on the front side of the CCD behaves Auto: Always closed expect during camera exposure Close: Closed Open: Opened Close on readout: Always open, except during readout of an image
Optical Black	On / Off	Optical black defines whether black lines and columns on the left and top side of the CCD (normally covered by black cover) are readout
Trigger number	Integer number	Number of triggers applied to the CCD camera before it readouts an image
Cycle Number	Integer number	Number of cycles for an EM-CCD camera
On Chip Amp	EM readout / CCD readout	Readout mode of an EM CCD camera.
Ext. Output Polarity	Positive / Negative	Polarity of the trigger output of an EM-CCD camera
Ext. Output Delay [s]	Floating point number	Delay of the trigger output of an EM-CCD camera in seconds
Ext.	Floating point	Width of the trigger output of an EM-CCD

Output Width [s]	number	camera in seconds
Photon imaging mode	No / middle / high	A mode to increase sensitivity under very weak light conditions (patented by Hamamatsu Photonics KK)
Mode	Internal Trig. / Ext.Trig./Time / Ext. Trig./Event / Ext. Trig./Level / Ext. Trig./Stop	Mode is comparable to a combination of Timing and trigger mode. These modes are special trigger modes of the C4880 standard camera. Please refer to the camera manual and the related timing diagrams for details.
BG sub	Checked/Unchecked	Starts or stops a background subtraction function inside the camera
Recursive filter	Checked/Unchecked	Starts or stops a recursive filter function inside the camera
Recursive filter number	Integer number	Defines the recursive filter number for the Recursive filter operation
High voltage	Checked / Unchecked	Applies or switches off high voltage
EB-Gain	Integer number	Sets the voltage on an electron bombarded CCD.
Sensitivity	Checked / Unchecked	Switches on or off the sensitivity parameter of a camera
Subtype	String	Defines which type of Flat panel is used. Normally this type is detected automatically.
New types	Checked/Unchecked	Sometimes Flat panels have older and newer types of the same device. Newer types normally have a Lemo connector at the camera front side.
Frame time	Any time in physical units	Interline CCD cameras: For interline cameras exposure and readout can be performed simultaneously. For exposure times \leq readout time the frame time is always the readout time. For exposure times $>$ readout time the frame time is always the exposure time. Full frame transfer CCDs: For full frame transfer CCD cameras exposure and readout is always done separated. This means frame time = exposure + readout time. CCIR/EIA video type cameras: The frame time for standard video readout cameras the frame time is always fixed (40 ms for CCIR and 33 ms for EIA).

The following parameters can be found on the camera acquisition dialog. They are however no parameters which are contributed by the cameras behavior, but are

parameter which are due to the acquisition process, therefore we call it acquisition parameters.

Start acquisition for sequence mode	"Next field", "Odd field", "Even field"	This parameter refers to images acquired by an analog CCD camera. It defines whether during sequence acquisition of full frames should start with an odd field an even field or the next field.
Wait for 2nd Frame in Acquire mode	Checked/Unchecked	This parameter defines whether for acquisition of flat panel images in acquire mode the system should ship the first frame and record only the second
Frame grabber trigger	Checked/Unchecked	This parameter defines whether acquisitions from an analog camera should start with an electrical signal applied to the FG trigger input (FG=frame grabber). This parameter applies only to standard video camera signals.
Frame trigger	Checked/Unchecked	This parameter defines whether acquisitions from a line sensor or TDI camera should start with an electrical signal applied to the frame trigger input.
Integrate after trigger	Integer number	This parameter defines whether the system should average a specified amount of frames after the reception of a trigger.
# of exposures	Integer number	In analog integration mode this parameter defines how many frames are summed up in memory
Clear on Start	Checked/Unchecked	This parameter defines whether the currently selected image should be cleared when starting image acquisition. If this parameter is not selected one can continue to acquire data starting with the data in the currently selected image
Photon counting threshold	Integer number	Threshold above parts in the image are regarded as photon spots
Use S/W trigger with single acquisitions	Checked/Unchecked	Cameras operated by DCAM module can have the option to be operated in Software Trigger mode. This is a mode where the software can control the start of an acquisition (Normally the exposure it either started by an internal mechanism or by external trigger). Together with single acquisitions this feature is used if the option "Use S/W trigger with single acquisitions" is selected. Otherwise software trigger mode is not used.

Rolling average	Checked/Unchecked	An image processing function which reduces the noise by combining the current image with the previous image (Live mode only).
Number (for Rolling average)	2- 255	Indicates the strength of the Rolling average function. The higher the number the stronger is the noise reduction but the more lag is introduced to the image.






Acquisition modes

Four image acquisition modes (Live, Acquire, Analog Integration and Photon Counting) can be controlled from the camera acquisition dialog.

There are two ways to open this control:

1. Choose the desired acquisition mode from the **Acquisition** menu.
2. Click on the desired acquisition mode button as described below




	Live mode	If the Live mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Live mode dialogue is open, image acquisition in Live mode will be started. The acquisition parameters as shown in the Live mode dialogue are used.
	Freeze Live mode	Live mode stops. The last image will be stored.
	Acquire	If the Acquire mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Acquire mode dialogue is open, an image acquisition in Acquire mode will be executed. The acquisition parameters as shown in the Acquisition mode dialogue are used.
	Analog integration	If the Analog integration mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Analog integration mode dialogue is open, an image acquisition in Analog integration mode will be executed. The acquisition parameters as shown in the Analog integration mode dialogue are used.
	Photon counting	If the Photon counting mode dialogue is not opened, it will be opened as you click this button for the first time. If you click it while the Photon counting mode dialogue is open, an image acquisition in Photon Counting mode will be executed. The acquisition parameters as shown in the Photon counting mode dialogue are used. See.

The acquisition buttons work different depending on the display status of the acquisition window:

If the window is closed, it will be opened by clicking one of the buttons.

If the window is already opened, the corresponding acquisition function will start after pressing an acquisition button. The acquisition will be made with the parameters set in the corresponding tab control.

The camera acquisition dialog has four tabs to control the different acquisition modes.



The camera and acquisition parameters of every acquisition mode is stored internally and reactivated if the acquisition mode is selected again.

The camera acquisition dialog can be changed in its size using the two pushbuttons on the right side of the dialogue.



If you want to see the whole window with all options, you have to choose the largest size.

However after you have finished initial settings, you may not need to see all setting features of the control. Then you can make it smaller and display just the controls you need for your work (e. g. start acquisition and set exposure).

If the Acquisition tab control is opened once, you can switch between the different modes either by selecting the desired mode from the main menu or by selecting the tab corresponding to the desired mode.

Each acquisition can be started by a hot-key and stopped by pressing the ESC key.

Hot-keys for start of an acquisition are:

- Ctrl +L Live
- Ctrl +F Freeze
- Ctrl +A Acquire
- Ctrl +I Analog integration
- Ctrl +C Photon counting

Live

Live mode is a mode in which images are acquired continuously. It is intended for monitoring images with maximum speed. It differs from other acquisition modes from display window handling. Live images are always displayed in the same window, while in other modes each new image can be displayed in a new display window.

The acquisition speed is displayed in the status bar. This feature can be used to control whether the camera is operated correctly.




Image quality can be improved by using real time background subtraction (RTBS), real time shading correction (RTSH) and Rolling average.

Acquire Mode

Acquire mode is very similar to Live mode Single image with two exceptions:

- 1.) It can contain different camera parameters as the Live mode. Thus it is easy to use two different camera parameter sets just by switching back and forth from Live to Acquire mode.
- 2.) Acquire mode does not use always the same window. Normally with every new acquisition the measurement data is written to a new image and displayed in an individual window. This prevents overwriting of already acquired data. Note: There is an option in the images options called “Acquire always to the same window”. If this option is checked all new acquisitions are directed to one image and are displayed in the same window which means that the previous image data is overwritten without warning.

Difference between Live and Acquire mode

Acquire mode is mainly intended to acquire measurement data whereas live mode is mainly intended for checking and adjustment purposes. In the case of a so called dual scan cameras like the C4880 and ORCA II there is an additional feature associated with the different acquisition modes:


Live mode and photon counting mode is operated with fast scan mode whereas Acquire mode and Analog integration mode is operated in slow scan mode. This enhances the meaning of Live mode as a fast operating mode for adjustment and checking purposes where Acquire mode is a measurement mode in which the camera is used in a high precision mode with maximum bit range.

Note: The feature that Live and Photon counting mode is operated in fast scan mode and Acquire and Analog integration is operated in slow scan mode is only available if the camera is used with internal drivers. If the camera is operated by a DCam module the setting of scan mode is done in the camera options and is valid for all acquisition modes.

Analog Integration

Analog Integration is a mode where a number of images from the camera are accumulated in the frame memory up to 16 or 32 bit depth. (16bit/32 bit).


The acquisition speed is displayed in the status bar. This feature can be used to control whether the camera is operated correctly.



Photon Counting

Photon Counting is a mode in which single photon events are added up in the frame memory. If the signal integration time is sufficiently long, a very high signal-to-noise ratio can be achieved. This mode also effectively suppresses certain crosstalk effects inherent to some image intensifiers. A dynamic photon counting mode allows temporal analysis of photon counting images.

The acquisition speed is displayed in the status bar. This feature can be used to control whether the camera is operated correctly.



See also the chapter Photon counting in this document for details.

Camera data types and data storage

Depending on the data delivered by the camera and the acquisition mode the result image will be stored in different types of images.

The data type of cameras range from 8 to 16 bit.

Data type of cameras

In practice the following data types can be delivered from the camera:

8, 10, 12, 14, 16 bit (all unsigned values)

Storage types

Image data is stored in three different types:

8, 16 and 32 bit

In other words: One pixel is either stored as 1 byte, 2 bytes or 4 bytes. More bytes per pixel means a higher range of pixel values but also means that an image needs more memory in RAM and on hard disk.

Depending on the data type of the camera and the acquisition mode the following storage types will appear.

Live and Acquire mode

Camera	Image	Value range
8	8	0 to 255
10	16	0 to 1023
12	16	0 to 4095

14	16	0 to 16393
16	32	-32768 to 32767

Analog Integration with option “32 bit images in Analog Integration” not selected

Camera	Image	Value range
8	16	-32768 to 32767
10	16	-32768 to 32767
12	16	-32768 to 32767
14	16	-32768 to 32767
16	32	Depending on integration number

Analog Integration with option “32 bit images in Analog Integration” selected

Camera	Image	Value range
8	32	Depending on integration number
10	32	Depending on integration number
12	32	Depending on integration number
14	32	Depending on integration number
16	32	Depending on integration number

Photon Counting with option “32 bit images in Photon Counting” not selected

Camera	Image	Value range
8	16	0 to 32767
10	16	0 to 32767
12	16	0 to 32767
14	16	0 to 32767
16	16	0 to 32767

Photon Counting with option “32 bit images in Photon Counting” selected

Camera	Image	Value range
8	32	Depending on integration number

10	32	Depending on integration number
12	32	Depending on integration number
14	32	Depending on integration number
16	32	Depending on integration number

Corrections

Corrections are used to compensate artifacts of CCD cameras, optics and other components used in the system. See also the chapter corrections in this document for more details.

Background subtraction correction

Most cameras delivers a nonzero signal even if they are not illuminated with light.

This has mainly two reasons:

- 1.) Image data is normally digitized by an A/D converter which outputs only positive numbers (including zero). Therefore it is an advantage to adjust the A/D offset in a way that already with no illumination there is a nonzero output. Otherwise small signals still leads zero as an output. As a result always a kind of noise or fixed pattern is visible in uncorrected images.
- 2.) Even with no illumination the sensor accumulates some charges which leads to a nonzero signal

As a consequence the most important correction method is a function which subtracts a dark data from the measurement data. See also the chapter Corrections - Background Subtraction in this document.

Standard background subtraction

Standard background subtraction is a function which subtracts a dark image from the given data. This is always executed after the image acquisition as a post process. The dark image can be either directly acquired from the camera or it can be read from an image file. There is an option to automatically call this post process after the image acquisition. This automated process is always called when an image acquisition finishes. This is:


In Live mode when the live mode has been freezed.

In Acquire mode when the image is completely acquired.


In Analog integration mode and photon counting mode when all images which should be summed up in memory are acquired

To operate background subtraction from camera

- 1.) Select Camera as the background source in the background options.



- 2.) Acquire an image
- 3.) Start background subtraction by either executing background subtraction from the menu:



- or -

Execute the Background subtraction icon on the toolbar:




RTBS (real time background subtraction)


Sometimes a background subtraction makes a so big difference on the image result that the user wishes that such background subtraction is also performed during live mode. This is a feature which can be called directly from the camera acquisition dialog. The data used as a background image can be either come from the camera directly or from an image file.

To execute real time background subtraction proceed as follows

- 1.) Start Live mode
- 2.) Select RTBS on the camera acquisition dialog.




If the background source is set to File the content of this file will be used



- 3.) Otherwise you may need to acquire new background data as follows:

Click Get BG Data.



This puts the next acquired image in the RTBS background buffer

Whenever you want to use a newer background image click Get BG Data.

Note: Real time background subtraction is also very useful to eliminate unwanted signal from an image.

Camera provided function

Some cameras have a built in background subtraction. In such cases the software only controls the camera feature but does not execute a background subtraction as an own function. Thus controlling the built in background subtraction with additionally executing the background subtraction functions from the software may lead to undesired results.

Cameras with built in background subtraction feature are

- C7190-21, C7190-23 and C7190-43 (video readout EB-CCD)
- C8000-20 (video readout BT CCD)
- X-Ray Line sensors like C7390, C8133, C8750, C9750 and C9133
- The C9100-13 (ImagEM)


Shading correction

The **Shading Correction** command corrects intensity non-uniformities of an image. In principle this is done by a pixel by pixel division of the intensity in the measurement image and the shading reference image. See also the chapter Corrections -Shading correction in the document.


To perform shading correction proceed as follows:

- 1.) Record a shading reference image. This image should represent the sensitivity distribution of the system.
- 2.) Perform background subtraction on this shading reference image.
- 3.) Save the image to hard disk.
- 4.) Specify the image as the shading image in the shading correction options. To do this click to the pushbutton Get at

the left side of the edit field for the shading file.



- 5.) Acquire your image data which you want to be shading corrected.
- 6.) Perform background subtraction on this image.
- 7.) Start Shading correction by either executing shading correction from the menu



- or -

Execute the Shading correction icon on the toolbar:



Note: Item 1) – 4) have to be performed only once if the system sensitivity does not change


Defect pixel correction

This function corrects defective pixel of a sensor by replacing them with not defective neighbors. Defects in this sense can be either single pixels complete lines or complete columns. A defect is regarded as either a hot pixel, which has too high intensity when not illuminated or a dead pixel which has too low intensity when illuminated. See also the chapter Corrections – Defect Pixel Correction in this document.

To activate the defect pixel correction proceed as follows:

- 1.) Create a defect pixel correction file by using the defect pixel correction tool and save it to a disk file. See the chapter Defect pixel correction tool.
- 2.) Specify the name of the defect pixel file in the defect pixel correction options and activate the checkbox Defect pixel


correction:



With every new acquisition the data is automatically corrected according to the information in the defect pixel file.

Rolling average

While all the a.m. corrections improve the image quality there is another operation which also improves image quality without being a correction in the strict sense: It is the **rolling average** function which is available in live mode only.



Rolling average is an image processing function which reduces the noise by combining the current image with the previous image (Live mode only).

The number which can be set from 2 to 255 indicates the strength of the **rolling average** function. The higher the number the stronger is the noise reduction but the more lag is introduced to the image.

Saving and loading images

This topic describes how to save or load images. Please see the Appendix Image File Formats for details.


Saving images

The **Save As..** command opens a dialogue where you can input file name, file type, drive and directory name for the image file to be saved. The complete image will be saved.

The **Save ROI As..** command saves the partial image within the currently selected rectangular ROI.

Choose **Save As...** from the **File** menu.

If you want to save a part of the image only, you have to select the ROI **before** you open the **Save ROI As...** dialogue!



Save As... dialogue


If you select **C**ancel the program cancels the operation.

The dialogue behaves like a standard Windows dialogue. It allows to select the directory where the file shall be saved, to define new directories, to display the file list in different formats etc.

After selecting the desired directory, you have to type the desired file name in the **File name** textbox.

Note: The correct file extension (e.g. “.tif” or “.img”) will be automatically appended to the file name. If you type a file extension which is not identical to the file extension as defined by the file type (e.g. you type “image.pix”), the correct file extension will be appended automatically (e.g. “image.pix.tif”).

Image data are saved in a format as selected in the **Save as type** list. Following formats are supported:



Save as type list

The file formats for image sequences will be available only, if the current image is a sequence image.

See the appendix of a description of all these file formats.

In the text box below the **Save As..** or **Save ROI As...**dialogue you can input a comment. This comment will be saved to disk together with the image file and it can be used to identify the image later on.

Loading images

This command allows loading and displaying an existing image file.

The **Open...** command pops up a dialogue similar to the **Save As...** dialogue. You have to select an image file type which you want to load and display. To select a file you can either input a file name in the file name text box or double-click on the filename in the file list.

If the selected file is part of an image sequence, the sequence, starting at the chosen file number will be loaded.

Following file types can be read by this program:

ITEX files and sequences, TIFF files and sequences (however not 24 bit color TIFF images), DPC files and DPC sequences.

File info

As soon as you have selected a file, detailed information about the actually selected file will be displayed in the **File info** table below the dialogue. The content of this dialogue is same as the **Image Status** display explained below.

Corrections

This topic describes the details of all image correction methods.

Overview

In the corrections menu you find commands which are related to image data corrections:

Background Subtraction	Strg+B
Shading Correction	Strg+S
Background/Shading assistant	

The menu items which can be executed from the correction menu are:

Background Subtraction, Shading correction and Background/Shading assistant

The Real time background subtraction and the real time shading correction can be invoked from the camera acquisition dialog.

The defect pixel correction is automatically performed once it has been activated in the defect pixel correction options or can be applied after image acquisition (offline defect correction).

From version 8.1 a background/shading assistant can help you to handle background and shading corrections easily.

Background Subtraction


The **Background subtraction** command subtracts a background image from the current image in the frame buffer. How this is done is determined by the settings in the background correction options.

If you want to automatically subtract background after each acquisition, you can select the **Auto Backsub** item in the background correction.

In Live Mode there is also the possibility to perform Real time background subtraction (RTBS), which displays the corrected image after every single acquire image.

The Background subtraction function allows you to subtract the camera dark current from an acquired image. The function can also be used to subtract offset signals of other origin such as undesired background light or stray light.

Background subtraction can be performed either "from camera" or "from file" depending on the setting in the background correction options.



Background Subtraction from Camera

When the option **Camera** is selected a new background image is acquired after each image acquisition and then subtracted from the acquired image. This means, the background image is taken after the data acquisition.

The background image is automatically acquired with the same parameters (such as exposure time and camera gain) as the acquired image. This method is easy to use and convenient if the exposure time is not too long. (In case of long exposure times you may wish to use the "from file" method instead.).

Following description refers to a configuration where a C4880 type camera (except C4880-8X and C4880-40) with a mechanical shutter is used.

The check box **Acquire image with open shutter** defines whether the C4880 shutter will be closed or opened during background acquisition. In most cases you will prefer that the shutter is closed automatically during background acquisition since this is the most convenient way. However, if you want to subtract also background light which has an external origin (such as stray light), you could disable this automatic shutter control by selecting the check box.

Following description refers to a configuration where an analog video camera or another camera without mechanical shutter (e.g. C3077, C4880-80, C4880-40 or C4742-95) is used.

If **Don't prompt the user before backgr. sub** is not selected the user is prompted to close any external shutter or switch off the illumination light or take a similar action to avoid that the light signal comes to the camera during background acquisition. If the option **Don't prompt the user before backgr. sub** is selected the user is **not** prompted to take any such action and Background subtraction immediately proceeds.

Background Subtraction from File

In case of long exposure times, the "from camera" method may become inconvenient since it would take the same time to acquire the background image after every data acquisition. In these cases, it is preferable to save a background image to hard disk once and load it every time it is needed. Of course, this makes sense only if you can assume that the same background image can be used for a series of data acquisitions (i.e. the acquisition conditions, like exposure time, will not change considerably.)


To perform background subtraction "from file" select the option **File**.

In some cases the readout parameters are different for different acquisition modes. In such cases you want to specify different files for the different acquisition modes (special case).


In many other cases (standard case), however, we do not want to specify several files for different acquisition modes and we assume that one file is suitable for every acquisition mode.

Note: When using the "from file" method the user must take care by himself that the background image contains suitable background data.

There is a checkbox where we can distinguish between these two cases. The default setting is that we do not distinguish between different acquisition modes and want to specify one single background file for every case.



Setting which uses one background file only (default)



Setting which uses different background files for different acquisition modes.

Press **Get** to open a file selection dialogue and select the desired file.

Save a background image

To obtain a valid background image you simply take an image under background light conditions with the desired image parameters which you want to use for your measurement, e.g. correct exposure time, scan area, binning factor etc. Save this image by using the **Save as...** command.

Warnings

If the file which you defined as background image file (see below) does not fit to your measurement conditions (e.g. if the binning factors are different) you get an warning message when executing **Background subtraction** and you are prompted to prepare a suitable file and try background subtraction again.

In some case, e.g. if temperature conditions or exposure conditions are very much different for the image in the background image file and in the frame buffer, you will get a warning message before the subtraction is performed and you can decide whether to proceed anyway or you want to prepare another background image. This should be just a hint from the program that the data may not be well suitable, but -as said above- the choice is yours.

Note: Background subtraction "from file" always subtracts pixels which have the same locations on the frame buffer one by one. Images which serve as background images should not have been shifted to another location (e.g. by loading it to a

different location and saving it), and conditions like binning factors and optical black parameters should not differ.

Details About Background Subtraction


Note: This chapter contains information about how the background subtraction process works.

If the original image or parts of it and the background image have similar intensities and/or if the original image is completely dark on some areas, negative intensity values may result in several pixels after background subtraction. There are, however, some reasons not to allow negative values in the image:

- Negative intensity values often do not have a physical meaning (there is no "negative light intensity").
- Image processing programs often cannot interpret negative values correctly.

Considering all these arguments HiPic handles the question of negative values as follows:

The user can select if the images should be clipped to zero or not during background subtraction by selecting or deselecting the checkbox **Clip to zero during background sub.** 8 bit images are always clipped to zero.



Considering Images which have acquired with a different integration count

This software allows integrating several camera images to one single image in Analog Integration mode. Software Versions prior to Version 8.1 did not allow to use images with different integration count for background subtraction, which leads to the necessity to acquire new background data for a different integration count. Version 8.1 and higher is now able to use any integration count number for the background image (as long as other parameters like image area or binning are suitable, as was described already earlier).


To increase precision and reduce noise it is recommended to use an integrated image as background image. The integration count should at least (or better higher) as high as in the data which we want to correct.

Constant

You can add a **constant** during background subtraction process.

This is especially useful when normal acquisition modes are used and negative values would be clipped to zero. By adding a constant you will see the image noise correctly even in dark areas. To add a constant during background

subtraction processing you simply enter the desired offset value in the edit box **Add constant after background subtraction.**



If you add a positive constant it may happen that the current bit range is exceeded. In this case the values are clipped to the maximum possible value within the bit range (e.g. 255 in 8 bit images, 1023 in 10 bit images and 4095 in 12 bit images). Therefore it is recommended to choose a constant not higher than necessary (the amount of noise is normally only a few counts). The constant is used for "from camera" type subtraction as well as for background subtraction "from file".

Real time corrections

There are two type of real time corrections: Real time background correction (RTBS) and Real time shading correction (RTSH).

Real time background correction

Background subtraction is a method of subtraction images directly from live images prior to display it. It can be invoked on the Live tap of a camera acquisition dialog. If "File" is selected as the background source in the background options, the specified file is used as the default data for real time background subtraction.

The "Get BG Data" command acquires new background data from the camera. The toggle switch "RT Backsub active" starts the background subtraction process. or stops it.

Note: The Get BG Data command is only available when Live mode is running.



Let's assume you have selected "File" as the background source in the background options and you have acquired new data during real time background subtraction. If you then close the program or start another background subtraction you will be asked whether you want to save this newly acquired background data permanently and use it for any subsequent background subtraction.

Real time shading correction

Once the real time shading correction (RTBS) is active, the real time shading correction (RTSH) becomes available. It uses the file specified in the shading correction options (see below for details). To start real time shading correction click to the toggle switch "RT shading corr. Active"

Shading Correction

The **Shading Correction** command corrects intensity non-uniformities of an image in the frame buffer. How this is done is determined by the settings in the shading correction options.

If you want to automatically correct shading after each acquisition, you can select the **Auto Shading** item in the shading correction options.

If you want to correct every image display in LIVE mode you can also use the real time shading correction (RTSH, see above).

Shading correction is a method to flatten a spatially non-uniform sensitivity of the complete imaging system.

Due to several reasons an imaging system can have a non-uniform sensitivity at different locations. The reasons may be uneven illumination, lens shading (vignetting) or different sensitivity of the CCD pixels. If an image intensifier is used, the image intensifier may have a spatially varying light amplification.

All these effects can be compensated by a multiplicative correction which we -by convention - call shading correction.

Get a Shading Reference Image


To be able to perform shading correction one must first acquire a shading reference image. A shading reference image is taken from a scene which is completely flat (homogenous illumination) within the field of view.

First, acquire this image by the imaging system, and then perform background subtraction on this image. Then you have to store the reference image using the **Save as...** command. Pay attention that the shading reference image should not have an intensity value of 0 at any pixel.

Sometimes, especially if one has to work under low light level conditions it may be very difficult or even impossible to obtain a shading reference image. In this case, shading correction can not be applied.

Declare an Image as Shading Reference Image

Specify the image as the shading image in the shading correction options. To do this click to the pushbutton Get at the left side of the edit field for the shading file.



Algorithm of Shading Correction

General

During shading correction the following calculation is performed:

$$C_{(x,y)} = \frac{D_{(x,y)} * K}{S_{(x,y)}}$$


C:	Corrected image data (output image)
D:	Uncorrected image data (input image)
K:	Constant
S:	Shading reference image

The calculation is performed with long integer arithmetic.

Constant K

The constant **K** can be specified by the user. At locations where $K = S_{(x,y)}$ the data will not be changed, at locations where $S_{(x,y)} < K$ the image intensity will be enhanced after shading correction, and at locations where $S_{(x,y)} > K$ the image intensity will be reduced after shading correction.

Depending on the setting in the shading options either the Upper LUT value of the shading reference image is used as the constant **K** or this constant will be calculated automatically according to a well defined criteria (Recommended).



When using the upper LUT value of the shading reference image the following is normally valid:

When setting the LUT values by the AUTO LUT function the upper LUT cursor will be set around the maximum value of the shading reference image. Normally shading data have a maximum intensity at the centre of the image and will slightly decrease at the image borders. Thus shading correction with a shading reference image where the LUT cursors have been set by the AUTO LUT function will keep the image intensity in the centre and enhance it at the borders. In almost every case this is the most convenient and easiest way of specifying the constant **K**.

When calculating the constant **K** automatically the average value of the center part (half size in both dimensions) of the shading image is used as the constant **K**.

Shading correction automatically accounts for chip type, frame-buffer type and binning factor, thus you can compensate an image which has been taken from the camera in normal mode with an image taken in 2 x 2 binning mode. Generally speaking the correction is always done with the pixels from the shading reference image which corresponds to the same location **on chip**.

The normal clipping is performed when the data exceed the current bit range. When zero is found in the shading reference image the corresponding pixel is not corrected (Shading reference images should not have zero values as a contents). The shading data should be background subtracted and should **not** be shading corrected itself. Otherwise an error message will be issued and the shading correction is not performed.

Defect Pixel Correction

The defect pixel correction corrects defective pixel of a sensor by replacing them with not defective neighbors

In a first step (by using the defect pixel correction tool), bad pixels, lines or columns of the camera are detected. The co-ordinates of the bad pixels are then stored in a special data file.

During normal operation with the sensor the defective pixels are then corrected if the pixel correction function is enabled.

The intensity value of a defective pixel, line or column is calculated by the intensity values of the neighboring pixel.

Defect pixel correction tool

The defect pixel correction tool is designed to generate coordinate data from image data showing defects. These coordinates are then used to correct images when they are acquired.

The procedure of getting coordinates of defects has to be done just once for a detector. In later sessions the stored pixel mask data will be used for correction.

Step 1

Acquire and store images with hot and/or dead pixel:

The recommended procedure is to acquire an image with a homogeneous illumination (only small grey level variation) at average pixel intensity of 50% of the maximum value. Whenever possible it is recommended to acquire the image in Image Accumulation mode with 8 or 16 images integration. Use the camera binning mode settings and exposure time settings which you will typically use for routine operation.

Save the image on your hard disk.

This image will be further used to calculate dead (pixel which have significantly reduced sensitivity) and hot pixel (which have a high dark current).

If you will use the detector with several binning modes, you have to record one image of any binning mode.

If you are planning to correct for hot pixel only, you may acquire an image without illuminating the detector (dark image).

You can also acquire a dark image and an image with illuminated detector and later use the dark image for hot pixel correction and the illuminated image for dead pixel correction.


Step 2

Defect Pixel localizing tool

Select **Defect pixels...** from the **Tools** menu to open the Defect pixel correction tool.



The Defect pixel correction tool will appear.



Step 3

At first you have to decide which method you want to apply. Select **Calculate hot pixels...** if you want to correct the pixels which have a higher dark signal than average pixel. Select **Calculate dead pixels...** if you want to correct the pixels which have a lower sensitivity than average pixel. Select **Both** if you want to correct hot and dead pixel (this is usually recommended).

Step 4

Select the defective pixel images in the defect pixel tool. If you have acquired different images for hot and dead pixel, please select the corresponding images in the sections.

Step 5

The system will suggest a reasonable threshold for discriminating normal pixel from bad pixel. The threshold is an intensity value which is related to the average intensity value and the standard discrimination of the reference image. The value suggested by the software can be manually changed if the correction map is not reasonable upon customer's considerations.

Press **ReCalculate** in the Hot pixel and/or Dark pixel section to calculate the suggested threshold value.

A separate value for singular pixel and lines/columns can be selected. These values may have to be chosen upon the features of individual sensors. If it turns out that the recommended values do not lead to satisfying results manual change is recommended. We can not give general recommendations for optimizing the threshold. User has to change and try.

An additional threshold is applied for lines/columns. Lines/columns are considered as defective if the percentage of defective pixel in a line or raw exceeds the value defined in Threshold lines/Columns.

Usually a value greater than 80% leads to reasonable results.

Step 6

After completing the threshold settings you have to click **Calculate defective pixel**. Then the defective pixel mask will be calculated. You can see the number of defective pixels and lines/columns indicated below the Calculate defective pixel button. The calculation will take a few seconds. As soon as „finished“ is displayed the calculation has ended.

Step 6a

Optionally you can shot the defects by clicking to the **Show** command. They are displayed as overflow values in either the hot or dead pixel file (A new image is created, the old image is not overwritten).

Step 7

Finally you have to save the defective pixel mask in a specific file (*.ini file).


Press **Save to *.ini file** and then select file name and directory.

If you additionally want to activate this defect pixel information you can use the **Save to *.ini file + Activate** command.

Note: When you generate different pixel masks with different binning factors of the same detector you have to choose the same file for all pixel masks. The data are then all collected in one *.ini file.

Enabling defect pixel correction


Specify the name of the defect pixel file in the defect pixel correction options and activate the checkbox Defect pixel correction:



With every new acquisition the data is automatically corrected according to the information in the defect pixel file. If you select **Save to *.ini file + Activate** on the Defect pixel tool this is automatically done.

Offline Defect Pixel Correction

There is also an offline defect pixel correction available for images which have not been corrected during acquisition. This offline defect pixel correction can be executed with the context sensitive menu.



If you execute this command your will be asked which defect pixel map you want to use for this correction. The file selected in the “Defect Pixel Options” will be selected as the default.

Background/Shading correction assistant

To make the process of acquiring data and defining the settings for background correction or background and shading correction more easy and getting familiar with these corrections the background/shading correction assistant can be used. It guides you through the process of acquiring background/shading data and makes

all necessary settings for you. It is strongly recommended for beginners to use the default/recommended settings initially. Once you got familiar and used to these settings you can easily use other more specific settings as well.

The Background shading correction assistant can be invoked with the Correction - Background/shading correction menu entry


Background Subtraction	Strg+B
Shading Correction	Strg+S
Background/Shading assistant	

A dialog will appear and asks you for several settings. We recommend to use the default/recommended settings. We also recommend to start with background subtraction only first and once you are familiar with it use Background/ and Shading correct.


Performing Background subtraction

To perform background subtraction with the background/shading correction assistant invoke the background/shading correction assistant with the menu command “Background/Shading assistant” from the corrections menu.


The following dialog will appear:




Select Background only (default) and click to Next. The next dialog will appear:



Select File (default) and click to Next. The next dialog will appear:



Select Manual (default) and click to Next. The next dialog will appear:




You should now acquire a background image which represents the background of your measurement. This image should be taken under the same conditions as the image you want to acquire for your measurement (i.e. binning, image size, analog gain settings etc.). Please use the camera acquisition dialog to acquire the image now. It is recommended to use the same acquisition mode as you want to use for your measurement (LIVE, ACQUIRE, ANALOGINTEGRATION).

Alternatively you could also use an image you have already saved to disk. Please load this image now or make it the current image (by clicking to it) if it has been loaded already.

If this is finished, click to next. The software will save this image to a default file name and select this file in the background options automatically. You will get an informative message similar to the following:





After clicking to OK, you will get the next dialog:




At this time you are finished with acquiring and setting up and you can acquire measurement data and correct them with the acquired background image.

To subtract background data from a given image you can use the Background subtraction menu item from the corrections menu



, the Background Toolbar Icon , the context sensitive menu of the image which you want to correct (Click with the right mouse button to the image), the Real time background subtraction feature or the automatic background subtraction mode (this has to be selected from within the background options).

Clicking with the right mouse button to the background/shading assistant dialog



will allow you to go directly to the background options. These options contain more specific settings related to background subtraction.

Performing Background subtraction and Shading correction

We assume that you are familiar with acquiring background data with the Background/Shading assistant. If not, please read the previous topic.


The Background shading correction assistant can be invoked with the Correction - Background/shading correction menu entry

Background Subtraction	Strg+B
Shading Correction	Strg+S
Background/Shading assistant	

A dialog will appear and asks you for several settings. We recommend to use the default/recommended settings.

To perform background subtraction and shading with the background/shading correction assistant invoke the background/shading correction assistant with the menu command “Background/Shading assistant” from the corrections menu.

The following dialog will appear:



Select Background & Shading and click to Next.

Select File (default) and click to Next. The next dialog will appear:

Select Manual (default) and click to Next. The next dialog will appear:

You should now acquire a background image which represents the background of your measurement..

If this is finished, click to next. The software will save this image to a default file name and select this file in the background options automatically. You will get an informative message.

After clicking to OK, you will get the next dialog:




Acquire a shading image now. This shading image should represent the sensitivity variation of your system, which you want to account for. You can use any acquisition mode for this acquisition, but you should select a mode which outputs the full camera image (e.g Subarray mode is not suitable).

Alternatively you could also use an image you have already saved to disk. Please load this image now or make it the current image (by clicking to it) if it has been loaded already.

If this is finished, click to next. The software will save this image to a default file name and select this file in the shading options automatically. You will get an informative message similar to the following:




After clicking to OK, you will get the next dialog:




At this time you are finished with acquiring and setting up and you can acquire measurement data and correct them with the acquired background and Shading image.


Subtract background data as described in the previous topic.

To perform shading correction from a given image you can use the Shading correction menu item from the corrections menu



, the Background Toolbar Icon , the context sensitive menu of the image which you want to correct (Click with the right mouse button to the image), the Real time Shading correction feature or the automatic Shading correction mode (this has to be selected from within the shading options)

Clicking with the right mouse button to the background/shading assistant dialog



will allow you to go directly to the shading options. These options contain more specific settings related to shading correction.

Reserved Filenames

As the Background/shading assistant saves and uses file to store background and Shading data there are several filenames you should not use by yourself. It is not recommended to replace or delete these files. The reserved file names are the following files within your application directory:

Reserved_BackgroundImage_General.img

Reserved_BackgroundImage_Live.img
Reserved_BackgroundImage_Acquire.img
Reserved_BackgroundImage_AnalogIntegration.img
Reserved_ShadingImage.img

Other Options

There are two more choices you can select if you run the background/Shading assistant.

In step 2 you can select “Camera” instead of “File”. In this case the background data is not saved to a file. Every time a background image is needed for a correction a new image is acquired from the camera. This setting does not influence the handling of shading data.

In step 3 you can select “Automatic” instead of “Manual”. In this case Automatic background subtraction and/or Automatic Shading correction is selected which means that after every acquisition immediately an background subtraction and/or Shading correction is performed. This options interferes with the real time background subtraction and the real time shading correction and should not be used if you intend to use the real time background subtraction and the real time shading correction.

There are several other options in the background and the shading options dialog. Please study these options if the standard behavior does not fit your needs.

Image Display & LUT


This topic contains details of image display and LUT and its operation.

Image display windows


Up to 20 images can be displayed at a time. Each image has its own status. All image processing functions refer to the currently selected image (called "**current image**" in this manual). If you click to an image it becomes the current image.

Each image has its own properties concerning the LUT and color selection. If several images are displayed on the screen and they have different color selection, not all images may be displayed correctly. Normally only the current image is displayed correctly. This behavior depends on the selected color depth of the Windows display.

The image is always displayed within a window with a specific zoom factor. Zoom factors from 1/8 to 8 are available. When the image with the current zoom factor does not fully fit into the window, scrollbars appear on the side and only a portion of the image is displayed (see).




The whole image displayed



Only a part of the image is displayed

When the window is larger than the area which is needed for the image a part of the window remains white.



When the image is smaller than the display window

You can freely select the size of the window. Sometimes one wants to set the window size to exactly match the image size. In such case it is sufficient to double-click to a grey area on the image window or to select the image and press the function key F2. When the window is smaller than the image there is a small grey area on the lower right corner where one can click.


The caption of the window contains information about the image (image mode, file name, zoom factor). Some functions resize the window automatically so that the image fits into the window automatically (e.g. Acquire, Open., Background Subtraction).

Caption

The image caption shows the image file name and the zooming factor.



If the image is not yet saved a default name is assigned derived from the acquisition with which it has been created and the internal image number.



Zoom

When the Zoom tool is selected the left mouse buttons can be used to change the zoom factor. The zoom factor is increased by a factor of approximately 1.4 when clicking with the left mouse button to the current image area. The zoom factor is decreased by a factor of two when clicking with the left mouse button while holding down the shift key.

The zoom factor ranges from 1/8 to 8.


The point where the user clicks to the image becomes the centre point of the image. Thus when you want to magnify a certain location it is sufficient to click to this location with the left mouse button.

If there are scroll bars on the right or bottom side of the window, just a part of the image is displayed. The scroll bars can be used to scroll to a location which is currently not displayed. It is however more convenient to zoom out with the left mouse button while holding down the shift key until the portion of the image


which should be displayed becomes visible. Then click to this location with the left mouse button to magnify it.

Rulers


Rulers can be used to display pixel numbers or physical units in the image. To show or hide rulers open the images options (right click in the image display and select Options...).



To select whether the rulers should display calibrated units or pixels select the calibrated checkbox.



After closing the Images Options with OK all images contains the desired rulers.




Quickprofile


The quickprofile is a profile which is directly displayed on the image window and updated with every new image display.

The quickprofile always shows the profile data extracted from the image based on the rectangular ROIs in this image.

Rectangular ROIs with full horizontal width (≡) are always displayed as a horizontal profile.

Rectangular ROIs with full vertical height () are always displayed as a vertical profile.


Arbitrary rectangle ROIs can be either displayed horizontally, vertically or not displayed. These three possibilities can be selected in the quick profile options:



To switch on the quick profile display toggle the status of the quick profile pushbutton on the toolbar or click with the right mouse button on the image display and select Quick Profile.




Normally the quick profile is displayed on the same area as the image.




Alternatively the quickprofile can also be displayed in a separate area outside the image display area. To select this alternative call the quick profile options and select Display quick profile outside of image.


Additionally the user must specify how large the display should be. A percentage has to be specified which determines how large the quick profile are should be relative to the image size.




The quick profile then appears outside of the image.




The size of the quick profile depends on the option settings for Adjust quick profile height:



If “Use LUT values” is selected the upper and lower LUT limits are used as the limiting values. If “Use Min/max values” is selected the minimum and maximum values are use as the limiting values. Portions of the quick profile which are beyond these limits are no displayed because they would lay outside the window which is used for the quick profile (This is either the image display window or a separate window neighboring the image display).



“Adjust quick profile height” set to “Use LUT values”



“Adjust quick profile height” set to “Use Min/max values”

Different from other functions not only the currently selected ROI contributes to the quickprofile display but all suitable rectangle ROIs. If the quick profile should be displayed outside of the image the outer areas appear and disappear automatically when needed.

Several quick profiles can be displayed in the same window. In this case the relation between ROI and Quickprofil can be easily recognized by the color.




FWHM


In some specific circumstances the user wants to get a Full width half maximum (FWHM) analysis of the currently selected ROI.

The result will be displayed on the left upper corner of the image.


This information can be used for checking and adjustment purposes. Its precision cannot be compared with quantitative analyses as data fitting, but it is a first hint of the situation.



The quick profile options contains some parameters to determine size, color and precision



The decision whether the FWHM should be displayed in pixels or in calibrated units is made in a parameter in the images options:








ROIs

Frequently you will have to select special areas on an image. For example you may sometimes want to save only a part of the image. In order to do this you will need a tool to select a special area of the image, a region of interest (ROI).


You can select up to 10 individual ROIs (Region Of Interest) for every image, but only one ROI is the currently active one.

Types of ROI


There are different types of ROIs which can be used for different purposes. The types of ROIs are:

-  Rectangle ROI
-  Arbitrary line ROI
-  Point ROI
-  Rectangle ROI with full size in horizontal direction
-  Rectangle ROI with full size in vertical direction


A Rectangle ROI is a ROI which starts at any arbitrary point and ends at another arbitrary point. The two points are combined to yield a rectangle ROI.




A Line is a ROI which starts at any arbitrary point and ends at another arbitrary point. The two points are combined to yield a rectangle.




A point ROI is a ROI defined at the point where the mouse is clicked.



A Rectangle ROI with full size in horizontal direction is a ROI which extends to the full horizontal width but can be arbitrary in its vertical starting and ending point.



A Rectangle ROI with full size in vertical direction is a ROI which extends to the full vertical height but can be arbitrary in its horizontal starting and ending point.



Create ROI

To create a ROI first of all select the ROI type on the toolbar:



Then draw the ROI with clicking to the starting point with the left mouse button, then draw the mouse to the end point of the ROI and release it. Once a ROI is created the ROI toolbar automatically switches to the pointer tool.

Move/Resize/Delete

To move the ROI select the Pointer on the ROI toolbar (after creation of a ROI this is automatically selected), click on the image display with the left mouse button and move the mouse.

To delete a ROI activate the image by clicking on its caption and type the Del Key or the d key.

To resize a ROI hold down the shift key and click and hold down the left mouse button and move the mouse until the ROI has the desired size.

Alternatively a ROI can be moved by using the cursor keys (left, right, up, down). Each time you press one of the keys the ROI is shifted by 1 pixel. If you press the CTRL key + one of the cursor keys, the ROI is shifted by 10 pixel, if you press CTRL+ALT + one of the cursor keys, the ROI is shifted by 100 pixel.


If you press the Shift key + one of the cursor keys, the lower right corner will be shifted (ROI size is changed).

If you want to expand an ROI over the entire horizontal or vertical image area, you can select an ROI and press H to extend the ROI over the maximum horizontal area, or V to extend it over the maximum vertical area, or F to extend it over the full image.

To precisely position or size an ROI the ROI Interface can be used.

Multiple ROIs

If the user wants to select several ROIs he has to select the desired ROI type again and create the next ROI as he did with the first one.



Select/Delete

Always one of the ROIs is the selected one and every image processing function (like the Save ROI as) is always using the selected ROI. The selected ROI is drawn with as dashed line.


A maximum of 10 ROIs can be selected at once and every ROI has its individual color.

To select a ROI different from the currently selected one just click with the mouse on or near the ROI or use the n key to change to the next ROI. Typing the n key (next) selects all ROIs one after the another.

Typing the Del-key or the d-key always deletes the currently selected ROI and automatically selects another one if there is still any ROI which is active.

Mouse Coordinates

While the mouse cursor is inside of an image window, the co-ordinates and the intensity value at the current mouse position are displayed.



Status display including actual pixel intensity information.

ROI I/F

The ROI Interface is a dialog which is designed to precisely define or read the coordinates of the currently selected ROI of the currently selected image.

To show the ROI I/F toggle the state of the ROI I/F button.




The ROI I/F will appear on screen and show the coordinates of the currently selected ROI. Changing the Selected ROI will automatically update the information in the ROI I/F.

ROI interface (rectangle ROI)				
	Start	End	Width	Unit
X	396	632	236	nm
Y	0	1023	1023	No unit
	Intensity	Area size	Diagonal	Slope
		241428		

Depending on the type of ROI either only the Start values (point ROI), or Start and End values (Line ROI) or Start, End and Width values (All three types of rectangle ROIs) can be inputted.

Image Status

Every image contains information about date and time, size, camera settings, calibration and many other parameters which are related to this image. To see this information select the Image Status menu from the file menu.




The image status will be displayed on screen.

The Image Status dialog box displays the following information:

- File info:** Name: Live [0], Type: 2 Bytes/Pixel, Size: 256 * 256
- Created with:** Softw.: HiPic 8.0.0 pf1, Camera: C4742-95-12NRK, FG: IC-PCI + AM-DIG, Date: 01-30-2006, Time: 16:51:51
- Exposure parameters:** Exp. Time: 31 ms, # Exp.: 1, # Trig.: 0
- Acquisition mode & LUT:** Acq. mode: Live, Data type: 12 bit (12 bit), LUT: 12 bit, B/W, Limits: 0, 4096, linear
- Corrections:** Background, Shading, Curvature
- Scaling:** X: Linear: 4 No unit/Pixel, Y: Linear: 4 No unit/Pixel
- Camera parameters:** [Camera]. AMD=N, NMD=F, EMD=E, SMD=S, ADS=12, SHT=1, FBL=1, EST=1, SHA=K, SFD=F, SPX=4, ATP=N, CEG=0, CEO=0, ESC=B, TimingMode="Internal timing", TriggerMode="Edge trigger", TriggerSource="BNC", VerticalBinning="1", TapNo="1", TriggerPolarity="neg.", CCDArea="1024 x 1024", Binning="4 x 4", ScanMode="Superpixel", NoLines=1024, CameraName="C4742-95-12NRK", Type=7, SubType=7
- Scan & other area properties:** Valid image [X,Y,DX,DY]: 0,0,256,256, Scan area [X,Y,DX,DY]: 0,0,256,256, Chip origin [X,Y]: 128,0, Frame buffer origin [X,Y]: 0,0, Binning factor [X,Y]: 4,4
- Comment:** (Empty text area)

Context sensitive menu


The context sensitive menu which appears if the user clicks with the right mouse button on the image display has the following entries:



Save As	Saves the image to file
Save ROI As	Saves the currently selected ROI of the image to file
Copy to new image	Copy the current image to a new image.
Copy ROI to new image	Copy the current selected ROI of the image to a new image.
Show Image status	Shows the image status
Zooming factor	Allows to set the zooming factor to the desired value
Do Auto LUT	Executes Auto LUT
LUT Color	Allows to set the LUT color to the desired value
Quick profile	Shows or hides the quick profile
Quick profile options	Shows the quick profile options
Do Background Subtraction	Executes background subtraction
Do Shading Subtraction	Executes shading correction
Options...	Shows the images options.

Print Images

The currently selected image can be printed by using the print menu command.




A printer setup dialog will appear. Make your selection and Click OK to print the image.

LUT

General

The LUT tool is used to control the image display by manipulating its brightness and contrast.



The LUT Tool

If you acquire images you will find that the images are frequently displayed with too low or too strong contrast. If you acquire images under low light level conditions the display may be too dark. Use the LUT tool to adapt the image display to the desired contrast. This is done by defining the lowest and highest grey value which shall be displayed.


The tool contains an area with two cursors which represents a lower (blue cursor) and an upper limit (red cursor) for the LUT. The meaning of the values of these cursors is described below. The cursors can be moved by clicking next to them and dragging them while the left mouse button is pressed. When you use the right mouse button instead, both cursors will move simultaneously, keeping their distance constant. An alternative way to move the cursors is provided by the two sets of slide bars at the left and right side of the LUT tool. Two display boxes (blue and red) display the exact numerical values of the two cursors.

You can easily see how the image contrast is changed when you change the upper or lower limit as the current image is updated after any change in the LUT control.

The bottom part of the LUT shows a scale representing the input intensity values. In 8 bit mode the scale ranges from 0 to 255, in 10 bit mode from 0 to 1023, in 12 bit mode from 0 to 4095 and in 16 bit mode from 0 to 32767.

At the centre part of the LUT tool (above the scale) a histogram is displayed where the user can see the frequency by which intensity values are present in the image. This histogram can be used as a guide for setting the LUT limits.

The LUT tool is re-sizeable horizontally. If the LUT tool is longer, the intensity scale can be seen and accessed with higher precision but the histogram display will need more time to update.



LUT tool when using 32 bit image data files

When 32 bit image data file format is used (e.g. with 16 bit digital cameras) the string “x16” in the lower left corner of the window indicates that the intensity figures shall be multiplied by 16 to give the real intensity data values.

Upper and lower limits


The software allows displaying images either in black and white or in pseudocolor mode (see). The meaning of the LUT limits is explained below for both modes.

Intensity values between the lower and the upper limit are displayed using the full range of the LUT (either in grey-scale or in color). The values above the upper limit (which exceed the current LUT) are displayed in yellow (in case of grey-scale LUT) or black (in case of color LUT).


The uppermost intensity value of the current acquisition mode is displayed in red (grey-scale LUT) or in white (color LUT). Areas in these colors indicate an overflow of the CCD camera. For 10 bit images the overflow value is 1023 and for 12 bit images it is 4095. Be sure to reduce your light intensity if you encounter this situation.

Display mode	<i>Smaller than blue cursor</i>	<i>Between cursors</i>	<i>Greater than red cursor</i>	<i>Saturation</i>
B/W	<i>Black</i>	<i>Black to white</i>	<i>Yellow</i>	<i>Red</i>
B/W inverted	<i>Yellow</i>	<i>White to black</i>	<i>white</i>	<i>Red</i>
Pseudocolor	<i>Black</i>	<i>Blue to red</i>	<i>Black</i>	<i>White</i>
Pseudocolor inverted	<i>Black</i>	<i>Red to blue</i>	<i>Black</i>	<i>White</i>
B/W without using colors	<i>Black</i>	<i>Black to white</i>	<i>White</i>	<i>White</i>
B/W without using colors inverted	<i>White</i>	<i>White to black</i>	<i>Black</i>	<i>White</i>

Table: Image color-coding in different display modes



Mapping of LUT values for black and white LUT



Mapping of LUT values for color LUT


LUT parameters

With this function you can select if images shall be displayed in black and white or pseudocolor.

Choose **LUT** from the **Display** menu to display the LUT parameters dialogue.



The LUT parameters dialog appears:



LUT parameter dialogue box

Within the dialogue box you can select the LUT type (either black-and-white or rainbow) by clicking to the **B/W** or **Rainbow** radio button in the **Color** frame.

B/W will produce a grey scale display where data overflow is marked by yellow and red color, while **Rainbow** will produce a pseudo-color display where different intensities are coded as different colors. **B/W without using colors for overflow** will produce a grey scale display where data overflow is not marked by yellow and red color.

Inverted look-up table

It is possible to apply an inverted look up table. This means that all intensity values are inverted (image looks like a negative image).

In order to switch from the normal to inverted mode, you have to open the dialogue **LUT** in the **Display** menu. Select the box **Inverted LUT**. The inverted look up table can be applied to linear and non-linear look-up table types.


Display negative values

In most cases negative images values, even though they may appear due to some operation like background subtraction, may not be of any interest for the display. In very special cases a user may want to display it as well. For this purpose the option “Display neg. values” is provided. If this option is selected, negative values will be displayed in the same brightness or color as the positive value with the same absolute value.

Non-linear contrast enhancement

The standard (default) enhancement mode offers a linear enhancement between the lower and upper limits.

Sometimes, however, a non-linear enhancement function may be preferred. This program supports two non-linear modes: Gamma and Sigmoid.



LUT dialogue with linear contrast enhancement

In order to switch from the linear to a non-linear mode, you have to open the dialogue **LUT** in the **Display** menu. Select the desired enhancement function in the **LUT Type** field. Then set the enhancement factor (gamma-factor or sigmoid

factor) by clicking and dragging the enhancement curve in the graphic box. Alternatively you can insert the desired factor into the textbox **Gamma**. Press the enter key to activate the new factor and recalculate the image (see also).

Choose LUT type according to your desire by selecting the related radio button in the LUT type section.

Set the Gamma or Sigmoid parameters either by clicking in the box where the transfer curve is drawn and dragging the curve (current image will be updated immediately) or by inserting the desired factor in the textbox Gamma followed by pressing the enter key to activate the new factor.

Optimize contrast enhancement

By clicking the small asterisk at the lower right side of the LUT tool or the corresponding button in the menu bar you can force LUT limits being calculated automatically according to the intensity values present in the image. In most cases this is an easy way to adapt the LUT settings to the current image and to obtain optimum digital contrast enhancement.


Selecting measurement parameters

This topic describes in detail how the user can influence the measurements and other preferences of the program.


Options

All options are accessible from one single options dialog. To call the options menu there are two ways:


- 1.) Call the options dialog from the menu:



- 2.) Click with the right mouse button to a dialog where you want to get the related options dialog. If for example click to the sequence dialog you will get the context sensitive menu for sequence which contains the Sequence options command.



General options




Restore window positions	If this option is checked, the program stores the position and size of every window when it is closed. The settings will be restored, when the window is opened for the next time
Background color: Set	A color selection window will be displayed where you can choose the background color of the program main window
Background color: Default	Sets the default background color of the program main window
USER function	When the UserFunctions option is checked, function calls to a DLL named CUSTOMER.DLL are enabled. With this option the user can implement his own control and data analysis functions into

CCD camera

The CCD camera options dialog shows all options which can be specified to define the CCD cameras operation. Depending on the type of CCD camera these options look different.

Acquisition



Wait till next display in LIVE mode [ms]:

Displays a new live image only if the specified time in ms has been elapsed since the last displayed image.

32 bit images in Analog Integration

If this option is set result images in Analog Integration mode will always be 32 bit.

Note: If the camera outputs 16 bit images the resulting images will be always 32 bit images (to allow storing signed values) independent of this option.

CCD Gain for Photon Counting

This option allows specifying the default CCD camera gain setting in photon counting mode. You can select between “Maximum gain” and “Minimum gain”. It is recommended to use “Maximum gain” unless the photon spots get saturated. The CCD gain can also be modified manually.

Write Dynamic Photon Counting file

If **Write dynamic photon counting file** is selected, the recording of photon counting images in the special DPC file format is enabled. In a DPC files the x-y co-ordinates of each photon and the time when it has been detected are recorded. This allows making a time dependent analysis of photon counting images.

32 bit images in Photon Counting

If this option is set result images in Photon counting mode will be 32 bit. Otherwise the result images will be 16 bit images.

Moiré reduction

This option specifies how strong Moiré reduction should take place. It is recommended to keep this parameters at its default value “Standard”.

Additional timeout [sec]:

The parameter **Additional timeout (sec)** can be used to define an additional timeout, after which an image

acquisition will be stopped, if no response comes from the camera.

Usually this time should be 0. Only in some special case when you get timeout errors before an image acquisition has been finished, you may increase this parameter.

Normally the timeout is automatically calculated:


Timeout = System defined timeout (derived from the exposure time) + Additional timeout (sec).

This parameter is very rarely needed in the current version of the software.

Deactivate grabber while not in use

Please keep this setting as it is (selected)

Calibration



Calibration method

Methods to tell the system which scaling it should apply for the X and Y axis. Existing methods are: No calibration. Square calibration, free calibration.

Horizontal

Calibration for the horizontal direction


Vertical

Calibration for the vertical direction

Current system calibration

The currently active calibration for the complete system

Sequence



Prompt before start

Normally a sequence acquisition is not started immediately after you click **Start Acquisition** in the sequence acquisition dialogue. Before the acquisition actually starts, an initialization procedure is executed. If **Prompt before start** is selected, the system will display a message box after the initialization is completed. Sequence acquisition will be started immediately after you click OK in this message box. This is useful if you want to precisely control the starting time of sequence acquisition.

Enable stop

If **Enable stop** is selected, it will be possible to stop sequence acquisition at any time. If it is not selected, sequence acquisition can be stopped only by the user function. Of course sequence acquisition will stop latest when the pre-selected number of images has been acquired. If **Enable stop** is selected, the sequence acquisition may be slightly slower on older computer.

Enable wrap around in Acquire to RAM mode


If **Enable wrap around** is selected, the sequence function will allow to acquire a sequence in the wrap around mode. This function is only available, when the **Enable stop** is selected.

Do not attempt to load a sequence to RAM

Normally loading a sequence will always try to load the sequence to computer RAM. If this succeeds operations and display can then be executed at maximum speed. However if after loading a part of the sequence there is not enough RAM all images are unloaded and only checked. If normally the sequences are too large to fit into the RAM this is just waste of time. If the option **Do not attempt to load a sequence to**

	RAM is active the system does not attempt to load the sequence.
Number of buffers	Number of buffers used for fast hard disk recording as an intermediate storage between acquisition and hard disk write.
Load HIS sequences after acquisition	HIS sequences are always written to a file. If this option is active the file is automatically selected as the current sequence.
Pack 10 and 12 bit images(increases performance)	10 and 12 bit images are packed internally, which increases the saving speed and the storage depth. Warning: If you want to perform background subtraction with the option “Clip to zero during background subtraction” not selected (in other words: if you want to create signed data) in a later step you should not use this option.
Always display image during acquisition	If this option is active the image is always displayed during the acquisition process.
Warning on	If warning on is active a display shows the remaining time in seconds until the next acquisition starts. Additionally the background color of the main window changes to red shortly before the next acquisition starts.
Play interval	Defines the time in ms between two images for the display in play mode.
Wrap during play	If during play of a sequence it encounters an end it restarts from the other end.
Do auto corrections after sequence	If the checkbox Do auto corrections after sequence is checked, all correction functions, as background or shading correction, are performed after a complete sequence has been acquired. A background image will be acquired only once at the end of the sequence. If this checkbox is unchecked, the corrections will be made after acquisition of each image. It may depend on your experiment, which mode to use: The mode Do auto correction after sequence allows a much faster processing of an image sequence. We recommend using it as default. However there may be circumstances which do not allow to use this mode, e.g. if the background image is changing after each acquisition. Then it is better to perform background subtraction immediately after an image has been acquired with a newly acquired background image.

Background



- Background source Defines the source of the background data. Available choices are Camera or File.

- Use individual background images for every acquisition mode (special case only) When we select to use File as a background source we can furthermore select whether we want to specify one single file for all cases (default) or whether we want to specify one background file for every acquisition mode. This can be useful if scanning options link scan speed are different in different acquisition modes.

- Take background image with open shutter Normally the CCD cameras shutter is closed during background subtraction from camera. If this options is selected the background image is taken with open shutter.

- Don't prompt the user before backgr. sub. If **Don't prompt the user before backgr. sub** is not selected the user is prompted to close any external shutter or switch off the illumination light or take a similar action to avoid that the light signal comes to the camera during background acquisition. If the option **Don't prompt the user before backgr. sub** is selected the user is **not** prompted to take any such action and Background subtraction immediately proceeds.

- Background file for: Live File for background subtraction in Live mode

- Background file for: Acquire File for background subtraction in Acquire mode

- Background file for: Analog Integration File for background subtraction in Analog Integration mode

- Add constant after background subtraction You can add a **constant** during background subtraction process. This is especially useful when normal

acquisition modes are used and negative values would be clipped to zero. By adding a constant you will see the image noise correctly even in dark areas. To add a constant during background subtraction processing you simply enter the desired offset value in the edit box **Add constant after background subtraction**.

If you add a positive constant it may happen that the current bit range is exceeded. In this case the values are clipped to the maximum possible value within the bit range (e.g. 255 in 8 bit images, 1023 in 10 bit images and 4095 in 12 bit images). Therefore it is recommended to choose a constant not higher than necessary (the amount of noise is normally only a few counts). The constant is used for "from camera" type subtraction as well as for background subtraction "from file".

Clip to zero during background subtraction

Negative intensity data may appear after background subtraction.

If this command is selected, negative intensity data are clipped to zero.

Background source for real time backsub

Defines the source for Real time back ground subtraction.


If camera is selected: Acquire a new background image in Live mode with your camera by clicking the button **Get BG data** in the **Live** menu while your camera is working in Live mode.

If file is selected: Use an image from file as background image (use this function if you have stored a background image on your hard disk before). After you press the **Get BG data** in the **Live** dialogue a file selection dialogue will be opened. Select the file which you want to use as background image for real time background subtraction.

Auto background subtraction

When the **Auto Backsub** option is checked, background subtraction is automatically performed after an image acquisition. Please refer to for details about Background subtraction.

Shading



Shading file

Image file used for shading correction

Shading constant

This option defines the method how to get the shading constant:

If “Calculate Automatically” is selected, the software calculates the constant according to the image data in the shading image.


If “Use upper LUT value of Shading file” is selected, the upper LUT value of shading file is used as the constant.

Note: In this case the upper LUT value within the shading file has to be selected correctly prior to saving the shading file.

Auto shading correction

If this options is selected a shading correction is executed automatically after image acquisition


Defect pixel



Defect pixel correction If this options is selected defect pixel correction is executed

Defect pixel correction file This file specifies the coordinates of the defects for the defect pixel correction

Images




Acquire always to the same window When the option **Acquire always to the same window** is selected the **Acquire** and **Analog Integration** functions always use the same image window for acquisition. This is

useful to avoid that the PC memory is quickly used up for storage of images. The default behavior concerning new images is as follows: The Live mode always acquires to the same window. The image load function opens always a new image. The Acquire and Analog Integration function opens up a new window every time they are invoked.

Default zoom factor	When a new window is opened a default zooming factor is used for the display.
Warn when unsaved images are closed	When a window is closed and the image is not or not fully saved the program issues a warning prior to closing the window (default). When the checkbox is not selected the program does not issue such a warning.
Calibrated (Rulers, ROIs, FWHM)	Defines whether the labeling of Rulers ROIs and FWHM are in calibrated units or in pixels.
LUT lower side = 0	When this option is checked the lower limit of the LUT will always be set to 0 when automatic LUT calculation is executed (e.g. by pressing the „*“ button of the LUT tool or by using AUTO LUT mode). If the option is unchecked the lower value will be calculated.
Auto LUT	When the Auto LUT option is checked the LUT is adapted automatically during image acquisition (e.g. during analog integration and photon counting), after acquisition, background subtraction and loading an image.
Auto LUT also during LIVE mode	When this option is checked the image contrast will be updated by the AUTO LUT function whenever a new image is displayed in Live mode.
Auto LUT optimized for selected ROI	When this option is checked the Auto LUT function will optimize the image contrast for grey values which are within the boundaries of the actual ROI. This allows optimizing the contrast based on the grey levels within a limited image area. If no ROI is selected, the whole image area will be considered as actual ROI. This function will be only active if Auto LUT is enabled.
Display horizontal ruler	Defines whether or not a horizontal ruler is displayed in the image
Display vertical ruler	Defines whether or not a vertical ruler is displayed in the image
Fixed header for ITEX images (10K)	The image status information is stored in a string in the image file header. Normally the length of this string varies with its content. However, in some cases it is desired that the header has a fixed length. If this option is checked the length of the header is always 10K.


Quick profile



Display FWHM on image	Defines whether or not a FWHM value should be displayed in the image
FWHM color	Defines the color in which the FWHM letters should be displayed
Font size	Defines the Font size in which the FWHM letters should be displayed
Number of digits	Defines the precision with which the FWHM should be displayed
Use minimum as zero for FWHM	Defines whether or not the Minimum of the profile should be used as the zero point for FWHM calculation
Display quick profile outside of image	Defines whether or not the Quick profile should be displayed outside of the image. If not it is displayed at the same area than the image.
Space for Quick Profile relative to image size [1%...100%]:	Defines how large the space for the quick profile should be relative to the image size. This option applies only to cases where the Quick Profile is displayed outside to image.
Display direction for general rectangle	This option defines whether and how an arbitrary rectangle should have a Quick Profile display or not. Possible values are: Don't display, Horizontal, Vertical
Adjust Quick profile height	Defines how to adjust the height of the displayed quick profile. If "Use LUT values" is selected the upper and lower LUT limits are used as the limiting values. If "Use Min/max values" is selected the minimum and

maximum values are use as the limiting values.

Analysis profiles



Show analysis cursors

When Show analysis cursors is selected, two cursors appear on the **Profile Display** Window and the **Profile Analysis** window appears on the screen.

Detailed analysis

If the **Detailed analysis** is selected under the **Options** menu several additional numerical values will be displayed.

Use minimum as zero

For some functions the "zero point" has to be determined. By default this point is exactly at intensity =0. However, if the option **Use Minimum as zero** under the **Options** menu is selected the minimum found in the region between the two cursors is regarded as the zero point. This is very useful if the profile contains an undesired offset.

Scaled display

When **Scaled Display** is active, the scaling data attached to the profile is used to get the correct physical value for every data point. If **Scaled display** is inactive corresponding pixel values are displayed. With this option you can even display profiles with different units. All analysis values are displayed in scaled values or in pixels according to this option.

Background black

Background black determines whether the profile background will be displayed black or white.

Auto zoom	If Auto zoom is selected, the system performs a zoom operation whenever a profile has been changed. This ensures that always the best display mode is selected to show all actual profiles within the display window with maximum size (like the Zoom all button).
Y-Axis	There are three different possibilities for display of the intensity axis. It can be either displayed linear or logarithmic with the base 10 or with the base e. This can be selected by choosing one of the options Linear , Log(base e) or Log(base 10) .

Saving and loading parameters

Permanent parameters

If the users closes a session and reopens it most of the parameters are remembered permanently and restored once the session is reopened. In other words most of the parameters used in the software are permanent.

There are a few parameters which may confuse the user if they are restored from one session to the next. These parameters are set to a default value when a new session is opened.

To demonstrate this we take the external trigger mode of a CCD camera. The external trigger mode switches the CCD camera to a mode where it only acquires images if a trigger is applied to the camera. We assume that a user which has operated the CCD camera in external trigger mode the day before reopens the session the next day and starts live mode. But due to the fact that the camera is switched to external trigger mode it does not acquire images. It may take long time for the user to remember that he applied triggers to the camera the day before and he may be confused.


Workfiles

In other cases the user may want to exactly reproduce the experimental conditions he had the day before. In this case he may use the Save/Load workfile feature of the software.

Saving a workfile saves all current parameters (even those which are not permanent in the previously explained sense). Loading a workfile restores all parameters stored previously in this file.

To save the current set of parameters to a workfile proceed as follows:


1) Call the Workfiles menu command.



- or -

Right click to the programs client area and select Workfiles...


The Workfile Dialog appears:



2) Select Save and enter a file name on the file dialog which appears

To Load a Workfile call the Workfile dialog and Select Load.

A warning will be issued saying that all settings will be overwritten by the new settings in the workfile.



Note:


If you are not sure do not load the new parameters

- or -

Save the current settings to another workfile first, and then load the new workfile.

Image positions

A special parameter in the general menu defines whether all window positions are restored once the window is opened again. This feature allows the user to tidy up his user interface and to place all dialogs at a location where he is used to search it.



If this option is not checked all windows are loaded to a default position.

Note: If it happens that dialogs are suddenly lost - because they are opened at locations outside of the main dialog -, switching off this feature, closing and reopening the software causes all dialogs to appear at default locations. Any “lost” dialogs can be found in this way.

Calibration & Analysis

These topics deals with the calibration of image and profile data and its analysis

Calibration

Calibration allows to assign physical units to image and profile data. First you can learn what calibration means in this program. Then you get an introduction to the different calibration methods (linear and table type) used in this software. In the topic "System, Image and Profile Calibration" you learn about the hierarchical ordering of calibration. The following chapters show how to make the set-up of calibration for the different methods.

Please see also the Appendix Calibration File Format for details.

What is Calibration?

When acquiring images with a two-dimensional camera, light intensities are measured at certain detector elements. These elements are called pixels. Of course there is some spatial relationship (correspondence) between the pixels and the real world. The camera may look into a room with an ordinary objective lens, or it may look at some microscope image. In all cases there is a specific spatial correspondence between the image on the camera and the real world. Calibration is a way to get quantitative information about the real world by extracting data from the camera image.

The geometrical transformations made during calibration process may be simple or more complicated depending on the type of transformation the physical measure undergoes in the complete system:

The simplest (and most often used) way of calibration is to attach a single calibration factor to the system. The pixel distance from one point to the next is then just multiplied by that factor.

The most general way of calibration is to allow a totally free mathematical co-ordinate transformation function like:


$$X_scaled = X_scaled (X_camera, Y_camera)$$

$$Y_scaled = Y_scaled (X_camera, Y_camera).$$

Such a very general calibration method, though it is the most versatile one, has three disadvantages:

- To perform such a calibration one would need two two-dimensional calibration arrays of the size of horizontal x vertical pixel number (x 4 byte at least).
- Lines with $X_scaled = const$ or $Y_scaled = const$ are no longer straight, making it very difficult to extract profiles along such lines (what the user typically wants to do).
- It is very complicated to handle this calibration scheme considering how the calibration data could be entered by the user.

Due to practical considerations a calibration scheme has been implemented which is easy to use but versatile enough for most cases. This calibration scheme has the following possibilities and restrictions:



Example image of calibration (non-linear in X and Y direction)

- Lines with $X_{scaled} = constant$ or $Y_{scaled} = constant$ always have to be parallel to the camera axes.
- The scales in X and Y can be linear or non-linear. In the latter case they have to be monotonous (increasing or decreasing).
- Non-linear calibration is realised by function tables (called Table Calibration) and any functional relationship can be modeled freely by the user, as long as the "monotony condition" is met.

Linear calibration in both dimensions enables the user to measure data along any arbitrary direction (also non-parallel to the axes).

Linear/Table

The system calibration contains two sets of calibration data: One for the X axis and one for the Y axis. Each calibration axis can either be of linear or of table type.

The **linear calibration** consists of a calibration factor and a unit. When the unit is the same for both directions and the calibration is linear in both axes, the calibration information can be attached to profiles generated in any arbitrary direction. When the dimension is not the same in both directions or when table calibration is used at least for one axis, only horizontal or vertical profiles can be scaled.

When linear calibration is used the origin for the calibration data of intensity profiles is always the starting point of the profile. The calibration of a profile therefore looks always like:

Pixel No.	Calibration value
0	0
1	1 * Factor

2 2 * Factor
n n * Factor

Thus, the absolute value will depend on the location of the starting point. If normal images are analyzed this is what the user typically wants. As a consequence, however, if you want to compare two profiles they should start at the same location. If you want to compensate an offset simply change the starting point.

The **table type calibration** consists of a table of **n** floating point values and a unit. **n** corresponds to the number of pixel in the axis. The **n** values reside in a calibration file with the extension **.scl** and are read automatically into memory when needed. The table may contain any values provided that they are strictly monotonous, either increasing or decreasing. (In other words, the function which the table resembles must be invertible.)

If you want to make non-linear calibration or attach absolute values to certain pixels or if you want to create calibration information which does not start at the value zero (e.g. for a spectrograph) you should use table calibration. Table calibration always yields the same scaled value at the same pixel location. Therefore absolute comparison is possible even if the profile windows do not start at the same location.

System / Image / Profile Calibration

The program distinguishes between the **system calibration**, the **image calibration**, and the calibration for every intensity **profile**.

System, image and all profiles individually can have different calibration information.

The calibration which is assigned to the system is always automatically applied to the current image at the moment the image is acquired. In case a profile is extracted from an image, the current image calibration is applied to the profile. The system passes its calibration to the image, and the image passes its calibration to the profiles.

For this transmission mechanism the following situations may appear:

- The system calibration attaches information to the pixels on the chip (camera). If the image origins on chip and on frame-buffer are not the same, the system calibration must be shifted accordingly.
- If binning is active during image acquisition, the system calibration is modified accordingly. This means, internally a table will be compressed accordingly, and a calibration factor will be multiplied by the binning factor. All these mechanisms are fully automatic and do not require special care from the user.
- When a profile is extracted from the frame-buffer the valid image area is regarded and automatically checked whether the profile sampling window is contained in the valid image area. If not, the user will be prompted to modify the profile


sampling window, extend the calibration data linearly or do not assign any calibration information to the profile.

The calibration is always automatically attached to the images and profiles and re-activated when images and profiles are reloaded from disk. In most cases the system calibration can also be reconstructed from the information in the image or profile.


Using Different Calibration Methods

There are several methods to tell the system which calibration it should apply for the X and Y axis.


To tell the system which calibration method he should use show the calibration options by calling the Options dialog from the menu



and select Calibration within the Tree dialog on the left side.



Calibration information can be different for the X and the Y Axis. The Calibration method deals with the way how to specify the calibration for both Axis.



There are three different methods:

- | | |
|----------------|---|
| No Calibration | No Calibration information is assigned to X and Y Axis |
| Square | A linear Calibration is assigned to X and Y Axis. Both calibration factors are identical |
| Free Scale | A different Calibration can be specified for the X and Y Axis. Both Axis can have linear or table calibration |

No Calibration


Select **No Calibration** to enable the square calibration method.

No Calibration is assigned to X and Y axis.

Note: Technically speaking this is not true. Instead a linear calibration with the factor 1 and the unit “No Calibration” is assigned to both X and Y Axis.

Square calibration


Select **Square** to enable the square calibration method.



You can specify a pixel distance a physical distance related to this pixel distance and the associated unit. If the currently selected image has a line ROI which is the currently selected one pushing the Get will calculate the pixel distance from the start and end point of this line ROI.

Free Calibration

Select **Free Scale** to enable the free calibration method.



The user can select the calibration type and the unit independently for both directions. For each direction five values can be specified:

Type, Pixel, Value, File and Unit.


When selecting the linear calibration, the meaning of **Pixel**, **Value** and **Unit** is the same as for square calibration. The only difference to square calibration is that you are specifying the calibration for each axis separately. Hence, for horizontal calibration you specify a horizontal pixel distance and for vertical calibration a vertical distance. For the ease of input you have to select a rectangular ROI where either its width or its height are used to calculate the pixel distance depending on the axis direction.

When selecting table calibration **Value** is no longer valid, the entry **File** is valid instead. By pushing the small button **Get** on the right side of the **File** entry you can choose a calibration file from a file list.

System Calibration

The lower part of the calibration dialog informs the user about the current system calibration (the calibration which was active before the calibration options dialog has been called).

After you have selected calibration values these are assigned to the system calibration. You additionally have the choice to assign these values to the currently selected image and to profiles.



Profiles

A profile is a one-dimensional data extracted from an image.

Quick profiles

The easiest way of extraction profiles is to use the quick profile feature where profiles extracted from horizontal or vertical Rios are displayed on the same window as the image.


Analysis profiles


Profile Analysis extracts one-dimensional intensity data from the image along user-defined positions, attaches scaling information and calculates and displays various information about these data.

In the following eight sections you learn how to use the profile analysis tools. First you get general information about the profile tools, followed by a detailed description of the profile control and profile display window. Then you learn how to acquire and display profiles, how to attach scaling to the profiles and how to make a more detailed analysis. Finally you get information about how to handle profile data and how to export them to other programs.

Please see also the appendix Profile File Format for details.


General Information on Profile Analysis

Choose **Profile** from the **Analysis** menu to display the Profile Control and the Profile Display dialogue boxes or toggle the state of the profile button on the toolbar .



Profile control dialogue box

The **Profile Control** dialogue box offers a variety of controls to acquire, display and handle profiles.



Profile display dialogue box

The **Profile Display** window is used to display the profiles. It can be freely resized by the user. A toolbar on top of the **Profile Display** dialogue box offers a variety of display and analysis options.

When the option **Analysis** on the **Profile Display** dialogue box is selected, the **Profile Analysis** dialogue box appears on the screen in its standard size.

Profile Analysis		
Cursor1:	7,730434 mm	
Cursor2:	69,57391 mm	
Difference:	61,84348 mm	
Memory:	M0	M1
Value at Cursor1:	1905,03	1952,65
Value at Cursor2:	1173,52	out

Standard type of the profile analysis dialogue box

If the option **Detailed analysis** is additionally selected the **Profile Analysis** dialogue box becomes larger and displays several additional analytical data.

Profile Analysis		
Cursor1:	7,730434 mm	
Cursor2:	69,57391 mm	
Difference:	61,84348 mm	
Memory:	M0	M1
Value at Cursor1:	1905,03	1952,65
Value at Cursor2:	1173,52	out
FWHM:	—	—
Lifetime (1/e):	—	—
Risetime(10%-90%):	—	—
Falltime(90%-10%):	—	—
Maximum:	1914,267	2017
Location of Max:	8,521739	8,180367
Integral:	96005,95	102604,8
Center of Gravity:	35,28615	35,03408

Large type of the profile analysis dialogue box

The Profile Control

The **Profile Control** dialogue box is designed in the sense that you first select the desired action (command button) and then the desired memory the action should work on. The last command remains active unless another command button is pressed. For example if you want to delete several profiles it is sufficient to click the **Del** push-button once and then the memory buttons whose profiles you want to delete.

Profile Control										
Get	Get H.	Get V.	Updt.		Load	Save				
Hide	Show	Del	Recall	Info	Zoom	Data				
Memory:	M0	M1	M2	M3	M4	M5	M6	M7	M8	M9
Color/Type:	Red	Green								
Auto update:	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Get profile with arbitrary direction										

Profile control dialogue box

The **Profile Control** window contains the following parts:

- A series of command buttons on the top of the window.
- 10 memory buttons
- 10 small fields below the memory buttons displaying the color and the type of the profile (these fields are invisible if the memory is empty)
- 10 check boxes to specify the **Auto update** function
- A text box displaying messages such as a brief description of the currently selected command or an error message

The command buttons activate the following commands:


<u>G</u>et	Acquires a profile with arbitrary starting and ending points (requires a line ROI)
<u>G</u>etH	Acquires a horizontal integrated profile (requires a rectangular ROI)
<u>G</u>etV	Acquires a vertical integrated profile (requires a rectangular ROI)
<u>U</u>pd	Updates the profile from the current image (at the same position and direction).
<u>L</u>oad	Loads a profile from disk into memory. If the profile is a part of a profile sequence the sequence will be loaded starting with the file you selected
<u>S</u>ave	Saves a profile from memory to disk. If the profile is acquired from an image sequence, profile data for all images of the sequence will be saved.
<u>H</u>ide	Hides a profile. It disappears from the Profile Display window.
<u>S</u>how	Shows a profile which previously has been hidden.
<u>D</u>el.	Deletes a profile from memory (<u>Caution</u> : this command is irreversible).
<u>R</u>ecall	Shows the profile sampling window on the current image and displays its co-ordinates in pixels and scaled units in the message field. The selected ROI becomes the current ROI. Its parameter can be changed now (e.g. using the ROI info tool).
<u>I</u>nf	Shows the Profile Info window, which contains the full information about the profile and the parent image
<u>Z</u>oom	Adjusts the scale of the Profile Display window so that the selected profile will best fit into it.
<u>D</u>ata	Displays the content of a profile memory in a text box. This can be used for Clipboard and DDE links.


Acquiring Profiles

There are three types of profiles: **Line** profiles, **Integrated Horizontal** profiles and **Integrated Vertical** profiles.

Before you acquire profiles, you may want to assign scales to the X-and Y-axis of the image.

To select any type of profile you first have to select an appropriate ROI. If more than one ROI is selected for an image all profile actions always refer to the currently selected ROI (this is the one which is displayed with a dashed line).

When you want to get a **Line** profile select a line-type ROI . Press **Get** and the memory button where you want to store the profile. The profile values are sampled along the line as described in the following drawing.






Sampling pixel data along a specified line

The profile is generated by sampling pixel data along the line with at a step width of one pixel. When square scaling (or free scaling with linear scaling in both directions with the same unit) is selected the measured distance will be always correct, independent of the orientation of the sampling line.

Note: As a side effect of this sampling technique the value of one image pixel may appear twice in two adjacent profile channels. The method of getting **Line** profiles is especially useful for imaging applications with image sensors having square or non-square pixels.

Integrated Horizontal profiles and **Integrated Vertical** profiles are sampled along the sensor axis either horizontally or vertically.


When you want to get an **integrated** profile select a rectangular ROI (Any of the three selections ,  or ). To get either of both profile types press the **GetH** or **GetV** command button and click to the memory button where to store the profile. Depending on the type of profile (horizontal or vertical) a profile is extracted from the image within the specified window. The pixel intensities are averaged (integrated and normalized by the sampling window width) along the direction perpendicular to the profile direction. Since the values are averaged they will be independent of the width of the profile sampling window and fractional pixel intensities may occur with integrated profiles.

When there is already a profile in the specified memory the previously selected profile sampling window is displayed as default. If you select the **Updt** command you automatically get the same profile type as it has been in the memory before with the previously selected profile sampling window, but the data are re-calculated from the current image (which may have changed).

When the **Auto update** checkbox is clicked for an already existing profile the profile is automatically updated when the image data have changed. This function is triggered for example by image acquisition, by dark subtraction, by loading an image from disk, and so on. This checkbox has to be clicked if the profile shall be saved during profile sequence acquisition.

Displaying Profiles

The **Profile Display** window is used for displaying the extracted profiles.










Profile display window


The full area of the **Profile Display** window is used to display the profiles. Several profiles can be displayed at a time, provided that the scaling unit is the same. (The scaling may be different, but the unit must be the same.) Hence it is possible to simultaneously display profiles which were derived from different images, even with different scaling.

*Note: The system is not able to compare profiles with the same physical dimensions but with different units like **mm** and **cm** or **m**.*

Above the profile display you will find a task bar with functions to save and load ROI sets, to print profiles and to set some preferences.


- | | | |
|---|----------------------|--|
|  | Print Profile | Prints the profile display window on your printer |
|  | Save ROI | Save the actual set of ROIs |
|  | Load ROI | Load a previously saved set of ROIs |
|  | Unzoom | Zooms out completely. (Then the intrinsic limits will determine the zoom ranges). |
|  | Zoom last | Restores the last zoom setting. This is useful if you have zoomed too deeply into a profile and you want to get back one step. |
| | | Note: By using this command repeatedly you can alternate between two defined zoom settings. |
|  | Zoom all | Optimizes the zoom factor to display all actual profiles in a maximum size |
|  | Options | Opens a dialogue box to set profile display options |

To Open the Profile options either call the options menu and select the Analysis Profile item on the tree control



or right click to any of the profile dialog windows.

The Analysis Profile dialog will appear.



Profile options

In the profile display options dialogue you can set following parameter:

When **Analysis** is selected, two cursors appear on the **Profile Display** Window and the **Profile Analysis** window appears on the screen.

The option **Detailed Analysis** and **Use Minimum as zero** are explained in the chapter "Profile Analysis".

When **Scaled Display** is active, the scaling data attached to the profile is used to get the correct physical value for every data point. If **Scaled display** is inactive corresponding pixel values are displayed. With this option you can even display

profiles with different units. All analysis values are displayed in scaled values or in pixels according to this option.

There are three different possibilities for display of the intensity axis. It can be either displayed linear or logarithmic with the base 10 or with the base e. This can be selected by choosing one of the options **Linear**, **Log(base e)** or **Log(base 10)**.

Background black determines whether the profile background will be displayed black or white.

If **Auto zoom** is selected, the system performs a zoom operation whenever a profile has been changed. This ensures that always the best display mode is selected to show all actual profiles within the display window with maximum size (like the Zoom all button).

By using the mouse you can interactively **zoom** into a portion of the profile as follows:

- Move the mouse cursor to one of the corners which define the area you want to zoom in.
- Press the **right** mouse down.
- Move the mouse to the opposite corner of the area you want to zoom in.
- Release the mouse button.

You can repeat this process as often as you like.


Every profile has its intrinsic **data limits**. These are defined by the first and last data point in the X direction and by the data type and the detailed circumstances of the exposure in the intensity direction. When the **Unzoom** function is executed these intrinsic limits are used to define the zoom ranges. If more than one profile is displayed the absolute maximum and minimum values of the intrinsic limits of all these profiles are determined. The intrinsic limits for different type of data are as follows:

If a memory contains a profile the "Mn" text on the push-button is underlined. In addition, the type of profile can be determined by the small symbol below the button.

- \ is displayed for line profiles,
- for integrated horizontal profiles and
- | for integrated vertical profiles.
- s profile sequence
- f profile file
- a averaged profile

If a memory contains a profile which is not hidden than the background color under the symbol is the same as the profile's color. If it is hidden, that color is white.

You can temporarily hide a profile by clicking the **Hide** button in the Profile control dialogue box. But the system may also hide a profile automatically. This happens if you display a scaled profile with a unit different to that of other profiles in the display. In this case, the other profiles will be hidden (until you **show** them again).

When selecting **Print Profile**  from the **Profile Display** window the systems sends a hard-copy of the current **Profile Display** window to the Windows printer currently installed.

Profile Scaling

In the moment a profile is generated the scaling information is transferred from the image to the profile. When an integrated horizontal or vertical profile is generated the scaling for the corresponding axis is transferred from the image to the profile. If a line profile (arbitrary direction) is generated scaling information is only transferred when the scaling for both axes are linear and have the same unit. If the scaling factor is not the same for both directions, the angle of the profile sampling line is taken into account and a new scaling factor which is corrected for the profile direction is calculated. If the starting point of the line is S_x, y , the ending point of the line is E_x, y and the scaling factors are SC_x and SC_y , then the new scaling factor SC_n is:

$$SC_n = \sqrt{\frac{((S_x - E_x) * SC_x)^2 + ((S_y - E_y) * SC_y)^2}{(S_x - E_x)^2 + (S_y - E_y)^2}}$$

The numerator is an expression for the scaled length of the line, while the denominator is an expression for its pixel length.

Scaling information cannot be attached to the profile when the profile sampling window is extended beyond the valid image area (because the image scaling is defined only inside the valid image area) and a warning message will appear.

The user can select among following possibilities:

- If the user selects **Clip profile to image area** the overlapping area between the valid image and the previously selected profile sampling window is used for the profile generation.
- If the user selects **Get new profile** he can define the profile sampling window again. The previously selected area will serve as the default.
- If he selects **Extend scaling data linearly** the scaling factor is calculated as if the profile sampling window had been within the valid image. This works only for linear scaling.
The user can select whether he wants to clip the profile sampling window to the valid image.

- If the user selects **Do not assign scaling data** the profile will obtain no scaling at all (will remain non-scaled, and its values will be pixels).

Profile Analysis

When **Analysis** from the **Option** menu of the Profile Display window is selected, the **Profile Analysis** window will appear on the screen in its standard size and two cursor will appear within the **Profile Display** window.

Profile Analysis		
Cursor1:	7,730434 mm	
Cursor2:	69,57391 mm	
Difference:	61,84348 mm	
Memory:	M0	M1
Value at Cursor1:	1905,03	1952,65
Value at Cursor2:	1173,52	out

Profile analysis window

You can move the cursors by clicking and dragging them with the left mouse button (like the cursors on the LUT tool). The profile analysis display shows the X-values (pixel numbers or scaled values) and the difference between the two cursors. This is useful for measuring distances. It also shows the intensity values of all profiles at the cursors locations. If the cursor is out of the range of a profile the strings "out>" or "<out" will appear instead of a number.

Note that the movement of the cursors is not bound to the profile channels. If you move them between two channels an interpolated value between the two neighbor channels will be used. The colors of the profiles are displayed in the text box showing the memory numbers. This Profile Analysis window can show the intensity values of all memories at one time (up to 10 memory buffers). If more memories are going to obtain profiles the width of the window will increase automatically.

When **Scaled Display** is selected under the **Options** menu all values are scaled values, otherwise all displayed values are pixel values. During movement of the cursors the values **Cursor1**, **Cursor2**, **Difference**, **Value at cursor1** and **Value at cursor2** are updated continuously.


If the **Detailed analysis** is selected under the **Options** menu several additional numerical values will be displayed.

Profile Analysis		
Cursor1:	7,730434 mm	
Cursor2:	69,57391 mm	
Difference:	61,84348 mm	
Memory:	M0	M1
Value at Cursor1:	1905,03	1952,65
Value at Cursor2:	1173,52	out
FWHM:	—	—
Lifetime (1/e):	—	—
Risetime(10%-90%):	—	—
Falltime(90%-10%):	—	—
Maximum:	1914,267	2017
Location of Max.:	8,521739	8,180367
Integral:	96005,95	102604,8
Center of Gravity:	35,28615	35,03408

Profile analysis window (extended type)

The profiles are analyzed in the region between the two cursors and the values for **FWHM** (full width at half maximum), **Lifetime**, **Risetime**, **Falltime**, **Maximum**, **Location of Max.**, **Integral** and **Center of Gravity**. These values are updated when the cursor is placed at a certain location and the mouse is released. They are not continuously updated during mouse movement.

All these analysis functions first search for the absolute maximum in the interval between the two cursors.



Data which are analyzed from an intensity profile

It does not matter whether cursor 1 is at the right side of cursor 2 or vice versa. The region for the analysis is always extended to the next measuring point. (This may be important to know when the selected area is very small.) After the absolute maximum has been found, the profile is scanned until its value becomes smaller or equal to a certain value. (For example, for the FWHM this value is just half the value of the maximum.) When the value lays between two data values a linear interpolation is made to find an estimation for the correct value.

For some functions the "zero point" has to be determined. By default this point is exactly at intensity =0. However, if the option **Use Minimum as zero** under the **Options** menu is selected the minimum found in the region between the two cursors is regarded as the zero point. This is very useful if the profile contains an undesired offset.

The values are calculated as follows:

Value	Search starts at	Direction	Search value	Search value if "Use Minimum as zero"
FWHM right	Maximum	to right	(Max)/2	(Max+Min)/2
FWHM left	Maximum	to left	(Max)/2	(Max+Min)/2
Lifetime top	=Maximum			
Lifetime right	Maximum	to right	Max/exp(1)	(Max-Min)/exp(1)+Min
Risetime	Maximum	to left	Max*0.1	(Max-Min)*0.1+Min

10%	Risetime	Maximum	to left	Max*0.9	(Max-Min)*0.9+Min
90%	Falltime	Maximum	to right	Max*0.9	(Max-Min)*0.9+Min
90%	Falltime	Maximum	to right	Max*0.1	(Max-Min)*0.1+Min
10%					

FWHM=FWHM right-FWHM left
 Lifetime=Lifetime right-Lifetime top
 Risetime=Risetime 90%-Risetime 10%
 Falltime=Falltime 10%-Falltime 90%


Note: The purpose of these simple analysis functions is to get fast information about the profile characteristics. The usefulness for precise analysis, however, is limited.

Integral calculates the integral (area under the profile) between the two cursors. (This can e.g. be used to calculate pulse energies).

Center of Gravity calculates the center of gravity of the profile between the two cursors. This can be used to calculate the position of a relatively broad pulse with greater precision than the Peak function.

This fast analysis cannot replace a thorough data analysis which is often necessary. Especially the fast analysis tends to be quite sensitive to noise. A thorough analysis would often require treatments like Fourier analysis, filtering, curve fitting, etc. which are out of scope of the functions built into the HiPic. We recommend to use specialized software for those higher requirements.

Displaying Profile Data



Profile display text box

When the **Data** command is selected on the **Profile Control** dialogue box the content of the profile data is displayed in a text box **Profile X Data** where X stands for the number of the profile (e.g. Profile 0 Data).

The data contains two entries for each data point (X and Y) separated by a TAB character. The data points are separated by CR+LF characters. X is the assigned scaling value and Y the intensity value.

Using the Clipboard and DDE Links

When the profile data is displayed in a text you can easily fill the clipboard or establish a DDE link to another application. To do so proceed as follows:


- Copy the content of the text box to the clipboard by selecting the **Edit Copy** (Ctrl C) command of the text box.
- Select the other application (like Microsoft Word for Windows, Microsoft Excel) and select either **Paste** (for pasting the clipboard contents) or the **Paste Link** function (some applications may use a different command name).
- Data are now transferred. In case of DDE, the link is now established. Whenever the data is changing, the other application gets the new data via DDE link.

This function can be used to make profile display or analysis within a separate program. Remember that the profile data are updated automatically after image acquisition if the auto update checkbox is checked within the profile control dialogue box.

Histogram

Histogram analysis is a statistical analysis of intensity data within an user-defined area of interest


Choose **Histogram** from the **Analysis** menu to display the Histogram window.



Select a rectangular Area of Interest) If no ROI is defined, the histogram will be calculated from the whole image area.

Click **Get Histogram** to get the intensity histogram.

The histogram window displays the calculated histogram. By clicking the right mouse button and moving the mouse into the histogram display a cursor will appear. You can see some statistical data for the selected intensity like how often this intensity value appears in the specified region (**Total Count**), average intensity value (**Average**) or standard deviation (**Stand. Dev.**). These data may sometimes be useful for signal characterization (e.g. noise characterization). Move the cursor to inspect the data of other intensity values.



Histogram window

A portion of the histogram can be displayed by zooming with the right mouse button. Clicking to zoom displays the whole histogram correctly. Zoom last restores the last zoom window.

Note: Due to technical reasons 14 and 16 bit image data can be analyzed by the Histogram function only with reduced accuracy.

3D Data


This function allows to display numerical intensity data from within a user-defined area of interest.

Extract 3D Data from an Image


The image data is displayed in a text box by this function. As the maximum number of characters in a text box is 32000 normally not all pixels values can be displayed. Therefore only every n^{th} Pixel value is displayed. n is determined by the size of the area to display or the area will be modified.

To extract 3D Data from an image proceed as follows:

1. Select a rectangle ROI within you image
2. Choose **3D Data** from the **Analysis** menu.




- If the Area is larger than the amount of data which can be displayed in the text box the user is prompted to select the area from where to extract the 3D data. According to the size of the area one has to select the factor n for every direction from the 3D data selection dialogue box. Alternatively the user can select whether he wants the software to automatically reduce the ROI while keeping the center point of this ROI.



3D data selection dialogue

- The minimum values for the horizontal and vertical reduction (means maximum amount of data) are predefined. If one wants to select less data points one has to increase the values for horizontal and vertical reduction. After you click OK the 3D Data will be extracted from the current image and its numerical values will be displayed in the 3D Data dialogue box.

The pixels values of one line are separated by TAB characters. Different lines are separated by CR+LF characters. The first line contains scaling information in the X direction. The first column contains scaling information in the Y direction.



3D Data dialogue box

Using the Clipboard or a 3D Data DDE Link

A DDE link can be established similar to the links which can be established with profile data. To do so please perform the following:

- Copy the content of the text box to the clipboard by selecting the **Edit Copy** (Ctrl C) command of the text box.
- Select the other application (like Microsoft Word for Windows, Microsoft Excel) and select either **Paste** (for pasting the clipboard contents) or the **Paste Link** function (some applications may use a different command name).
- Data are now transferred. In case of DDE, the link is now established. Whenever the data is changing, the other application gets the new data via DDE link.

This feature can be used to display and analyze 3D Data.

Note: An alternative method to get image data in Text-Format is to save images in ASCII format. With this format the full image can be saved but the amount of data will be rather high.

Special acquisition modes & special setups

This topic describes the operation of special configurations or special setups in the sense that these setups behave somewhat different from what is described in other topics of this document. The setups described here differ especially in the way image data is generated or displayed.

This is true for line sensor data and TDI mode camera where – even though we still can think of image data delivered by the camera – not individual images are delivered from a camera but a continuous stream of line data.

In this sense also a PIV mode camera behaves a little bit different from ordinary cameras because the TDI mode camera always delivers a pair of images with every trigger.

Finally the way images are created by Hamamatsu X-ray Flat panel sensors are a little bit different because these sensors require two triggers to trigger a single frame and also these sensors operate best if they are triggered continuously. Therefore a special acquisition mode is provided if synchronization with a pulsed X-Ray source is needed.

Line sensors

This topic describes how to acquire image data from a X-ray Line sensor and how to set camera and acquisition parameters.

Introduction

The line sensor differs somewhat in the way how image data is acquired. Instead of images of a well defined size the line sensor delivers continuously line data.

Scrolling display

To cope with this feature a different type of image display is used to visualize the image data. Image data is still displayed in a two dimensional window; however the data is updated as a new bundle of lines has arrived from the grabber to show always the most recently acquired lines. The size of this window is identical with the number of pixel in horizontal direction. The size of this window in vertical direction is – however – arbitrary and can be defined by the user as long as the memory limitations of the frame grabber board is not exceeded.

Corrections

Corrections as background and shading corrections are performed inside the line sensor. The way how to perform the image corrections thus differs a little bit from the way corrections are done with ordinary area scan cameras. The software built in features available from the toolbar and the menu should no be used. This may lead to undesired results. Instead correction should be performed by using the commands on the camera acquisition dialog as described in the related topic.


These functions allow executing functions inside the line sensor to acquire and store correction data.

As a consequence of the acquisition process also the defect pixel correction cannot be used.

Line sensor Options

In the X-ray line sensor options dialogue several parameters of the sensors can be set.

Click the File – options menu and select camera on the tree control of the left side.




The Line Sensor Options dialog will appear on screen.

The Line Sensor Options dialog contains several tabs for groups of camera and acquisition parameters.

These groups are:

- | | |
|-------------|--|
| Measurement | Defines how basic measurements are performed. |
| Correction | Defines how corrections are performed |
| Areas | Allows specifying areas along the line sensor where correction or acquisition is treated in a special way. |
| Info | Display information about the currently select line sensor |

Line sensor measurement parameters



Line Sensor Options dialog – Measurement tab

Several basic settings of the Sensor can be selected in measurement tab of the line sensor options dialog. The user normally wants to select these settings once and keep them for some time. Therefore these settings are stored and automatically used for your next session with. Press **OK** to memorize and activate the options parameter after you have made any changes in the options dialogue.

Trigger Mode

Three different trigger modes are offered. In **Internal** mode the lines are read out as defined by the speed parameters as explained in the section below. In **External** mode the start of the readout of a single line is triggered with an external trigger signal applied to the sensor. The exposure time is defined by the external trigger interval. In **External+constant exp.time** mode the start of the readout of a single line is also triggered by an external trigger signal applied to the sensor. The exposure time is defined by the Line Sensor Speed parameter as explained below.

Timing in Internal mode

The sensor timing in internal mode can be defined in two different ways:

- by using the magnifying geometry parameters or
- by using the Line sensor speed parameter

Choose the method **Magnifying Geometry Parameters** if the line sensor is not in proximity to the samples. In this case a magnification factor has to be considered.

To apply this method you have to click at the **Magnifying geometry** box. Then you can insert three parameters which define the geometrical magnification on your set-up. You have to insert the distance between X-ray source and object in meter (**Distance X-Ray tube <->object**), the distance between X-ray source and sensor in meter (**Distance X-Ray tube <-> sensor**) and the speed in which your samples are moved in meter per minute (**Speed of conveyor belt**). The most suitable sensor timing will be automatically calculated based on these data. The timing will be made to obtain a square size pixel format.

Choose the method **Line sensor speed** if the sensor is very close to the sample plane. In this case no geometrical magnification has to be considered.

This mode is active if the box **Magnifying geometry** is unchecked.

You have to define only a single parameter: the speed with which the samples are moved in meter per minute (**Line sensor speed**). The timing parameter to obtain a square size pixel format will be automatically calculated. Make sure that the inserted parameter is within the permitted range which is displayed below the Line sensor speed input box.

The calculated line frequency and exposure time will be displayed on the lower left side of the options dialogue.


The parameter **Line sensor speed** will be used also to calculate the exposure time when the sensor is used in the trigger mode **External+const exp. time**.

The parameter **No. of acquired lines per image** defines how many lines will be displayed in an image. A typical range is 500 to 1000 lines. The range may be limited by the frame grabber used in your configuration.

The **Frame Trigger** option defines whether an image is acquired starting from an electrical trigger. This can be used to start scanning whenever the object of interest comes close to the line sensor.


The **Scrolling Live Display** option defines whether the system acquires all lines of an image prior to display it (non scrolling mode) or whether the image is automatically updated as soon as a (small) portion of lines appear (scrolling mode).

To set no scrolling display set the following:



A screenshot of a software settings dialog box. It features a light gray background. On the left, the text 'Lines per image' is followed by a text input field containing the number '1024'. To the right of this field is a checkbox labeled 'Scrolling live display', which is currently unchecked.

To set the Scrolling display set:



A screenshot of a software settings dialog box. It features a light gray background. On the left, the text 'Lines per image' is followed by a text input field containing the number '2048'. Below it, the text 'Bundle height:' is followed by a text input field containing the number '256'. To the right of these fields is a checkbox labeled 'Scrolling live display', which is currently checked.

The value of “Bundle height” can be set according to the personal taste, however there are two recommendations:

If the bundle height is very small this can lead to data corruption inside the Grabber module.


For the eye it is convenient if the bundle height is set to a value which leads to frame update rate of a few Hz, which means the bundle height should be around the order of 1/10 of the image height.

Important Note: Depending on computers speed, line rate, bundle height and overall image size data corruption may occur in Scrolling live display mode. This cannot be avoided for all possible combination of parameters. In such cases please switch to non scrolling display mode.

The parameter **Displayed intensity data during measurement uses:** defines the type of data displayed in the acquisition dialogue. The data may help to optimize the setup by choosing an adequate intensity level. Three data types can be displayed.

- Average intensity within a line
- Minimum intensity within a line or
- Maximum intensity within a line

Line sensor correction parameters



Sensor Options dialog – Correction tab

In the Correction dialogue several set-up parameter for background and shading correction can be defined.

Please make sure that these settings are correctly made before you execute background or shading correction.

Before background or shading correction will be executed the software checks if reasonable conditions are set to perform these corrections. For background

correction the X-ray source has to be switched off, for sensitivity correction the X-ray source has to be switched on and no sample shall be in the image area.


The background correction will be executed only if the maximum intensity level within a line does not exceed the intensity value defined in the **Background data acquisition:Maximum level in any pixel** parameter. Typically this may be set to a value 200 to 300.

The sensitivity correction will not be executed if the minimum intensity level in a line reaches an intensity value higher than defined **under Sensitivity data acquisition:Minimum level in any pixel**.

The parameter **Number of times to check** defines how many lines will be checked with regards to above threshold parameters before background correction or shading correction will be started. Only if the condition is true for the defined number of lines correction will start.

The parameter **Number of integration for correction data acquisition** determines how many lines will be averaged in order to calculate the background and shading correction parameters. A typical value is 32.

Line Sensor area specification



Sensor Options dialog – Area tab


In the **Area** dialogue you can exclude certain parts of the image from processing.

The section **Correction confirmation range** allows excluding certain image areas from the background and shading correction. The section **Output mask** allows masking image areas in the image display. The determined areas will be masked displayed in black.

Set **No. of areas**: to define how many areas you want to define (maximum is 4).

Then define the start and end pixel of each range.

Line sensor information display



Line Sensor Options dialog – Info tab

In this dialog the camera version and some important features of the sensor are displayed.

Image Acquisition

Execute Live and Acquire


Choose **Live** from the **Acquisition** menu to display the acquisition dialogue.

In this dialogue you can control image acquisition and acquisition of correction data.

Press **Live** to display the actual live image. When you click **Freeze** the image acquisition stops and the last image acquired is displayed.

Press **Acquire** to acquire a single image.

The number of lines displayed in one image is defined in the options dialogue.



Live mode dialogue

Display a quick profile

A **QuickProfile** can be displayed during **Live** mode and after image acquisition has been terminated.

To show a quick profile proceed as follows:

- 1.) Select a rectangle ROI with full horizontal width on the toolbar:



- 2.) Define a ROI on the Live mode image with clicking and dragging with the left mouse button.


- 3.) Switch on the Quickprofile by clicking on the quickprofile toolbar button.



- 4.) Adjust the LUT to display the desired range of intensities. A simple way of adjusting the LUT in an optimum way is to click to the AutoLUT toolbar pushbutton.



This is how an uncorrected line sensor image with maximum LUT contrast enhancement normally looks like.



Execute Corrections

You can also start acquisition of correction data from this dialogue.

Normally you start by acquiring a background data set.


To do this you have to switch off your X-ray source and press **Get Background**. Now the system will check if the conditions as defined in the options dialogue are true. If so, background data will be acquired.

Second step is to set the A/D converter offset.

To do this you have to switch off your X-ray source and press **Get A/D offs**. Now the system will check if the conditions as defined in the options dialogue are true. If so, A/D offset data will be acquired.

You can also do the first and second step in one sequence by pressing **the Get A/D offset + background** button.

After getting background data successfully the image of the Line sensor – again with maximum LUT contrast enhancement - will look like this:



Third step is acquiring the shading data set.

To do this you have to switch on your X-ray source and make sure that there is no object in the image area. Then press **Get sensitivity**.

Now the system will check if the conditions as defined in the options dialogue are true. If so, shading data will be acquired.

Correction data will be automatically used during **Live** or **Acquire** mode.

You can save a set of correction data by pressing the **Save correction** button and following the normal file save procedure.

You can recall a correction data set which you have saved earlier by pressing the **Recall corr.** Button.

Press **Clear correction** in order to clear the actual correction data set.

After the correction data are applied to the raw image data the intensity values may be drastically changed. In order to correct this you can choose one of 5 **Correction methods** in order to decide which mean intensity level will be used for correction. This defines a correction factor used for image display.

-None will not apply any additional correction factor.

-Average will correct the intensity values to reach the mean intensity value of the uncorrected image.

-Max. will correct the intensity values to the maximum intensity value of the uncorrected image.

Min. will correct the intensity values to the minimum intensity value of the uncorrected image.

User will correct the intensity values to the intensity value defined in the **Value** box.

PIV mode

Introduction

PIV (particle image velocitometry) mode is a special acquisition mode to acquire a pair of image data acquired within short time starting from an external trigger.


Note: PIV is not a feature of standard cameras. Only special cameras offer this feature. If you cannot find the PIV feature in the camera options this may be an indication that your camera does not support this feature.

Normal image acquisition algorithms like LIVE or Acquire mode acquire only one image at a time so we have to acquire the data in a special way if we want to get both images.


Setting the camera to PIV mode

PIV mode can be selected in the camera options dialog. Once this mode is selected every trigger will cause the camera to output two images.

To Select PIV mode open the camera options dialog by clicking the File – options menu and select camera on the tree control of the left side.



Select PIV in the camera options dialog and close this dialog by clicking to OK.




To check the functionality you can use the LIVE mode. With every trigger you should see both images appearing shortly one after another.

Recording two PIV images in sequence mode


To record both images coming from the camera with one trigger we use the sequence mode.

To call the sequence dialog execute the Sequence menu.




The sequence dialog will appear.

Select Live for the acquisition mode, 2 for the number of loops (or a multiple of 2 if you want to record several pairs of PIV images) and Full speed for the Acquisition timing.




Then select Store to RAM as the data storage option.



Then click to Acquire on the acquisition TAB.

After execution of the first trigger the sequence dialog will show 2 acquired images




Please refer to the topic sequence acquisition for details of sequence mode.

X-ray Flat Panels


Options

In the Flat Panel options dialogue several parameters of the sensors can be set.

Click the File – options menu and select camera on the tree control of the left side.



The Flat Panel options dialog will appear



Sensor

options dialogue (left internal grab mode, right external mode)

Several basic settings of the Sensor can be selected here. The user normally wants to select these settings once and keep them for some time. Therefore these settings are stored and automatically used for your next session with HiPic.

Subtype

Select the sensor type you use. Normally the subtype is automatically detected according to the frame time in internal mode.

Detect Subtype at Startup

This option defines whether the subtype should be detected at startup or not. If the system cannot detect the subtype due to some reason it is an advantage to select this option and to manually select the subtype. These settings are remembered permanently and activated in future sessions.

New Types

Click the check box **New Types** if you use a sensor with „-02“ type number. New types have a little bit different number of active pixels and a Lemo connector for external trigger input. If this option is checked the choice “External lemo” appears in the timing mode selection box.

Timing mode

If **Internal timing** is used the exposure time is controlled by the sensor. The exposure time is then determined by the sensor model and binning mode. If **External Grb.** is used, exposure time and/or start of exposure is controlled by a trigger signal generated by the frame grabber. The user can set the exposure time in the acquisition dialogue

If **External Lemo** is selected, exposure time and start of exposure is controlled by an electrical signal applied to the Trigger signal connector at the Flat Panel (“Lemo connector”). For details please refer also to the sensor manual.

Binning


Several binning formats can be selected depending on the selected subtype. There are also sensors which do not support binning. Frame rate depends on the binning mode selection.

Wait for 2nd Frame in Acquire mode

This mode can be used only if an External timing mode was selected.

If the user clicks to Acquire the following happens.

- 1.) The HiPic starts image acquisition and waits until the first frame has been readout completely.
- 2.) Then it displays a text in red color indicating the user that he now can apply a pulsed X-Ray source. The X-Ray pulse can be applied as long as the red color is displayed.



This feature can be used to manually synchronize the acquisition with an external X-Ray source. The X-Ray Source must be synchronized with the flat panel and the


acquisition process (The software must be the master). Due to the design of the flat panel there is currently no way to synchronize the Flat panel with the X-ray source as a master.

Camera information

Here the camera version and some important features of the sensor are displayed.

Live Mode


Select Live from the Acquisition menu.



The Flat Panel camera acquisition dialog will appear and show the Live Tab.

In this dialogue you can select **Live**, **Freeze** and **Single Exp.**. Press **Live** in order to start image acquisition. When you click **Freeze** the image acquisition stops and the last image acquired is displayed.

Press **Single Exp.** to acquire a single image.



Flat Panel camera acquisition dialog.

Display a quick profile

A **QuickProfile** can be displayed during **Live** mode and after image acquisition has been terminated.

To show a quick profile proceed as follows:

- 5.) Select a rectangle ROI with full horizontal width on the toolbar:



- 6.) Define a ROI on the Live mode image with clicking and dragging with the left mouse button.

- 7.) Switch on the Quickprofile by clicking on the quickprofile toolbar button.



Adjust the LUT to display the desired range of intensities. A simple way of adjusting the LUT in an optimum way is to click to the AutoLUT toolbar pushbutton.



Exposure time

Exposure time can be set if trigger mode is set to External Grb. in the options dialogue. In Internal mode the exposure time is defined by the sensor, in External Lemo mode it is defined by the electrical signal.

Real-time background subtraction

Real-time background subtraction is an operation where a background image is continuously subtracted from the acquired live images. This function can be used to subtract a background like mottle image, intensity offset signals etc. from images.

To execute background subtraction you have to decide if you want to acquire a new background image or if you want to use any saved image as background: Open the dialogue **Background-Corrections** options and select or deselect the checkbox **Get real-time backsub data from camera**. The way how to get the background image depends on this setting:

- **Get real-time backsub data from camera** is checked:
Acquire a new background image in Live mode with your camera by clicking the button **Get BG data** in the **Live** menu while your camera is working in Live mode.
- **Get real-time backsub data from camera** is unchecked:
Use an image from file as background image (use this function if you have stored a background image on your hard disk before). After you press the **Get BG data** in the **Live** dialogue a file selection dialogue will be opened. Select the file which you want to use as background image for real-time background subtraction.

To start real-time background subtraction press **Do RT Backsub** while the camera is working in **Live** mode.

The subtracted image is now displayed. In many cases you may have to change the contrast settings with the LUT tool now in order display the resulting images in good contrast.

Press **Do RT Backsub** again to stop real-time subtraction.

Acquire mode

Choose **Acquire** from the **Acquisition** menu to display the Flat Panel camera acquisition dialog.


The Flat Panel camera acquisition dialog will appear and show the Live Tab.

The commands in this box are similar as explained under the chapter with the exception that only single images can be acquired by pressing the **Acquire** button. The next incoming image will be stored.

As the X-ray Flat Panel Sensors offer limited control of acquisition timing (also depending on the Timing mode chosen in the options dialogue), you may sometimes get unexpected image results, especially if you expect to synchronize the image acquisition to external events (e.g. a pulsed X-ray source). In order to offer a better timing control the option **Wait for 2nd image in Acquire mode** has been implemented.

If this has been activated the following will happen:

- 1) The HiPic starts image acquisition and waits until the first frame has been readout completely.
- 2) Then it displays a text in red color indicating the user that he now can apply a pulsed X-Ray source. The X-Ray pulse can be applied as long as the red color is displayed.



Acquire mode dialogue (left standard mode, right “Wait for 2nd Frame” mode)

Analog Integration

Analog Integration with X-ray Flat Panels adds up a number of images in the frame.

Choose **Analog Integration** from the **Acquisition** menu to display the Analog integration dialogue.

Its parameters are similar to those of the **Live** window. Additionally you can specify the parameter **# of Exposures** in the section **Integration**, which determines how many images will be accumulated.

If **Clear on start** is checked, the last image which has been acquired will be cleared. If it is not checked, new images will be added to the last image.

Press **Integrate** to start acquisition and integration of a number of images (as defined under **# of Exposures**).

TDI mode for ORCA

Introduction

TDI (time delayed integration) is a mode where an area sensor is operated similar than a line sensor. The object is moved along the camera and trigger pulses have to be applied to the camera to shift the charges in the same speed than the moving object.


Note: TDI mode is not a feature of the standard ORCA camera. Only special cameras offer this feature. If you cannot find the TDI mode feature in the camera options this may be an indication that your camera does not support this feature.

An Orca operated in TDI mode differs somewhat in the way how image data is acquired. Instead of images of a well defined size the line sensor delivers continuously line data.

To cope with this feature a different type of image display is used to visualize the image data. Image data is still displayed in a two dimensional window; however the data is updated as often as possible to show always the most recently acquired lines. The size of this window is identical with the number of pixel in horizontal direction. The size of this window in vertical direction is – however – arbitrary and can be defined by the user as long as the memory limitations of the frame grabber board is not exceeded.


Selecting TDI mode

When starting the HiPic the user has to select ORCA TDI as a camera.



Note: The Orca TDI camera is only supported by IC-PCI + Am-DIG and the PCDig frame grabber.

Select the TDI mode in the camera options.



Specify the number of lines in the vertical direction which is used to display the image and the binning factor in vertical direction.

Live and Acquire mode

If you start Live mode an image is displayed which is updated as fast as possible and shows always the most recent lines acquired by the frame grabber.

Acquire mode acquires the specified number of lines and outputs them in an image.

Display a quick profile

A **QuickProfile** can be displayed during **Live** mode and after image acquisition has been terminated.

To show a quick profile proceed as follows:

- 8.) Select a rectangle ROI with full horizontal width on the toolbar:



- 9.) Define a ROI on the Live mode image with clicking and dragging with the left mouse button.

- 10.) Switch on the Quickprofile by clicking on the quickprofile toolbar button.



- 11.) Adjust the LUT to display the desired range of intensities. A simple way of adjusting the LUT in an optimum way is to click to the AutoLUT toolbar pushbutton.



TDI Mode for C10000

TDI mode cameras from Hamamatsu are designed to scan images with high speed, high sensitivity. Special models are even sensitive in the UV. Even though their operation principle is very similar to a full frame CCD camera from the users point of view they behave very similar to a line scan camera.

The Hamamatsu TDI sensor C10000 has two readout modes: Area and TDI mode. In area mode they operate like a CCD camera and output a two dimensional image. In TDI mode they scan the image synchronized to a moving sample (e.g. a conveyer belt). Besides others the sensor has the following unique parameters:


Sensor Mode	Defines in which mode the sensor is operated. Can be area or TDI mode.
Sensor Bundle Height	Defines how many lines are transferred from the DCAM module to the application at once (This property is controlled

completely by software and is therefore not visible in any of the dialogs)


Readout Direction Defines which readout registers are used to readout the sensor. Can be Forward or Reverse
Exposure Time Defines the readout speed of the sensor

Area Mode

In Area Mode the sensor is readout as an ordinary area sensor. This is useful for focusing the optics. Set the Sensor Mode to Area and the exposure time to an appropriate value. For focusing purposes this is sufficient. The sensor has a specific pattern which is visible in area mode. If needed, background and shading correction can be applied to improve image quality. To setup area mode set the property Sensor mode to Area in the camera options.



Setup of the C10000 for area mode.



Operating the C10000 in area mode

TDI mode

TDI mode is the actual measurement mode of this sensor. Even though the TDI mode is a kind of line sensor mode, which means that line data is continuously output without vertical blank the outputted data is organized in larger quantities which define one image. The image height which should be used for display is arbitrary but should be large enough to be able to view items which should be on screen for the same time. The maximum height of such an image is set to 4096.


To setup TDI mode the following properties have to be set:

Set sensor mode to TDI


Set readout direction to forward or reverse

Set the number of lines which should be readout for one image

Set the exposure time to an appropriate value



Setup of the C10000 in TDI mode.




Operating the C10000 in TDi mode.

Normal and scrolling LIVE display.


The HiPic provides two ways of live display:

The normal Live display In this mode a new image is displayed when the specified number of lines per image is reached. This is identical to conventional area sensors

The scrolling Live display In this mode the displayed image contains the specified number of lines per image. However if a defined number of lines are reached the complete image is shifted to the top and the new arrived lines are displayed on the bottom part.



No scrolling display image 1 and 2



Scrolling display image 1 and image 2

To set no scrolling display set the following:



The property “sensor bundle height” is automatically set to “Lines per image”

To set the Scrolling display set:



The property “sensor bundle height” is automatically set to “Bundle height”

The value of “Bundle height” can be set according to the personal taste, however there are two recommendations:

If the bundle height is very small this can lead to data corruption inside the DCDam module.

For the eye it is convenient if the bundle height is set to a value which leads to frame update rate of a few Hz, which means the bundle height should be around the order of 1/10 of the image height.

Important Note: Depending on computers speed, line rate, bundle height and overall image size data corruption may occur in Scrolling live display mode. This cannot be avoided for all possible combination of parameters. In such cases please switch to non scrolling display mode.


Exposure time in TDI mode

The exposure time in TDI mode is to be set so that the scanning is synchronized to the object movement.

The Exposure time (Exp) can be calculated using the number of Lines used for scanning (N) and the line frequency (F).

$$\text{Exp} = N / F$$

Using a magnification not equal to 1 we get the following relations:



Ls: Distance Sensor – Objective lens

Lo: Distance Object – Objective lens

$$M = l_s / l_o \text{ (Magnification)}$$

The moving object is imaged on the Sensor thus the speed of the object (S_o) translates to the so called line speed S_s (a virtual movement) on the sensor.

$$S_s = S_o * M$$

The line speed translates to the line frequency (F) using the sensor cell size (C). In other words if the sensor is readout with this line frequency the object image gets square pixels

$$F = S_s / C$$

We then get the exposure time:

$$Exp = N / F = N * C / S_s = (N * C) / (S_o * M)$$

$$Exp = (N * C) / (S_o * M)$$

where

N: Number of lines used for TDI scanning [1]

C: Cell size [m]

S_o : Speed of the moving object [m/s]

M: Magnification on the object [1]

In the case of the currently used sensor $N = 128$ and $C = 12 * e-6$ m (=12 μ m).


Using Sequences

Both in scrolling and non scrolling mode Live streaming to HIS file (Hamamatsu Image sequence file) can be used. Also traditional sequence recording to RAM and hard disk is available. When loading such sequence into the HiPic a scrolling display for this kind of line data is available. Such sequences have the Tag $IsLineData=-1$ in the section [Acquisition]. If this tag is found the HiPic offers some addition controls to play and switch images in the sequence dialog.


With normal images the sequence dialog looks like this:



With such kind of line data the sequence dialog get some more controls:



If the checkbox “Move # of lines” is activated not only the saved images are displayed as they are saved, but positions between the individual images are displayed. If the Line Offset is nonzero starting from this line number the current selected image is displayed and a part of the next image is added on the bottom part of the image (very similar to the scrolling live display). The Backward, Play backward, Forward and Play Forward pushbuttons then act no longer on full images but “move” image just the specified number of lines.



Original images saved on disk and displayed image with “Line Offset” (Red marks are only used for illustrating purposes they are not displayed in the real image).

Important notes using image sizes

The displayed image size is always the same as acquired during sequence acquisition and cannot be enlarged during display. Therefore it is important to select the value “Lines per image” large enough to see all details which the user wants to see at the same time.

The “lines per image” value defines the number of buffers used internally for the scrolling display. To avoid image data loss it is important to make this value sufficiently large.


Forward/Reverse scanning

The HiPic automatically flips image data in the case reverse scanning is selected.

Generic camera


Selecting Generic Camera

By selecting this camera model one can input image data from a third party camera by a frame grabber.



Generic Camera Options Dialog

To adapt the frame grabber to the parameters of the camera the generic camera options dialog can be used.



The generic camera options dialog has the following options:


Image Width	Width of the image in pixels
Image Height	Height of the image in pixels
X-Offset	Horizontal offset in pixels
Y-Offset	Vertical offset in pixels
X-Max-Area	Number of pixels outputted by the camera in X-direction.
Y-Max-Area	Number of pixels outputted by the camera in Y-direction.
Bits per Pixel	Number of bits per Pixel
Output mode	Can either be line or area mode
Data Valid	Sets the data valid option (Camera Link only)
Tap Config.	Tap configuration. Describes the number of TAPs and how they are combined to the final image
Camera Name	Name of the camera. Can have an arbitrary value. Is only used for documentation purposes
Exposure time	Value which is used to calculate timeouts and defines the acquisition strategy. Most frame grabbers does not need this value.
Readout time	Value which is used to calculate timeouts and defines the acquisition strategy. Most frame grabbers does not need this value.
Mode control: Mode0-2	TTL value outputted at a freely programmable output. Possible values are low or high
Trigger	TTL output which can be Low, High or set to a defined frequency (Pattern).
Interval time	Time which defines the frequency of the outputted pattern ($f=1/T$).

There is the possibility to communicate with the camera by a serial port which has been opened for this camera by using the RS232 tool.

Note: If a camera requires grabber parameters which are not mentioned here, they could eventually be set by using a configuration file. Most frame grabbers support such configuration files. Please the configuration program of the frame grabber how to set parameters and its meaning.

Generic camera acquisition dialog

Camera image can be acquired with the generic camera acquisition dialog.




Sequence acquisition

In sequence acquisition mode a series of images and/or profiles can be acquired. Dependent on the settings the sequence is recorded to the computer RAM (fast method, but length of recording is limited by RAM size) or to hard disk. The sequence function can save and load sequences and replay the sequence. Also sequences of intensity profiles can be recorded.


Note: At any time only **one** sequence can be active (either image or profile sequence).

There are three different ways (RAM, HD one file per image, HD one single streaming file) of storing sequence images. Depending on the complexity of image acquisition not all three ways of storing are available.

Sequence operation is controlled from the sequence dialog. To call the sequence dialog execute the Sequence menu.



The sequence dialog will appear.




The sequence acquisition dialog has three tabs where different types of parameters can be specified

Acquisition	Defines which and how images acquire and in which timing this should be done.
Data Storage	Defines which types of data should be saved (images and/or profiles) and how and where these should be stored
Processing	Allows to process the collected sequence data (average)


Selecting acquisition parameters

To specify the acquisition parameters select the Acquisition tab.




Acquisition mode and camera parameters

The first thing to specify is in which acquisition mode the image data is acquired.




Images are acquired with the acquisition parameters specified in the camera acquisition dialog with the parameters specified in the corresponding tab:



Also the camera parameters specified in the camera options are valid for the selected acquisition.

Cycle definition

The next things to specify are the properties of the cycle. The user has to specify how many loops he wants to perform and whether this should take place with full speed or predefined fixed interval.



No. of loops specifies how many images are acquired and stored.

Choices of Acquisition timing are **Full speed** or **Fixed intervals**. If you select **Full speed** the camera is running freely and the sequence is recorded with this timing. Immediately after the acquisition and storage of an image the next acquisition will be executed. The time interval is just defined by the exposure time, the readout time, the response time of the mechanical shutter and the data storage time.

Note: Fast hard disk recording works with **Live** mode and **Full speed** setting only.

If you select **Fixed intervals** you can define the interval from one acquisition to the next (examples of allowed formats: 500ms, 2.5s, 1m). This interval should of course be larger than the time needed for acquisition, readout and data storage.



Other information on the acquisition tab

Depending on the data storage the program may display other information on the acquisition tab.

These may be information about remaining disk space if saving to HD

Info	
Disk space (bytes):	8468303872
Bytes/image	2098240
Free images:	4035

In the case fast hard disk recording is used the available buffers during acquisition is displayed and if there were lost frames due to limited hard disk performance.

Avail. buffers:	5
Lost Frames:	0

In the case of long interval time and if the options Warning on is selected in the sequence options the remaining time to the next acquisition is displayed.

Time left (s):	17
----------------	----

Start sequence acquisition


If all parameter are selected clicking to start starts the sequence acquisition process.




The precise behavior of this process is defined by several parameters in the sequence options.

Prompt before start

One option defines whether the sequence option should start immediately after the Start button is pressed or not.



If this option is selected the software shows a dialog and does not start until the User Clicks to OK



This feature can be used to precisely start the sequence.

Enable stop

Set **Enable stop** in order to enable the possibility to stop the image acquisition at a desired time.



On older computers this function may sometimes slow down the sequence acquisition a little bit. This option should be selected normally. Only in cases where cameras with very high frame rate are used it may be an advantage to deselect this option.

Wrap around

A sequence acquisition to RAM offers the possibility to use the specified number of images as a ring buffer. Once the last image has been reached it restarts from the first image overwriting older images by newer one. To activate this option the user has to select Enable wrap around in Acquire to Ram mode in the sequence options dialog

Enable wrap around in Acquire to RAM mode

Then a checkbox at the acquisition tab will appear.

Wrap around

If the user selects this option the sequence acquisition will restart automatically once it has reached the end of the buffer. In this measurement mode the sequence acquisition does not end until the user stops the acquisition with the Stop button. The most recent acquired images are kept after stopping the acquisition.

Note: This option is only available if the images are acquired to RAM.

Display during acquisition

Normally all sequence images are displayed during acquisition. Sometimes the frame rate is so high that the image acquisition is slowed down just because of the image display. In such case the image option “Always display image during acquisition” should be switched off.

Always display image during acquisition

Sometimes the image display speed can be influenced by the image zoom factor. Many graphics boards have the feature that image display with zoom 1 is faster than image display with other zooming factors.

Note: In Fast hard disk recording mode (Live streaming) you can also use the option “Wait till next display in LIVE mode [ms]” to avoid slowing down the acquisition only due to highly repeating image display.

Fast hard disk recording

Fast hard disk recording starts the LIVE mode and writes all images to the specified frame.


Important! For the optimum performance of the hard disk write operation the routine to check the hard disk performance must have been executed before.

The fast hard disk recording functions assumes that there is no acquisition in progress. If there is an acquisition in progress - even a LIVE mode acquisition – the fast hard disk recording issues an error message and does not start.

If the camera has a very high frame it may happen that image display speed slows down the overall acquisition speed. There is a way to completely switch off the


image display. However this is not what a user normally wants. There is another way to acquire images with high speed without displaying all images in full speed.

There is an option in the acquisition options which allows specifying an interval time. Even if new live images are acquired they are not displayed unless the specified time has elapsed.



If for example we take a camera which outputs 150 frames / sec, we may want to save all images to hard disk (and this is normally possible because the image size is small in such case) but we do not need to display 150 frames / sec. If we specify 50 ms as a waiting time, the image display is limited to 20 frames / sec which is still sufficient for observing the image but does not slow down hard disk recording.

If the grabber model which you use supports image buffering (PCDig, PCCamLink, National Instruments boards, DCAM drivers), the buffers can be used and the buffer size can be defined by the user.



The number of buffers must be large enough to buffer timing interval where the hard disk is busy. Depending on the hard disk this is typically in the order of $\frac{1}{2}$ or 1 second. The number of buffers has to be selected in a way that it can buffer images which arrive in this interval. If we acquire images from an ORCA in 2x2 binning mode we have a frame rate of approximately 16 frame / sec. If we assume a maximum busy time of the hard disk of $\frac{1}{2}$ sec we have to provide at least 8 buffers. This number should not be larger than necessary because it consumes RAM memory but it should also be large enough to prevent lost frames.

During acquisition a display of the currently free buffers is visible and an information whether and how many frames have been lost already.



Optimize acquisition speed in fast hard disk recording

If lost frames appear during fast hard disk recording there are the following countermeasures which could be eventually done:

- Execute Check perform. HD
- Switch off the image display during sequence or lower the display rate
- Increase the number of buffers
- Use a different (faster) hard disk
- Use display data instead of the full camera data
- Acquire only an ROI


- Use a camera mode which produces less data like binning
- Increase frame time by increasing the exposure time

If there are still lost frame it is recommended to acquire image data to RAM.

Data storage options


To specify the data storage parameters select the data storage tab.

On this tab the user can define whether the software should store image and/or profiles and how and where to these data are stored.



Storing images

In most cases the user wants to store images



To do so select “Acquire images”. As an additional feature the user can specify to store only part of the image. For this purpose he can specify a ROI from where image data is stored. To do so select a Rectangle ROI on the image click the pushbutton ROI and select the selection ROI only.

If you select **Store to RAM** the sequence will be recorded to the computer RAM. Use this in order to acquire a sequence with the highest speed. However you may also have to consider that the RAM is limited so the number of images you can record is limited. If you specify a loop number larger than the amount of physical memory, the software only acquires a smaller number of images and stops then.

For storing to hard disk there are two possibilities:

- HD (streaming, LIVE only)
- HD (individual files, all modes)

HD (streaming, LIVE only) is also called fast hard disk recording or just streaming is the fastest way to save image data to file. It is however limited to LIVE mode and full speed acquisition.

In other cases like Analog Integration mode or fixed interval you have to use HD (individual files, all modes). In such cases the Store to selection does not offer the possibility of HD (streaming, LIVE only). Please see also the appendices ITEX Sequence, TIFF Sequence, ASCII Sequence and HIS Image sequence file.

Differences between storage modes

There are several differences between these three modes. The following table shows the most important differences.

Topic	Store to RAM	HD (streaming, LIVE only)	HD (individual files, all modes)
Space	Limited to free physical RAM space	Limited only to HD space	Limited only to HD space
Permanent Storage	Has to be saved to individual files after acquisition	One single file	Many files linked only by naming
Acquisition mode	All acquisition modes	LIVE mode and full speed acquisition only	All acquisition modes
Speed	Fast, can acquire in full speed from every camera	Fast, but limited to hard disk performance. A sophisticated buffer management tries to avoid lost frames	Slow
ROI	Can acquire ROI	Can acquire ROI. Allows to save only display data (data reduced to 8 bit)	Can acquire ROI
File type	*.img, *.tif,*.dat,*.his	*.his	*.img, *.tif,*.dat

Corrections	All corrections available	Images have to be saved to individual files first, then all corrections available	All corrections available
File management	The Windows Explorer may take long time for file management operations if the number of files is larger than 10000.	No problems with Windows Explorer. File size may be eventually too large to move (Several GB can be acquired within short time)	The Windows Explorer may take long time for file management operations if the number of files is larger than 10000.
Security against hangup	In case of hangup data is lost.	In case of hangup data is lost.	Data acquired until hangup is available.

Additional features and limitations of fast hard disc recording

Please consider following limitations of the **Fast hard disc** recording mode: It applies to images acquired in Live mode only. It applies only to sequences acquired at full speed. If the grabber model which you use supports image buffering (PCDig, PCCamLink, National Instruments boards, DCAM drivers), the buffers can be used and the buffer size can be defined by the user. A precise time stamp can be expected only if PCCamlink or PCDig frame grabber boards are used. With other configurations the time stamp may not be accurate.


You have the additional possibility to save data as 8 bit only if you additionally check **Display data only**. The images will be saved in a contrast setting as displayed on the screen.

Please make sure that the LUT setting is reasonable before you acquire an image with this option selected. If you are using a camera with more than 8 bit per pixel output, the **Display data only** function will reduce the amount of data to be saved. Even a computer with limited performance may be fast enough to save such sequences.

The data saved in fast hard disc recording mode is written to one single file. This avoids delays due to file management but may be difficult to handle once the file is very large.

Check performance of hard disk recording

You can check the performance of your recording system for the fast hard disk recording mode.



If you click the **Check perform. Video** button, the time needed to display an image on your monitor will be measured and displayed.

If you click on **Check perform. HD** the performance (byte per second) of your hard disk will be checked with a file size as determined by the **HD space (MB)** parameter. Before you click on the **HD** button you have to insert a reasonable figure in the **HD space (MB)** text box. The measurement may take several minutes depending on your computer performance and file size selected.

The result of the analysis is shown in a table. You can use this information to check if your configuration is fast enough to record images in LIVE mode on your hard disk.


Important Note: The information gained in the **Check perform. HD** will be used to optimize data access to the hard disk. Executing this test is essential for the optimum performance of the hard disk recording function. It is also important that the **HD space (MB)** parameter is set sufficiently large. 100 MB is the minimum reasonable value, any value larger than this increases the precision and significance of this test.

Specifying file names

If a sequence should be saved to a file or a set of files additional controls for file specification will appear.


In the case of fast hard disk recording one has to specify a single file with the extension *.his. To browse the hard disk for a directory and file name click to the pushbutton Get.

Additionally the user can specify whether only 8 bit data (as displayed on screen) should be saved. To do so select the option Display data only.



In the case of storing images to individual files the images are stored to a set of images linked by the same base name added by a number indicating the individual sequence image.


In this case the user has to specify a single file with the extension *.img or *.tif or *.dat but with a number for the start of the sequence. To browse the hard disk for a directory and file name click to the pushbutton Get.



Note: Please make sure the number of digits is large enough to number all images.
Example: If you want to acquire 5000 images you need to specify at least four digits in the file name.

Storing profiles

In other cases the user wants to store profiles (or images and profiles).




If you select **Acquire and store defined profiles**, all intensity profiles with defined and displayed ROIs will be stored. You have to select the file name in the file text box. Profile sequences are always stored to disk.





Image Sequence display (normal images)



To control the display and replay of a sequence the controls on the lower part of the sequence dialogue are used. If you run the program in demo mode you still can display and replay existing sequences.


Unless a sequence is already opened, you have to open a sequence using the File – Open dialogue or acquire a new sequence.

Select **Images** in the Sequence control – Acquisition dialogue to display a sequence of images (only necessary if an image sequence and a profile sequence are currently loaded).




The controls in the center are used to specify the image within the sequence. The text box **Current** shows the index of the image currently shown. You can enter the desired image number, switch to the first (with ) or last image (with ) or switch to the next (with ) or previous image (with ). You can

automatically replay the sequence backward or forward with the controls  and . The parameter **Interval** defines the time for which each image will be displayed.

If the parameter **Wrap** is selected, the sequence display will be continued until the user stops it by pressing the  (=Stop) button.

If you want to extract a single image from the sequence you can use the push-button **Copy to new image** to create a single image with the same data as the current image. This is displayed in a separate window.

The currently displayed image within the sequence contains some useful information in the window title bar: The word "Sequence", the file name (or the image number if the sequence is not yet saved), the time and the zooming factor. If the image is the first image within a sequence the time parameter shows the time when the image was acquired (11:08:09). If it is another image it shows the offset to the first image in milliseconds (e.g. 147 ms). Therefore an exact timing can be reconstructed. If **Full speed** was selected during acquisition the time when the image has been read out completely is registered, if **Fixed Timing** was selected the start time of the acquisition command is registered. This is normally more precise.




Title bar of a sequence image


Image Sequence display (Linesensor and TDI data)

In the case of line sensors sequence data is saved similar to image data in portions of individual images with fixed width and height however the meaning of the data is more like a continuous stream of lines. As a consequence the user may want to see not only the (more or less arbitrary) packages of saved data but any position between these packages. The software provides a method to select not only data which starts at full image borders but also between these image borders. The image header contains information whether the data has the character of line data.

With such kind of line data the sequence dialog get some more controls:




If the checkbox “Move # of lines” is activated not only the saved images are displayed as they are saved, but positions between the individual images are displayed. If the Line Offset is nonzero starting from this line number the current selected image is displayed and a part of the next image is added on the bottom part of the image (very similar to the scrolling live display). The Backward, Play backward, Forward and Play Forward pushbuttons then act no longer on full images but “move” image just the specified number of lines.



Original images saved on disk and displayed image with “Line Offset” (Red marks are only used for illustrating purposes they are not displayed in the real image).

Sequence image file naming conventions

When you select the **Save as..** function within the **File** menu you will find several formats to save image sequences:



Save as.. file type list

1. ITEX sequence

All images of a sequence are saved to separate ITEX-files. These files are only linked together by the naming convention. There are two possible naming conventions:

*.001 format

Example1: Starting name: TEST.001

Sequence: TEST.001, TEST.002, TEST.003 ...

Example2: Starting name: TEST.031

Sequence: TEST.031, TEST.032, TEST.033 ...

Note: This format should not be used any more, since it may not be supported in future versions of this program. It is only provided due to compatibility reasons with previous versions.

*001.img format

Example: Starting name: TEST145.IMG

Sequence: TEST145.IMG, TEST146.IMG, TEST147.IMG etc.

The default extension is 001 and the user only has to select a file base name. In this case the files are named NAME.001, NAME.002 etc.

When saving a sequence one has to select the starting name. All other file names are derived starting from this file name:

Note: The numbering can have up to 7 digits

2. TIFF sequence format.

There are two ways to save image sequences in TIFF format:

2.1. **Data to TIFF** sequence

Images are saved in 8 or 16 bit TIFF format, depending on the depth of the image buffer. If the image buffer is 8 bit, TIFF will be 8 bit deep, otherwise it will be 16 bit deep.

"Data to TIFF" preserves the measurement data. The saved image data do not depend on the LUT settings, although the LUT settings are saved in the TIFF header and restored from it when re-loading such an image.

2.2. **Display to TIFF** sequence

Images are saved in 8 bit TIFF format. This is a color-palette format. If the LUT is grayscale when saving the TIFF file, the palette contains gray scales, otherwise it contains hues colors.

"Display to TIFF" does not preserve the measurement data. Its purpose is to store pictures which look like the image on the monitor screen. So, the contents of image data depend on the current LUT settings. Use "Display to TIFF" if you want to save images for presentation purposes, but never try to store your measurement data this way!

3. **ASCII** sequence

Images are saved in ASCII format. Note that ASCII format does not preserve any calibration information.

The files created can be read with the **Open** function (however ASCII files can not be read by this program). Depending on the file format you choose, you can either display single images or a whole sequence.

If you execute image processing operations on a sequence all sequence images are treated in parallel. This is true for **Background subtraction** (the background image is taken and is subtracted from all images), **Shading correction**, **Arithmetic operations** and **Map values by LUT**.


4. **HIS** sequence

All images of a sequence are saved in one single file. This allows using full performance of your computer system. It is the fastest way to record images on hard disk. You have to use this format if you use the **Fast hard disk** recording mode.

Profile sequences

Profile sequence acquisition

Select **Acquire and store defined profiles** in the Sequence – Data storage dialogue to enable recording of profile sequences.



Also a base file name or the name of the first profile within the sequence has to be defined. During sequence acquisition the system adds four digits to the base name: The first digit specifies the profile memory number and the last three or four digits the number of the series. For example a profile name test6001.prf indicates the profile memory number 6 and the first exposure


Before you can start to record a sequence of profiles, you have to get a sample images and define the profiles.

Profile sequences are recorded for all profiles for which ROIs are defined.

Then you have to define the timing parameters of a sequence as described in the chapter above. After pressing **Start Acquisition** a series of profile files will be stored.

Profile sequence display

Unless you have just recorded a sequence, you need to load one from disk before you can display it . Use the **Load** function of the profile dialogue in order to open a profile sequence. You have to select Profile sequence in the file type selection, otherwise only a single profile is loads. If there have been several profiles where a ROI was specified than several profiles are loaded at the same time. All these profiles together are called one sample (In the case of images as sample is always on image).




Select **Profiles** in the Sequence control – Acquisition dialogue to display a sequence of profiles.

The way how to use the profile sequence dialogue is analog to the image display.

Processing sequences

It is possible to do some post processing operations with image and profile sequences. Choose the **Processing** tab from the **Sequence control** to display the dialogue for sequence processing. It is possible to

- Average a number of images or profiles



Sequence processing dialogue

Averaging

If you press **Average**, all data (images or profiles) of the current sequence will be averaged. This can be used for creating noise reduced data.

It is possible to exclude some data of the sequence from averaging. This is e.g. useful, if one or several data of the sequence show unwanted contents.

Display the data you want to exclude from averaging and click **exclude this sample** in order to eliminate this data from averaging.

Finally, hit the **Proceed** button.

In case of image data, the result will be displayed in a new image window. In case of profiles, the result will be displayed in the **Profile Display** window.

The result image is always a 16 bit type. All images are averaged, and in the case that the original images are no 16 bit images they are multiplied by an appropriate factor to use the full 16 bit depth without exceeding it. In the case of photon counting images the images are only summed up without averaging. The following formulas are applied to normal and photon counting images:

$$R_N = \frac{\sum I_n * f}{n}, R_{PC} = \sum I_n$$

Where R_N is the result pixel value for normal images, R_{PC} is the result pixel value for photon counting images, I_n is the pixel value of the image n (at the shifted location), f is a factor indicated below and n is the number of samples.

Data type of the original images	f
16 bit	1
14 bit	2
12 bit	8
10 bit	32
8 bit	128

Photon Counting

Introduction

Photon Counting is a mode in which single photon events are added up in the frame memory. If the signal integration time is sufficiently long, a very high signal-to-noise ratio can be achieved. This mode also effectively suppresses certain crosstalk effects inherent to some image intensifiers. A dynamic photon counting mode allows temporal analysis of photon counting images.


Version 8.2 of this software has enhanced the features of photon counting adding real time background subtraction, Moiré reduction, intelligent border treatment, support of fast two tapped cameras (like the C9300-201 with 150 fps) and offline photon counting.

Required Hardware


The following drawing shows the typical hardware required for photon counting.

It shows the vacuum tube including Photocathode MCP and phosphor screen as well as the output optics and CCD camera.

Starting from a photon an electron is created at the photocathode. This electron is multiplied inside the MCP and finally hits the phosphor screen. A small light spot is emitted at the phosphor screen and imaged to a CCD camera by the output optics.



The CCD camera finally registers a small light spot for every photon (see the following screenshot).



Conditions for Photon Counting Mode

Sensitivity to single photons

Photon Counting requires that the system sensitivity is high enough to see single photon events. This is not the case if a CCD camera is used without any other light amplification device. Such a light amplification device could be an image intensifier whose gain is set above a certain threshold.

Low light level

Photon counting can only be performed properly if the light intensity is low enough.

In case that the photon spot areas of two photons overlap at the same image on the CCD they will be recognized as only one photon event. To avoid such a counting error the probability of a "two-photon event" should be low.

Basic principle of photon counting

Threshold

The first principle of photon counting is that only the portions of the image are regarded where the intensity exceeds a certain **threshold**. Only these portions contribute to the result. The setting of this threshold is depending on features of the camera and of the photon spots which are recorded by the camera. This feature results in high processing speed which allows us to use photon counting in real time.

Reduce the spot information to a single count

Another important feature of the used algorithm is that the whole information of a single photon spot is reduced to a single count.

Coordinate determination


The third principle is that for every photon spot we determine as precise as possible the coordinates of the photon. This allows us to perform a two dimensional photon counting.

Photon counting algorithm and implementation


Background subtraction and threshold calculation

In a first step a background subtraction is performed and every pixel below a certain threshold is set to zero. The background subtraction reduces the effect coming from uneven CCD camera background or hot pixels.

The raw data of photon counting processing looks like:




All processing steps are done with two dimensional data but for the ease of illustration we just use a one dimensional view of the data:




The green one is the raw photon data whereas the yellow one is the background data.

After background subtraction we get:



If we apply a threshold of 100 we get the following:



Coordinate calculation

Further processing steps uses only data which is non-zero, which increases processing speed drastically. For every photon we find its coordinates x and y . In our one dimensional example the location of the photon would be pixel 20.

Advanced processing

Advanced processing steps reduces moiré pattern which may appear in systems with discrete pattern of different size (in our case MCP and CCD) cells and handles bordering effect.

Data Storage

The data is stored in an image which means that the frame buffer is increased for every photon at the photon location or it is written to a so called DPC file which allows performing dynamic photon counting analysis.

Counting errors and maximum count rate


Correct threshold

An important factor for precise photon counting is the correct selection of the threshold. If the threshold is too low the system will count events which come from noise of the CCD system and which are no real photons. Be sure that the photon counting acquisition really counts no photons if there is no light.

As the number of electrons arriving at the output has a certain distribution (we call this pulse height distribution) the brightness of the photons spots vary strongly. On other words: There are photons which produce a higher signal then others. If the threshold is too high, not all of the photons are counted.

Overlapping photon spots

The photon counting algorithm can count photons as individual photons only if they are separated from each other. The following screenshot will illustrate this. The red and green lines represent intensity distributions originating from two photon whereas the yellow line can be judged as a single photon only (even though there may have been two photons)



The probability of double photons which are counted as a single one because they overlap is a measure of the counting error.

A good indication of this probability is the percentage of the image which is "covered" by photons, in other words the percentage of the image which exceeds the threshold. The smaller this percentage is the smaller is the probability that a "new" photon hits an already recorded one. To keep this value small the "single

photoelectron" probability should be made small by reducing the intensity of the light signal or minimizing the single photoelectron spot size. Increasing the frame rate reduces the number of photons detected on a single CCD frame thus decreasing the counting error.

Practical experiments have showed that the counting error have approximately the same value than the “above threshold” value.

This program displays the current percentage of pixels above the threshold in a specified ROI. This value should not exceed a few percent (around 5%).

Maximum count rate

The maximum count rate we can achieve with a give system depends on the size of the photon spots, the CCD camera frame rate and – of course – on the required precision.

Operating Photon Counting

Photon Counting Setup

Before you start photon counting image acquisition **background data has to be acquired** and a suitable **threshold has to be set**. A sophisticated automatic setup routine has been implemented to acquire background data and automatically searches for the most suitable threshold.


Open the Photon Counting dialogue and press **Setup** to start this routine. The program automatically acquires background data and performs histogram analysis. Then it starts photon counting acquisition in test mode to get an optimal threshold for the given system . Finally it proposes certain values necessary for the photon counting, as described in the chapter. Depending on the detector settings and type of detector, the procedure may take several minutes.

Note 1 During this setup the amplification device (MCP) should not be active otherwise wrong thresholds are calculated.

Note 2 If you did not perform Photon Counting Setup before you start photon counting acquisition, the system will prompt you to do so first. If you prefer to enter the threshold manually, you can ignore this message. This prompt is re-activated after you changed some parameters (like CCD exposure time) that can affect the proper threshold determination.

Note 3 In some rare cases the threshold has to be modified after photon counting setup to yield good results. In any case of doubt please contact your local Hamamatsu representative.


In case of a dual tapped camera like the C9300-201 two thresholds can be specified separately. The photon counting setup determines these two thresholds automatically.



Photon counting, Proposed threshold after automatic threshold detection

Photon counting acquisition


After you have executed Photon counting setup you can start to acquire photon counting data. Set your amplification device (MCP) to maximum gain, select the number of integrations and start photon counting.



The exposure time can be set freely but it is recommended to set it close or equal to the camera frame time. If the exposure time is set smaller then the frame time signal is wasted (this may be desired if the signal is strong), if the exposure time is set larger then the frame time the maximum possible count rate is decreased.

The photon counting process is continuously monitored and two types of messages are output.


The first is information about the number of pixels which exceed the threshold. As explained earlier this value should not exceed 5%.



The system also monitors whether the CCD camera becomes saturated, while looking at photon spots (Saturation always means loss of information). In such cases a warning message appears at the top part of the acquisition dialog and it is recommended to reduce the CCD camera gain (see also the option “CCD Gain for photon counting”). This message is only output in case of a problem.

Options related to photon counting

There are the following options related to photon counting:



The Selection box for **CCD Gain for photon counting** allows specifying the default CCD camera gain setting in photon counting mode. You can select between “Maximum gain” and “Minimum gain”. It is recommended to use “Maximum gain” unless the photon spots gets saturated.

If **Write dynamic photon counting file** is selected, the recording of photon counting images in the special DPC file format is enabled

If **32 bit images in Photon Counting** is set result images in Analog Integration mode will always be 32 bit.

The selection of **Moiré Reduction** specifies how strong Moiré reduction should take place. It is recommended to keep this parameters at its default value “Standard”

Dynamic Photon Counting (time resolved 2-D photon counting)

The standard photon counting mode just accumulates the photons detected during the acquisition period. If you want to record also the time when a photon has been detected, you can use the Dynamic Photon Counting acquisition mode. If this mode is enabled, a DPC-type file will be generated where the x-y co-ordinate of each photon is registered as well as the time when it has been detected. In a later analysis step this file can be used for temporal analysis of photon counting images.

Please see also the Appendix DPC File Format for details.

Acquisition

To activate the dynamic photon counting mode you have to select the item **Write dynamic photon counting file** in the menu Setup - Options- Display.

Then start photon counting image acquisition as described above.

Before the acquisition will start, a file name selection dialogue will be shown. There you have to define the file name of the file where the image data will be stored. For these files a special file extension is used: *.dpc. In this file the x-y coordinates of each photon and the time (CCD frame number) when it was detected are stored. The time resolution is limited by the frame rate of the camera you use for image acquisition.

In parallel to the data recording in the DPC file, an accumulated photon counting image is generated as described above.


Analysis

There are two ways to use dynamic photon counting image files of the DPC format:

If you want to see the accumulated photon counting image without using the time information you can load the file by selecting the “dpc file” format in the File – Open dialogue.

If you want to analyze the temporal information of a DPC file you have to open the file by selecting the “dpc sequence” format in the File – Open dialogue. Then an image sequence will be generated from the file.

After you opened the file a dialogue is displayed where you can set several parameters of the image sequence which will be generated:



Dynamic photon counting load dialogue

In the upper part of this dialogue several basic parameters of the DPC file are displayed, like file size, exposure time of one frame, number of exposures and total acquisition time.

In the lower part you can set parameters of the sequence which will be generated:

Start frame sets the frame of the photon counting image at which the sequence will start. **Frames /image** define the number of frames which will be accumulated for one image of the sequence and **No images** defines the number of images of the sequence.

Example (shown in above dialogue): If data were recorded with 15000 camera frames and you intend to split the data into 10 images, each image will be calculated from 1500 frames.

Set the **No images** parameter to define the number of sequence images you want to generate. If the number exceeds the value (No.exposures / Frames/image), the **No images** parameter will be automatically corrected. The header of the sequence images will show the time of the first image of the time frame.


Click **OK** to generate the sequence after you defined the sequence parameters. The generated sequence is a normal image sequence and can be treated in the standard fashion.

For information about handling and analysis of image sequences please refer to the chapter.


You can find a description of the DPC file format in Appendix F. You can use this if you want to use your own program for analysis of DPC files.

Offline photon counting

Offline photon counting can be performed on the basis of sequence data. To perform Offline photon counting right click with the mouse to a sequence and select:



Then please select the photon counting parameters:



Select the appropriate parameters and click to OK. A message box will appear telling you the progress of offline photon counting analysis and the value of the **Above Threshold** parameter.

Reserved Filename

Please note that the photon counting setup saves the background data to a special file inside the application directory. Do not delete or modify the following reserved file:

`Reserved_PhotonCountingBackground.img`

Otherwise photon counting may no longer work as expected.


Image processing

This topic describes how to apply image processing function.


Arithmetic

Arithmetic operations on an image or among two images can be made by using the Arithmetic commands.

Choose **A**rithmetic from the **P**rocessing menu to display the Arithmetic dialogue.



The Arithmetic dialog will appear.



Arithmetic dialogue

The arithmetic functions can be used for tasks like subtracting or multiplying two images. These functions will always work on the current image and will change that image.

Caution: Save your original image first if you need to keep it.

Both, unary and binary operations are possible. In case of binary operations the first operand is always the current image and the second operand is always an image in a file on a disk.

Note: You can select a partial area of the image and the function will be performed only inside that area.

In case of binary operations the system always automatically calculates the overlapping area between the two images and will perform the function only on that area. After the operation is finished the systems valid image area will be set to that area.

In addition, for all functions (also the unary ones) a constant may be defined as a calculation parameter. This constant has to be input into the edit box named **Const.**

For the binary operations the name of the file which contains the second operand image must be defined. Press the **Get** button to open a dialogue and select an image file.

You may also specify which type of arithmetic is used for the calculations, (long) integer or floating-point arithmetic. If the given constant is an integer value by default the arithmetic type is integer but you can force the program to use floating-point arithmetic, if you like. Choose the desired radio button in the frame **Type of arithmetic used.**

Note: This function can also be used to convert the bit range of an image without altering it otherwise. Just select the function $F=F+C$ with $C=0$, for example.

After you pushed the **Proceed** button the program calculation will start.

Note: The time required for a calculation will depend on the speed of your computer.

The following unary operations are available:

- | | |
|--------------------------------|--|
| $F=F+C$ | Adds a constant (integer only) |
| $F=-F+C$ | Inverts the frame-buffer contents (integer only) |
| $F=F*C$ | Multiplies the frame-buffer by a constant (integer or float) |
| $F=F/C$ | Divides the frame-buffer by a constant (integer or float) |
| $F=\ln(F)*C$ | Calculates the natural logarithm (float only) |

The following binary operations are available:

- | | |
|-----------------------------|--|
| $F=F+D+C$ | Adds a constant and a disk file to the frame-buffer (integer only) |
| $F=F-D+C$ | Adds a constant and subtracts a disk file from the frame-buffer (integer only) |
| $F=F*D/C$ | Multiplies the image in frame-buffer with a disk file (integer or float) |
| $F=F/D*C$ | Divides the image in frame-buffer by a disk file (integer or float) |

Based on these Functions one predefined operation is provided to invert the image content.

All integer operations are performed with long integers (4 byte) and are done in a way to minimize rounding errors (e.g. multiplication is done before division etc.).

After the calculation clipping is performed to keep the values within the range of the current data type. The lower and upper clipping values for the different data types are as follows:

Data type	Lower	Upper
8 bit	0	255
10 bit	0	1023
12 bit	0	4095
14	0	16383
16 bit	-32768	+32767
32 bit	-2.147483.648	2.147483.647

Predefined operations are operations on the image content which have well defined operators and constants. Currently there is only one predefined operation, it inverts the image data. If for example the image is a 12 bit image the inverted image is calculated by

$$I_{\text{new}} = (4095 - I_{\text{old}})$$

which leads zero for the highest intensity and 4095 for the lowest intensity. A similar effect can be achieved if the LUT option “Invert” is selected, with the difference that the LUT operation does not change the data whereas the processing functions changes the data.

Superimpose

Superimpose is a function which can overlay two images. This chapter is split into two parts. In the first part you can learn how Superimpose technically works. In the second part you learn how to operate the Superimpose function.

Technical Explanation of Superimpose

Two layers within the image memory are used to display two images simultaneously: The **bottom layer** and the **top layer** (see figure below).


The bottom layer is displayed in the color scale currently selected (either B/W or color). The display of the bottom layer can be modified with the LUT tool like at any other image.

The top layer overlays the bottom layer. The top layer is displayed in colors when the bottom layer is B/W and it is displayed in B/W when the bottom layer is colored. Thus, you can easily distinguish between the two layers.

The bottom layer is displayed with 8 bit grey level resolution and the top layer is displayed with 4 bit grey level resolution. The image containing both layers is 12 bit deep and can additionally contain an overlay.

When images are superimposed, the source images must be reduced in its dynamic range to 8 bit for the bottom layer and 4 bit for the top layer. This is done in a similar way as the output LUT displays images: The full dynamic range of data (e.g. 12 bit) is reduced to 8 bit which is also output on the monitor. The LUT

cursors saved with the source images are used to reduce the dynamic. This means if the LUT settings were appropriate to display the image correctly the images with reduced dynamic range also yield a correct image.




Scheme: How two source images are superimposed

Sometimes it happens that the top image hides wide areas of the bottom image. Therefore a possibility has been provided to make the display of intensity values representing low intensity values of the top image transparent. If more and more of the lower intensity values are transparent only the high intensity parts of the top image remain visible and you can more and more see all parts of the bottom image.

The Superimpose Control

Choose **Superimpose** from the **File** menu to display the Superimpose dialogue.



Superimpose dialogue

The **Superimpose Control** dialogue contains several controls to handle the superimpose function.


With the two push-buttons on the top of the Superimpose Control dialogue **FB → Top + File → Bottom** and **FB → Bottom + File → Top** a superimposed image can be created in one step. The current content of the frame buffer (current image) is transferred to one layer whereas the other layer is created using the data from a file.

With the push-button **FB → Top + File → Bottom** the current image is transferred to the top layer and the bottom layer is loaded from a specified file.

With the second push-button **FB → Bottom + File → Top** the current image is transferred to the bottom layer and the top layer is loaded from a specified file.

To specify the file which should be used as bottom or top layer image a file name must be entered into the text box below the corresponding push-buttons. The **Get** push-button on the right side can be used to select a file name out of a file list.

Use the LUT tool to change the brightness of the bottom image and **the 8/12 bit Top layer color** slider of the LUT parameter dialog to change the top image color settings.



LUT parameter dialogue

Example: A typical application of the superimpose function is present in the field of fluorescence microscopy. One wants to display an image of the illuminated sample showing the sample structure in detail together with an image of the fluorescence labeled parts.

The normal procedure for to get a superimposed image is as follows:

One first acquires the image with the illuminated sample and saves it to a file, using the **Save Image** command.

This file name has to be selected as file reference for the bottom layer by inserting the file name into the text box below the button **FB → Top + File → Bottom**.

Then acquire an image of the fluorescence light.

By pressing the push-button **FB → Top + File → Bottom** the superimposed image is created.

Average

Average is a function which is based on an image or profile sequence. Please see the chapter Averaging for details.

Creating data for special purposes / Tools

Tools are dialogs which are used to create special data which is used inside the program for different purposes.


Calibration file editor

The calibration file editor tool is used to create, view and edit calibration files (tables).

Calibration files are files which contain floating point entries (old type calibration files contain 1024 or 1280 entries) in either ascending or descending order. They must be strictly monotonous. They are used to assign calibration information to the system. Please see also the Appendix Calibration File Format for details.

If you want to add new calibration tables, be sure to have the necessary data (polynomial coefficients) prepared.

Choose **Calibration File Editor** from the **Tools** menu to display the Calibration File Editor.



Calibration File Editor

It displays the values of the floating point numbers as a graphical curve. The display can be either integral or differential. If the display is integral the values themselves are displayed. If the display is differential the differences of neighboring values are displayed. One can look at values by moving a cursor to a special location. Additionally the dialog shows the filename (if any) in the

caption, the pixel number of the actually displayed value, the value itself in numerical form, and an indication whether the display is integral or differential.

The dialog contains the following pushbuttons

Title	Command	Function
File	Load	Loads a calibration file into the memory and displays it.
	Save	Saves the current set of data to a calibration file on disk.
Create	Diff.	Opens a dialogue where a set of calibration data starting from a differential polynomial can be calculated.
	Integr.	Opens a dialogue where a set of calibration data starting from an integral polynomial can be calculated.
Arith	Add	Adds a value to the polynomial.
	Multiply	Multiplies the polynomial with a value.
Various	Print	Print calibration file data.
	Zoom	Adjusts the zoom in a way that the data is fully visible within the Calibration Files editor's window.

File Load loads a calibration file into memory and displays it. This function automatically checks whether the specified file is a valid calibration file or not.

File Save saves the current set of data to a calibration file on disk. It is only enabled when valid calibration data are defined.

Create Differential Polynomial allows calculating a set of calibration data starting from a polynomial specifying the differences of data values.

$$\text{Diff} = A_0 + A_1 \cdot i + A_2 \cdot i^2 + \dots$$


A₀: Coefficient 0

A₁: Coefficient 1 etc.

A_n: Coefficient n (if the order of the polynomial is n)

If you push the **Create Diff.** pushbutton the dialogue **Create Differential Polynomial** is opened

To create a polynomial select the order of polynomial first. Then insert the coefficients. Finally press **OK** to confirm the data input.



Create differential polynomial

Create Integral Polynomial allows calculating a set of calibration data starting from a polynomial specifying the data values.

$$Value = \sum_0^x A0 + A1 \times i + A2 \times i^2 + \dots$$


A0: Coefficient 0

A1: Coefficient 1 etc.

An: Coefficient n (if the order of the polynomial is n)

If you push the **Create Integr.** pushbutton the dialogue **Create Integral Polynomial** is opened

To create a polynomial select the order of polynomial first. Then insert the coefficients. Finally press **OK** to confirm the data input.



Create integral polynomial

The maximum polynomial order allowed is five for the integral and the differential polynomials.

Artih. Add Adds a value to the polynomial. If you click to the pushbutton Artih. A input box opens where the user can input a number. The number is added to the current values.

Artih. Multiply multiplies the polynomial entries with a value. If you click to the pushbutton Artih. Multiply an input box opens where the user can input a number. The polynomial entries are multiplied with this number.

A combo box labeled **Display** defines whether the calibration data is displayed in an integral (All items are displayed as they are) or a differential way (The difference of two adjacent entries is displayed).

By zooming with the right mouse button a part of the calibration data can be viewed

Defect pixel correction tool

The defect pixel correction function corrects defective pixel of a sensor by replacing them with not defective neighbors. The Defect pixel correction tool allows creating defect pixel data. See also the chapter Defect pixel correction tool for details.

Programming techniques

This topic deals with sophisticated programming techniques which can be used by expert users to extend the features of this software.

User function


User function is a feature that allows the user to implement any own functionality at specific locations within HiPic E.g. control of other devices or online analysis of data.

General

When the **UserFunctions** option is checked, function calls to a DLL named CUSTOMER.DLL are enabled. With this option the user can implement his own control and data analysis functions into HiPic

Choose **User Function** from the **Processing** menu to display the User function dialogue.

This function is kept for compatibility reasons with elder versions.



User function dialogue

User function is a feature of HiPic which allows the user to implement his own routines into the system. When the option **User Functions** from the **Options** menu is selected, the system calls the function `UserFunction` within a DLL with the name `CUSTOMER.DLL` from various locations. The definitions of the Function `UserFunction` is given below:

```
int FAR PASCAL UserFunction(int iIndex, char far *
sStatusString, int iStringLength, char far * sBackString,
int nBackStringMaxLen, int far * nBackStringLength, float
ftPar1, float ftPar2, float ftPar3, float ftPar4, int
iCycleIndex);
```

A dummy DLL is already delivered with the distribution disk and can be used as a prototype to built user-specific functions. The function **UserFunction** is called from various locations from within HiPic e.g. before and after an acquisition. The dummy DLL contains a list of all locations from where the Function is called. The parameter `iIndex` indicates from which location within HiPic the call came (e.g. Start of Acquire, Live etc.)

Function Parameters

The meaning of the parameters is as follows:

Parameter	Meaning
<code>iIndex</code>	Location from where the call came
<code>sStatusString</code>	Status string indication information about the image
<code>iStringLength</code>	Length of the status string in bytes
<code>sBackString</code>	String which can be handled back to HiPic Under certain circumstances this string is appended to the comment area of the status string (see below).
<code>nBackStringMaxLen</code>	Maximum number of bytes for the string of <code>sBackString</code>
<code>nBackStringLength</code>	Length of the string handled back (must not be bigger than <code>nBackStringMaxLen</code>)
<code>ftPar1</code>	Parameter which can be handled to the Function*
<code>ftPar2</code>	Parameter which can be handled to the Function*
<code>ftPar3</code>	Parameter which can be handled to the Function*
<code>ftPar4</code>	Parameter which can be handled to the Function*
<code>iCycleIndex</code>	Cycle number (when the Function called from within the cycle, Index numbers 25, 26, 27 only)

*: The user has to input these parameters in the user function dialogue

Memory and Display parameters

- handled to the `UserFunction` with the `sBackString` parameter -

2,-1,0,0,1000,1018,0,0,1000,1018,0.25,3288,1688,
0 1 2 3 4 5 6 7 8 9 10 11 12

4,140,297,485,309,148176908,2097152,150274064,1048576,
13 14 15 16 17 18 19 20 21

123896328,1318,96622516,10340,LIVE (0),0,3
22 23 24 25 26 27 28

0: iBytesPerPixel
1: fDisplayOnVGA
2: areSource.iX
3: areSource.iY
4: areSource.iDX
5: areSource.iDY
6: areImgToDisplay.iX
7: areImgToDisplay.iY
8: areImgToDisplay.iDX
9: areImgToDisplay.iDY
10: ftZoom
11: FormImageDisplay.hwnd
12: FormImageDisplay.PicImage.hwnd
13: iROIType
14: areROI.iX
15: areROI.iY
16: areROI.iDX
17: areROI.iDY
18: IDataHandle
19: IDataBytes
20: IDisplayHandle
21: IDisplayBytes
22: IStatusHandle
23: IStatusBytes
24: IScalingHandle
25: IScalingBytes
26: sFileName
27: fDataSaved
28: iPRFDirection

Notes:

areSource:	Coordinates of the stored image.
fDisplayOnVGA:	0=Image is not displayed, -1 is displayed on the VGA screen
areImgToDisplay:	Coordinates of the image which is display within the window FormImageDisplay.PicImage
ftZoom:	Zooming factor range: 0.25, 0.5, 1, 2, 4
FormImageDisplay.hwnd:	Handle of the Window containing the displayed image including caption, white areas and scroll bar.

FormImageDisplay.PicImage.hwnd:	Handle of the Window containing the image area (areImgToDisplay) only.
iROIType:	NONE =0, ZOOM = 1, ROIPOINT = 2, ROILINE = 3, ROIRECTANGLE = 4, ROIMODIFY = 5
areROI:	Area of the user specified ROI.
IDataHandle:	Memory is allocated with the GlobalAlloc function. To use the memory call lpvData=GlobalLock(IDataHandle) to unlock it call lpvData=GlobalUnlock(IDataHandle).
IDataBytes:	Number of Bytes allocated for the data.
LDisplayHandle, IDisplayBytes:	Same for Display data (8 bit per pixel)
LStatusHandle, IStatusBytes:	Memory area used to store the status string in the case another image is activated.
IScalingHandle, IScalingBytes:	1280 single (=float, 4 byte) values containing the scaling table, in the case table scaling is used
sFileName:	Name which appears in the image caption. In the case the image is saved it contains the full path of the file name.
fDataSaved:	0 if the image data is created or modified, -1 if the image is saved.
iPRFDirection:	Direction of the quick profile (HORINTEGRPROFILE = 2, VERINTEGRPROFILE = 3)

If the image does not exist (e.g. if the UserFunction is called prior to the acquisition) the sBackString is empty and iLenBackString=0.

If the string exceeds 255 characters it will be truncated to 255 characters.

User Function Dialog

The **User function** dialogue can be used to check the DLL function. It allows to enter the function index and the 4 parameters and returns the Cycle number, the index number and the returned string. The 4 free parameters which are entered at the corresponding locations are always passed to the function independent whether the function is called from the User function dialogue or from within HiPic. These parameters can be used to rule specific behaviour or the user function.

Prototype

The following is a printout of the prototype of the `UserFunction`

```
#include "windows.h"

short FAR PASCAL UserFunction(short iIndex, char * sStatusString, short iStringLen, char * sBackString, short
nBackStringMaxLen, short * nBackStringLen, float ftPar1, float ftPar2, float ftPar3, float ftPar4, short iCycleIndex);

BOOL WINAPI DllMain( HANDLE hModule, ULONG ul_reason_for_call, LPVOID lpReserved)
{
    switch(ul_reason_for_call)
    {
        case DLL_PROCESS_ATTACH:
            return(TRUE);
            break;
        case DLL_THREAD_ATTACH:
            return(TRUE);
        case DLL_THREAD_DETACH:
            return(TRUE);
        case DLL_PROCESS_DETACH:
            return(TRUE);
    }
    return 1;
}

short FAR PASCAL UserFunction(short iIndex, char * sStatusString, short iStringLen, char * sBackString, short
nBackStringMaxLen, short * nBackStringLen, float ftPar1, float ftPar2, float ftPar3, float ftPar4, short iCycleIndex)
{
    short iLen;

    iLen=0;

    switch (iIndex)
    {
        case 0: /* Called from the Function panel */
        {
            break;
        }
        case 1: /* Start of C4880 LIVE mode */
        {
            break;
        }
        case 2: /* Start of C4880 ACQUIRE mode */
        {
            break;
        }
        case 3: /* Start of C4880 ANALOG INTEGRATION mode */
        {
```

```

        break;
    }
case 4: /* Start of C4880 PHOTON COUNTING mode */
{
    break;
}
case 5: /* Start of C4880 BackSub LIVE mode */ Note:obsolete
{
    break;
}
case 6: /* Start of C4742 LIVE mode */
{
    break;
}
case 7: /* Start of C4742 acquisition (except Live mode) */
{
    break;
}
case 8: /* Start of analog camera LIVE and Acquire mode */
{
    break;
}
case 9: /* Start of analog camera ANALOG INTEGRATION mode */
{
    break;
}
case 10: /* Start of analog camera PHOTON COUNTING mode */
{
    break;
}
case 11: /* Start of analog camera Sync LIVE mode */
{
    break;
}
case 12: /* Start of analog camera Single Exposure Sync LIVE mode */
{
    break;
}
case 13: /* Background subtraction for C4880 */
{
    break;
}
case 14: /* Background subtraction for C4742 */
{
    break;
}
case 15: /* Background subtraction for analog cameras */
{

```

```

        break;
    }
    case 16: /* Start background subtraction from file */
    {
        break;
    }
    case 17: /* End background subtraction from file */
    {
        break;
    }
    case 18: /* End Acquisition for C4880 (all modes) */
        /* The return string is pasted to the comment of the image! */
    {
        break;
    }
    case 19: /* End Acquisition for C4742 and analog cameras (all modes) */
        /* The return string is pasted to the comment of the image! */
    {
        break;
    }
    case 20: /* Start shading correction */
    {
        break;
    }
    case 21: /* End shading correction */
    {
        break;
    }
    case 22: /* Start curvature correction */
    {
        break;
    }
    case 23: /* End curvature correction */
    {
        break;
    }
    case 24: /* Start cycle loop */
        /* iCycleIndex only valid with 24, 25 and 26 */
    {
        MessageBox (NULL, "Start of Cycle", "Cycle", 0);
        break;
    }
    case 25: /* End acquisition within cycle loop */
        /* The return string is pasted to the comment of the image! */
        /* iCycleIndex only valid with 24, 25 and 26 */
        strcpy (sBackString,"User Function index 25 executed");
        iLen=strlen(sBackString);
    {

```

```

        break;
    }
    case 26: /* End cycle loop */
        /* iCycleIndex only valid with 24, 25 and 26 */
        {
            MessageBox (NULL, "End of Cycle", "Cycle", 0);
            break;
        }
    case 27: /* Start of a sequence */
        {
            break;
        }
    case 28: /* A single acquisition during sequence mode has ended */
        /* If the return string is "stop" the sequence ends! */
        {
            /* This ends the sequence: */
            /* strcpy (sBackString,"stop"); */
            /* iLen=strlen(sBackString); */
            break;
        }
    case 29: /* After each new image in Live mode*/

        }
    default:
        {
            strcpy (sBackString,"No valid index number specified");
            iLen=strlen(sBackString);
            break;
        }
    }

    * nBackStringLen=iLen;
    return 1;
}
}

```

User function in Sequence mode

Depending on the submode used for sequence acquisition different cycles are executed (please refer also to “Details of sequence mode” in the chapter “Sequence”):

Mode S (fastest mode)

UserFunction 1 (Live mode) or 2 (Acquire mode)

```
UserFunction 27
Do
    Single acquisition
    UserFunction 28
Loop for all images
UserFunction 18 (C4880) or 19 (all other cameras)
```

Mode F (most flexible mode)

```
UserFunction 27
Do
    UserFunction 24
    GetImage (within this function the user function is
    additionally called as if the normal
    corresponding acquisition mode is called)
    UserFunction 25
Loop for all images
UserFunction 26
```

Remote control programming

To perform similar tasks and even more complex tasks the powerful **Remote Control SDK** is now available.

This Remote Control is based on Microsoft ActiveX Interface and allows to control all features from another program.

The Remote Control SDK is sold together with a training to users who are familiar with High Level language Programming like Microsoft Visual Basic or Visual C++.

Please contact your supplier for more detailed information.

RemoteEx

RemoteEx is a small program and allows to control the HiPic from another application via a text based communication which is exchanged by TCP-IP. This allows to control the software from another application on the same computer or an application on another computer. The other application can also be based on a different operating system like UNIX or Linux provided this operating system is able to communicate via TCP-IP. See the “Remote Ex programmers Handbook” in your application directory for details.


Script programming

The RemoteEx program also allows running script files to generate easily automated measurement cycles. Please see the “Remote Ex programmers Handbook” for details.

Accessing the CCD camera directly

This command is not available if the DCam driver configuration is used. The RS232 command can be used to watch or send control commands to a camera which is controlled by a serial command I/F. This feature is operable if the serial commands are transferred to the camera by a Windows COM port or by a dedicated com port of the frame grabber (like the PCI 1428 or PCCamLink frame grabbers) It is mainly intended for diagnostics purposes and should not be used under normal circumstances.

If you choose **RS232** from the **Info** menu the Communication dialogue appears.



It shows the commands sent to the camera by the program and the response strings received from the camera.




Also, the user can send commands to the camera manually (see the camera manual for details on the commands and their syntax). To send a command, write your command into the text box **Command to send** and push the **Send** button. The received response string is displayed in the display box **Received string**.

As long as the input focus is at the RS232 dialogue the user no longer can see commands sent and strings received in the meantime (Otherwise the user's commands would be permanently overwritten by the system's commands and could not be entered correctly). Once the user sets the focus to another window the RS232 dialogue continues to display all system messages.

RS232 communication dialogue

If, at any time during program operation, there should be a time-out problem with the serial communication to the camera, the program will display the following **RS232 Timeout** dialogue.



RS232 Timeout dialogue

You can either select **R**etry to retry to communicate via RS232, **C**ontinue to continue the program without sending the command, or **C**lose to close the communication link. When selecting **R**etry you can change to another COM port.

Appendix

Image file formats

The following describes the image file formats used in the program

ITEX image (*.img format).

This is one of the standard image formats (the other standard format is “Data to TIFF image”) which maintains the full information of all images and is compatible with all data processing functions of this program. If it is not intended to export image data to other programs, we recommend using this format.

Notes: If a camera with 16 bit data output is used, the image files will be saved in 16bit/pixel format. Images saved in the 32 bit data format can not be loaded in software versions lower than 6.1.

The .IMG file format used in HiPic is compatible to the ITEX format.

Bytes	Content
0-1	Characters IM
2-3	Comment length (byte)
4-5	Width of the image in pixels
6-7	Height of the image in lines
8-9	X-Offset
10-11	Y-Offset
12-13	File type: 0=8 Bit, 1=Compressed (Not used), 2=16 Bit, 3=32 Bit
14-64	Reserved
64- <i>nnn</i>	Comment area containing the status string
<i>nnn</i> +1- End	Data Area (one or two bytes per pixel stored in row order from the top to the bottom of the image)

Note: If Table Calibration is used, the Calibration tables are stored after the end of this section. Please see the Appendix Status String Format for details.

TIFF format

TIFF (*Tagged Image File Format*) is a widely used image format and supported by most image-processing and word-processing software packages.

The software also supports this file format. Two types of TIFF files can be generated.

TIFF images are stored using TIFF format version 6.0. When Calibration tables are used, these are appended to the normal data according to the status string which is saved in the comment tag.

Note: If Table Calibration is used, the Calibration tables are stored after the end of this section.

Display to TIFF image

If **Display to TIFF image is selected**, a (palletized)8 bit TIFF image will be saved.

The function of this command is to create a TIFF file which **looks** exactly like the image which you see on the display screen, including contrast enhancement by the LUT. It does, however, **not** contain the full data bit depth of the image, and hence it should not be used to store measurement results for other than display purposes.

ROI overlay and QuickProfile will not be saved.

The TIFF image looks exactly like the image which you see on the monitor.

Note: If you want to load TIFF images within this program later on, always use the **Data to TIFF image** format.

Note: Display to TIFF will reduce the bit depth of the images!

If you want to store image data in another file format than IMG or TIFF, save them as TIFF first and use a graphic program (PhotoShop, Paint Shop Pro, Corel PhotoPaint) to convert the format.

Note: If the current image is part of a sequence, you have the choice, if you want to save the full sequence or only the current image.

Data to TIFF image

Data to TIFF image allows to save image data with full data depth in TIFF format. The file will contain **8 or 16 bit** data without the LUT transformation.

Note: Though 16 TIFF files are true TIFF files according to their definitions, not many other programs support this relatively new data format yet. If you are not sure whether your image analysis program can handle these files try it and use standard TIFF, as generated with the **Display to TIFF image** in the case your program cannot handle them.

Images stored in 32 bit format will be saved in 16bit/pixel format. Most upper 16 bits are used. Lower bits are skipped.

Note: If the current image is part of a sequence, you have the choice to save the full sequence or only the current image.

ITEX Sequence

*This command allows saving a sequence of images which has been previously recorded with the **Sequence** function.*

You have to choose the name for the first image in the sequence. The system will **automatically** name all other images. Sub-areas of images (ROIs) can be selected.

Note: There are two naming conventions for naming a sequence: NAMEXXX.IMG and NAME.**XXX** where XXX is a numeric expression (e.g. 001).

Select the file format ITEX Sequence (*001.img) to save a sequence in the format NAMEXXX.img.

Select the file format ITEX Sequence (*.001) to save a sequence in the format NAME.XXX. **Note:** This format may not be supported in future versions of this program. It is only provided for compatibility reasons with older program versions!

TIFF Sequence

also supports TIFF file format for image sequences. Two types of TIFF file sequences can be generated:

Display to TIFF Sequence

If **Display to TIFF Sequence** is selected, 8 bit TIFF files will be saved.

The function of this command is to create a sequence of TIFF files which **look** exactly like the images which you see on the display screen, including contrast enhancement by the LUT. They do, however, **not** contain the full data bit depth of the original images, and hence this format should not be used to store measurement results for other than display purposes.

Note: If you want to load TIFF images within this program later on, always use the **Data to TIFF image** format.

The sequence will be recorded in the naming convention NAMEXXX.tif.

Note: Display to TIFF will reduce the bit depth of the images!

If you want to store image data in another file format than IMG or TIFF, save them as TIFF first and use a graphic program (PhotoShop, Paint Shop Pro, Corel PhotoPaint) to convert the format.

Data to TIFF Sequence

Data to TIFF Sequence allows to save image sequence data with full data depth in TIFF format. The files will contain 8 or 16 bit data without the LUT transformation.

Note: Though 16 bit TIFF files are true TIFF files according to their definitions, not many other programs support this relatively new data format yet. If you are not sure whether your image analysis program can handle these files try it and use standard TIFF, as generated with the **Display to TIFF Sequence** in the case your program cannot handle them.

Images stored in 32 bit format will be saved in 16bit/pixel format. Most upper 16 bits are used. Lower bits are skipped.

The sequence will be recorded in the naming convention NAMEXXX.tif.

ASCII image

Image data can be saved in ASCII format. The data are stored without any header information. Each pixel data is separated by a TAB character. The end of each line is indicated by carriage return + line feed.

The data can not be read by this program. ASCII files can be read by many spreadsheet programs like Excel and data analysis programs like Origin.

The extension *.dat is used for indicating ASCII image files.

ASCII Sequence

Image sequences can be saved in ASCII format using this format. The data are stored without any header information. Each pixel data is separated by a TAB character. The end of each line is indicated by carriage return + line feed.

The data can not be read by this program. ASCII files can be read by many spreadsheet programs like Excel and data analysis programs like Origin.

The naming convention NAMEXXX.dat is used for this file format.

HIS image sequence file

The HIS format is used for image sequences acquired with the fast hard disk recording mode. A full image sequence is stored in one single file. Please refer also to.

Following is a detailed description of the HIS file format:

A) General

There are two types of headers (with almost identical format)

First header of every channel

Header describing all following images within a sequence

The content is:

1. First header of every channel

- Image format (width etc.)
- Channel number
- Number of images in the file
- Number of additional cannels
- Time stamp
- Marker
- Info about channel x (Comment string, variable length)

This header has variable length size (Length is determined by comment)

2. Header describing all following images within a sequence

- Image format (width etc.)
- Channel number
- Time stamp
- Marker

This header has variable length size (Length is determined by comment)

The header is always immediately followed by the data of one image.

B) Universal case

1. First header of every channel

Bytes	Content
0-1	Character IM
2-3	Comment length in Bytes (ComLen)
4-5	Width of the image in pixels (iDX)
6-7	Height of the image in lines (iDY)
8-9	X-Offset (iX)
10-11	Y-Offset (iY)
12-13	File Type 1=8bit, 2=16bit, 3=32bit, 6=12bit ¹⁾ , 11=24bit RGB, 12=48bit RGB, 13=96bit RGB (not used now), 14=36bit RGB ¹⁾
14-17	Number of images in this file for this channel (-1=unknown)
18-19	Number of additional channels in the file
20-21	Channel number
22-29	Time stamp of image 0 (double)
30-33	Marker
35-64	Additional information (can be used freely by the application)
64 to 64+ComLen-1	Comment area can contain any information. There is no restriction in length (except that a two byte variable is used for specifying the length). The HiPic/HPDTA stores a string in the format of an INI file and scaling tables here (if scaling is table type).

64+ComLen to 64+ComLen+DatLen-1	Data of image
---------------------------------	---------------

2. Header describing all following images within a sequence

Bytes	Content
0-1	Character IM
2-3	Comment length in Bytes (ComLen)
4-5	Width of the image in pixels (iDX)
6-7	Height of the image in lines (iDY)
8-9	X-Offset (iX)
10-11	Y-Offset (iY)
12-13	File Type 1=8bit, 2=16bit, 3=32bit, 6=12bit ¹⁾ , 11=24bit RGB, 12=48bit RGB, 13=96bit RGB (not used now), 14=36bit RGB ¹⁾
14-17	=0
18-19	=0
20-21	Channel number
22-29	Time stamp of image 0 (double)
30-33	Marker

34-64	Additional information (can be used freely by the application)
64 to 64+ComLen-1	Comment area can contain any information. There is no restriction in length (except that a two byte variable is used for specifying the length). The HiPic/HPDTA stores a string in the format of an INI file and scaling tables here (if scaling is table type).

64+ComLen to 64+ComLen+DatLen-1	Data of image
---------------------------------	---------------

Remarks:

- ¹⁾ 12 bit grayscale images and 36 bit RGB images are stored bit-packed.

C) Standard case (one channel)

1. First header of every channel

Bytes	Content
0-1	Character IM
2-3	Comment length in Bytes (ComLen)
4-5	Width of the image in pixels (iDX)
6-7	Height of the image in lines (iDY)
8-9	X-Offset (iX)
10-11	Y-Offset (iY)
12-13	File Type 1=8bit, 2=16bit, 3=32bit, 11=24bit RGB, 12=48bit RGB, 13=96bit RGB (not used now)
14-17	Number of images in this file (-1=unknown)
18-19	=0
20-21	=0
22-29	Time stamp of image 0 (double)
30-33	Marker
34-64	Additional information (can be used freely by the application)
64 to 64+ComLen-1	Comment area can contain any information. There is no restriction in length (except that a two byte variable is used for specifying the length). The HiPic/HPDTA stores a string in the format of an INI file and scaling tables here (if scaling is table type).

64+ComLen to 64+ComLen+DatLen-1	Data of image
---------------------------------	---------------

2. Header describing all following images within a sequence

Bytes	Content
0-1	Character IM
2-3	Comment length in Bytes (ComLen)
4-5	Width of the image in pixels (iDX)
6-7	Height of the image in lines (iDY)
8-9	X-Offset (iX)
10-11	Y-Offset (iY)
12-13	File Type 1=8bit, 2=16bit, 3=32bit (not used now), 11=24bit RGB, 12=48bit RGB, 13=96bit RGB (not used now)
14-17	=0
18-19	=0
20-21	=0
22-29	Time stamp of image 0 (double)
30-33	Marker
34-64	Additional information (can be used freely by the application)

64 to 64+ComLen-1	Comment area can contain any information. There is no restriction in length (except that a two byte variable is used for specifying the length). The HiPic/HPDTA stores a string in the format of an INI file and scaling tables here (if scaling is table type).
----------------------	---

64+ComLen to 64+ComLen+DatLen-1	Data of image
------------------------------------	---------------

Profile File Format

The format of the profiles files written by **Fehler! AutoText-Eintrag nicht definiert**.is as follows:

```
;"HiPic 6.2 Profile"                                1.)
;"HiPic,5.0,1,4.0,3,6,3,3,373,3868,1,01-28-1994..." 2.)
;585,70,799,411                                    3.)
;342,0,3                                            4.)
;1,1.515152,"ps      ", "scal1  "                 5.)
0,357.4605
1.515152,360.8232
3.030303,354.1535
4.545455,352.8047
```

- 1) Identification line (contains ;"HiPic 6.2 Profile").
- 2) Status string in "", preceded by semicolon.
- 3) Start and end position of the profiles in the form: StartX, StartY, EndX, EndY. For Integrated profiles it indicates the position of two opposite corners of the rectangle used for integration, preceded by semicolon.
- 4) Number of data points, X-Offset (always 0) and Profile Type:
1=line, 2=Integrated horizontal, 3=Integrated vertical, preceded by semicolon.
- 5) Calibration Type (1=Linear, 2=Table), Calibration Factor (linear Calibration only), Unit in "", Calibration file without extension (table Calibration only), preceded by semicolon.
- 6) Data in subsequent lines, as many lines as number of data points. Format: X, Y value

Calibration File Format

Calibration files are used to provide nonlinear or special Calibration. They contain a list of floating point numbers (4 byte type, called float or single). Each number corresponds to one pixel on the chip. The Calibration file does not contain a unit nor the information for which direction the Calibration is applied. Calibration files must always have the extension **.SCL**. The floating point numbers must be strictly monotonous (ascending or descending), otherwise the file is not accepted by the system as a valid Calibration file. The format of the file is as follows:

Bytes 0-3	Value0
Bytes 4-7	Value1
Etc.	Etc,
Bytes n*4	Value
- Bytes n*4+3	n

Note: The Calibration data is directly written to the image file if the Calibration type is table. The tokens `ScalingXFile` and `ScalingYFile` contains an address where the Calibration table is written in the file. An asterik „*“ or a plus „+“ indicates the address. The asterik indicates that the Calibration has 1024 entries, the plus indicates 1280 entries. If it contains e.g. the entry `*473533` the Calibration data is written in the image file at an offset of 473533 bytes. The format `#xxxxxxxxyyyy` means the table with yyyy entries and address xxxxxxxx.

DPC File Format

Dynamic photon counting images are saved in the DPC file format. In a DPC file the x-y coordinate of each photon and the time when it has been detected are recorded.

Bytes	Content
0-1	Characters IM
2-3	Comment length in bytes (ComLen)
4-5	Width of the image in pixels (iDX)
6-7	Height of the image in lines (iDY)
8-9	X-Offset (iX)
10-11	Y-Offset (iY)
12-13	File type: 2=16 Bit
14-64	Reserved
64-nnn	Comment area can contain any information. It is used by the program to store the status string and the Calibration tables (if any)

nnn+1-End	Data Area (Starts at address 64+ ComLen)
-----------	--

The **Data Area** looks like n times the following diagram, where there is one set of such data for every frame. The first entry is the time stamp relative to the origin in ms, followed by a the coordinates of all photons counted within this frame. The end of the data for the first frame is indicated by a delimiter (0xFFFFFFFF). This is repeated for all frames recorded.

	Byte	Coordinates	Remark
03	00-	Timestamp	Time in ms from first image (long int)
05	04-	X0	Photon 0, X-Coord.
07	06-	Y0	Photon 0, Y-Coord.
09	08-	X1	Photon 1, X-Coord.
11	10-	Y1	Photon 1, Y-Coord.
	Etc.	Etc.	Etc.
	Etc.	Etc.	Etc.
24	20-	0xFFFFFFFF	Delimiter (Long int, all bits set to 1)

Status String Format

The status string is a string which is attached to an image and to profiles derived from an image. It contains all information about the image. The following is a sample string and the description of the different information. The status string contains only ASCII strings separated by comma.

The image status in this version is formatted differently from former versions.

The image status is saved as one string and it is organized like a *.INI file. It contains different sections where every section can contain tokens with assigned values. Other programs like remote control client programs or just any different image processing programs can add their own sections to save special data without disturbing the HiPic. It is, however, important to know that several of the entries in the status string are mandatory. Otherwise the file cannot read by the programs correctly. So be careful if modifying any of the tokens mentioned here. This is especially important for the Sections [Acquisition] and [DisplayLUT].

As an example we take the following status string (It was created when acquiring an image from the C7342-95 camera):

[Application],Date="06-03-2004",Time="16:03:46",Software="HiPic",Application=1,ApplicationTitle="High Performance Image Control System",SoftwareVersion="7.0.0",SoftwareDate="19.05.2004"

[Camera],AMD=N,NMD=S,EMD=E,SMD=N,ADS=12,SHT=150,FBL=1,EST=1,SHA=F,SFD=F,SPX=2,ATP=N,CEG=0,CEO=0,ESC=B,TimingMode="Internal timing",TriggerMode="Edge trigger",TriggerSource="BNC",TriggerPolarity="neg.",CCDArea="1280 x 1024",Binning="2 x 2",ScanMode="Normal",CameraName="C4742-95-12NRB",Type=7,SubType=7

[Acquisition],NrExposure=1,NrTrigger=0,ExposureTime=16ms,AcqMode=1,DataType=3,DataTypeOfSingleImage=3,ShadingCorr=0,CurveCorr=0,areSource="0,0,1280,1024",areGRBScan="0,0,1280,1024",pntOrigCh="0,0",pntOrigFB="0,0",pntBinning="1,1",BytesPerPixel=2,BacksubCorr=0

[Grabber],ConfigFile="C:\32BitApp\HiPic32_700PCDig.txt",Type=3,SubType=1

[DisplayLUT],EntrySize=3,LowerValue=0,UpperValue=4096,BitRange="12 bit",Color=3,LUTType=0,Gamma=1,First812OvlCol=1,Lut16xShift=0,Lut16xOvlVal=32767

[Scaling],ScalingXType=1,ScalingXScale=1,ScalingXUnit="No unit",ScalingXScalingFile="No scaling",ScalingYType=1,ScalingYScale=1,ScalingYUnit="No unit",ScalingYScalingFile="No scaling"

The name of every section is enclosed in brackets []. Every token is separated from it's value with a "=" character. The tokens are separated from each other with comma. After one section there can be a Carriage Return – Line Feed combination, but this is optional.

As an example we take the "SoftwareVersion" token with the "Application" section. It's value is "7.0.0". The value can be enclosed by quotes if necessary, but this is optional. Lets take the "EntrySize" token from the "DisplayLUT" section. It's value is "3" (But no quotes are used). When quotes are used the value can even contain commas like the "pntOrigCh" token within the "Acquisition" section. Section and token names are case sensitive which means "camera" and "Camera" are different.

Images created from the HiPic can contain the following sections with the following tokens

Application section

[Application]

Token	Type	Description	Example
Date	String	Image creation date	"06-03-2004"
Time	String	Image creation time	15:41:51
Software	String	Creating software	HiPic
Application	Integer	Creating software (see ApplicationType)	1
ApplicationTitle	String	Title of creating software	"High Performance Image Control System"
SoftwareVersion	String	Version of creating software	"7.0.0"
SoftwareDate	String	Version date of creating software	"19.05.2004"

Camera section

[Camera]

This section can have different values depending on the camera used for acquisition:

Analog Camera

In the case of an Analog Camera it has the following tokens:

Token	Type	Description	Example
CameraName	String	Name of the camera	Analog camera
Type	Integer	Camera type	3
SubType	Integer	Camera subtype	0
VideoInput	Integer	Video input	0
FrameTime	String	Video frame time	40ms
ExternalSync	Integer	External Sync (flag)	0
StartAcquisition Field	Integer	Starting acquisition field in sequence mode	ICP_NEXT

C4742-95

In the case of a C4742-95 it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (N=Normal, E=External)	N
NMD	String	Normal mode (N=Normal, S=Electrical Shutter, F=Frame Blanking)	F
EMD	String	External mode (E=Edge Trigger, L=Level Trigger)	L
SMD	String	Scan mode	S
ADS	String	AD select (possible values are 8 10 or 12)	12
SHT	String	Shutter (number of lines, possible values are depening on superpixel settings)	1
FBL	String	Frame blanking number (number of frames, possible values are depening on superpixel settings)	1
EST	String	External shutter (number of lines, possible values are depening on superpixel settings)	1
SHA	String	Scan horizontal area (F=Full, K=Kilo)	K
SFD	String	Scan front dummy (flag)	F
SPX	String	Superpixel (possile values are depending on camera version)	2
ATP	String	Acquire trigger polarity (P=Positiv, NB=Negativ)	N
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
ESC	String	External source (B=BNC, D=D-Sub, I=Digital I/F)	I
Timing Mode	String	Timing mode	Internal timing
Trigger Mode	String	Trigger mode	Level trigger
Trigger Source	String	Trigger source	Digital I/F
Trigger Polarity	String	Trigger polarity	neg.

CCDArea	String	CCD area	1024 x 1024
Binning	String	Binning mode	2 x 2
Scanmode	String	Scan Mode (Normal/External)	Normal
Camera Name	String	Camera name	C4742-95
Type	integer	Camera type	7
SubType	integer	Camera subtype	7

C4880

In the case of a C4880 it has the following tokens:

Token	Type	Description	Example
SSP	String	Scan speed (H=High, S=Slow)	H
SAG	String	Scan amplitude gain (L=Low, H=High, SH=Super High)	H
SMD	String	Scan mode (N=Normal, A=Subarray, S=Superpixel, B=Binning)	N
SHA	String	Scan horizontal area (F=Full, HC=Half-center, HL=Half-left, HR=Half-right, QC=Quarter-center, QL=Quarter-left, QR=Quarter-Right, EC=Eight-center)	F
SVO	String	Scan vertical offset	0
SVW	String	Scan vertical width	1018
SVB	String	Scan vertical binning	2
SHB	String	Scan horizontal offset	2
SPX	String	Superpixel	2
SOP	String	Scan optical black (flag)	I
AET	String	Acquire exposure time	20ms
ASH	String	Acquire shutter (A=Auto, C=Close, O=Open)	A
AMD	String	Acquire mode (I=Internal, E=External Event, T=External Time, S=External Stop, L=External Level)	I
ATP	String	Acquire trigger polarity (N=Negativ, P=Positiv)	N
ATN	String	Acquire number	1
ACN	String	Acquire cycle number	1
CSW	String	Cooler switch (O=On, F=Off)	F
PSW	String	Panel switch (O=On, F=Off)	O
SHC	String	Shading control (O=On, F=Off)	F
TST	String	Temperature set	-30
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
CEC	String	Contrast enhancement control	F

		(V=Volume, E=through CEG and CEO commands, F=Off)	
IGC	String	Image intensifier gain control	0
IIG	String	Image intensifier gain	0
TDY	String	Trigger delay	0
AGT	String	Acquire gate time	0
Temperature	String	Temperature	-27.7
CVG	String	Contrast volume gain 0	
CVO	String	Contrast volume offset 31	
Camera Name	String	Camera name	C4880
Type	Integer	Camera type	1
SubType	Integer	Camera subtype	1

C4880-8X

In the case of a C4880-8X it has the following tokens:

Token	Type	Description	Example
ACN	String	Acquire cycle number	1
AET	String	Acquire exposure time	0
AMD	String	Acquire mode (I=Internal, E=External)	I
ATP	String	Acquire trigger polarity (P=Positiv, N=Negativ)	N
CEC	String	Contrast enhancement control (V=Volume, E=through CEG and CEO commands, F=Off)	V
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
SAG	String	Scan amplitude gain	H
SAR	String	Scan area	0,0,659,493,1
SHT	String	Shutter (number of lines)	500
SMD	String	Scan mode (N=Normal, E=Extended)	N
SSP	String	Scan speed (H=High, S=Slow)	H
Chip	String	Chip	ICX074
VolGain	String	Volume gain	0
VolOffset	String	Volume offset	0
ROM	String	ROM version	2.1D
Version	String	Firmware version	2.5C
Gain	String	Gain 0	
Offset	String	Offset	0
Contrast Enhancement	String	Contrast enhancement	Potentiometer
Camera	String	Camera name	C4880-81

Name			
Type	integer	Camera type	8
SubType	integer	Camera subtype	3

C4742-98

In the case of a C4742-98 it has the following tokens:

Token	Type	Description	Example
SSP	String	Scan speed (H=High, S=Slow)	H
AET	String	Acquire exposure time	20ms
AMD	String	Acquire mode (N=Normal, E=External)	N
SMD	String	Scan mode (N=Normal, S=Super-Pixel, A=Subarray)	N
SPX	String	Superpixel (1, 2, 4, 8)	2
SHA	String	Scan horizontal area (F=Full, HC=Half-center, HL=Half-left, HR=Half-right, QC=Quarter-center, QL=Quarter-left, QR=Quarter-Right, EC-Eight-center)	F
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
EMD	String	External mode (E=Edge, L=Level)	L
SFD	String	Scan front dummy (O=On, F=Off)	F
ATP	String	Acquire trigger polarity (P=Positiv, N=Negativ)	N
SVO	String	Scan vertical offset	0
SVW	String	Scan vertical width	1024
SHO	String	Scan horizontal offset	0
SHW	String	Scan horizontal width	1024
Temperature	String	Temperature	50.0
Timing Mode	String	Timing mode	Internal timing
Trigger Mode	String	Trigger mode	Level trigger
Trigger Polarity	String	Trigger polarity	neg.
Camera Name	String	Camera name	C4742-98
Type	integer	Camera type	10
SubType	integer	Camera subtype	11

C7300

In the case of a C7300 it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (N=Normal,	N

		E=External)	
NMD	String	Normal mode (N=Normal, S=Shutter, F=Frame blanking)	F
EMD	String	External mode (E=Edge, L=Level)	L
SMD	String	Scan mode (N=Normal, S=Superpixel, A=Subarray)	N
ADS	String	AD select (8, 10 or 12)	12
SHT	String	Shutter (number of lines)	1
FBL	String	Frame blanking mode (number of frames)	1
EST	String	External shutter (number of lines)	1
SHA	String	Scan horizontal area (F=Full, K=Kilo)	K
SFD	String	Scan front dummy (O=On, F=Off)	F
ATP	String	Acquire trigger polarity (P=Positiv, N=Negativ)	N
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
ESC	String	External source (B=BNC, D=D-Sub, I=Digital I/F)	I
SVW	String	Scan vertical width	1024
SVO	String	Scan vertical offset	0
Timing Mode	String	Timing mode	Internal timing
Trigger Mode	String	Trigger mode	Level trigger
Trigger Source	String	Trigger source	Digital I/F
Trigger Polarity	String	Trigger polarity	neg.
Camera Name	String	Camera name	C7300-10-12NRP
Type	Integer	Camera type	19
SubType	Integer	Camera subtype	12

C8000-10

In the case of a C8000-10 it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (N=Normal, E=External)	N
EMD	String	External mode (E=Edge Trigger, L=Level Trigger)	E
SMD	String	Scan mode	N
SHA	String	Scan horizontal area (F=Full, K=Kilo)	F
SFD	String	Scan front dummy (flag)	F
SPX	String	Superpixel (1,2 ,4 ,8)	2
ATP	String	Acquire trigger polarity	N

		(P=Positiv, N=Negativ)	
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
SSP	String	Scan speed (H=High, M=Medium, S=Slow)	H
SHO	String	Scan horizontal offset	0
SHW	String	Scan horizontal width	640
SVO	String	Scan vertical offset	0
SVW	String	Scan vertical width	480
Timing Mode	String	Timing mode	"Internal timing"
Trigger Mode	String	Trigger mode	"Edge trigger"
Trigger Polarity	String	Trigger polarity	"neg."
CCDArea	String	CCD Area	"High"
Binning	String	Binning	"2 x 2"
ScanMode	String	Scan Mode	"Normal"
SubarrayHOffs	Integer	Subarray horizontal offset	0
SubarrayHWidth	Integer	Subarray horizontal width	640
SubarrayVOffs	Integer	Subarray vertical offset	0
SubarrayVWidth	Integer	Subarray horizontal width	480
Camera Name	String	Camera name	„C8000-10“
Type	integer	Camera type	21
SubType	integer	Camera subtype	23

C8000-20

In the case of a C8000-20 it has the following tokens:

Token	Type	Description	Example
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
REC	String	Recursive filter	F
RCN	Integer	Recursive filter number	2
BGO	Integer	Background offset	0
HVC	String	High voltage control	F
HVN	Integer	High voltage number	0
Camera Name	String	Camera name	„C8000-10“
Type	integer	Camera type	20
SubType	integer	Camera subtype	19

e			
---	--	--	--

X-Ray Flat panel

In the case of a Flat panel module it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (I=Internal, E=External)	I
SPX	String	Superpixel (1, 2, 4)	2
Timing Mode	String	Timing mode	"Internal timing"
Binning	String	Binning	"2 x 2"
Camera Name	String	Camera name	„Flat panel module“
Type	Integer	Camera type	22
SubType	Integer	Camera subtype	24

C8800

In the case of a C8800 it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (N=Normal, E=External)	N
EMD	String	External mode (E=Edge Trigger, L=Level Trigger)	L
SMD	String	Scan mode	S
SPX	String	Superpixel (1, 2, 4, 8)	2
ATP	String	Acquire trigger polarity (P=Positiv, N=Negativ)	N
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
ESC	String	External source (M=Multiport, I=Digital I/F)	I
SHO	String	Scan horizontal offset	0
SHW	String	Scan horizontal width	1024
SVO	String	Scan vertical offset	0
SVW	String	Scan vertical width	1024
Timing Mode	String	Timing mode	"Internal timing"
Trigger Mode	String	Trigger mode	"Edge trigger"
Trigger Source	String	Trigger Source	"Multiport"
Trigger Polarity	String	Trigger polarity	"neg."
Binning	String	Binning	"2 x 2"
ScanMode	String	Scan Mode	"Normal"
SubarrayHOffs	Integer	Subarray horizontal offset	0
Subarray	Integer	Subarray horizontal width	640

yHWidth			
SubarrayVOffs	Integer	Subarray vertical offset	0
SubarrayVWidth	Integer	Subarray horizontal width	480
Camera Name	String	Camera name	„C8800“
Type	integer	Camera type	27
SubType	integer	Camera subtype	38

C9100

In the case of a C9100 it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (N=Normal, E=External)	N
EMD	String	External mode (E=Edge, L=Level)	L
SMD	String	Scan mode (N=Normal, S=Superpixel, A=Subarray)	N
SPX	String	Superpixel	2
ATP	String	Acquire trigger polarity (N=Negativ, P=Positiv)	N
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
EMG	Integer	Electron multiplying gain	20
ESC	String	External source (M-Multiport, I=Digital I/F)	I
SHO	String	Scan horizontal offset	0
SHW	String	Scan horizontal width	1024
SVO	String	Scan vertical offset	0
SVW	String	Scan vertical width	1024
TimingMode	String	Timing mode	"Internal timing"
TriggerMode	String	Trigger mode	"Edge trigger"
TriggerSource	String	Trigger Source	"Multiport"
TriggerPolarity	String	Trigger polarity	"neg."
Binning	String	Binning	"2 x 2"
ScanMode	String	Scan Mode	"Normal"
SubarrayHOffs	Integer	Subarray horizontal offset	0
SubarrayHWidth	Integer	Subarray horizontal width	640
SubarrayVOffs	Integer	Subarray vertical offset	0

VOffs			
Subarray VWidth	Integer	Subarray horizontal width	480
CameraName	String	Camera name	„C9100-01“
Type	Integer	Camera type	26
SubType	Integer	Camera subtype	34

C9266

In the case of a C9100 it has the following tokens:

Token	Type	Description	Example
Type	Integer	Camera type	“C9266_02”
SubType	Integer	Camera subtype	59
TriggerMode	String	Internal / External Edge 2 / External Level 1 / External Level 2 / External Edge 1	Internal
Polarity	String	Positive / Negative	Negative
Binning	String	1 x 1 / 2 x 2	1 x 1

C9300

In the case of a C9300 it has the following tokens:

Token	Type	Description	Example
AMD	String	Acquire mode (N=Normal, E=External)	N
AMD	String	Acquire mode (N=Normal, E=External)	N
NMD	String	Normal mode (N=Normal, S=Electrical Shutter, F=Frame Blanking)	F
EMD	String	External mode (E=Edge Trigger, L=Level Trigger)	L
SMD	String	Scan mode	S
ADS	String	AD select (possible values are 8 10 or 12)	12
SHT	String	Shutter (number of lines, possible values are depending on superpixel settings)	1
FBL	String	Frame blanking number (number of frames, possible values are depending on superpixel settings)	1
EST	String	External shutter (number of lines, possible values are depending on superpixel settings)	1
SHA	String	Scan horizontal area (F=Full, K=Kilo)	K
SFD	String	Scan front dummy (flag)	F
SPX	String	Superpixel (possible values are depending on camera version)	2
TNS	String	Tap number select	1
ATP	String	Acquire trigger polarity (P=Positive, NB=Negative)	N
CEG	String	Contrast enhancement gain	0
CEO	String	Contrast enhancement offset	0
ESC	String	External source (B=BNC, D=D-Sub,	I

		I=Digital I/F)	
SHO	String	Scan horizontal offset	0
SHW	String	Scan horizontal width	1024
SVO	String	Scan vertical offset	0
SVW	String	Scan vertical width	1024
Timing Mode	String	Timing mode	Internal timing
Trigger Mode	String	Trigger mode	Level trigger
Trigger Source	String	Trigger source	Digital I/F
Trigger Polarity	String	Trigger polarity	neg.
CCDArea	String	CCD area	1024 x 1024
Binning	String	Binning mode	2 x 2
Scanmode	String	Scan Mode (Normal/External)	Normal
SubarrayHOffs	String	Scan horizontal offset	0
SubarrayHWidth	String	Scan horizontal width	1024
SubarrayVOffs	String	Scan vertical offset	0
SubarrayVWidth	String	Scan vertical width	1024
Camera Name	String	Camera name	C9300
Type	integer	Camera type	29
SubType	integer	Camera subtype	44

C10000

In the case of a C10000 it has the following tokens:

Token	Type	Description	Example
ADS	integer	8 / 12	8
AMD	String	N / T	N
ATP	String	P / N	P
CEG	integer	0 ... 15	1
FBL	Integer	2 ... 16	2
HLN	Long	600 ... 65535	10000
TMD	String	I / E	I
TRD	String	F / R	F
Timing Mode	String	Area / TDI	TDI
Readout Direction	String	Forward / Reverse	Forward

n			
BitsPerPixel	integer	8 / 12	12
Trigger Polarity	String	Pos. / neg.	neg.
Trigger Mode	String	Internal / External	Internal
LinesPerImage	integer	Line number in case of TDI scanning	1024
Camera Name	String	Camera name	„C10000“
Type	integer	Camera type	57

DCam

In the case of a camera supported by DCam it has the following tokens:

Token	Type	Description	Example
CameraName	String	Camera name	„Hamamatsu 1394 C8484“
Type	integer	Camera type	23
Binning	Integer	Binning	"1"
BinningKeyVal	String	Binning	"2 x 2"
LightMode	Integer	Light mode (0/1)	0
Scan Speed	Integer	Scan Speed	1
Scan Mode	Integer	Scan Mode	2
ScanModeKeyVal	String	Scan Mode	"Subarray"
Hoffs	Integer	Horizontal Offset	0
HWidth	Integer	Horizontal Width	512
VOffs	Integer	Vertical Offset	0
VWidth	Integer	Vertical Width	512
Trigger Mode	Integer	Trigger mode	3
TriggerModeKeyVal	String	Trigger mode	"Internal"
TriggerPolarity	Integer	Trigger Polarity	"1"
TriggerPolarityKeyVal	String	Trigger Polarity	"neg."
Gain	Integer	Gain	0

DCAM 3.0

In the case of a camera supported by DCAM 3.0 module it has the following tokens:

Token	Type	Description	Example
CameraName	String	Camera name	„Hamamatsu 1394 C8484“
Type	integer	Camera type	23
Prop_XXX	Double	All available properties	Prop_SubarrayHpos=256

		following the name from the DCAM 3.0 module (Without blanks).	
--	--	---	--

Acquisition section

[Acquisition]

Token	Type	Description	Example
NrExposure	Integer	Number of exposures (integrated in memory)	1
NrTrigger	Integer	Number of trigger	1
ExposureTime	String	Exposure time	20ms
AcqMode	Integer	Acquisition mode (see AcqMode)	1
DataType	Integer	Data type (see DataType)	2
DataTypeOfSingleImage	Integer	Data type of a single exposure (see DataType)	2
BacksubCorr	Integer	Background corrected (flag)	0
ShadingCorr	Integer	Shading corrected (flag)	0
CurveCorr	Integer	Curvature corrected (flag)	0
areSource	String	Source area (valid image, X, Y, DX, DY)	0,0,1000,1018
areGRBScan	String	Scan area (for grabber, X, Y, DX, DY)	8,0,1000,1018
pntOrigCh	String	Origin on chip (X, Y)	0,0
pntOrigFB	String	Origin on frame buffer, (X, Y)	0,0
pntBinning	String	Binning (X, Y)	1,1
BytesPerPixel	Integer	Bytes per pixel	1

Grabber Section

[Grabber]

Token	Type	Description	Example
ConfigFile	String	Configuration file	D:\Program Files\HiPic32\digital.cnf
Type	Integer	Grabber type (see FrameGrabber)	2
SubType	Integer	Acquisition module type (See AcquisitionModule)	1

Display LUT Section

[DisplayLUT]

Token	Type	Description	Example
EntrySize	Integer	Size of LUT (see LUTSize)	2

LowerValue	String	Lower value (lower cursor)	0
UpperValue	String	Lower value (upper cursor)	1024
BitRange	String	Bit range	10 bit
Color	Integer	Color (See LUTColor)	1
LUTType	Integer	LUT Type (See LUTtype)	0
Gamma	String	Gamma value	1
First812OvlColor	Integer	First overlay color in 812 bit mode (overlaid images only)	1
Lut16xShift	Integer	Shift value in case of Lut16x mode	1
Lut16xOvlVal	Long	Overlay vlaue for Lut16x mode	32767

Calibration Section

[Scaling]

Token	Type	Description	Example
ScalingXType	Integer	Calibration type in X direction (see ScalingType)	SCALING_LINEAR
ScalingYType	Integer	Calibration type in Y direction (see ScalingType)	SCALING_LINEAR
ScalingXScale	String	Calibration factor in X direction	1.57
ScalingYScale	String	Calibration factor in Y direction	1.57
ScalingXUnit	String	Calibration unit in X direction	mm
ScalingYUnit	String	Calibration unit in Y direction	mm
ScalingXScalingFile	String	Calibration file in X direction	scaling1
ScalingYScalingFile	String	Calibration file in Y direction	scaling2

Comment Section

[Comment]

Token	Type	Description	Example
UserComment	String	Comment string inputted by the user	This is my comment

Constants

Used constants:

ApplicationType:
APPLICATIONHIPIC = 1
APPLICATIONONTA = 2

APPLICATIONEM = 3

CameraType:

NOCAMERA = 0
C4880 = 1
C4742 = 2
ANALOG0 = 3
ANALOG1 = 4
ANALOG2 = 5
ANALOG3 = 6
C474295 = 7
C488080 = 8
C474298 = 10
C7300 = 19
C800020 = 20
C800010 = 21
FlatPanel = 22
DCam = 23
OrcaHR = 24
C8484 = 25
C9100 = 26
C8800 = 27
LineSensor = 28
C9300 = 29
C474295TDI = 30
USBFlatPanel = 31
GenericCam = 32
DCam30 = 33
DCam213Plus = 34
C10000 = 35
C9266 = 36

CameraSubType:

NoSubType = 0
C4880_00 = 1
C4880_60 = 2
C4880_80 = 3
C4880_91 = 4
C4880_92 = 5
C4880_93 = 6
C4742_95 = 7
C4880_60OU = 8
C4880_1K2K = 9
C4742_98 = 11
C7300_10 = 12
C4880_10 = 13
C4880_20 = 14
C4880_21 = 15
C4880_30 = 16
C4880_40 = 17
C7190_10 = 18
C8000_20 = 19
C4742_95ER = 20
C4880_31 = 21
C4880_50 = 22
C8000_10 = 23
FlatPanel_C7942 = 24
FlatPanel_C7943 = 25
C4742_95HR = 26
C8484_01 = 27
FlatPanel_C7921 = 28
C4742_98ER = 29

C4742_98BT_K = 30
C4742_98BT_L = 31
C7300_10_NRK = 32
FlatPanel_C7930DP = 33
C9100_01 = 34
C9100_02 = 35
C9100_11 = 36
C9100_12 = 37
C8800_01 = 38
C8800_21 = 39
C7390 = 40
C8133 = 41
C8750 = 42
C9133 = 43
C9300_001 = 44
C9300_101 = 45
C9300_201 = 46
C9300_501 = 47
FlatPanel_USB = 48
FlatPanel_C9250DP = 49
FlatPanel_C9312 = 50
FlatPanel_C9321 = 51
FlatPanel_C9311 = 52
GenericCam_01 = 53
FlatPanel_C10013 = 54
C9100_13 = 55
C10000_01 = 56
C9300_502 = 57
C9266_01 = 58
C9266_02 = 59
FlatPanel_C9252DK = 60
FlatPanel_C10322D = 61
C9750 = 62
FlatPanel_C10500D = 63
FlatPanel_C10502D = 64
C10000_401 = 65
FlatPanel_Dummy = 67

AcqMode:

LIVE = 1
ACQUIRE = 2
PHOTONCOUNTING = 3
ANALOGINTEGRATION = 4

DatType:

DAT8 = 1
DAT10 = 2
DAT12 = 3
DAT16 = 4
DAT812 = 5
DAT14 = 6
DAT16u = 7
DAT32 = 8

FrameGrabber:

grbNone = 0
grbAFG = 1
grbICP = 2
grbPC = 3
grbNI = 4
grbDCam = 5
grbUSBFlatPanel = 6

GrbPhx = 7
GrbSap = 8
grbC9266 = 9

AcquisitionModule:

AMDIG = 1
AMVS = 2
CamLink = 3

LUTSize:

LUTSize8 = 1
LutSize10 = 2
LutSize12 = 3
LUTSize16 = 4
LUTSize812 = 5
LUTSize14 = 6
LUTSize16x = 9

LUTColor:

LUTColorBW = 1
LUTColorRainbow = 2
LUTColorBWwithoutColor = 3

LUTType:

LUTTypeLinear
LUTTypeGamma
LUTTypeSigmoid

CalibrationType:

SCALING_LINEAR = 1
SCALING_TABLE = 2

Details of sequence mode

Following information may be useful to use the sequence mode most effectively:

The program for sequence mode has been designed to allow fast and flexible image acquisition. However these two targets can not be reached under all operating conditions at the same time. Therefore two submodes of the sequence mode can be distinguished:

A mode optimized for highest speed (Mode S) and a mode optimized for best flexibility (Mode F).

The Mode S is the default mode and is used when **all** of following conditions are true:

- The images are saved to RAM
- No profiles are saved
- Live or Acquire mode is used
- Interval is “full speed” or < 1 second
- No auto corrections are used or “Do auto correction after sequence” is selected

The mode F is used when **one** of following conditions is given:

- Analog or Photon counting mode is used
- No images are saved (neither to RAM nor to hard disk)
- Images are saved to hard disk

- Profiles are saved to hard disk
- Fixed interval is used with interval time > 1 second
- Auto corrections are executed directly after image acquisition

Additional remarks:

Wrap during acquisition is only available when Live mode is used.

The image is displayed and the Auto LUT is executed after every sample if either **Always display image during acquisition** is selected or fixed interval with interval time >2 seconds is used or if no images or profiles are saved (neither to RAM nor to hard disk).

Histogram and Profile update is done if fixed interval with interval time > 2 seconds is used.

If fixed interval with interval time > 10 seconds is used the timing information in the image caption (and in the header of the image files) is formatted like hh:mm:ss, using the system standard time. Otherwise it is formatted like xxxxms showing the time passed after acquisition of the first image in milliseconds.

Only if images are stored to hard disk every image status can have its individual information. Only in these cases the sBackString of the UserFunction will be appended to the status string and saved to the image header.


LUT

There are three types of LUT:

Linear LUT (default)

This LUT makes a linear grayscale transformation between the input image data and the displayed data.


$$Int_{out} = \frac{Int_{in} - L}{H - L} * 253$$



Gamma LUT

This LUT makes a nonlinear grayscale transformation according to following algorithm:

$$Int_{out} = \left[\frac{Int_{in} - L}{H - L} \right]^k * 253$$




Sigmoid LUT

This LUT makes a nonlinear grayscale transformation according to following algorithm:

$$Int_{out} = \left[\frac{Int_{in} - L}{\frac{H - L}{2}} \right]^k * 126 \quad \text{for } Int_{in} < \frac{H - L}{2} + L$$

and

$$Int_{out} = 253 - \left[\frac{H - Int_{in}}{\frac{H - L}{2}} \right]^k * 126 \quad \text{for } Int_{in} \geq \frac{H - L}{2} + L$$



Error handler

Though the program is carefully designed, planned, coded, tested and debugged, software which exceeds a certain amount of code cannot be absolutely bug free in every circumstances and with all available parameter settings. As a consequence this program also may have certain bugs which are not detected during test but which may appear at customers site. To allow easy debugging and ensure program quality a powerful error handler has been established within this program. When a fatal error occurs the chain of function calls and other important data is written to a file before the program is ending. This file can be used to locate the origin of the error as quickly as possible.

When such an error occurs a message-box appears informing the user that such an error has appeared and all information has been written to the file **ERRORS.TXT** within the directory in which the programs EXE file is located. If such bug has occurred at your side please send this file to your local **Hamamatsu** dealer.

Glossary of Terms

Acquire Mode

A scanning mode of C4880 type cameras. In the C4880 terminology it is called "Slow Scan Mode". This mode allows to acquire single images with the highest possible quality.

Analog Contrast Enhancement

C4880 cameras and C4742 cameras have an analog contrast enhancement circuit built in. It allows to increase the analog amplifier gain (before the signal is digitized) and to subtract a analog offset signal. In case of C4880 cameras the enhancement function works in Live Mode only.

This function can be used to spread the contrast of an image within a desired intensity range.

Analog Video Camera

In this manual standard video cameras, which offer an analog video output signal according to RS 170 or CCIR (System B) standard are called "Analog video cameras". Hamamatsu offers for example the types C3077 (2/3"), C5405 (1/2") and C 5403 (1/3").

Background Subtraction

Subtraction of a background image from an acquired image. The background image is typically a image acquired under complete darkness. This image contains the camera readout noise, dark current and offset signals. Sometimes this function is called "Dark subtraction" or "Dark current subtraction".

Binning Mode

A special readout mode of C4880 and C4742 cameras. Already in the camera the charges from a predefined number of pixel is accumulated. The accumulated charge signal is then digitized and read out. Binning allows to get a higher camera sensitivity and reduce the amount of data. Of course this function works in expense of spatial resolution.

Current image

The image window which is the active window, or which was the last active image window before another (non-image) window was opened.

DCAM, DCAM driver

DCAM driver are unified software drivers supplied for most of the cameras. The Hamamatsu DCAM driver concept allows to change easily between different camera and interface versions.

DDE

Dynamic data exchange. This function allows to send measurement data automatically to another program (like Excel), where they can be further processed. DDE is supported for export of profile and 3-D analysis data.

Fast Scan Mode

Same as Live mode.

High Precision Scan Mode

Same as Acquire mode.

IC-PCI

Frame grabber board with PCI bus. Different versions are available: 2M/DIG with 2 MByte frame memory for digital input signals, 2M/VS with 2MByte frame memory for analog video signals and 4M/DIG with 4 MByte frame memory for digital input signals.

Interlaced

A way of image data readout from standard analog video cameras. One image (frame) is read out as two fields.

Live Mode

A mode to acquire and display images continuously. This mode is available for all cameras. In case of C4880 type cameras, the camera is switched to the fast readout mode.

LUT (Look up table)

A conversion table which is internally calculated to convert intensity values. It can be used to convert black and white images into pseudocolor images or to make contrast enhancement.

Menu Bar

Main application control window. From this window you can control all sub-menus and commands.

Photon Counting

A special image processing mode which can be used in combination of image intensifier with C4880, C4742 or Analog video cameras. If the gain of the intensifier is high, single photon events can be detected by the camera as a relatively bright spot. The photon counting function can process these images and record each spot a single photon event. This mode allows to accumulate low light level images with best signal to noise ratio and highest sensitivity.

Realtime background subtraction

A background image is continuously subtracted from images acquired in Live mode. This function is useful to subtract mottle images or dark current images.

ROI

A Region Of Interest (ROI) is a part of the image which was selected by the user using the ROI definition tools

Calibration

If physical dimensions are assigned to images or profiles, we call it scaling (for example μm or m, wavelength etc.). The scales may be of linear or non-linear type.

Scanning mode

The C4880 allows to work with two different scanning modes. Fast scanning mode is mainly applicable for image adjustment (focusing), while the high precision mode is used for precise image acquisition. In HiPic the fast scanning mode is called "Live mode" and high precision mode is called "Acquire mode"

Sequence

The Sequence function allows to automatically record large numbers of image and/or profile data, stored either in RAM or on hard disk.

Shading Correction

A uneven image brightness is called shading. It can be corrected by multiplying each pixel of an image by a correction factor. Typically the correction factors are obtained from an reference image.

Subarray Readout Mode

A readout mode of C4880 cameras. Only a user-predefined area of the CCD is read out. This function can be used to reduce the amount of data and speed up the readout, if only a part of the image is necessary for analysis.

Superimpose

A image processing function to overlay two images.

Super Pixel Mode

A special binning mode where a square number of pixel is binned (e.g. 2 x 2, 4 x 4). This function can be used to increase camera sensitivity, increase the readout rate (frame rate) and reduce the amount of data.

TIFF

Image file format. TIFF is widely used as image format. Most application programs can import TIFF files.

User Functions

In order to extend the functionality of the program to users needs, it is possible to call user-programmed DLL functions from many locations within the program, especially within the acquisition procedures. These function calls are called user function calls.

Index

3

3D Data 155

A

About 3

Above Threshold in Photon counting 214

Accessing the CCD by RS232 238

Acquire mode 69

 Difference between Live and Acquire mode 69

Acquiring images 21

Acquisition and Options dialog 38

Acquisition modes 67

Analog integration mode 70

Analysis Profiles 143

Arithmetic 219

ASCII images 244

Auto LUT 120

Averaging 205

B

Background correction 72

Background subtraction 79, 126

Background Subtraction 127

Basic operations 19

C

Calibration 124, 137

 File Editor 225

Calibration File Editor 225

Calibration file format 251

Camera and Frame grabber selection 7

Camera data types 70

Camera dialogs 40

CCD camera control 38

Context sensitive menus 33, 115

Corrections

 Background correction 72

 Background subtraction 79

 Defect pixel correction 75

 Defect pixel correction 86

 Defect pixel correction 86

 Realtime background subtraction 73, 83

 Shading correction 74, 84

D

Data storage 70

Data to TIFF images 242

DCam API 15

DCAM API 15

DDE Links 153

Defect pixel correction 75

Defect pixel correction 86

Defect pixel correction 86

Digital contrast enhancement 116

Display negative values 119

Display to TIFF images 242

DPC file format 251

Dynamic file format 251

Dynamic Photon Counting 215

 Acquisition 216

 Analysis 216

 Enable 122, 123, 215

E

Error handler 272

F

Falltime 151

Fast hard disk recording 194

 Optimize speed 195

Features of HiPic 4

File comment 78

File info 78

FWHM 108, 151, 152

G

Getting started 19

 Acquiring images 21

 Acquisition and options dialog 25

 Live mode 21

 LUT 23

 Options 27

 Quick profile 28

 saving and loading images 26

 Setting camera parameters 24

 Zooming 24

H

Help 3

HIS image sequence file format 245

HIS sequence format 202, 203

Histogram 154

I

Image display

 Displaying mouse coordinates 113

 FWHM 108

 Quickprofile 105

 ROIs 109

 Rulers 105

 The caption 104

 Zoom 104

Image Display 103

Image display windows 32

Image file formats 241

Image processing 219

Image sequences 196

Image status 114

Installation

- Internal driver 11
- Installation
 - FireWire cameras 11
 - General 9
 - Setup launcher 9
 - Software 9
- Installation
 - FireWire cameras 11
- Installation
 - Hardware 14
- Installation
 - Hardware 14
- Installation
 - Changing hardware after installation 15
- Installation
 - Changing hardware after installation 15
- Internal drivers 7
- Inverted LUT 119
- ITEX Images 241
- L**
- Lifetime 151
- Line sensors 159
- Linear and table scaling 138
- List of camera parameters 60
- Live mode 21, 69
 - Difference between Live and Acquire mode 69
- Loading images 77
- Look up table 116
- LUT 23, 24, 32, 116
 - Gamma 270
 - Linear 270
 - Sigmoid 270
- M**
- Measurement parameters 121
- N**
- Nonlinear LUT 119
- O**
- Options 27, 121
 - Auto Backsub 126, 128
 - Auto LUT 130, 131
 - User Function 229
- P**
- Permanent parameters 134
- Photon counting 207
 - Condition 209
 - Conditions 208
 - Setup 213
- Photon Counting
 - Conditions 208
 - Setup 213
- Photon counting mode 70
- PIV mode 169
- Print images 115
- Profile analysis
 - Displaying Profile Data 153

- Profile Analysis 151
 - Acquiring Profiles 145
 - Displaying Profile Data 153
 - Displaying Profiles 147
 - Profile Scaling 150
- Profile file format 250
- Profile sequences 200
- Profiles 143
- Q**
- Quick profile 28
- Quickprofile 105
- R**
- Realtime background subtraction 73, 83
- Remote control programming 237
- Restore window positions 136
- Risk of overexposure 2
- ROI 109
- ROIs 32, 109, 113
- RS 232 238
- RTBS 73
- Rulers 105
- S**
- Saturation in Photon counting 214, 215
- Save and restore measurement parameters 134
- Safety intructions 2
- Saving images 76
- Selecting Regions of Interest 109
- Sequence 189
 - Averaging 205
 - Details 269
 - Image display 200
 - Image file naming conventions 202
 - Processing 205
 - Profile acquisition 203
 - Profile display 204
- Sequence acquisition 189
- Sequence acquisition parameters 190
- Sequence data storage 196
- Serial communication
 - Problems during startup 20
 - RS232 timeout 20, 21
- Setting camera parameters 24
- Shading correction 74, 84
- Specifying ROI coordinates 113
- Start the program 19
- Status string format 252
- Superimpose 221
 - Bottom layer 222
 - Top layer 221
- Supported CCD cameras 6
- Supported frame grabbers 6
- System, image and profile calibration 139

T

TDI mode 176

TDI mode camera 177

The LUT 116

The menu

- Acquisition menu 35

- Analysis menu 36

- Correction menu 36

- Display menu 36

- File menu 35

- General 34

- Info menu 37

- Main menu 34

- Processing menu 36

- Tools menu 35

- Window menu 37

The ROI I/F 113

The statusbar 30

The toolbar 30, 37

The user I/F

- Child windows and system modal windows 30

- Context sensitive menus 33

- Handling Windows 31

- Options, acquisition and tool dialogs 30

- ROIs 32

- The image display 32

- The LUT 32

- The menu 29

- The statusbar 30

- The toolbar 30

The user interface 29

threshold 209

TIFF Images 242

Tools

- Calibration

 - File Editor 225

- Defect pixel correction tool 227

U

User Function 229

User function in sequence mode 236

W

Workfiles 134

X

X-Ray Flat panels 171

X-ray Line Sensors

- Acquisition 165

X-Ray line sensors 159

Z

Zoom 104