

FLAT PANEL SENSORS



Digital X-ray image sensors developed as key devices
for real-time X-ray imaging applications requiring high sensitivity and high image quality

Contents

1. X-ray Image Sensor Technology	1
1-1. CMOS Process Technology	1
1-2. Scintillator Technology	1
1-3. Assembly Technology	1
1-4. IC Design Technology	2
2. Line-up of Flat Panel Sensors	3
2-1. Radiology	3
Rotational radiography	3
Biochemistry	7
Diffraction	9
2-2. Non-destructive Inspection	11
General X-ray application	11
High tube voltage application	13
3. Description of Terms	15
4. Operating Principle and Structure	16
4-1. Features	16
4-2. Operating Principle	16
4-3. Active Pixel Construction	18
4-4. Coupling Method of Scintillator	18
4-5. Scintillator	19
4-6. Window Material of Top Cover	19
4-7. X-ray Radiation Damage	19
5. Features	20
5-1. Spatial Resolution	20
5-2. Dynamic Range	20
5-3. Detective Quantum Efficiency	20
6. System Set-up	21
6-1. Set-up Example	21
6-2. Precaution for Power Source Configuration	21
6-3. Frame grabber cables	21
7. Image Samples	22
8. Related Products	28



HAMAMATSU Flat Panel Sensor



HAMAMATSU has been a leading supplier of opto-semiconductor sensors of diverse types and categories including sensors for high energy physics, UV sensors, visible light sensors, infrared sensors, position sensors, and linear or image sensors since 1957. Over the last 20 years, intense research and development have been made into various techniques for capturing X-ray images with solid state image sensors and some remarkable results have been achieved. Besides being ideal for digital imaging, solid state image sensors are expected to replace conventional image acquisition methods based on ortho film or imaging tubes. After years of laboratory development, HAMAMATSU has deemed digital X-ray imaging is ready. A unique flat panel sensor based on the single crystal CMOS technologies was developed in 1987. Then a prototype of CMOS flat panel sensor was fabricated in 1993. This sensor had a large pixel size and still retained an analog interface, but it was a milestone in the digital radiography field. The dedicated world's first monolithic CMOS flat panel sensor with 12 cm photodiode array area was developed in 1999. The latest line-up of flat panel sensor still inherits the advantages of single crystal Si yet has expanded analog and digital performances achieved through CMOS technology and outstanding scintillator.

Line-up of Flat Panel Sensors

Application	Type No.	Scan mode	A/D (bits)	Pixel size (μm)	Photodiode area [(H) · (V) mm]	Frame rate *1 (frames/s)	Scintillator	Energy range	Top cover material (thickness)	CE marking *2	Page No.	
Radiology	C10900D	Fast	13	200	124.8 · 124.8	35	Direct deposition CsI	20 kVp to 90 kVp	—	—	3	
		Partial				70						
		Fine	12	100		17						
		Panoramic				280						
	C9252DK-14	Whole	13	200	249.6 · 124.8	30	CsI FSP	20 kVp to 80 kVp	Carbon fiber (1.0 mm)	—	5	
		Partial	12	100		146						
	Biochemistry	C7942CK-22	—	12	50	120 · 120	2	Direct deposition CsI	17 keV (Mo source)	Carbon fiber (1.0 mm)	○	7
		C9730DK-10	—	14		52.8 · 52.8	4					
		C9732DK-11	—			120 · 120	1					
	Diffraction	C9728DK-10	—	14	50	52.8 · 52.8	3	Direct deposition CsI	18 keV Max.	Carbon fiber (1.0 mm)	○	9
C10158DK		Whole	14	50	118.8 · 118.8	3						
	Partial	12										
Non-destructive inspection (offline) *3	General X-ray	C7921CA-29	12	50	52.8 · 52.8	4	CsI FSP	20 kVp to 100 kVp	Aluminum (1.0 mm)	○	11	
		C7942CA-22			120 · 120	2						
		C7943CA-22			124.8 · 124.8	7						
	High tube voltage	C7942SK-25	—	12	50	120 · 120	2	GOS-deposited FOP	20 kVp to 150 kVp	Carbon fiber (1.0 mm)	○	13
		C9312SK-06	—			124.8 · 115.2	8		20 kVp to 110 kVp			

*1: Single operation

*2: Types marked with a circle (○) conform to the European EMC directives: EN61326-1 Class A.

*3: Not recommended for use in inline applications where large radiation damage may occur.

Note:

Please consult us for other energy ranges not listed above.

The HAMAMATSU flat panel sensor has not been legally approved for medical applications. This means that the flat panel sensor alone cannot be used as medical equipment. When incorporating the flat panel sensor in systems for medical applications, be sure to obtain any required legal approvals.

1. X-ray Image Sensor Technology

1-1. CMOS Process Technology

HAMAMATSU develops and manufactures unique flat panel sensors with large area, low noise and high resolution for X-ray imaging field. The latest Si process line realizes the world largest monolithic CMOS chip where optimized high performance electronics are obtained. HAMAMATSU can provide high-quality X-ray images for various X-ray imaging applications.



1-2. Scintillator Technology

HAMAMATSU has succeeded in developing an optimal and high-sensitivity X-ray device ideal for indirect X-ray detection. A CsI:TI crystal plate with needle structure mounted on the largest formatted photodiode array allows the scintillation to propagate through the fiber-like-crystals. This structure offers advantages in light propagation over other scintillators of compound. Scintillator thickness has been optimized according to the application.



1-3. Assembly Technology

Assembling and packaging of flat panel sensor is of critical importance. It is considerably more challenging than for integrated circuits since flat panel sensors must interact with X-ray irradiation and other harsh environments in higher reliability. The latest robotics techniques were utilized to attain precision alignment, higher reliability assembling and stable quality in large size active area ranging from 5 · 5 cm to 12 · 24 cm.

Figure 1: Wire bonding





Figure 2: Sensor board assembly



1-4. IC Design Technology

In recent years, large formatted CMOS solid state image sensors of monocrystalline material have become available. These offer many options compared to amorphous Si and conventional imaging devices. We have developed many types of flat panel sensors that show the tremendous potential offered by large size image sensor devices made in a 0.6 μm standard CMOS process. These devices have an anti-blooming (overflow drain), a correlated double sampling (CDS) circuit, external frame start and binning functions.

Figure 3: CAD image of circuit design






Figure 4: CAD image of pattern design



Rotational radiography

2-1. Radiology

The C10900D was developed for CT and panoramic imaging. It operates in 4 selectable scan modes to capture X-ray images: "Fast mode" and "Partial mode" with a pixel size of $200 \cdot 200 \mu\text{m}$, and "Fine mode" and "Panoramic mode" with a pixel size of $100 \cdot 100 \mu\text{m}$.



Features

- Four selectable scan mode
- High sensitivity
- High-speed frame rate: 280 frames/s (panoramic mode)
- Flat panel structure without image distortion
- Scintillator: direct deposition CsI
- For assembly into equipment: supplied without case

Applications

- Cone beam CT, etc.
- Digital radiography, etc.

(Ta=25 °C, Typ. unless otherwise noted)

Parameter	C10900D				Unit
	Fast mode	Partial mode	Fine mode	Panoramic mode	
Pixel size	200		100		μm
Photodiode area	124.8 (H) · 124.8 (V)				mm
Number of pixels	624 (H) · 624 (V)		1248 (H) · 1248 (V)		pixels
Number of active pixels	608 (H) · 616 (V)	608 (H) · 310 (V)	1216 (H) · 1232 (V)	1216 (H) · 72 (V)	pixels
Frame rate	35	70	17	280	frames/s
Frame rate external	Sf(int)*1 to 10		Sf(int)*1 to 10	Sf(int)*1 to 50	frames/s
Noise *2	2900		1300		electrons
Sensitivity *3	6000		3500		LSB/mR
Saturation charge	10.5		2.3		M electrons
Resolution *4	2.5		4.5		line pairs/mm
Dynamic range	3600		1700		-
Defect line *5	8 Max.				lines
X-ray tube voltage	20 to 90				kVp
Digital output (interface)	13 (LVDS differential)		12 (LVDS differential)		bits
Scintillator	Direct deposition CsI				-

*1: Frame rate for internal mode

*2: Internal trigger mode


*3: At 80 kVp, acrylic filter 170 mm

*4: Spatial frequency at CTF=5 %

*5: Adjacent defect lines are not allowed. For the definition of defect line, see Chapter 3, "Description of Terms".

Note: For detailed data on the products listed in this catalog, see their datasheets that are available from our website www.hamamatsu.com or contact our sales office.

Dimensional outlines (unit: mm, tolerance: ±1 mm)



KACCA0268EA

Rotational radiography

2-1. Radiology

This flat panel sensor was developed for CT and panoramic imaging application up to 90 kVp. C9252DK-14 has a tiling construction for large image area capability. C9252DK-14 is whole scan and partial scan selectable for CT and panorama application respectively.



Features

- High sensitivity
- High-speed frame rate
- Flat panel structure without image distortion
- Scintillator: direct deposition CsI
- Large area size with tiling construction
- Partial scan and whole scan selectable

Applications

- Panoramic imaging
- CT imaging, etc.

2. Line-up of Flat Panel Sensors

(Ta=25 °C, Typ. unless otherwise noted)

Parameter	C9252DK-14	Unit
Pixel size	200	μm
Photodiode area	249.6 (H) · 124.8 (V)	mm
Number of pixels	1248 (H) · 624 (V)	pixels
Number of active pixels	1216 (H) · 616 (V)	pixels
Frame rate	30	frames/s
Frame rate external	Sf (int)*1 to 10	frames/s
Noise *2	2600	electrons
Sensitivity *3	6000	LSB/mR
Saturation charge	8.3	M electrons
Resolution *4	2.5	line pairs/mm
Dynamic range	3200	-
Defect line *5	20	lines
X-ray tube voltage	20 to 90	kVp
Digital output (interface)	13 (LVDS differential)	bits
Scintillator	Direct deposition CsI	-
Top cover material (thickness)	Carbon fiber (1.0 mm)	-
Suitable frame grabber cable	A8406-6* series	-

*1: Frame rate for internal mode *2: Internal trigger mode *3: At 80 kVp, acrylic filter 170 mm *4: Spatial frequency at CTF=5 %
 *5: Adjacent defect lines are not allowed. For the definition of defect line, see Chapter 3, "Description of Terms".


■ Partial scan

Pixel size	100	μm
Number of active pixels *6	2432 (H) · 100 (V)	pixels
Partial frame speed	146	frames/s
Noise	1300	electrons
Saturation charge	4.2	M electrons
Sensitivity *7	1500	LSB/mR
Dynamic range	3200	-
Resolution	4.5	line pairs/mm
Digital output (Interface)	12 (LVDS differential)	bits

*6: The position of the active area under the partial scan mode is described in dimensional outline *7: 80 kVp, acrylic filter 170 mm

Note: For detailed data on the products listed in this catalog, see their datasheets that are available from our website www.hamamatsu.com or contact our sales office.

Dimensional outline (unit: mm, tolerance: ±1 mm)



KACCA0221EA

Biochemistry

2-1. Radiology

The biochemical application series C9730DK-10 and C9732DK-11 have fine 50 μm pixel pitch and directly deposited CsI scintillator onto a large formatted photodiode array optimized to 17 keV radiation. They have a cassette shape with a carbon fiber top cover. C7942CK-22 installs CsI flipped scintillator plate (FSP) which is optimized for low absorption material radiography of less than 80 kVp energy range.



Features

- High resolution
- No dead area (insensitive area) due to seamless structure
- Scintillator: CsI flipped scintillator plate (C7942CK-22)
direct deposition CsI (C9730DK-10, C9732DK-11)
- Top cover material: carbon fiber

Applications

- Soft X-ray radiography
- X-ray imaging, etc.
(C9730DK-10, C9732DK-11: optimized to 17 keV)

2. Line-up of Flat Panel Sensors

(Ta=25 °C, Typ. unless otherwise noted)

Parameter	C7942CK-22	C9730DK-10	C9732DK-11	Unit
Pixel size	50	50	50	μm
Photodiode area	120 (H) · 120 (V)	52.8 (H) · 52.8 (V)	120 (H) · 120 (V)	mm
Number of pixels	2400 (H) · 2400 (V)	1056 (H) · 1056 (V)	2400 (H) · 2400 (V)	pixels
Number of active pixels	2240 (H) · 2344 (V)	1032 (H) · 1032 (V)	2368 (H) · 2340 (V)	pixels
Frame rate (single operation)	2	4	1	frames/s
Frame rate (2 · 2 binning)	4	-	-	frames/s
Frame rate (4 · 4 binning)	9	-	-	frames/s
Frame rate external (single operation)	Sf (int) ^{*1} to 0.1	Sf (int) ^{*1} to 0.5	Sf (int) ^{*1} to 0.5	frames/s
Noise ^{*2}	1100	1250	1250	electrons
Sensitivity	35 ^{*3}	65 ^{*4}	65 ^{*4}	LSB/mR
Saturation charge	2.2	6.4	6.4	M electrons
Resolution ^{*5}	8	10	10	line pairs/mm
Dynamic range	2000	5100	5100	-
Defect line	20 Max. ^{*6}	10 Max. ^{*6}	16 Max. ^{*7}	lines
X-ray tube voltage	20 kVp to 80 kVp	Less than 35 kVp (17 keV Mo source)		-
Digital output (interface)	12 (RS-422 differential)	14 (USB 2.0)		bits
Scintillator	CsI flipped scintillator plate	Direct deposition CsI		-
Top cover material (thickness)	Carbon fiber (1.0 mm)			-
Suitable frame grabber cable	A8406-5* series	-	-	-

*1: Frame rate for internal mode

*2: Internal trigger mode and single operation

*3: At 80 kVp, without filter

*4: At Mo target 30 kVp, without filter


*5: Spatial frequency at CTF=5 %

*6: Adjacent defect lines are not allowed. For the definition of defect line, see Chapter 3, "Description of Terms".

*7: Only one pair of adjacent defect lines is allowed. For the definition of defect lines, see Chapter 3, "Description of Terms".

Note: For detailed data on the products listed in this catalog, see their datasheets that are available from our website www.hamamatsu.com or contact our sales office.

Dimensional outlines (unit: mm, tolerance: ±1 mm)



KACCA0259EA

KACCA0222EA

KACCA0276EA

Diffraction

2-1. Radiology

These flat panel sensors feature a low noise and direct deposition CsI scintillator for applications where diffraction is critical.



Features

- Active pixel sensor
- Low noise: 80 electrons (C9728DK-10)
- High quality image
- With USB 2.0 interface (C9728DK-10)
- Scintillator: direct deposition CsI
- Top cover material: carbon fiber
- Partial scan and whole scan selectable (C10158DK)

Applications

- Diffraction
- Radiography, etc.

2. Line-up of Flat Panel Sensors

(Ta=25 °C, Typ. unless otherwise noted)

Parameter	C9728DK-10	C10158DK	Unit
Pixel size	50	50	μm
Photodiode area	52.8 (H) · 52.8 (V)	118.8 (H) · 118.8 (V)	mm
Number of pixels	1056 (H) · 1056 (V)	2376 (H) · 2376 (V)	pixels
Number of active pixels	1032 (H) · 1032 (V)	Whole scan	2352 (H) · 2352 (V)
		Partial scan	2352 (H) · 528 (V)
Frame rate	3	Whole scan	3
		Partial scan	12
Frame rate external	Sf (int)* ¹ to 0.1	Sf (int)* ¹ to 1.0	frames/s
Noise * ²	80	180	electrons
Sensitivity * ³	450	220	LSB/mR
Saturation charge	0.45	1.1	M electrons
Resolution * ⁴	10	10	line pairs/mm
Dynamic range	5600	6000	-
Defect line * ⁵	10 Max.	30 Max.	lines
Adjacent defect line	0 Max.	2 Max.	lines
X-ray energy range	18 Max.	18 Max.	keV
Digital output (interface)	14 (USB 2.0)	14 (LVDS)	bits
Scintillator	Direct deposition CsI		-
Top cover material (thickness)	Carbon fiber (1.0 mm)		-
Suitable frame grabber cable	-	A8406-5* series	-

*1: Frame rate for internal mode

*2: Internal trigger mode

*3: At Mo target 25 kVp, without filter


*4: Spatial frequency at CTF=5 %

*5: For the definition of defect line, see Chapter 3, "Description of Terms".

Note: For detailed data on the products listed in this catalog, see their datasheets that are available from our website www.hamamatsu.com or contact our sales office.


Dimensional outlines (unit: mm, tolerance: ±1 mm)

C9728DK-10



KACCA0224EA

C10158DK



KACCA0225EA

General X-ray application

2-2. Non-destructive Inspection

These flat panel sensors install high quality CsI flipped scintillator plate (FSP) and aluminum top cover. These flat panel sensors offer 50 or 100 μm pixel size and high resolution of 8 or 5 line pairs/mm for general X-ray applications. Maximum radiation energy is less than 100 kVp and guaranteed dose is 10^6 Roentgens. These flat panel sensors are available to use with microfocus X-ray sources having energy range up to 100 kV.



Features

- High quality image
- Wide dynamic range
- High-speed imaging
- Flat panel structure without image distortion
- Scintillator: high quality CsI flipped scintillator plate
- Top cover material: aluminum

Applications

- Non-destructive inspection
- Biochemical imaging
- X-ray microscopy
- Digital X-ray photography, etc.

2. Line-up of Flat Panel Sensors

(Ta=25 °C, Typ. unless otherwise noted)

Parameter	C7921CA-29	C7942CA-22	C7943CA-22	Unit
Pixel size	50	50	100	μm
Photodiode area	52.8 (H) · 52.8 (V)	120 (H) · 120 (V)	124.8 (H) · 124.8 (V)	mm
Number of pixels	1056 (H) · 1056 (V)	2400 (H) · 2400 (V)	1248 (H) · 1248 (V)	pixels
Number of active pixels	1032 (H) · 1012 (V)	2240 (H) · 2344 (V)	1216 (H) · 1220 (V)	pixels
Frame rate (single operation)	4	2	7	frames/s
Frame rate (2 · 2 binning)	8	4	15	frames/s
Frame rate (4 · 4 binning)	16	9	30	frames/s
Frame rate external (single operation)	Sf (int)*1 to 0.1	Sf (int)*1 to 0.1	Sf (int)*1 to 0.1	frames/s
Noise *2	1000	1100	2300	electrons
Sensitivity *3	18	25	28	LSB/mR
Saturation charge	2.9	2.2	10	M electrons
Resolution *4	8	8	5	line pairs/mm
Dynamic range	2900	2000	4300	-
Defect line *5	10 Max.	20 Max.	10 Max.	lines
Power source	A.vdd=+5.0 V, D.vdd=+5.0 V	A.vdd=+5.0 V, D.vdd=+5.0 V, V (±7.5)=±7.5 V		-
X-ray tube voltage	20 to 100			kVp
Digital output (Interface)	12 (RS-422 differential)			bits
Scintillator	CsI flipped scintillator plate			-
Top cover material (thickness)	Aluminum (1.0 mm)			-
Suitable frame grabber cable	A8406-5* series			-

*1: Frame rate for internal mode

*2: Internal trigger mode, single operation


*3: At 80 kVp, without filter

*4: Spatial frequency at CTF=5 %

*5: Adjacent defect lines are not allowed. For the definition of defect line, see Chapter 3, "Description of Terms".

Note: For detailed data on the products listed in this catalog, see their datasheets that are available from our website www.hamamatsu.com or contact our sales office.

Dimensional outlines (unit: mm, tolerance: ±1 mm)



	C7942CA-22	C7943CA-22
A	112	121.6
B	117.2	122
C	33.5	31.1
D	22.9	20.5

KACCA0261EA

High tube voltage application

2-2. Non-destructive Inspection

These flat panel sensors install shielded FOP with GOS scintillator to protect CMOS chip against higher energy radiation. These sensors have 50 μ m pixel pitch. C7942SK-25 has maximum radiation energy up to 150 kVp and C9312SK-06 has 110 kVp Max. Guaranteed dose is $6 \cdot 10^6$ Roentgens.



Features

- High quality image
- Wide dynamic range
- High-speed imaging
- Energy range:
150 kVp Max. (C7942SK-25)
110 kVp Max. (C9312SK-06)
- Scintillator: GOS-deposited FOP
- Top cover material: carbon fiber

Applications

- Non-destructive inspection (offline)
- Radiography, etc.

The C9312SK-06 is warranted for a period of 12 months after the date of shipment. Even within the warranty period, if the total radiation dose exceeds 6 million Roentgen against 110 kVp X-ray irradiation, it is not warranted. In case of 8-hour operation in a day, it is equivalent to 92 days on the condition that the 110 kVp continuous X-ray is exposed to the detector as the intensity of 4096 LSB output can be taken at 8 frames per second (maximum frame rate for single operation). In a real operation, the dark and light field images may be taken frequently for offset and sensitivity calibration. The flat panel sensor may be exposed to X-ray before or after the valid image is taken. It is needed to count up all these shots for proper estimation. Please calculate the life time considering the operation condition in your system.

2. Line-up of Flat Panel Sensors

(Ta=25 °C, Typ, unless otherwise noted)

Parameter	C7942SK-25	C9312SK-06	Unit
Pixel size	50	50	μm
Photodiode area	120 (H) · 120 (V)	124.8 (H) · 115.2 (V)	mm
Number of pixels	2400 (H) · 2400 (V)	2496 (H) · 2304 (V)	pixels
Number of active pixels	2316 (H) · 2316 (V)	2472 (H) · 2184 (V)	pixels
Frame rate (single operation)	2	8	frames/s
Frame rate (2 · 2 binning)	4	-	frames/s
Frame rate (4 · 4 binning)	9	-	frames/s
Frame rate external (single operation)	Sf (int)* ¹ to 0.2	Sf (int) * ¹ to 1.0	frames/s
Noise * ²	1300	700	electrons
Sensitivity * ³	58	20	LSB/mR
Saturation charge	1.1	2.25	M electrons
Resolution * ⁴	8	9	line pairs/mm
Dynamic range	850	3200	-
Defect line	30 Max. * ⁵	20 Max. * ⁶	lines
X-ray tube voltage	20 to 150	20 to 110	kVp
Digital output (Interface)	12 (RS-422 differential)	12 (LVDS differential)	bits
Scintillator	GOS-deposited FOP		-
Top cover material (thickness)	Carbon fiber (1.0 mm)		-
Suitable frame grabber cable	A8406-5* series	A8406-4* series	-

*1: Frame rate for internal mode

*2: Internal trigger mode and single operation

*3: At 150 kVp, without filter (C7942SK-25)
At 110 kVp, without filter (C9312SK-06)


*4: Spatial frequency at CTF=5 %

*5: Only one pair of adjacent defect lines is allowed. For the definition of defect line, see Chapter 3 "Description of terms"


*6: Adjacent defect lines are not allowed. For the definition of defect line, see Chapter 3 "Description of terms"

Note: For detailed data on the products listed in this catalog, see their datasheets that are available from our website www.hamamatsu.com or contact our sales office.

Dimensional outlines (unit: mm, tolerance: ±1 mm)



KACCA0262EA



KACCA0179EB

3. Description of Terms

Active pixels

Pixels in the active area used for imaging excluding the shielded portion surrounding the photodiode area.

Binning Mode

Binning mode is a method for taking images as superpixels that are the sum of several pixels, for example, 2 · 2 or 4 · 4 pixels. See 4-2 for details.

Defect line

A defect line is a horizontal or vertical line containing 4 or more consecutive pixels located at the opposite side of an amplifier array or a shift register, that produce 1/8 * of the sensitivity average of the surrounding pixels.

Note 1: Generally, the opposite side of an amplifier array and a shift register array corresponds to the 1st row and 1st column of the image data. On two-port devices, the image data format is different so refer to the individual datasheet.

Note 2: The port is the number of output digital data sets that are sent in parallel to a frame grabber board. The C9312SK-06 and C9252DK-14 are two-port flat panel sensors.

* The definition for some products differs from this value. Please consult the individual datasheets.

Dynamic Range

The ratio of saturation charge to noise.

Fill Factor

The ratio of active area (no-shielded area) to total pixel area. A larger fill factor can reduce light detection loss.

Frame rate

The vertical sync (Vsync) frequency determines the frame rate. In internal trigger mode, video and sync signals always flow from the flat panel sensor at the highest frame rate. External trigger mode can drive the flat panel sensor at a lower frame rate. See 4-2 for details.

Note: The actual frame rate in individual software depends on the CPU processing performance as well as the display board or software load. The software might skip a few frames if the processing performance is lower than system requirements.

Noise

The rms noise swing width in units of electrons, measured at 25 °C in internal trigger mode and single operation.

Photodiode area

Area on the chip where a photodiode array is formed.

Resolution

Spatial frequency at CTF=5 %, expressed in LP/mm (LP: line pairs). Resolution is measured using a 60 kV microfocus X-ray source.

Saturation charge

The maximum output bit in units of electrons.

Scintillator

A material that emits a small flash of light (or scintillation) when struck by nuclear particles or radiation. HAMAMATSU offers various scintillators to cover a wide energy range of incident X-rays.

Sensitivity

Output bit per 1 mR irradiation, expressed in units of LSB/mR (LSB: least significant bit). The full scale of a 12-bit video output is 4095 LSB, 8191 LSB for 13-bit type and 16383 LSB for the 14-bit type. See tube voltages or filter conditions described in each data sheet.

X-ray Tube Voltage


The voltage applied across the anode and cathode of an X-ray tube to accelerate electrons emitted from the cathode filament in order to hit the X-ray tube target. When a voltage V_p is applied to the X-ray tube, X-rays are generated at an energy corresponding to V_p eV. The X-ray tube voltage range guaranteed for flat panel sensors will differ depending on the sensor type.

4. Operating Principle and Structure

4-1. Features

Flat panel sensor is a lightweight and compact X-ray imaging sensor module consisting of a sensor board and a control board. The sensor board has plural blocks of charge-sensitive amplifier arrays with CDS (Correlated Double Sampling) circuit and horizontal shift register. Amplified analog video signals sequentially flow from each block and transferred to the control board. The analog signal is converted into a digital signal and outputs to an external frame grabber through the parallel port.

Figure 5: Example for block diagram



4-2. Operating Principle

The CMOS flat panel sensor chip consists of large formatted photodiode matrix, low noise charge integration amplifier and vertical shift register.


A trans-impedance amplifier with a feedback capacitor converts photocurrent into an output voltage proportional to the input photon intensity. The output voltage of charge integration amplifier is proportional to the product of the photon intensity and integration time. The long integration time allows detecting extremely faint radiation even with a tiny pixel size. The integration time in internal mode is specified in each data sheets. The integration time in external mode is determined by interval of the rising edges of a couple of successive external trigger pulses. The flat panel sensor cannot detect a change of radiation within the integration time. The vertical shift pulse is applied to every row of photodiode matrix sequentially. The start timing of the integration time of each row has a time lag that is equal to the cycle time of a horizontal synchronous pulse.

The capacitor of each photodiode element has limited value, so the integration charge is also limited. This is defined as the saturation charge in the data sheet. The linear output cannot be obtained if over radiation dose corresponds to the value is irradiated.

Each vertical shift register pulse, charge amplifier drive pulse, horizontal shift register pulse, A/D converters and memory chip drive pulse are generated from associated electronic circuits within the sensor case. The user does not have to design and build an intricate timing generator but needs only connect it to a digital grabber card to acquire an X-ray image.


The integration time is the reciprocal of the frame rate. The frame rate specified in internal mode is the maximum speed. It corresponds to the shortest integration time. The longer integration time is attained in external trigger mode, which supplies a longer cycle of external trigger pulse to the flat panel sensor.

Figure 6: Timing chart



In the flat panel sensor drive sequence, the high-sensitivity photodiode matrix generates carriers and accumulates them in the junction capacitance of each pixel when they receive scintillation light from the scintillator material. All pixels have an overflow drain to protect against blooming in adjacent pixels even when a portion of the pixels is saturated. The standard CMOS process implemented with a photodiode matrix achieves a high fill factor ratio of 79 % for 50 · 50 μm pixels and 87 % for 100 μm pixels. (See Figure7.)

Figure 7: Photodiode array



Shift pulses are supplied by the vertical shift register for sequential scanning and a row of photodiodes flows the carrier into each data lines. The on-chip extremely-low noise charge amplifier array with CDS circuit is somewhat complicated in terms of structure and operation, however it can cancel out the offset components by finding the signal difference between the accumulated charge and zero level. This makes a huge improvement in image output uniformity. Usually, some kinds of corrections are made to an acquired image before it is displayed, but the flat panel sensor was designed to have a high level of image quality even before making those corrections. The output noise is mainly determined by the charge amplifier itself and the data line capacitance. The ENC (total Equivalent Noise Charge) is given by the following equation.

$$ENC = 8KT/3 gm \cdot (Ct)^2 \dots\dots\dots (1)$$

$$Ct = Cd + Cp + Cf \dots\dots\dots (2)$$

- K : Boltzman constant $1.3806 \cdot 10^{-23}$
- T : Absolute temperature [K]
- gm : Transconductance [S]
- Cd : Junction capacitance of photodiode [F]
- Cp : Data line capacitance [F]
- Cf : Feedback capacitance of the charge amplifier [F]

● **Binning mode**

The flat panel sensor can deliver high-resolution images in the single pixel drive mode. On the other hand, some types of the flat panel sensors have a binning drive mode for high-speed and high-sensitivity operation. A user can select 2 · 2 or 4 · 4 binning mode. In the 2 · 2 binning mode, the neighboring 2 · 2 pixels are read out together. In the 4 · 4 binning mode, the neighboring 4 · 4 pixels are read out together. So the user can drive the module with higher frame rate controlling the binning mode. Moreover when it is driven at the same frame rate, the amount of signals becomes 4 times larger by 2 · 2 binning mode and 16 times by 4 · 4 binning mode.

● **Trigger Mode**

The flat panel sensors have two trigger modes (internal/external trigger mode). Under the internal trigger mode, video signal and synchronous signals always flow out from the flat panel sensor at the highest frame speed of Sf (int), which is described in the electrical specification in its data sheets. On the other hand, the user can choose the external trigger mode to synchronize the flat panel sensor with the X-ray source or to expand the accumulation period to make the sensitivity higher. The period of this external trigger signal is equal to that of the integration time. The sensor starts to forward the video data to the frame grabber board at the rising edge of this signal. The duty of this trigger signal should be settled between 1 % and 99 % (50 % is recommended). The voltage of this signal should be compatible with TTL-level.

When the user selects internal trigger mode, the signal of “IntExt” or “IntExtTrg” should be set to Low. If “IntExt” or “IntExtTrg” is set to High, the flat panel sensor works as the external trigger mode.

In case of the external trigger mode, there are two methods to control integration period. One is supplying the synchronous signal through “ExtTrgGrb” from PC. The other is supplying the synchronous signal through “ExtTrgLemo” or “ExtTrgIO” from X-ray source. Please refer to the data sheet for the flat panel sensor.

Under the external trigger mode, the first image, which follows the first external trigger signal, is bright and useless because the integration time is not defined. It means that charge from dark current continuously accumulated until the trigger pulse is supplied. Therefore, the second image or later ones are valid.


4-3. Active Pixel Construction

HAMAMATSU CMOS flat panel sensors adopt two kinds of amplifier circuit constructions ; one is passive pixel type for ordinary application, and the other is active pixel type for high-end use.

The passive pixel type incorporates amplifier arrays existing at one side of photodiode array and each amplifier is connected to one column of photodiode array via address switches. This simple structure allows high fill factor, high productivity and high radiation tolerance although there is limitation to reduce thermal noise of amplifier dramatically due to its large input capacitance caused by long data line.

The active pixel type solves the drawback of the passive pixel type. Its low noise or high S/N feature is appealing on low energy or low exposure application like an high quality imaging of low contrast objects or usage of an open-type microfocus X-ray source. This type has an amplifier built in each pixel, and accumulated charges are converted into voltage. The noise level is less than 180 electrons rms, which is less than 1/6 of the passive pixel type.


Figure 8: Internal circuit of CMOS chip



4-4. Coupling Method of Scintillator

HAMAMATSU flat panel sensor utilizes a needle structure CsI as an X-ray conversion layer. There are three methods for coupling a scintillator structure with a CMOS device. The first one is direct deposition of CsI onto the photodiode matrix, the second one is a CsI flipped scintillator plate (FSP). Both of these approaches utilize the needle structure CsI scintillator. A combination of CsI scintillator and CMOS large formatted photodiode array is essential for high-resolution devices, and the direct deposition CsI flat panel sensor in particular shows outstanding resolution. On the other hand, the FSP offers the advantage of relatively low cost. Another scintillator component is GOS-deposited FOP. It can be used higher tube voltages application.

Figure 9: Cross sectional figure of indirect and direct conversion devices



4-5. Scintillator

Achieving a spectral response characteristic that matches the peak wavelength and spectral range of the scintillator emission is a critical factor in photodiode matrix design. By controlling the impurity profile and passivation, HAMAMATSU has succeeded in developing an optimal high-sensitivity X-ray device ideal for indirect X-ray detection.

The second to the last letter of each product type number (e.g. C7942CA-22) indicates which type of scintillator is used.

Type C utilizes a cost effective and high resolution scintillator plate. Type D is an direct deposition type flat panel sensor which realizes high resolution and high sensitivity. Type S has GOS deposited X-ray shielded fiber optics plate against higher input irradiation.

C: CsI:TI [CSI Flipped Scintillator Plate (FSP)], D: CsI:TI (Directly deposited on photodiode array)
S: GOS (GOS is deposited onto an X-ray shielded FOP)

Figure 10: Scintillation spectrum of CsI and photo sensitivity of a pixel




Figure 11: The needle-like crystal structure of CsI



4-6. Window Material of Top Cover

The last letter of each product number (e.g. C7942CA-22) indicates which type of window is used.

A: Aluminum (1 mmt), K: Carbon fiber (1 mmt)

The top cover protects the device against physical damage, extends the service life of the CsI and photodiode array, and acts as a filter for the scattered soft X-rays.

The aluminum top cover is adequate for radiation of the X-ray tube voltage up to 100 kVp. The carbon fiber top cover, on the other hand, have high transmittance (more permeable to X-rays) and are the right option for detecting lower energy X-ray of less than 80 kVp.

The carbon fiber is also applied for the flat panel sensor on which X-ray shielded FOP is mounted.

4-7. X-ray Radiation Damage

The flat panel sensor shows the degradation of characteristics by X-ray irradiation. In general X-ray sensitivity decreases because of the browning of scintillator and the dark output increases by the radiation damage of the photodiode. In some cases, partial dark output increase is also observed. In the case of C7921CA-29, for example, sensitivity decrease around 25 % and the dark output increase of around 200 LSB/frame at internal frame rate are observed after the exposure of 1 million Roentgen at the tube voltage of 100 kVp. "SK" type is relatively tolerant against X-ray than the "CA" type because of the high absorption efficiency of X-ray shielded FOP. The X-ray life time of the "C9312SK-06", for example, is 6 million Roentgen at 110 kVp. In case of 8-hour operation in a day, it is equivalent to 92 days on the condition that the 110 kVp continuous X-ray is exposed to the detector as the intensity of 4096 LSB output can be taken at 8 frames per second (maximum frame rate for single operation).

$$\frac{6 \text{ [million Roentgen]} \cdot 16^{*1} \text{ [LSB/mR]} \cdot 0.9^{*2}}{4096 \text{ [LSB/frame]} \cdot 8 \text{ [frames/sec]} \cdot 3600 \text{ [sec/hour]} \cdot 8 \text{ [hour/day]}} = 92 \text{ [day]}$$

*1: Minimum X-ray sensitivity value. Refer to the individual datasheets for the values.

*2: Calculated assuming that the sensitivity decrease rate by X-ray irradiation is 0.9.

In a real operation, the dark and light field images may be taken frequently for offset and sensitivity calibration. The flat panel sensor may be exposed to X-ray before or after the valid image is taken. It is needed to count up all these shots for proper estimation. The upper limit of X-ray tube voltage and the life time of the sensor are determined by scintillator type, window material and the structure of sensor ship. Refer to the datasheet for each type of the flat panel sensor and application manual for more details.

Not recommended for use in inline applications where large radiation damage may occur.

5. Features


5-1. Spatial Resolution

The flat panel sensor spatial resolution is determined by the frequency function of photodiode array and the applied scintillator layer. The contrast transfer function (CTF) curve is an index of the total resolution of the flat panel sensor (See Figure 12). It appears from 0 line pair/mm to the Nyquist frequency (Fn). The frequency is defined by the pixel pitch.

$$F_n = \frac{1}{2 \cdot d} \dots\dots\dots (3)$$

d: Pixel pitch [mm]

Figure 12: Resolution




5-2. Dynamic Range

An advanced design is essential for attaining a wide dynamic range with the 12-bit A/D converter which provides 4096 levels of gray or the 14-bit A/D converter which provides 16384 levels of gray. The value of dynamic range shown in each sensor type specification corresponds to the maximum bit depth of each type. Amplifier noise is minimized by lowering noise voltages in initial stage of the amplifiers. A wide dynamic range is obtained by keeping the swing width of the amplifier within the range of the 5 V supply used in the standard CMOS process and by optimizing the saturation voltage for the photodiode array. However, improving the gain in the initial stage of amplifier for the tiny pixels requires reducing the charge amplifier feedback capacitance (Cf). The capacitance of a column gate switch drain is also added on the data line capacitance, along a large active area making up large input capacitance (Ct).

X-ray image sensors usually require large active area. This causes the Ct value to increase significantly so that a large format active area and very fine pixels are drawbacks in terms of noise. However, by optimizing the transconductance (gm) of the initial amplifier stage, and using patterning that miniaturizes the data lines and reduces stray capacitance, a 100 μm type flat panel sensor was developed having a wide dynamic ranges of 4300 as shown above.

Figure 13: Linearity



5-3. Detective Quantum Efficiency

Detective Quantum Efficiency (DQE) is one of the physical parameters used for the evaluation of image quality and it was derived from the idea of S/N. The DQE is defined as shown in equation (4).

$$DQE = (S / N_{out})^2 / (S / N_{in})^2 \dots\dots\dots (4)$$

Generally a detector is never ideal because the information fed into the detector is degraded due to the incomplete absorption efficiency, noise and contrast loss in the detection system. The DQE analysis is important to determine the effectiveness of the imaging system. The DQE essentially involves the conversion gain of the system, resolution characteristics (MTF: Modulation Transfer Function) and noise characteristics (NPS: Noise Power Spectrum). Since the MTF and NPS are the function of the spatial frequency, DQE is also represented as a function of the spatial frequency. The DQE(f) is calculated by the equation below,

$$DQE(f) = G^2 \cdot q \cdot MTF(f)^2 / NPS(f) \dots\dots\dots (5)$$

where G is a constant which is defined as the average conversion gain [mV/quantum] of the system, q is the number of X-ray quanta incident per pixel upon the detector during one acquisition period, MTF(f) is the Modulation Transfer Function which is another representation of the spatial resolution obtained from the Line Spread Function (LSF) by applying the Fourier transform method, and NPS(f) is the noise power spectrum obtained from the image data by applying the Fourier transform method.


6. System Set-up

6-1. Set-up Example

Connect the flat panel sensor to the grabber board that was installed in a PC according to the manufacturer's instructions. A highly flexible grabber board can control binning bit, IntExt, ExtTrg as binning or triggers operation though its digital I/O lines. A power cable* which is terminated with the applicable plug at one end and open at another end is shipped as an accessory with the flat panel sensor.

* Power cable is not supplied with some products. Please refer to the individual datasheets.

Figure 14: Connection example



KACCC0269EB

6-2. Precaution for Power Source Configuration


In case of using a power source that has large ripple noise, horizontal line noise may appear in an image. A power source whose ripple noise is less than 1 mVp-p should be used.

For the high quality image acquisition, analog +5 V and digital +5 V sources should be separated, because all power source ground are connected to the frame grounds at the flat panel sensor side, they should be disconnected from the frame ground or another power source ground at the power source side (refer to Figure 15).

The voltages described above are specified at the flat panel sensor side. The impedance of the power cable attached with the flat panel sensor is low enough, but it causes 0.1 V drop. Therefore the voltage at the power source side should be set 0.1 V higher than the voltage specified above. Avoid to extend the power cable because it causes voltage drop.

An AC line filter should be placed at the AC input line in the power source to prevent the surge from the AC line. A power source with CE mark is recommended.

Figure 15: Example for power cable connection



KACCC0312EA

6-3. Frame grabber cables

HAMAMATSU provides the A8406 series as frame grabber board cables. The A8406 series is available in several types with different cable lengths and terminations. Check the cable specifications that meet your needs. For detailed information, refer to the A8406 series datasheet.

Power Supply Unit

Please consult with our sales office.

7. Image Samples

Dog paw

Taken with radiology type flat panel sensor

Photodiode area: 124.8 · 124.8 mm

Pixel size: 200 μm

Scintillator: direct deposition CsI

Flame rate: 0.5 frames/s

X-ray source: HAMAMATSU L6622-02

(Focus size: 40 μm , window Be: 200 μm) 50 kV, 300 μA



Hornet

Taken with non-destructive inspection type flat panel sensor

Photodiode area: 52.8 · 52.8 mm

Pixel size: 50 µm

Scintillator: CsI flipped scintillator plate

Frame rate: 0.35 frames/s

X-ray source: RIGAKU 150CPT

(Focus size: 300 µm, window Be: 1 mm) 35 kV, 1.7 mA



Fish

Taken with radiology type flat panel sensor

Photodiode area: 120 · 120 mm

Pixel size: 50 μm

Scintillator: direct deposition CsI

Frame rate: 0.1 frames/s

X-ray source: HAMAMATSU L6622-02

(Focus size: 10 μm , window Be: 200 μm) 35 kV, 200 μA



Chicken wing

Taken with radiology type flat panel sensor

Photodiode area: 124.8 · 124.8 mm

Pixel size: 200 μm

Scintillator: direct deposition CsI

Frame rate: 0.1 frames/s

X-ray source: HAMAMATSU L6622-02

(Focus size: 10 μm , window Be: 200 μm) 30 kV, 150 μA



Leaf

Taken with radiology type flat panel sensor

Photodiode area: 120 · 120 mm


Pixel size: 50 μ m

Scintillator: direct deposition CsI

Frame rate: 0.1 frames/s

X-ray source: HAMAMATSU L6622-02

(Focus size: 10 μ m, window Be: 200 μ m) 30 kV, 150 μ A



Cell phone

Taken with non-destructive inspection type flat panel sensor

Photodiode area: 120 · 120 mm


Pixel size: 50 μm

Scintillator: CsI flipped scintillator plate

Frame rate: 0.5 frames/s

X-ray source: HAMAMATSU L7901

(Focus size: 7 μm , window Be: 150 μm) 90 kV, 100 μA



8. Related Products

HAMAMATSU has a full line of X-ray sources and cameras for X-ray inspection of electronic components, industrial products as well as for a wide range of fields to meet needs in food processing and security, etc. For detailed information and PDF catalogs on these products, please access our homepage at www.hamamatsu.com

Microfocus X-ray source
(sealed type, 100 kV)



Microfocus X-ray source
(sealed type, 150 kV)



X-ray I. I. digital camera unit



Berillium window
X-ray I. I. camera unit




Compact X-ray CCD camera
"X-CUBE™"



X-ray line scan camera



Dual-energy
X-ray line scan camera



X-ray TDI camera



High resolution
X-ray imaging system



Scintillator plate
"FOS, ACS, ALS"



Energy differentiating
64 ch CdTe radiation line sensor



Line-up of Opto-semiconductors

Photosensor

 <p>Si photodiodes, APD, MPPC</p>	 <p>Photo IC and Related Products</p>	 <p>Image Sensors</p>
 <p>PSD (Position Sensitive Detectors)</p>	 <p>Infrared Detectors</p>	 <p>X-ray Flat Panel Sensors</p>

LED



Devices for Optical Communication



Devices for Car Applications



Mini-spectrometers



Related Products



We welcome your access to our web site



www.hamamatsu.com

Please access our web site to check various information about our latest product catalogues, news, technology introduction and corporate outline.
Some of the new/developmental products in this catalogue may not be available on our web site. Please consult your local sales office for more information.

Notice

- Copies of the full warranty can be obtained prior to the purchase of products by contacting your local HAMAMATSU sales office.

HAMAMATSU makes no other warranties, and any and all implied warranties of merchantability, or fitness for a particular purpose, are hereby disclaimed. The customer is responsible for use of the product in accordance with HAMAMATSU's instructions and within the operating specifications and ratings listed in this catalogue. HAMAMATSU shall not be responsible for the customer's improper selection of a product for a particular application or otherwise. No warranty will apply if the products are in any way altered or modified after delivery by HAMAMATSU or for any intentional misuse or abuse of the products. Proper design safety rules should be followed when incorporating these products into devices that could potentially cause bodily injury.

HAMAMATSU's liability on any claim for loss or damage arising out of the supplying of any products, whether based on contract, warranty, tort (including negligence and for property damage or death and bodily injury) or other grounds, shall not in any event exceed the price allocable to such products or a part thereof involved in the claim, regardless of cause or fault. In no event shall HAMAMATSU be responsible to the customer or any third party for any consequential, incidental or indirect damages, including but not limited to loss of profits, revenues, sales, data, business, goodwill or use, even if the company has been advised of the possibility of such loss or damage. The limitation of liability set forth herein applies both to products and services purchased or otherwise provided hereunder. This warranty is limited to repair or replacement, at the sole option of HAMAMATSU, of any product which is defective in workmanship or materials used in manufacture. All warranty claims must be made within 1 year from the date of purchase or provision of the products or services.

Products that are amenable to repair shall be done so either under warranty or pursuant to a separate repair agreement. Some products cannot be repaired either because of the nature or age of the product, the unavailability of spare parts, or the extent of the damage is too great. Please contact your local HAMAMATSU office for more details.

- The products described in this catalogue should be used by persons who are accustomed to the properties of photoelectronics devices, and have expertise in handling and operating them. They should not be used by persons who are not experienced or trained in the necessary precautions surrounding their use.
- The information in this catalogue is subject to change without prior notice
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
- No patent rights are granted to any of the circuits described herein.

HAMAMATSU

HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1, Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558, Japan
Telephone: (81)53-434-3311, Fax: (81)53-434-5184

www.hamamatsu.com

Main Products

Si photodiodes
APD
MPPC
Photo IC
Image sensors
X-ray flat panel sensors
PSD
Infrared detectors
LED
Optical communication devices
Automotive devices
Mini-spectrometers
High energy particle/X-ray detectors
Opto-semiconductor modules

Hamamatsu also supplies:

Photoelectric tubes
Imaging tubes
Light sources
Imaging and processing systems

Information in this catalogue is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein.

© 2010 Hamamatsu Photonics K.K.

Quality, technology, and service are part of every product.

Sales Offices

JAPAN:
HAMAMATSU PHOTONICS K.K.
325-6, Sunayama-cho, Naka-ku
Hamamatsu City, 430-8587, Japan
Telephone: (81)53-452-2141, Fax: (81)53-456-7889
E-mail: intl-div@hq.hp.k.co.jp

U.S.A.:
HAMAMATSU CORPORATION
Main Office
360 Foothill Road, P.O. BOX 6910,
Bridgewater, N.J. 08807-0910, U.S.A.
Telephone: (1)908-231-0960, Fax: (1)908-231-1218
E-mail: usa@hamamatsu.com

Western U.S.A. Office:
Suite 200, 2875 Moorpark Avenue
San Jose, CA 95128, U.S.A.
Telephone: (1)408-261-2022, Fax: (1)408-261-2522
E-mail: usa@hamamatsu.com

United Kingdom, South Africa:
HAMAMATSU PHOTONICS UK LIMITED
Main Office
2 Howard Court, 10 Tewin Road, Welwyn Garden City,
Hertfordshire AL7 1BW, United Kingdom
Telephone: (44)1707-294888, Fax: (44)1707-325777
E-mail: info@hamamatsu.co.uk

South Africa office:
PO Box 1112
Buccleuch 2066
Johannesburg, South Africa
Telephone/Fax: (27)11-802-5505

France, Portugal, Belgium, Switzerland, Spain:
HAMAMATSU PHOTONICS FRANCE S.A.R.L.
19, Rue du Saule Trapu, Parc du Moulin de Massy,
91882 Massy Cedex, France
Telephone: (33)1 69 53 71 00
Fax: (33)1 69 53 71 10
E-mail: infos@hamamatsu.fr

Swiss Office:
Dornacherplatz 7
4500 Solothurn, Switzerland
Telephone: (41)32/625 60 60,
Fax: (41)32/625 60 61
E-mail: swiss@hamamatsu.ch

Belgian Office:
Scientific Park,
7, Rue du Bosquet
B-1348 Louvain-La-Neuve, Belgium
Telephone: (32)10 45 63 34
Fax: (32)10 45 63 67
E-mail: epirson@hamamatsu.com

Spanish Office:
C. Argenters, 4 edif 2
Parque Tecnológico del Valles
E-08290 CERDANYOLA, (Barcelona) Spain
Telephone: (34)93 582 44 30
Fax: (34)93 582 44 31
E-mail: infospain@hamamatsu.es

Germany, Denmark, Netherlands, Poland:
HAMAMATSU PHOTONICS DEUTSCHLAND GmbH
Arzbergerstr. 10,
D-82211 Herrsching am Ammersee, Germany
Telephone: (49)8152-375-0, Fax: (49)8152-265-8
E-mail: info@hamamatsu.de

Danish Office:
Lautruphøj 1-3
DK-2750 Ballerup, Denmark
Telephone: (45)70 20 93 69, Fax: (45)44 20 99 10
E-mail: sfeldskov@hamamatsu.de

Netherlands Office:
PO BOX 50.075, 1305 AB ALMERE, The Netherlands
Telephone: (31)36-5382123, Fax: (31)36-5382124
E-mail: info@hamamatsu.nl

Poland Office:
02-525 Warsaw,
8 St. A. Boboli Str., Poland
Telephone: (48)22-646-0016, Fax: (48)22-646-0018
E-mail: jbaszak@hamamatsu.se

North Europe and CIS:
HAMAMATSU PHOTONICS NORDEN AB
Smidesvägen 12
SE-171 41 Solna, Sweden
Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01
E-mail: info@hamamatsu.se

Russian Office:
Vyatskaya St. 27, bld. 15
Kosmodamianskaya nab. 52/1, 14th floor
RU-127015 Moscow, Russia
Telephone: (7) 495 258 85 18, Fax: (7) 495 258 85 19
E-mail: info@hamamatsu.ru

Italy:
HAMAMATSU PHOTONICS ITALIA S.R.L.
Strada della Moia, 1/E
20020 Arese, (Milano), Italy
Telephone: (39)02-935 81 733
Fax: (39)02-935 81 741
E-mail: info@hamamatsu.it

Rome Office:
Viale Cesare Pavese, 435
00144 Roma, Italy
Telephone: (39)06-50513454, Fax: (39)06-50513460
E-mail: inforoma@hamamatsu.it

Hong Kong:
HAKUTO ENTERPRISES LTD.
Units 709-712, 7/F, Miramar Tower,
132 Nathan Road, Tsimshatsui, Kowloon, Hong Kong
Telephone: (852)2578-4921, Fax: (852)2807-2498

Taiwan:
HAKUTO TAIWAN LTD.
6F, No.308, Pa teh Road, Sec. 2,
Taipei, Taiwan R.O.C.
Telephone: (886)2-8772-8910
Fax: (886)2-8772-8918

KORYO ELECTRONICS CO., LTD.
9F-7, No.79, Hsin Tai Wu Road
Sec.1, Hsi-Chih, Taipei, Taiwan, R.O.C.
Telephone: (886)2-2698-1143, Fax: (886)2-2698-1147

Republic of Korea:
SANGKI TRADING CO., LTD.
Suite 431, World Vision BLDG.
24-2 Yoido-Dong
Youngdeungpo-Ku
Seoul, 150-877
Telephone: (82)2-780-8515
Fax: (82)2-784-6062

Singapore:
HAKUTO SINGAPORE PTE LTD.
Block 2, Kaki Bukit Avenue 1, #04-01 to #04-04
Kaki Bukit Industrial Estate, Singapore 417938
Telephone: (65)67458910, Fax: (65)67418200



Cat. No. KOTH0008E08
Aug. 2010 DN
Printed in Japan (2,000)