

INSTRUCTIONS

MP-65180(SCL)

Scaler

No. IMP65180-1

Notice

- The information in this manual, which is based on specifications believed correct at the time of publications, is subject to change without notice due to improvements made in the instrument.
- In order to assist us in preparing future documentation, please advise your nearest JEOL service office if you find any errors in this manual.
- This manual has been prepared using a desktop publishing system in which photographs and drawings are read by the use of an image scanner.
Therefore, details of some of the photos may not be as clear as those of the originals.
- In no event will JEOL Ltd. be liable for, indirect, incidental or consequential damages resulting from the use of software of this instrument must not be modified without leave.
JEOL Ltd. assumes no responsibility for effect of such modifications.
- This instrument must not be modified, and products other than those manufactured by JEOL Ltd. must not be attached to this instrument, without prior written permission.
If any such modification or attachment is made, all the stipulated warranties and services contacted by JEOL Ltd. or its affiliated company will be void.
- This manual is copyrighted. All rights are reserved. This document may not, in whole or part, be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine readable form without prior consent, in writing, from JEOL Ltd.
- When this manual or the software described in this manual is furnished under a license agreement, it may only be used or copied in accordance with the terms of such license agreement.

Copyright ©1997 JEOL Ltd.

SELLING

JEOL Ltd.

1-2 Musashino 3-chome, Akishima-Shi, Tokyo, 196-0021 Japan

Telephone: 81-42-542-2187

Facsimile: 81-42-546-5757

MANUFACTURE

JEOL technics Ltd.

6-38 Musashino 2-Chome, Akishima-Shi, Tokyo, 196-0021 Japan

For servicing or inquires, contact your local service office.

Contents

1. General	1
2. Specifications	1
3. Composition	1
4. Functions of the Scaler and Multi Point Measurement Menus	2
4.1 Scaler Menu	2
4.2 Multi Point Measurement Menu	3
5. Measurement	4
5.1 X-Mode	4
5.2 Y-Mode	5
5.3 R-Mode	6
6. Multi Point Measurement	7
6.1 Line Mode	7
6.2 Circle Mode	8

Please note that the drawings given in this instruction manual are somewhat different from actual components.

1. General

The scaler is capable of setting markers at any positions on an image, and measuring and displaying the length between two markers or the diameter of a circle. Linewidths (in the X-, Y-, and diagonal directions), distances between any two points, and diameters of circles can be measured.

2. Specifications

Linewidth Measurement


- Measured items: X(horizontal)•Y(vertical)•R(diagonal)
- Measuring range: The whole magnification range
- Measurement accuracy: Equivalent to the magnification accuracy of SEM.
- Operation: Carried out by using the Measurement menu opened by GUI.
- Display of measured values: Digitally displays X, Y, and R in three digits on the monitor.
- Recording of measured values: The frozen image on CRT can be recorded by taking a photograph or printing on a commercially available printer (recommended one by JEOL).

Multi Point Measurement

- Measured items: Distances between any two points (Line), diameters of circles (Circle)
- Measuring range: The whole magnification range
- Measurement accuracy: Equivalent to the magnification accuracy of SEM.
- Operation: Carried out by using the Measurement menu opened by GUI.
- Display of measured values: Digitally displays Line/Circle in three digits on the monitor.
- Recording of measured values: The frozen image on CRT can be recorded by taking a photograph or printing on a commercially available printer (recommended one by JEOL).

3. Composition


- Set-up floppy disk



- Measuring modes:
- X-mode: Measures a distance in the horizontal direction.
(Two vertical line-cursors are displayed.)
 - Y-mode: Measures a distance in the vertical direction.
(Two horizontal line-cursors are displayed.)
 - R-mode: Measures a distance in the diagonal direction.
(Two pairs of vertical and horizontal line-cursors are displayed.)

Write: By clicking [Write], the cursors and measured value(s) are pasted on the image and the Measurement menu is closed.

Cancel: By clicking [Cancel], the state before activating this function returns.



- Measuring modes:
- Line mode: Measures the length of any line (the distance between two points).
 - Circle mode: Measures the diameter of any circle.
- Up to 10 measurements in total can be carried out in the line and circle modes.
- Write: By clicking [Write], the line or circle and the measured value are pasted on the image and the Multi Point Measurement menu is closed.
- Cancel: By clicking [Cancel], the state before activating this function returns.

5. Measurement

5.1 X-Mode

- 1 Freeze the image or open a BMP file.

To open a BMP file, select [Open...] of [File] in the menu bar.

- 2 Click [Image] in the menu bar and select [Scaler...] from the pull-down menu.

The Measurement menu opens.

- 3 Click [X] and drag two cursors to enclose the line to be measured.



The measured value is displayed in real time at the right top of the screen.

- 4 Click [Write].

The cursors and measured value are pasted on the original image and the Measurement menu is closed. The image on which the cursors and measured value are pasted can be

photographed, printed, and saved. (Refer to "Chapter 3" of the main instruction manual.)

By clicking [Cancel], the display of the original image (the one before activating this function) returns.



- 1 Freeze the image or open a BMP file.

To open a BMP file, select [Open...] of [File] in the menu bar.

- 2 Click [Image] in the menu bar and select [Scaler...] from the pull-down menu.

The Measurement menu opens.


- 3 Click [Y] and drag two cursors to enclose the line to be measured.

The measured value is displayed in real time at the right top of the screen.

- 4 Click [Write].

The cursors and measured value are pasted on the original image and the Measurement menu is closed. The image on which the cursors and measured value are pasted can be photographed, printed, and saved. (Refer to "Chapter 3" of the main instruction manual.)

By clicking [Cancel], the display of the original image (the one before activating this function) returns.



- 1 Freeze the image or open a BMP file.

To open a BMP file, select [Open...] of [File] in the menu bar.

- 2 Click [Image] in the menu bar and select [Scaler...] from the pull-down menu.

The Measurement menu opens.

- 3 Click [R] and drag the cursors to enclose the line to be measured.


A pair of vertical and horizontal line-cursors can be moved together by grabbing at their cross point. The measured value is displayed in real time at the right top of the screen.

(The X- and Y- values are also displayed simultaneously.)

- 4 Click [Write].

The cursors and measured values are pasted on the image and the Measurement menu is closed. The image on which the cursors and measured values are pasted can be photographed, printed, and saved. (Refer to "Chapter 3" of the main instruction manual.)

By clicking [Cancel], the display of the original image (the one before activating this function) returns.



6. Multi Point Measurement

6.1 Line Mode

- 1 Freeze the image or open a BMP file.

To open a BMP file, select [Open...] of [File] in the menu bar.

- 2 Click [Image] in the menu bar and select [Multi Point Measurement...] from the pull-down menu.
The Multi Point Measurement menu opens.


- 3 Click [Line] and drag with the mouse.

The position where the left mouse button is pressed becomes the start point and the position where the mouse button is released the end point. When the mouse button is released, the measured value (the distance between the two points) is displayed.

- 4 Click [Write].

The line and the measured value are pasted on the image and the Multi Point Measurement menu is closed. The image on which the line and the measured value are pasted can be photographed, printed, and saved. (Refer to "Chapter 3" of the main instruction manual.)

By clicking [Cancel], the display of the original image (the one before activating this function) returns.



The measured value or length of the displayed line can be changed by pointing the cursor to the measured value and dragging the mouse with the left button held.
The line itself can be deleted by clicking the right mouse button while

6. Multi Point Measurement

6.2 Circle Mode

- 1 Freeze the image or open a BMP file.

To open a BMP file, select [Open...] of [File] in the menu bar.

- 2 Click [Image] in the menu bar and select [Multi Point Measurement...] from the pull-down menu.

The Multi Point Measurement menu opens.

- 3 [Click [Circle] and drag with the mouse.


A circle having the diameter corresponding to the distance you drag the mouse is displayed.

When the mouse button is released, the measured value (the diameter of the circle) is displayed.

- 4 Click [Write].

The cursors and the measured value are pasted on the image and the Multi Point Measurement menu is closed. The image on which the circle and the measured value are pasted can be photographed, printed, and saved. (Refer to "Chapter" 3 of the main instruction manual.)

By clicking [Cancel], the display of the original image (the one before activating this function) returns.



The measured value or diameter of the displayed circle can be changed by pointing the cursor to the measured value and dragging the mouse with the left button held.

The measured value or displayed circle itself can be deleted by clicking the right mouse button while pointing the cursor to the measured value or the circle itself.