



MKS RGA Communications Library

**Software Development Kit
For Application Programmers**

Contents

Introduction.....	5
Installing The CD-ROM Software.....	7
Hardware Management	9
<i>Hardware Programmable Features of RGA Control Units</i>	<i>10</i>
<i>The MKS RGA Device Manager Software for Microsoft Windows® Computers</i>	<i>15</i>
<i>Working with a Firewall.....</i>	<i>17</i>
A 'Quick Start' Introduction to The RGA Communications Library	18
<i>The Significant Events</i>	<i>19</i>
<i>Significant Properties and Methods.....</i>	<i>22</i>
RGAComms ActiveX Object Model Reference.....	23
<i>Object Hierarchy</i>	<i>23</i>
<i>RGAConnection</i>	<i>24</i>
<i>RGASensor.....</i>	<i>28</i>
<i>ProtocolInfo</i>	<i>31</i>
<i>Filaments.....</i>	<i>31</i>
<i>Measurements.....</i>	<i>32</i>
<i>Measurement.....</i>	<i>36</i>
<i>BarchartMeasurement.....</i>	<i>37</i>
<i>AnalogMeasurement.....</i>	<i>37</i>
<i>SinglePeakMeasurement.....</i>	<i>38</i>
<i>PeakJumpMeasurement.....</i>	<i>39</i>
<i>PeakJumpPeak.....</i>	<i>40</i>
<i>PeakJumpPeaks.....</i>	<i>40</i>
<i>Scan.....</i>	<i>41</i>
<i>ScanMeasurements.....</i>	<i>42</i>
<i>Readings.....</i>	<i>42</i>
<i>SourceSettings</i>	<i>43</i>
<i>SourceSettingsCollection</i>	<i>44</i>
<i>DetectorSettings.....</i>	<i>44</i>
<i>Inlet.....</i>	<i>45</i>
<i>Inlets.....</i>	<i>45</i>

<i>AnalogInput</i>	46
<i>AnalogInputs</i>	47
<i>AnalogOutput</i>	47
<i>AnalogOutputs</i>	47
<i>Cirrus</i>	48
<i>RVC</i>	48
<i>Rollover</i>	49
<i>RF</i>	50
<i>Multiplier</i>	50
<i>DigitalPort</i>	51
<i>DigitalPorts</i>	51
<i>TotalPressure</i>	51
<i>Enum rgaFilamentTrips</i>	52
<i>Enum rgaFilamentSummary</i>	52
<i>Enum rgaFilamentXTripMode</i>	52
<i>Enum rgaMeasurementTypes</i>	53
<i>Enum rgaFilterModes</i>	53
<i>Enum rgaInfoAvailableStates</i>	53
<i>Enum rgaPressureUnits</i>	53
<i>Enum rgaLinkDownReasons</i>	53
<i>Enum rgaCirrusPumpStates</i>	54
<i>Enum rgaCirrusHeaterStates</i>	54
<i>Enum rgaRVCPumpStates</i>	54
<i>Enum rgaRVCHeaterStates</i>	54
<i>Enum rgaRVCValveModes</i>	54
The Java Libraries	55
Redistributing the MKS RGA Libraries	56
An Introduction to the Examples	57
The VB Example	58
The C++ Example	63
<i>Using the sample</i>	63
<i>The CRGACONNECTION class</i>	64
<i>CPPSample.cpp</i>	65
The Java Example	69

The LabView Example	72
Using Virtual RGA Systems	75
<i>What is a Virtual RGA System?</i>	<i>75</i>
<i>How are Virtual RGA's Installed?</i>	<i>75</i>
<i>How do Virtual RGA's Work?</i>	<i>76</i>
<i>Are there any limitations to Virtual RGA's?</i>	<i>76</i>
<i>Using a Virtual System</i>	<i>77</i>
Vacuum Profile 'Replay' File Format	96
Background Information	100

Chapter 1

Introduction

This manual describes the ***RGA Communications Library Software Development Kit (SDK)*** provided by MKS Instruments that enables the integration of MKS RGA Control Units into other software applications.

Who should read this manual?

If you are a software developer working in one of the popular programming languages such as Java, Visual Basic, Visual C++ or LabView, and you wish to control an MKS RGA control unit and use the data that it acquires directly in your own software, you should read this manual.

What does the manual include?

The manual includes

- A description of the RGA control unit hardware.
- What you need to know about installing and configuring the control unit hardware
- The principles of how to use the software tools provided by the SDK
- A full reference guide to the RGA Communications library object model
- A description of the example programs and how you can adapt them for your own use.

What does the manual not include?

The manual does not include any information about MKS' own software products, Process Eye[®] and EasyView[®], which also fulfil the function of controlling MKS mass spectrometers.

What is in the SDK for Windows[®] computers?

In addition to the Java and ActiveX RGA communications libraries themselves, the SDK includes the complete RGA Device Management software for installing and configuring both MKS network control units – e-Vision, e-Vision Plus, Microvision-IP and HPQ-IP – and serial RS232 control units.

Finally the SDK contains examples of the use of the libraries, which you can use as the basis for your own development.

What is in the SDK for non-Windows[®] computers?

If your operating system is not Microsoft Windows[®], for example Linux, Unix or Mac, you should copy the entire contents of the 'Linux' folder to your computer and read the documentation therein. You will find the Java JAR library file and the Java example files. MKS Instruments does not provide any hardware management software for non-Windows operating systems.

However e-Vision and e-Vision Plus control units require no external management because browsing to the control unit with your favourite browser software performs all user-configurable control unit options. All network instruments have browser interfaces and support the applet that can control the RGA, acquire data, tune the instrument and run diagnostic tests.


You cannot use the SDK on non-Windows computers to communicate with serial RS232 control units such as Microvision Plus and HPQ2. You can only use the SDK with Microvision-IP, HPQ2-IP, e-Vision and e-Vision Plus control units.

Chapter 2

Installing The CD-ROM Software

Installing on a Windows® Computer

Insert the CD-ROM in a suitable drive. It should auto-run, but if it does not, double-click 'setup.exe' in the root folder.



Read the options on the CD-browser window. If you do not already have Acrobat Reader software installed on your computer, please select the option to install this first, since all SDK documentation is in PDF format.

The default option is to install the SDK. Selecting this option will launch the install package, which offers a complete choice of the features to be installed. Several components can be installed to a folder location of your choice. If you are new to writing software for RGA's, be sure to include installation of the examples in your chosen language(s).

If, after installing the software, you want to install additional components that were not installed initially, you can do so by launching Add/Remove Programs from Control Panel.

Installing on a non-Windows® Computer

The software contained in the 'Linux' folder of the CD-ROM should be self-explanatory and simply requires copying to an appropriate location on an individual developer's computer.

Chapter 3

Hardware Management

The tools provided with this SDK can be used with any of the following control units

Micovision Plus and HPQ2 with serial RS232 connections.

You can work with these control units only if they are connected to Windows computers. These control units should have a minimum of version 5.0 embedded core and application software. You must access them on a computer that is running at least version 5.10 of the 'RGAServer.exe' application software. This application is the data server to which you connect. (RGAServer.exe can be installed as part of the SDK.)

Micovision-IP and HPQ-IP with network connections.

All computer operating systems can control both these control units provided they have version 5.0 embedded core and application software and version 5.00 of the internal data server software as a minimum.

e-Vision® and e-Vision Plus® control units with network connections.

You can work with all versions of these control units. Once again, the data server to which your software connects is embedded in the control unit.

The Network Control Units

All the control units listed above that have network connections provide a web interface accessible from any browser software. This web interface permits control of the RGA and provides a graphical view on the data being acquired in real time. Typically you would use this as a method for quickly determining the correct behaviour of the RGA and also for 'tuning' its mass alignment and the resolution of the individual peaks. The web interface provided by these control units is an example of a Java applet written by the same MKS software developer that has developed the software for this SDK.

The IP addressing options can also be set by browsing to the control unit. All the network instruments support DHCP, but it is likely that you will be assigning static IP addresses because the instruments will be on a small local network without a DHCP server. All network instruments are shipped with the default address of 192.168.0.250. The client computer must be on the same subnet as the control unit for it to serve up the pages correctly.

If it becomes necessary to update the embedded software in a network control unit, it is normally set into FTP mode and an FTP client application is used to transfer the new files. MKS often provides 'download packages' that contain several files and which perform the operation by means of a 'Wizard' – but these 'Package Downloads' are only available for Microsoft Windows® computers.

The Serial RS232 Control Units

You will be controlling these control units by connecting to the network address of the computer to which the control unit(s) are connected. This computer must be running Microsoft Windows® and at least version 5.10 of the RGA Server software. One instance of this server application can support multiple control units.

Obviously, none of the serial control units have a web interface. All aspects of hardware configuration and set-up are undertaken using the supplied Windows applications. Unless you wish to include mass alignment and tuning algorithms in your own software, you will probably need to have at least one instance of the full EasyView or Process Eye suite on a computer that can be used for diagnostic and RGA maintenance purposes. This software is available separately from the MKS Spectra RGA Group.

Software Protocol Version

All the different control units share the same software protocol versioning strategy. When you make a connection to any data server, you get two important pieces of information.

- The installed version
- The minimum compatible version.

You can only connect to a server if the protocol version of the libraries that are installed on the computer (the required version) is greater than or equal to the minimum compatible version reported by the sensor. This checks that the computer software is not too old for the control unit. It is up to the computer software itself to check that the control unit is not too old to work with the current version of the client application.

Hardware Programmable Features of RGA Control Units

You will almost certainly find it necessary to program the ion source settings of the analyser connected to the control unit that you are working with. The ion source has eight adjustable properties that affect the behaviour and signal response of the analyser – and hence its calibration. These properties are

Low mass alignment	0-65535
High mass alignment	0-65535
Low mass resolution	0-65535
High mass resolution	0-65535
Ion energy	0 to
Electron energy	0 to 130 volts
Extract potential	0 to 130 volts
Filament emission current	0 to 5 mA

The ion source settings determine the ion amps that arrive at the detector and, therefore, directly determine the magnitude of the peaks measured. All analysers have a Faraday detector – a simple detector that generates a current by collecting the ions that fall onto it. Some analysers have an additional ‘multiplier’ detector. A multiplier detector generates an output current that is a multiple of the incident current. The multiple is referred to as the Detector Gain. This gain is not a fixed number but is itself adjustable by increasing or decreasing the voltage applied to the detector. The voltage-gain response is non-linear. Typically the gain is about 50 at –700V and 1000 at –900V. Unfortunately multiplier detectors have a finite life in a direct relationship to the amount they are used. As time goes by increasingly negative voltages must be applied to achieve the same gain that they had when new.

In addition to any gain introduced by the detector (the Faraday detector, by definition, always has a gain of 1), all MKS control units feature at least two electronic gain ranges, the smallest of which also has, by definition, a value of 1.

The calibration of an analyser is essentially the transfer function that converts ion amps injected into the electronics by the detector to the partial or total pressure value that created the signal when the electronic gain is 1.

It should be apparent that with all these adjustable parameters, it is necessary to create some ordered arrangement whereby one of a number of commonly used settings can be called up at a particular time. E-Vision control units store one set of ion source settings. All other control units store a table of six. Each of these source settings has one (Faraday only detector) or four (dual Faraday and multiplier detector) detector settings. Each detector setting has two calibration factors – one for each of the two filaments.

The default settings for all the control units are shown in Table 1 overleaf.

Please note the following important points that affect you as a developer. All the ion source settings are positive numbers even though some of them in hardware terms are negative voltages. However, the multiplier detector voltage is negative both in reality and in software.

The only fixed item in the ion source settings table, which in software is an array (0 to 5) of six items, is the **Name**. All other values are completely flexible within the allowable range of the property concerned. So you can change the electron energy of the “Standard electron energy” row to be any number you like, but it obviously makes sense that if you want to work at a low electron energy, you should be adjusting the second row in the array in order to keep the meaning correct.

The **ion energy** value can affect the sensitivity and peak shape. If you are particularly interested in low masses you may benefit from using a slightly higher value. Conversely high masses may benefit from a lower value. The sensitivity of the analyser will increase with increased **emission current** – but only up to a point. Beyond a certain emission current the signal will probably fall. Lower **electron energies** may start to result in lower signals generally, but the main advantage is that they will usually result in much lower ‘doubly ionised’ fragments. The commonest occasion that this is apparent is when Argon is present. At 70eV electron energy a significant portion of the argon will appear doubly ionised at an m/e value of 20, whereas at 40eV all the argon will appear at an m/e value of 40. The **extractor potential** on an open source analyser should always be numerically greater than the electron energy. Beyond that, its impact is secondary compared with the other three properties.

The table of six source settings (only one in an e-Vision) is stored in non-volatile RAM within the control unit. Therefore it persists as control units are moved around or de-powered. The 48 calibration factors are stored at the data server, which is effectively the control unit itself on e-Vision and –IP control units. However RS232 serial instruments have their data server on the computer to which they are connected. This means that when you move an RS232 control unit from one computer to another, or perform a new computer installation, the calibration is ‘lost’. However, the Windows applications that are associated with the Windows server provide ways of backing up the settings from one PC and restoring the same settings file to another PC.

Ion Source Settings Table for Standard and High Resolution Open Source Analysers

ID	Name	Electron Energy	Ion Energy	Extract Potential	Emission Current	Low mass alignment	High mass alignment	Low mass resolution	High mass resolution
0	Standard Electron Energy	70	5.5	112	1.0	32767	32767	32767	32767
1	Low electron energy	40	5.5	112	1.0	32767	32767	32767	32767
2	User defined 1	70	5.5	112	1.0	32767	32767	32767	32767
3	User defined 2	70	5.5	112	1.0	32767	32767	32767	32767
4	User defined 3	70	5.5	112	1.0	32767	32767	32767	32767
5	EasyView	70	5.5	112	1.0	32767	32767	32767	32767

Ion Source Settings Table for Standard and High Resolution Enclosed Source Analysers

ID	Name	Electron Energy	Ion Energy	Extract Potential	Emission Current	Low mass alignment	High mass alignment	Low mass resolution	High mass resolution
0	Standard Electron Energy	40	7	20	1.0	32767	32767	32767	32767
1	High electron energy	70	5	20	1.0	32767	32767	32767	32767
2	User defined 1	40	7	20	1.0	32767	32767	32767	32767
3	User defined 2	40	7	20	1.0	32767	32767	32767	32767
4	User defined 3	40	7	20	1.0	32767	32767	32767	32767
5	EasyView	40	7	20	1.0	32767	32767	32767	32767

Ion Source Settings Table for HPQ Analysers

ID	Name	Electron Energy	Ion Energy	Extract Potential	Emission Current	Low mass alignment	High mass alignment	Low mass resolution	High mass resolution
0	Base pressure RGA	70	7	110	0.39	32767	32767	32767	32767
1	Mid pressure RGA	40	7	61	0.7	32767	32767	32767	32767
2	Leak hunt (helium0)	85	10	130	1.0	32767	32767	32767	32767
3	HP Analytical (HPQ2S)	35	5	58	0.1	32767	32767	32767	32767
4	HP General (HPQ2S)	40	8	61	0.16	32767	32767	32767	32767
5	EasyView	70	7	110	0.39	32767	32767	32767	32767

Ion Source Settings Table for Analysers with e-Vision Control Units

ID	Name	Electron Energy	Ion Energy	Extract Potential	Emission Current	Low mass alignment	High mass alignment	Low mass resolution	High mass resolution
0	Standard Electron Energy	40 / 70	7.5	125	1.0	32767	32767	32767	32767

Detector Settings Table for Standard and High Resolution Open Source Analysers and e-Vision Control Units

ID	Name	Factor
0	Faraday	1.5e-6
1	Multiplier1	4.5e-5
2	Multiplier2	4.5e-4
3	Multiplier3	4.5e-3

Each Source setting table entry has two of these tables, one for each filament. By default the factors for all rows of the source settings table are the same.

Detector Settings Table for Standard and High Resolution Enclosed Source Analysers

ID	Name	Factor
0	Faraday	1.5e-8
1	Multiplier1	4.5e-7
2	Multiplier2	4.5e-6
3	Multiplier3	4.5e-5

Each Source setting table entry has two of these tables, one for each filament. By default the factors for all rows of the source settings table are the same.

Detector Settings for HPQ analysers on Source Settings ID=0 and ID=1 and ID=5

ID	Name	Factor
0	Faraday	4.5e-7

Detector Settings for HPQ analysers on Source Settings ID=2 and ID=4


ID	Name	Factor
0	Faraday	2.25e-7

Detector Settings for HPQ analysers on Source Settings ID=3

ID	Name	Factor
0	Faraday	7.5e-8

The MKS RGA Device Manager Software for Microsoft Windows® Computers

The '**RGA Device Manager**' application, which is included as part of the SDK on Windows computers, provides complete installation and configuration capabilities for the entire range of control units, whether connected to this computer or on the network. Specifically this includes serial RS232 instruments connected to this computer or a different computer on the network, as well as all network instruments on the local subnet. The main screen of the application is shown below.



All the control units referred to above can be 'auto-discovered' by the '**RGA Device Manager**' application. If you are installing control units for the first time, switch them on and, if necessary, start the RGA Server application on any remote computer that has serial instruments. The local instance of the RGA Server will start up automatically as required. When you launch the application it will 'discover' all these control units. Network instruments are delivered with an IP address of 192.168.0.250. Even if this does not match the subnet of the computer, **RGA Device Manager** will still find the control unit, but the first thing you will have to do is to use the application to re-program the IP address so that it matches your subnet. You will only be able to connect to instruments using the SDK software libraries if the control unit is on your subnet. For example, if your computer address is 192.168.1.100, **RGA Device Manager** will find it and alert you to the fact that the address of the control unit should be changed to a 192.168.1.xxx address.


The **RGA Device Manager** application has full on-line help by pressing the F1 key or clicking the '**Help**' button.

Depending on the particular control unit, you can use **RGA Device Manager** to perform the following tasks.

- Change the 'friendly name' of a control unit
- Change an IP address
- Troubleshoot an RGA to which you cannot connect
- Modify the ion source settings, detector settings and calibration directly
- Download embedded software – core, application or server
- Download 'packages' of upgrade software supplied by MKS from time to time

- Configure a control unit's digital or analog inputs
- Change the trip settings (external trip, emission trip, RF trip etc)
- Configure any additional inlets

If the control unit of interest is not shown with a green tick, double click the item in the list. This will display more detailed information. If the unit requires some kind of configuration, click the 'Configure Now' button.



You can right mouse click on an item in the main window and a context-dependent menu will offer you the relevant options for that control unit.

Passwords

You will need a password to make changes to the settings of a control unit. For Microvisions and HPQ instruments the password is **mksrga1** (case-insensitive). For e-Vision and IP control units the browser page password is **MKS** (all upper case). The web browser applet password for changing the tuning settings is **peaks** (all lower case).

What to do if your unit is not discovered

If the **RGA Device Manager** does not discover your control unit try the following.

- Check the power cables and serial or network connections.
- Read the auto-discovery log on the Details button. This will give you information about the serial ports opened and all messages received from on-line servers.
- Press the control unit Reset button. If it was in ftp mode, this will ensure it reboots in normal mode.
- Click the 'Troubleshoot ...' button to run the troubleshooting wizard. This can diagnose the most common reasons why a particular control unit does not

- appear in the list. These include the action of the Windows XP firewall in blocking communications between networked items.
- Make sure that you can 'ping' the IP address, using the ping command in a DOS window.

Working with a Firewall

If your computer is protected by a firewall, you will need to program the firewall to allow access to certain applications and/or ports. Windows XP service pack 2 introduced a firewall to Windows by default. The table below shows the port settings that are required in order to communicate with the RGA hardware.

If you have the Windows XP firewall, these settings are automatically installed when you install the software on your computer. The installer sets the scope to 'unrestricted' by default. If you use a different firewall, you will need to unblock these ports.

If you are using RS232 serial instruments

Name	Port	Protocol	Scope (where restricted)
RGA Server UDP	10013	UDP	At least your local subnet, if you want other computers to have access to the control units connected to the server.
RGA Server 1	10014	TCP	Ditto
RGA Server 2	10015	TCP	Ditto

All instruments

Name	Port	Protocol	Scope (where restricted)
RGA Network Listener	10014	UDP	In addition to your local subnet you need to ensure that the 192.168.0.250 address is accessible in order to auto-detect new or un-configured instruments.

Network instruments only

In order to enable package downloads supplied by MKS from time to time as a means of upgrading control units, the [ftp.exe](#) application, which is located in the folder referenced by the environment variable 'comspec' needs to be unblocked.

A 'Quick Start' Introduction to The RGA Communications Library

The RGA Comms library features an extensive set of methods and properties that deal with all aspects of making a connection to a control unit, taking control, adding measurements, constructing a scan and acquiring data.

The object model is essentially the same whether you use the ActiveX DLL to program in, for example, Visual Basic or the Java libraries to target a non-Windows computer. In all these cases the software installation includes an Object Browser that you can use to explore all the features of the object model.

For this reason this chapter will give you a general introduction to RGA programming fundamentals, whatever your programming language.

Please note the following important point. Different programming languages have their own conventions for differentiating properties and methods. Java and LabView are rather strict – properties must return one of the basic data types. Anything that returns another object, for example, is a method. VB is more lax in that properties can return objects. This example shows the difference.

```
VB:  
MyConnection.SerialNumber(0)
```

This item is an indexed property

```
Java:  
MyConnection.getSerialNumber(0)
```

This is a method

```
LabView:  
MyConnection.SerialNumber(0)
```

This is a method

This manual adopts the VB approach since it is one of the commonest languages, but you need to remember that, when programming in your chosen language, some functions listed as properties may in fact be methods.

Having made that point clear, let us now take a 'Quick Start' tour of the main elements of the object model.

The two main high level objects are the **RGASensor** (which is accessed from the **RGAConnection** object) and the **Scan** (which is accessed from the **RGASensor** object).

The **RGAConnection** object is the server to which you connect. If it is a server embedded in a network instrument there is obviously only one **RGASensor** object associated with it and by default it is selected for use. However, a server running on a Windows computer could have multiple sensors connected to RS232 ports and therefore multiple **RGASensor** objects. Thus the sensor that a program wishes to

query or control must first be selected. Once connected to a server and having selected a sensor the **RGASensor** object has many properties and methods for controlling it.

The chances are that your program will want to use the RGA to acquire partial pressure data. There are a several alternative ways that data can be acquired by the RGA and the object model uses the concept of different measurement types to manage each of these acquisition methods. One or more individual measurements can be grouped together and run as a Scan. The Scan object allows an application to control which measurements have data acquired, in what order and when.

The four different measurement types are:

Barchart	Consecutive masses are measured from a starting mass through to an ending mass and only one reading is taken per mass.
Analog	Consecutive masses are measured starting from startmass-0.5 through to end mass+0.5. 32, 16, 8 or 4 readings can be returned per mass
Peakjump	Individual masses are grouped together and measured with 1 reading per mass being returned. The data acquisition is essentially the same as the barchart mode but arbitrary masses can be chosen without the linear sweep through the mass range.
Singlepeak	A single reading can be taken from any position in the mass range of the sensor usually within 1/32 AMU. This mode uses a default zeroing mode which must be re-triggered. Typically this mode is used for leak checking.

When creating measurements the methods return a reference to a measurement object that is not fully initialised. When the measurement is fully initialised (or failed to initialise) the MeasurementCreated event will fire. While it is possible to create as many measurements as are needed in quick succession it is not possible to add measurements to the Scan until they have been successfully created.

Once the measurements have been added to the scan the Scan's Start method initiates data acquisition; data can be polled by accessing the measurements **Data** property to obtain a reference to the **Readings** object but normally it is more useful to process the entire scans data in the **EndOfScan** event which is fired when all measurements in the scan have completed their data acquisition.

The coding requirements for this are explained in more detail below.

The Significant Events

The Connected Event

This event is fired after calling the `Connect` method on the Connection object. So to connect to the server at IP address 10.10.100.25, you would call the following method.

```
MyConnection.Connect "10.10.100.25"
```

You would then wait for the Connected event to fire

```
MyConnection_Connected (bSuccess As Boolean)
```

The `bSuccess` parameter indicates whether the connect attempt succeeded. If `bSuccess` is false, there was a problem connecting to the server at the address specified.

The Controlling Event

This event is fired after invoking the `Control` method on your selected sensor.

```
MyConnection.SelectedSensor.Control "Vbrga", "version 1.0"
```

Code execution would then continue with the Controlling event

```
MyConnection_Controlling (bSuccess As Boolean)
```

The `bSuccess` parameter indicates whether the control attempt succeeded. If `bSuccess` is False there could be several reasons the most likely of course being that the control unit is in use by another user. This would have been apparent before attempting the `Control` method by examining the `State` property of the `SelectedSensor`

```
If MyConnection.SelectedSensor.State = Available Then ...
```

The MeasurementCreated Event

This event is fired after invoking one of the methods that add a measurement to the control unit. These methods are

```
AddBarchart  
AddPeakJump  
AddAnalog  
AddSinglePeak
```

Code execution would then continue with the `MeasurementCreated` event. Note that it is not uncommon to add several measurements to the control unit even though you may not wish to use them all immediately. A typical code fragment might look like this

```
With MyConnection.SelectedSensor  
  .AddBarchart "Barchart", 1, 20, PeakCenter, 5, 1, 0, 0  
  With .AddPeakJump "PeakJump", PeakCenter, 5, 1, 0, 0  
    .Masses.Add 28  
    .Masses.Add 32  
    .Masses.Add 40  
  End With  
End With  
WaitForMeasurementCount = 2
```

The code execution would continue with the `MeasurementCreated` event

```
MyConnection_MeasurementCreated (bSuccess As Boolean, theMeasurement As  
Measurement)
```

Typically this code would check for the `bSuccess` parameter and decrement the global `WaitForMeasurementCount` variable. Complete success would require

`bSuccess` to be true at each firing of the event and when `WaitForMeasurementCount` reached zero the code would continue with the next task.

The StartScan and EndOfScan Events

These events are fired once the control unit is scanning and acquiring data. You start the scanning process by invoking the `Start` method on the `Scan` object. You stop scanning by calling the `Stop` method.

```
Scan.Start 2
```

At the end of each scan the `EndOfScan` event is fired and you can use this event to access the data for the current scan

```
MyConnection_EndOfScan (theScan As Scan)
```

The `theScan` reference to the `Scan` object provides access to many properties, but in particular the data acquired during the scan, as follows

```
Mass28Value = theScan.Measurements("Barchart").Data.Readings("Mass 28")
```

The FilamentChange Event

This event is the prime means of knowing when the filament state changes. It is not necessary to poll for the filament status. After initially discovering the filament status by reading the property value on the `Filament` object the `FilamentChange` event will inform you of any changes in state as they happen.

A code fragment for the `FilamentChange` event might look like

```
Sub MyConnection_FilamentChange (theFilaments As RGAFilaments)
    IsOn = theFilaments.filamentOn
    Xtripped = (theFilaments.SummaryState And Xtripped) <> 0
End Sub
```

The LinkDown Event

This event is the prime means of handling a broken connection. Some control units (even some network ones) have both a TCP/IP connection and a real or virtual serial connection to the source of the data. Either of these connections can fail and the result is effectively the same – there will be no more data! You do not need to take any action to stop scanning. It is important, however, that if you have kept any references to the `MyConnection.SelectedSensor` object that these be released since it is no longer available.

A code fragment for the `LinkDown` event might look like

```
Sub MyConnection_LinkDown (Reason as LinkDownReason)
    Set MySensorReference = Nothing
End Sub
```

Significant Properties and Methods

You will find that you make extensive use of these properties and methods when you code your event handlers.

For the Connection Object

The *Connect* method is used to initiate a new connection, and the *Select* method is required if there is more than one sensor at that address. The *SelectedSensor* property gives you a reference to the specific sensor as an *RGASensor* object.

For the RGASensor Object

You will need to check properties such as *MaximumMass* and *EgainCount* to ensure that you work within the range of the instrument. You can obtain a considerable amount of information about the sensor even when another user is controlling it. However you will need to invoke the *Control* method in order to change settings or acquire data. You will use the *SourceSettings* property to set the working condition of the ion source. You will access the scan data by obtaining a reference to the *Scan* object and control the filaments by obtaining a reference to the *Filaments* object.

You will certainly want to download measurement definition(s) to the sensor from which you can construct the *measurements* that make up the scan.

For the Scan Object

The scan object has methods to add sensor *measurements* into the scan. It has methods to *start* and *stop* scanning. The scan object has properties that lead directly to the partial pressure data *readings*.

This completes the Quick Tour of the object model. You can browse both the Java and ActiveX implementations in their entirety by clicking the relevant shortcut on the Start | Programs menu. The following chapter provides a complete reference to all the features.


Chapter 5

RGAComms ActiveX Object Model Reference

The RGAComms ActiveX component provides an extensive list of objects that can be used to manipulate the RGA. This reference documents all the objects, their properties and methods. The following diagram outlines the structure of the object model and allows you to link to the appropriate objects reference text. The diagram only lists properties that reference other objects, each object may have many more properties and methods but they tend to deal with base types such as integers, floats and strings and have been left out to make it easier to see how the objects fit together in the system. As you can see the hierarchy of objects starts from the RGAConnection class. This is the only object in the component that can be created by your program, all others are 'sub-objects' that are created dynamically as the connection finds out about the RGA it is communicating with. Having made a tcp/ip connection to a server and implicitly or explicitly selected a sensor to use the program will spend most of it's time dealing with the RGASensor object referenced by the SelectedSensor property of the RGAConnection. For most data acquisition tasks the program will then concentrate on the Filaments, Measurements, Scan and Data objects, while there are many other objects in the system they tend to be for information, tuning and calibration of the sensor which is done less frequently.

Object Hierarchy

The **Object Hierarchy** looks like this



RGACConnection

This is the starting point for all interaction between an application and the RGA. In order to make the component easy to use in as many languages as possible all events for all of the RGA objects are fired from this object. The CLASSID of the object is RGACComms.RGACConnection, the following code snippets show how to create the object in various languages.

VB (ensure the DLL is added to the references for the project)

```
Dim WithEvents myRGA As RGACConnection
Set myRGA = New RGACConnection
myRGA.Connect("127.0.0.1")
```

VBScript (running under WScript)

```
Dim myRGA
Set myRGA = WScript.CreateObject("RGACComms.RGACConnection","myRGA_")
myRGA.Connect("127.0.0.1")

Sub myRGA_Connected(bConnected)
    If bConnected Then
        MsgBox "Connected successfully to sensor at 127.0.0.1"
    Else
        MsgBox "Failed to connect to sensor 127.0.0.1"
    End If
End Sub
```

C++ (using Visual C++ with COM smart pointers)

```
#import "RGACComms.dll"

using namespace RGACCommsLib;
IRGACConnectionPtr pRGACConnection;
HRESULT hRes = pRGACConnection.CreateInstance(__uuidof(RGACConnection));
if (hRes==S_OK)
{
    printf("Created RGACConnection object OK\n");
}
else
{
    printf("Failed to create RGACConnection [0x%x]\n",hRes);
}
```

[Back to Object tree](#)

Properties

SelectedSensor As [RGASensor](#) [Read-Only]

The selected sensor if the connection is to a sensor server, or the actual sensor if the connection is to a dedicated sensor like eVision or MicroVisionIP. Upon connecting to a sensor server this property will be NULL and the program must decide which sensor to select by checking the SensorCount, SensorSerialNo and SensorName properties. Having decided which sensor to select the Select method should be used to select the sensor for this connection.

Protocol As [ProtocolInfo](#) [Read-Only]

The sensor protocol information. If an attempt to connect to a sensor fails then it might be due to the sensor not supporting the appropriate protocol version. This object allows an application to see the version required by this version of the RGACOMMS library and if the sensor did accept a TCP/IP connection then you can also find out the sensors protocol version details.

SensorCount As Long [Read-Only]

The number of sensors available on this server, for eVision and MicroVision IP sensors this will be 1 but for a windows server supporting MicroVision+ units it can be any number.

SensorSerialNo As String [Read-Only]

This property allows access to the serial number of a sensor connected to the server, it is accessed like an array where valid indexes are between 0 and SensorCount-1. For MicroVision+ sensors where the exact sensor must be selected this property can be used to verify that the server currently knows about the serial number you would like to select.

```
Dim n As Long
For n=0 To myRGA.SensorCount-1
    Debug.Print "Serial# " & myRGA.SensorSerialNo(n)
Next
```

SensorName As String [Read-Only]

This is similar to the SensorSerialNo property except that it returns the friendly name assigned to the sensor.

[Back to Object tree](#)

Methods

Connect(Address As String)

Connects over TCP/IP to the sensor or sensor server at the address specified. This method returns immediately and the Connected event will fire some time later indicating success or failure of the connection operation.

Refresh

This method forces the server to re-read its connections and to re-acquire the settings from the control unit. Normally the server automatically detects if instruments come back on-line, for example. Sometimes it may be beneficial to force the server to perform a Refresh, which is effectively the same as having it start over again.

Select(Index As Variant)

When a connection has been made to a windows server that is in control of MicroVision+ sensors connected via the RS232 interface it is necessary to select the appropriate sensor. Index can be the numeric index in the SensorSerialNo or SensorName properties or it can be the name or serial number text. While this method will work with all sensor types it is redundant for network enabled sensors such as the eVision or MicroVision IP. See Connected event description for more details of use.

[Back to Object tree](#)

Events

Connected(bSuccess As Boolean)

Following a call to the Connect method this event will eventually fire. The bSuccess parameter indicates whether a successful connection was made or not. In the event of bSuccess being False then the ProtocolInfo property can be checked to see if the connection failed as a result of a protocol version issue or because of a lack of network connectivity.

If bSuccess is True and the application is using an eVision or MicroVision IP sensor then the SelectedSensor property will be valid and allow access to the information about the sensor. If the application is communicating with a windows server that is managing MicroVision+ sensors then it will need to use the Select method to select the appropriate sensor. Having selected a sensor the Connected event will fire once again, indicating either a failure to select and download information about the sensor or success in which case SelectedSensor will be valid.

Controlling(bInControl As Boolean)

Following a call to the selected sensors Control method this event will fire indicating whether control was granted or not. If control was not granted then the sensor is either in use by another connection or is not in a state where it can be controlled, for example it might require configuration.

MeasurementCreated(bSuccess As Boolean, theMeasurement As [Measurement](#))

This event fires following a call to one of the RGASensor.Measurements methods to create and add a measurement. If a measurement fails to be created at the sensor then bSuccess will be False. The theMeasurement parameter always references the measurement object that was returned from the call to AddBarchart, AddAnalog, AddPeakJump or AddSinglePeak methods.

Note that it is possible to create many measurements in one go by calling the AddXXX methods successively but until the measurements are fully created they cannot be added to a scan. Similarly for PeakJump measurements, masses cannot be added/removed or changed until the underlying measurement has been fully created.

StartingScan(theScan As [Scan](#))

This event fires when a scan is just starting. The theScan parameter references the RGASensor.Scan property and provides access to the current scan number and how many scans are still to be done before a call to Scan.Resume will be needed. The purpose of this event is to allow any UI to be updated to show scan progress and to keep the scan running by calling Scan.Resume as required.

FilamentChange(theFilaments As [Filaments](#))

This event is fired whenever the filaments state changes. This allows an application to update any filament UI or abort what it is doing if the filaments trip off for some reason.

EndOfScan(theScan As [Scan](#))

This event is fired when the last piece of data for a scan has been acquired. Access to the data is through the Data property of the individual measurements in the scan. Note that no more data will be processed until this event returns, instead it will be buffered and processed upon completion of this events processing.

FilamentOnTime(theFilaments As [Filaments](#))

If the filament is configured to have a maximum on time and it is on this event will fire approximately every 2 seconds. The Filaments object allows access to the remaining time

which can be reset. The purpose of this event is to allow update of any UI that might allow a user to reset the filament on time or to automatically reset it if the application is in a state where it must ensure filaments don't accidentally switch off.

DigitalInputChange(thePort As [DigitalPort](#))

Whenever digital input lines change, this event will fire. The thePort parameter references the DigitalPort object whose input lines have changed.

InletChange(ActiveInlet As [Inlet](#))

Whenever the active inlet changes as a result of a valve change from an RVC or VSC device this event will fire. The ActiveInlet parameter is a reference to the newly selected inlet. Note that all data will have the appropriate inlet factor applied automatically. This event allows for UI update and re-scaling of any charts to reflect the new pressure ranges that might be acquired with the new inlet setting.

LinkDown(Reason As [rgaLinkDownReasons](#))

This event fires whenever the link to the sensor or server is lost. The Reason parameter identifies where the problem arose. For example if the sensor is a MicroVision+ and the RS232 cable was unplugged then the Reason value would indicate a problem with the Serial link.

This is a pretty catastrophic problem and the connection will need to be re-made when the sensor is available again.

RefreshComplete(bSuccess As Boolean)

This event is fired on completion of the server refresh method. If bSuccess is True the information about all the connected instrument(s) will be up to date.

AnalogInputReading(theInput As [AnalogInput](#))

Whenever an analog input reading is sent from the sensor this event will fire allowing easy access to the new value of the input.

TotalPressureReading(theTotalPressure As [TotalPressure](#))

Whenever a total pressure reading is sent from the sensor this event will fire allowing easy access to the new total pressure reading value.

RFTripState(theRF As [RF](#))

Fires when the RF Trip state changes

MultiplierState(theMultiplier As [Multiplier](#))

If the sensor has a multiplier fitted then this event fires when the Multiplier status changes

CirrusState(theCirrus As [Cirrus](#))

If the sensor has Cirrus hardware then this event fires when the Cirrus status changes

RVCPumpState(theRVC As [RVC](#))

If the sensor has an RVC fitted then this event fires when the RVC pump state changes

RVCH HeaterState (theRVC As [RVC](#))

If the sensor has an RVC fitted then this event fires when the RVC heater state changes

RVCValveState (theRVC As [RVC](#))

If the sensor has an RVC fitted then this event fires when the RVC valve state changes

RVCInterlocksState (theRVC As [RVC](#))

If the sensor has an RVC fitted then this event fires when the RVC interlocks state changes

RVCStatus (theRVC As [RVC](#))

If the sensor has an RVC fitted then this event fires when the RVC status state changes

RVCDigitalInputState (theRVC As [RVC](#))

If the sensor has an RVC fitted then this event fires when the RVC digital input state changes

[Back to Object tree](#)

RGASensor

This object represents a particular sensor and is accessed through the RGAConnection.SelectedSensor property. Initially the RGAConnection selects the sensor and downloads information about the sensors configuration. All this information can be accessed but no data can be acquired until control of the sensor is taken. Only one application can control a sensor at any time.

[Back to Object tree](#)

Properties**SerialNumber As String [Read-Only]**

Returns the serial number of the sensor.

State As rgainfoAvailableStates [Read-Only]

The state of the sensor when the connection was made and information downloaded.

Name As String [Read-Only]

Returns the friendly name assigned to the sensor

Version As String [Read-Only]

Returns the sensor software version, this can be used to ensure an application is working with the appropriate version of software if necessary.

SourceSettings As [SourceSettingsCollection](#) [Read-Only]

Gets at the RGA's source settings collection.

Inlets As [Inlets](#) [Read-Only]

Gets at the RGA's inlets collection.

EGainCount As Integer [Read-Only]

Number of electronic gains that the sensor can apply when acquiring data. See the Hardware Management section for more details on how the RGA electronics acquire data.

EGain As Single [Read-Only]

Gets the actual electronic gain value for a given index from 0 to EGainCount-1.

MaximumMass As Integer [Read-Only]

Maximum mass that the sensor can monitor. The complete range of the sensor will be 0.5 to MaximumMass+0.5 for Analog or Single Peak measurements and 1 to MaximumMass for the Barchart and PeakJump measurement types.

DetectorType As String [Read-Only]

A descriptive string for the type of detector that the sensor has. Possible values are:

Faraday only
Dual Faraday and multiplier
Dual Faraday and SCEM
Dual Faraday and advanced SCEM

PeakResolution As Integer [Read-Only]

The peak resolution of the sensor. This is the number of points across one AMU that can be sampled. Current sensors from MKS use a value of 32.

ConfigurableIonSource As Boolean [Read-Only]

Whether the ion source settings can be modified from software or if they are fixed in hardware. For MicroVision+ and MicroVision IP units it is possible to modify parameters for Ion Energy, Electron Energy, Extract Volts and Emission Current but for the eVision and eVision+ units they have just one set of fixed values.

FilamentType As String [Read-Only]

Type of filaments fitted:

Tungsten
Thoria
Rhenium
Yttrium

Filaments As [Filaments](#) [Read-Only]

Accesses the Filaments object which supports reading back of the current state and configuration of the filaments when not in control of the sensor. When in control of the sensor the filaments can be manipulated through this object.

InControl As Boolean

Gets whether this connection is in control of the sensor. Control can be released by setting this property to False. The Controlling event will fire when the sensor acknowledges that control has been released.

Measurements As [Measurements](#) [Read-Only]

The sensors measurements collection. Measurements can be created and removed from this collection. Measurements from this collection can be added to a scan. When a scan is stopped the scan's list of measurements is cleared but the measurements are still available in this collection.

Scan As [Scan](#) [Read-Only]

The sensors scan, it controls data acquisition of 1 or more measurements from the Measurements collection.

PressureUnits As [rgaPressureUnits](#)

Gets/Sets the pressure units used to report all pressure readings and factors. The default is Pascal but this can be changed depending on the units you are used to working in or would like to display.

AnalogInputs As [AnalogInputs](#) [Read-Only]

The sensors analog inputs collection. If the sensor has no analog inputs then the collection will be empty.

AnalogOutputs As [AnalogOutputs](#) [Read-Only]

The sensors analog outputs collection. If the sensor has no analog outputs then the collection will be empty.

RVC As [RVC](#) [Read-Only]

If the sensor has an RVC fitted then this property will return an object that allows access to the RVC's state. If no RVC is fitted then this property will return nothing.

Cirrus As [Cirrus](#) [Read-Only]

If the sensor has Cirrus hardware then this property will return an object that allows access to the Cirrus state. If the sensor does not have the Cirrus hardware then this property will return nothing..

Rollover As [Rollover](#) [Read-Only]

If the sensor supports rollover compensation (e.g. HPQ2s and HPQ-IP) then this property returns an object that allows reading and writing of the rollover compensation algorithm variables. If the sensor does not support rollover compensation then this property will return nothing.

DigitalPorts As [DigitalPorts](#) [Read-Only]

The sensors DigitalPorts. If the sensor has no digital I/O capability then this collection will be empty.

TotalPressure As [TotalPressure](#) [Read-Only]

Gets the sensors total pressure object if the sensor is configured to have an external total pressure gauge. If the sensor is not configured with a gauge this property will return nothing.

[Back to Object tree](#)

RF As [RF](#) [Read-Only]

Gets at the sensors RF configuration and state

[Back to Object tree](#)

Multiplier As [Multiplier](#) [Read-Only]

Gets the Multiplier object if the sensor has a multiplier fitted

[Back to Object tree](#)

Methods

Control(AppName As String, AppVersion As String)

Attempts to take control of the sensor. The application should pass it's name and a version string which will be returned by the sensor to other connections allowing them to see what is controlling the sensor and what IP address they are at. This method will return immediately but eventually the Connected event will fire indicating success or failure to take control.

[Back to Object tree](#)

ProtocolInfo

The ProtocolInfo object is always available to indicate the required version of the protocol that a sensor must implement in order for the control to be able to connect. In the Connected event the fields of SensorVersion and SensorMinCompatibleVersion will be filled in if the bConnected parameter is True. If the bConnected parameter is False then they may be filled in which indicates that the control was able to make a tcp/ip connection to the sensor but it's protocol version is either incompatible with the control or some other network error occurred after making the initial connection.

[Back to Object tree](#)

Properties

RequiredVersion As String [Read-Only]

The required protocol version, this is currently set at 1.2.

SensorVersion As String [Read-Only]

The sensors protocol version

SensorMinCompatibleVersion As String [Read-Only]

The sensors minimum compatible version

[Back to Object tree](#)

Filaments

The Filaments object provides information on the configuration and state of the sensors filaments as well as allowing control of them when the application is in control of the sensor.

[Back to Object tree](#)

Properties

FilamentOn As Boolean

Gets/Sets the filaments on/off state. Setting this property will attempt to put the filaments into the requested on or off state, this will result in one or more FilamentChange events.

SelectedFilament As Integer

Gets/Sets the selected filament 1 or 2. Switching between filaments 1 and 2 will result in one or more FilamentChange events.

TripState As [rgaFilamentTrips](#) [Read-Only]

Indicates if the filament is in tripped state or not and if so, the reason for the trip.

SummaryState As [rgaFilamentSummary](#) [Read-Only]

Indicates the summary state of the selected filament. Only on and off are guaranteed states but some sensors such as the MicroVision+ and IP will provide transient warming up and cooling down states.

XTripMode As [rgaFilamentXTripMode](#) [Read-Only]

Indicates the external trip configuration.

XTripEnabled As Boolean [Read-Only]

Indicates whether the external trip input is enabled or not.

EmissionTripEnabled As Boolean [Read-Only]

Indicates if the emission trip is enabled or not.

MaxOnTime As Long [Read-Only]

The maximum time in milliseconds that the filaments will stay on without being re-set. 0 indicates the filaments will stay on indefinitely.

RemainingOnTime As Long

Gets/Sets the time remaining in milliseconds for filaments to stay on if MaxOnTime is not 0.

[Back to Object tree](#)

Measurements

This collection allows an application to access the measurements of a sensor as well as create new ones or delete already existing ones. Measurements should not be created or deleted while the sensor is scanning.

[Back to Object tree](#)

Properties**Count As Long [Read-Only]**

Gets the number of measurements currently set up in the sensor

Item(Index As Variant) As [Measurement](#)

Gets a specific measurement from the collection. Index can be a numeric index from 0 to Count-1 or the name of the measurement. Note that accessing measurements from this collection returns the base Measurement type which all measurements are derived from. If you need to get at the properties of a specific measurement type then you must query for the appropriate interface, for late bound languages such as script you can access methods and properties of the specific measurement type as they are queried for at runtime anyway. If code needs to check the type of a measurement it can use the Measurement.Type property. The following example shows how one might access a measurement that is known to be a PeakJumpMeasurement in both VB and C++.

VB

```
Dim peakjump as PeakJumpMeasurement
Set peakjump = myRGA.SelectedSensor.Measurements("peaks")
peakjump.Masses(0).Mass = 18           ` Change first peak mass to be 18
```

C++ (Without smart pointers)

```
_variant_t vIndex("peak");
IMeasurement* pIMeasurement=NULL;
if (pIMeasurements->get_Item(vIndex, &pIMeasurement) == S_OK)
{
    IPeakJumpMeasurement* pIPeakJump=NULL;
    pIMeasurements->
    >QueryInterface(__uuidof(IPeakJumpMeasurement), (void**) &pIPeakJump);
    // Use the peakjump measurement
    ...
    pIPeakJump->Release();
}
```

[Back to Object tree](#)

Methods

AddPeakJump(Name As String, FilterMode As [rgaFilterModes](#), Accuracy As Integer, EGainIndex As Integer, SourceIndex As Integer, DetectorIndex As Integer) As [PeakJumpMeasurement](#)

Adds a peakjump measurement to the collection of measurements. The measurement object returned is in an uninitialised state until it is acknowledged by the sensor. At this point the MeasurementCreated event will fire on the RGAConnection object indicating success or failure.

Name	Each measurement must be given a unique name by the application.
FilterMode	Determines how the data for the peaks is acquired.
Accuracy	A number between 0 and 8 where 0 is the fastest but least accurate and 8 is slower but more precise, the exact speeds and data improvements vary from model to model of sensor.
EGainIndex	The index in the EGain array for the electronic gain that should be used when acquiring data. See the DetectorSettings EGainUsable and MaxPressure properties for help on selecting the correct gain. Not all electronic gains will be allowed depending upon the configuration and type of sensor.
SourceIndex	The index into the SourceSettings collection for the settings that should be used while acquiring data for this measurement.
DetectorIndex	The index into the chosen SourceSettings Detectors array. This is always between 0 and 3 where 0 is the Faraday detector and 1,2 and 3 use a multiplier with a configurable gain.

AddAnalog(Name As String, StartMass As Integer, EndMass As Integer, PointsPerPeak As Integer, Accuracy As Integer, EGainIndex As Integer, SourceIndex As Integer, DetectorIndex As Integer) As [AnalogMeasurement](#)

Adds an analog measurement to the collection of measurements. The measurement object returned is in an uninitialised state until it is acknowledged by the sensor. At this point the MeasurementCreated event will fire on the RGAConnection object indicating success or failure.

Name	Each measurement must be given a unique name by the application.
StartMass, EndMass	The mass span that will be acquired by this measurement. Whole AMU's are specified but the actual range spanned will be StartMass-0.5 to EndMass+0.5
PointsPerPeak	The number of points to take across each AMU. The maximum for this is specified by the sensors PeakResolution property. The value must be a power of 2 so typical values are 2, 4, 8, 16 and 32.
Accuracy	A number between 0 and 8 where 0 is the fastest but least accurate and 8 is slower but more precise, the exact speeds and data improvements vary from model to model of sensor.
EGainIndex	The index in the EGain array for the electronic gain that should be used when acquiring data. See the DetectorSettings EGainUsable and MaxPressure properties for help on selecting the correct gain. Not all electronic gains will be allowed depending upon the configuration and type of sensor.
SourceIndex	The index into the SourceSettings collection for the settings that should be used while acquiring data for this measurement.
DetectorIndex	The index into the chosen SourceSettings Detectors array. This is always between 0 and 3 where 0 is the Faraday detector and 1,2 and 3 use a multiplier with a configurable gain.

AddBarchart(Name As String, StartMass As Integer, EndMass As Integer, FilterMode As [rgaFilterModes](#), Accuracy As Integer, EGainIndex As Integer, SourceIndex As Integer, DetectorIndex As Integer) As [BarchartMeasurement](#)

Adds a barchart measurement to the collection of measurements. The measurement object returned is in an uninitialised state until it is acknowledged by the sensor. At this point the MeasurementCreated event will fire on the RGAConnection object indicating success or failure.

Name	Each measurement must be given a unique name by the application.
StartMass, EndMass	The mass span that will be acquired by this measurement. Whole AMU's are specified but the actual range spanned will be StartMass-0.5 to EndMass+0.5
FilterMode	Determines how the data for the peaks is acquired.
Accuracy	A number between 0 and 8 where 0 is the fastest but least accurate and 8 is slower but more precise, the exact speeds and data improvements vary from model to model of sensor.
EGainIndex	The index in the EGain array for the electronic gain that should be used when acquiring data. See the DetectorSettings EGainUsable and MaxPressure properties for help on selecting the correct gain. Not all electronic gains will be allowed depending upon the configuration and type of sensor.
SourceIndex	The index into the SourceSettings collection for the settings that should be used while acquiring data for this measurement.
DetectorIndex	The index into the chosen SourceSettings Detectors array. This is always between 0 and 3 where 0 is the Faraday detector and 1,2 and 3 use a multiplier with a configurable gain.

AddSinglePeak(Name As String, mass As Single, Accuracy As Integer, EGainIndex As Integer, SourceIndex As Integer, DetectorIndex As Integer) As [SinglePeakMeasurement](#)

Adds a single peak measurement to the collection of measurements. The measurement object returned is in an uninitialised state until it is acknowledged by the sensor. At this point the MeasurementCreated event will fire on the RGACollection object indicating success or failure.

Name	Each measurement must be given a unique name by the application.
Mass	The mass that will be measured. This can be fractional.
Accuracy	A number between 0 and 8 where 0 is the fastest but least accurate and 8 is slower but more precise, the exact speeds and data improvements vary from model to model of sensor.
EGainIndex	The index in the EGain array for the electronic gain that should be used when acquiring data. See the DetectorSettings EGainUsable and MaxPressure properties for help on selecting the correct gain. Not all electronic gains will be allowed depending upon the configuration and type of sensor.
SourceIndex	The index into the SourceSettings collection for the settings that should be used while acquiring data for this measurement.
DetectorIndex	The index into the chosen SourceSettings Detectors array. This is always between 0 and 3 where 0 is the Faraday detector and 1,2 and 3 use a multiplier with a configurable gain.

Remove(Index As Variant)

Removes a single measurement from the collection. Index can be the numeric index in the collection of the measurement or the name of the measurement.

RemoveAll()

Removes all measurements from the collection

[Back to Object tree](#)

Measurement

All measurements are derived from this base object which holds common properties that are shared by all measurements. It is possible to access data when the measurement is part of a scan, the specific type of the measurement to aid in casting to the appropriate specific measurement type in type safe early bound languages.

[Back to Object tree](#)

Properties

Name As String [Read-Only]

The name of the measurement

Data As [Readings](#) [Read-Only]

Accesses the data of the measurement if it is currently involved in a scan

Type As [rgaMeasurementTypes](#) [Read-Only]

The real type of this measurement

Accuracy As Integer

Gets/Sets the accuracy of the measurement

EGainIndex As Integer

Gets/Sets the selected electronic gain for the measurement

SourceIndex As Integer

Gets/Sets the selected source settings for the measurement

DetectorIndex As Integer

Gets/Sets the selected detector for the measurement

UseRolloverCorrection As Boolean

Determines if the measurement uses total pressure readings as a way to compensate for quadrupole rollover effects. This property can only be set True for sensors that support rollover correction, e.g. the HPQ2s or HPQ-IP.

Initialised As Boolean

Gets whether the measurement is fully initialised or not. For a measurement to be added to a scan it must have been created and initialised successfully. See the MeasurementCreated event for more details.

[Back to Object tree](#)

BarchartMeasurement

The BarchartMeasurement object allows acquisition of a series of consecutive masses.

[Back to Object tree](#)

Properties

Name As String [Read-Only]

The name of the measurement

Data As [Readings](#) [Read-Only]

Accesses the data of the measurement if it is currently involved in a scan

Type As [rgaMeasurementTypes](#) [Read-Only]

The real type of this measurement

Accuracy As Integer

Gets/Sets the accuracy of the measurement

EGainIndex As Integer

Gets/Sets the selected electronic gain for the measurement

SourceIndex As Integer

Gets/Sets the selected source settings for the measurement

DetectorIndex As Integer

Gets/Sets the selected detector for the measurement

FirstMass As Integer

Gets/Sets the first mass that will be acquired. When setting the FirstMass property the LastMass value will remain as it is unless FirstMass is set higher, in which case LastMass will be moved to match FirstMass.

LastMass As Integer

Gets/Sets the last mass that will be acquired. When setting the LastMass property the FirstMass value will remain the same unless LastMass is set lower, in which case FirstMass will be moved to match LastMass.

FilterMode As [rgaFilterModes](#)

Gets/Sets the way in which the measurements data is acquired and reported.

[Back to Object tree](#)

AnalogMeasurement

The AnalogMeasurement object allows acquisition of a series of consecutive masses where readings are taken across each mass. The actual mass range acquired by the measurement is from FirstMass-0.5 to LastMass+0.5.

[Back to Object tree](#)

Properties

Name As String [Read-Only]

The name of the measurement

Data As [Readings](#) [Read-Only]

Accesses the data of the measurement if it is currently involved in a scan

Type As [rgaMeasurementTypes](#) [Read-Only]

The real type of this measurement

Accuracy As Integer

Gets/Sets the accuracy of the measurement

EGainIndex As Integer

Gets/Sets the selected electronic gain for the measurement

SourceIndex As Integer

Gets/Sets the selected source settings for the measurement

DetectorIndex As Integer

Gets/Sets the selected detector for the measurement

FirstMass As Integer

Gets/Sets the first mass that will be acquired. When setting the FirstMass property the LastMass value will remain as it is unless FirstMass is set higher, in which case LastMass will be moved to match FirstMass.

LastMass As Integer

Gets/Sets the last mass that will be acquired. When setting the LastMass property the FirstMass value will remain the same unless LastMass is set lower, in which case FirstMass will be moved to match LastMass.

PointsPerPeak As Integer

Gets/Sets the number of points to be acquired for each mass. This must be a power of 2 with the maximum value being defined by the RGASensor.PeakResolution property. Usual values are 2, 4, 8, 16 and 32.

[Back to Object tree](#)

SinglePeakMeasurement

The SinglePeakMeasurement acquires a single reading from anywhere on the sensors mass range. Valid range is 0.5 to MaximumMass+0.5. The peak resolution of the sensor determines the smallest valid increment 1/PeakResolution which is normally 1/32.

[Back to Object tree](#)

Properties

Name As String [Read-Only]

The name of the measurement

Data As [Readings](#) [Read-Only]

Accesses the data of the measurement if it is currently involved in a scan

Type As [rgaMeasurementTypes](#) [Read-Only]

The real type of this measurement

Accuracy As Integer

Gets/Sets the accuracy of the measurement

EGainIndex As Integer

Gets/Sets the selected electronic gain for the measurement

SourceIndex As Integer

Gets/Sets the selected source settings for the measurement

DetectorIndex As Integer

Gets/Sets the selected detector for the measurement

Mass As Single

Gets/Sets the mass value that will be acquired by the measurement.

[Back to Object tree](#)

PeakJumpMeasurement

The PeakJumpMeasurement object manages a set of masses that will be acquired in turn which all share the same basic acquisition settings. When the measurement is initially created it contains no list of peaks to measure. The Masses property accesses the PeakJumpPeaks collection which allows masses to be added or removed.

[Back to Object tree](#)

Properties

Name As String [Read-Only]

The name of the measurement

Data As [Readings](#) [Read-Only]

Accesses the data of the measurement if it is currently involved in a scan

Type As [rgaMeasurementTypes](#) [Read-Only]

The real type of this measurement

Accuracy As Integer

Gets/Sets the accuracy of the measurement

EGainIndex As Integer

Gets/Sets the selected electronic gain for the measurement

SourceIndex As Integer

Gets/Sets the selected source settings for the measurement.

DetectorIndex As Integer

Gets/Sets the selected detector for the measurement.

Masses As [PeakJumpPeaks](#) [Read-Only]

Gives access to the masses that are being measured by this measurement.

FilterMode As [rgaFilterModes](#)

Gets/Sets the way that individual mass readings are acquired and reported.

[Back to Object tree](#)

PeakJumpPeak

The PeakJumpPeak object represents an individual mass that is being acquired by a PeakJumpMeasurement. Through this object it is possible to alter the mass value.

[Back to Object tree](#)

Properties**Mass As Integer**

Gets/Sets the mass for a given peak

Index As Long [Read-Only]

The index of this peak in all the peaks of the measurement

[Back to Object tree](#)

PeakJumpPeaks***Properties*****Count As Long [Read-Only]**

Gets the number of peaks that will be scanned by this PeakJumpMeasurement.

Item(Index As Long) As PeakJumpPeak

Accesses the individual peaks of the measurement by their 0 based index.

[Back to Object tree](#)

Methods**AddMass(Mass As Integer)**

Adds a mass to the peakjump measurement.

Scan

The data acquisition of the RGA Sensor is controlled by the Scan object, you can build up the data that you want to be measured by adding measurements to the scan. Calling Scan.Start starts the data acquisition process where all measurements are taken in the order that they were added. When each new scan starts the StartingScan event will fire which gives an opportunity to re-arm the scan so that it continues acquiring data indefinitely, alternatively the event might be ignored if individual scans need to be triggered from some asynchronous external input. After each measurement has had it's data acquired the EndOfScan event will fire which provides a convenient place to access all of the data from each of the measurements.

Properties

Measurements As [ScanMeasurements](#) [Read-Only]

The measurements that are part of this scan.

IsScanning As Boolean [Read-Only]

Indicates if the scan is actually taking readings or is in a paused but ready to take readings state.

Number As Long [Read-Only]

The current scan number.

Remaining As Long [Read-Only]

The number of scans remaining before the scan becomes paused.

Methods

Start(NumScans As Integer)

Starts the scan running. NumScans indicates how many scans the sensor will run before pausing.

Stop()

Stops the scan running and removes all measurements from it's list.

Resume(NumScans As Integer)

Resumes a paused scan or keeps an already running scan running. This is commonly used from the StartingScan event in order to keep a surplus of scans remaining so there is no delay at the sensor in acquiring data. For example if an application wants to free-run the scan so that data is coming in as quickly as it is acquired it can use Scan.Start(x) where x is some number of scans, say 10. Then in the StartingScan event the application can look at the Scan.Remaining count to see when it falls below some threshold, say 5 and call Scan.Resume(10). This way the sensor should always be busy acquiring data even if the application experiences some delays while processing data.

If on the other hand an application wants to start each individual scan when some other external event occurs so that acquired data is synchronised with the external event then it can use `Scan.Start(1)` and `Scan.Resume(1)` at the appropriate times.

[Back to Object tree](#)

ScanMeasurements

The collection of measurements being managed by the Scan object.

[Back to Object tree](#)

Properties

Count As Long [Read-Only]

Number of measurements that are part of the current scan

Item(Index) As [Measurement](#)

Provides access to the individual scan measurements by name or numeric index.

[Back to Object tree](#)

Methods

Add(Measurement As [Measurement](#))

Adds a measurement to the scan, the measurement must have been created and fully initialised. See the `MeasurementCreated` event and `Measurement.Initialised` property.

[Back to Object tree](#)

Readings

The Readings object provides access to the data acquired by a measurement that is part of a scan. It is possible to poll the data to see whether a measurement has completed its data acquisition although this is only recommended for updating progress information as there is a possibility of missing data. For safe access to all of the data an application should use the `EndOfScan` event as during this event no more data will be written to the Readings buffer even if it is being received from the sensor.

[Back to Object tree](#)

Properties

NumReadingsTotal As Long [Read-Only]

The number of readings that the measurement collects in a complete scan.

NumReadingsAvailable As Long [Read-Only]

The number of valid readings the measurement has currently collected in the current scan

Measurement As [Measurement](#) [Read-Only]

The measurement that these readings belong to

Reading(Index) As Single [Read-Only]

Gets an individual data point from the measurements data. Index can be a numeric index from 0 to NumReadingsTotal-1.

[Back to Object tree](#)

SourceSettings

The configuration for a given source entry. If the sensor has configurable ion source settings then it will typically have 6 SourceSettings entries that can be used when acquiring data depending upon the environment where acquisition is taking place. For sensors without configurable parameters there will be just one SourceSettings entry that allows control over the instruments mass resolution and alignment settings. See the Hardware Management chapter for more details.

[Back to Object tree](#)

Properties

DetectorCount As Integer [Read-Only]

The number of detector options available with these source settings. This is typically 1 or 4 depending on whether the sensor has a multiplier fitted or not.

DetectorSettings(Index) As [DetectorSettings](#) [Read-Only]

Gets the DetectorSettings object. Index 0 is always the faraday detector while 1,2 and 3 if present are different configurations of the multiplier.

Name As String [Read-Only]

The name of these source settings

ElectronEnergy As Single

The Electron Energy setting. 0 – 100eV.

IonEnergy As Single

The Ion Energy setting. 0 – 10eV.

ExtractVolts As Single

The Extract Voltage setting. -130 – 0V.

ElectronEmission As Single

The Electron Emission setting. 0 – 5mA

LowMassAlignment As Long

The Low Mass Alignment setting

HighMassAlignment As Long

The High Mass Alignment setting

LowMassResolution As Long

The Low Mass Resolution setting

HighMassResolution As Long

The High Mass Resolution setting

MaxRecommendedPressure As Single [Read-Only]

The maximum pressure recommended by MKS that can be measured using these source settings.

[Back to Object tree](#)

SourceSettingsCollection

The collection of SourceSettings objects available on the sensor. The index of each SourceSettings object is used when creating measurements to determine the settings used when the measurements data is acquired, the settings can be changed for each measurement through the SourceIndex property common to all measurements.

[Back to Object tree](#)

Properties

Count As Long [Read-Only]

The number of source settings available for the sensor. Typically 1 for a sensor with fixed settings such as the eVision and 6 for configurable sensors such as MicroVision+ or IP.

Item(Index) As [SourceSettings](#)

Gets a specific SourceSettings object from the collection. Index is 0 based numeric index.

[Back to Object tree](#)

DetectorSettings

The DetectorSettings object maintains calibration information and settings for faraday or multiplier detector using the specific SourceSettings. It also has some helper properties that indicate valid settings for data acquisition and the pressure range that those settings will work to.

[Back to Object tree](#)

Properties

Name As String [Read-Only]

Name of the detector (e.g. Faraday).

DefaultFactor As Single [Read-Only]

The default factor used to convert current to pressure.

DefaultVoltage As Single [Read-Only]

The default detector voltage. Unused for faraday detectors.

Factor(Filament As Integer) As Single

The factor presently being used to convert current to pressure.

Voltage(Filament As Integer) As Single

The detector voltage being used, only for Multiplier settings.

CalDate(Filament As Integer) As Date

Gets/Sets the calibration date associated with the detector factor/voltage. A date value of 0.0 indicates that a calibration date has never been applied. All date values should be read and written as UTC times.

MaxPressure(EGainIndex As Integer,[InletIndex As Integer]) As Single [Read-Only]

Calculates the max pressure that can be measured using this detector with the given electronic gain setting. If InletIndex is not specified or -1 passed then the currently active inlet will be used for the calculation.

EGainUsable(EGainIndex As Integer) As Boolean [Read-Only]

Determines whether an electronic gain setting is valid for this detector.

[Back to Object tree](#)

Inlet

Each sensor is considered to have at least one inlet but can be configured to have more if an external RVC or VSC is fitted. Each inlet has a pressure reduction factor that is applied to the data that is acquired when the inlet is the active one.

[Back to Object tree](#)

Properties

Fixed As Boolean [Read-Only]

Indicates if the inlet is fixed or variable.

CanCalibrate As Boolean [Read-Only]

Indicates if the pressure reduction factor can be changed.

TypeName As String [Read-Only]

Indicates the type of inlet

DefaultFactor As Single [Read-Only]

The default pressure reduction factor for the inlet.

Factor As Single

The pressure reduction factor being used by the inlet.

[Back to Object tree](#)

Inlets

Collection of Inlet objects, each Inlet maintains the calibration values for each inlet. The active inlet is maintained by this collection and if an inlet changes the InletChange event will fire to notify your app.

[Back to Object tree](#)

Properties

Count As Long [Read-Only]

The number of inlets for the sensor.

ActiveInlet As Long [Read-Only]

Gets the currently active inlet index.

Item(Index) As [Inlet](#)

Gets a specific inlet settings object.

[Back to Object tree](#)

AnalogInput

Settings of a specific analog input. When the analog input is enabled, values will be sent periodically from the sensor and the AnalogInputReading event will fire with the new value.

[Back to Object tree](#)

Properties

Enabled As Boolean

Enables/Disables the analog input reading

MinVolts As Single [Read-Only]

Returns the minimum voltage value that can be measured by this input

MaxVolts As Single [Read-Only]

Returns the maximum voltage value that can be measured by this input

Resolution As Integer [Read-Only]

Returns the resolution in bits of the ADC of this input

AverageCount As Byte

Gets/Sets the number of readings that are taken and averaged before being sent from the sensor

Interval As Long

Gets/Sets the time interval in microseconds between readings of the input. Time between readings is Interval x AverageCount

Value As Single [Read-Only]

Returns the most recent value of the input

Index As Long [Read-Only]

Returns the index of the input in the AnalogInputs collection

[Back to Object tree](#)

AnalogInputs

Collection of AnalogInput objects. If the sensor has no analog input capability then the collection will be empty.

[Back to Object tree](#)

Properties

Count As Long [Read-Only]

Gets the number of AnalogInput objects in the collection

Item(Index As Long) As [AnalogInput](#) [Read-Only]

Gets a specific AnalogInput instance from the collection

[Back to Object tree](#)

AnalogOutput

Settings of a specific analog output. The voltage of the output can be set via the Value property.

[Back to Object tree](#)

Properties

MinVolts As Single [Read-Only]

Returns the minimum voltage that this output can be set to

MaxVolts As Single [Read-Only]

Returns the maximum voltage that this output can be set to

Resolution As Integer [Read-Only]

Returns the resolution in bits of the DAC used by the output

Value As Single

Gets/Sets the value of the output voltage

Index As Long [Read-Only]

Returns the index of the output object in the AnalogOutputs collection

[Back to Object tree](#)

AnalogOutputs

Collection of AnalogOutput objects. If the sensor has no analog output capability then the collection will be empty.

[Back to Object tree](#)

Properties

Count As Long [Read-Only]

Gets the number of AnalogOutput objects in the collection

Item(Index As Long) As [AnalogOutput](#) [Read-Only]

Gets a specific instance of an AnalogOutput object from the collection

[Back to Object tree](#)

Cirrus

If the sensor has Cirrus hardware then this object allows access to the configuration and current state.

Properties

PumpState As [rgaCirrusPumpStates](#)

Gets/Sets the state of the pumps

HeaterState As [rgaCirrusHeaterStates](#)

Gets/Sets the state of the heater

CapillaryHeaterOn As Boolean

Gets/Sets the state of the capillary heater

ValveCount As Integer [Read-Only]

Gets the number of valves

ValvePosition As Integer

Gets/Sets the current valve position if ValveCount is not 0

ChamberPressure As Single [Read-Only]

Gets the chamber pressure

[Back to Object tree](#)

RVC

If the sensor has an RVC fitted then this object provides access to the RVC's configuration and state.

Properties

PumpState As [rgaRVC PumpStates](#)

Gets/Sets the state of the pumps

HeaterState As [rgaRVCH HeaterStates](#)

Gets/Sets the state of the heater

ValveOpen(ValveIndex As Integer) As Boolean

Gets/Sets the state of the valves

InterlocksOn As Boolean [Read-Only]

Gets whether the RVC is using its interlock logic or not

ValveMode As rgaRVCValveModes

Gets/Sets the current valve mode

Status(Index As Integer) As Boolean [Read-Only]

Gets the status (Hi Vac (0) / Lo Vac (1)) False = Pressure OK, True = Pressure too high

DigitalInput(Index As Integer) As Boolean [Read-Only]

Gets the current state of one of the digital inputs

AlarmOutput As Boolean

Gets/Sets the current alarm output state of the RVC

Methods**CloseAllValves**

Closes all the valves controlled by the RVC

[Back to Object tree](#)

Rollover

If the sensor supports quad rollover compensation then this object provides access to the variables used in the rollover compensation algorithm.

Properties**M1 As Single**

Gets/Sets the M1 value

M2 As Single

Gets/Sets the M2 value

B1 As Single

Gets/Sets the B1 value

B2 As Single

Gets/Sets the B2 value

BP1 As Single

Gets/Sets the BP1 value

PSF(Mass As Integer) As Single

Gets/Sets the peak scale factor for a given mass

Methods

UpdateVariables

Updates the sensors copy of M1, M2, B1, B2 and BP1 variables

[Back to Object tree](#)

RF

Holds configuration and current status of the RF.

Properties

TripEnabled As Boolean [Read-Only]

Gets whether the RF trip is enabled

Tripped As Boolean [Read-Only]

Gets whether the RF is currently tripped or not

[Back to Object tree](#)

Multiplier

If the sensor has a multiplier fitted then this object will be accessible. It maintains the configuration settings and current status of the multiplier. The LockedBySoftware property also allows custom control of blocking the multiplier for applications that need that control.

Properties

InhibitWhenFilamentOff As Boolean [Read-Only]

Gets whether the multiplier is configured to not be used if the filaments are off

InhibitWhenRVCh HeaterOn As Boolean [Read-Only]

Gets whether the multiplier is configured not to be used if the RVC heater is on

MultiplierOn As Boolean [Read-Only]

Gets whether the multiplier is currently on

Locked As Boolean [Read-Only]

Gets whether the multiplier has been locked for some reason

LockedBySoftware As Boolean

Gets/Sets the locked status of the multiplier under software control

LockedByFilament As Boolean [Read-Only]

Gets whether the multiplier has been locked due to the filament being off

LockedByRVC As Boolean [Read-Only]

Gets whether the multiplier has been locked because the RVC heater is on

[Back to Object tree](#)

DigitalPort

Holds the settings for a specific digital port and allows output bits to be set. When input bits change the DigitalInputChange event will fire.

Properties

Name As String [Read-Only]

Gets the name of this port

Index As Long [Read-Only]

Gets the index of this port in the DigitalPorts collection

ConnectedMask As Byte [Read-Only]

Gets the mask of bits that are configured to be connected

OutputMask As Byte [Read-Only]

Gets the mask of bits that are configured as outputs

Value As Byte

Gets/Sets the current value for a range of bits

[Back to Object tree](#)

DigitalPorts

Collection of DigitalPort objects. If the sensor has no digital I/O capability then the collection will be empty.

Properties

Count As Long [Read-Only]

Gets the number of digital ports.

Item(Index) As [DigitalPort](#) [Read-Only]

Gets a specific digital port by index or name.

[Back to Object tree](#)

TotalPressure

If the sensor is configured to have an external total pressure gauge then this object is accessible from the SelectedSensor object. It allows access to the latest reading for total pressure, calibration factor for the readings and the settings for how often the pressure is measured. Also see the TotalPressureReading event which is fired when readings are received from the sensor.

Properties

Value As Single [Read-Only]

Gets the latest total pressure reading sent from the sensor.

CalibrationFactor As Single

The calibration factor applied to the readings from the gauge. Usually this is left at the default value of 1.0 but in cases where a more accurate pressure value is available it may be desirable to calibrate against that value.

CalibrationDate As Date

Gets/Sets the calibration date associated with the total pressure factor (UTC). A value of 0.0 indicates that no date has been applied for a calibration.

AverageCount As Byte [Read-Only]

The number of readings that are measured and averaged before a reading is sent back from the sensor, this value is set at configuration time.

ReadingInterval As Long [Read-Only]

The interval in microseconds between successive readings of the ADC. Note that AverageCount readings are taken and averaged before the value is returned from the sensor.

[Back to Object tree](#)

Enum rgaFilamentTrips

Filament Trip States. These indicate if the filaments are tripped and if so, why.

- tripNone=0 Filaments are not tripped
- tripEmission=1 Filaments are tripped due to bad emission
- tripExternal=2 Filaments are tripped by external trip input
- tripRVC=3 Filaments are tripped by RVC

[Back to Object tree](#)

Enum rgaFilamentSummary

Filament summary state.

- filamentOff=1 Filament is off
- filamentWarmUp=2 Filament is warming up
- FilamentOn=3 Filament is on
- filamentCoolDown=4 Filament is cooling down
- filamentBadEmission=5 This state occurs when emission trip is disabled but the emission is bad - filaments are still on

[Back to Object tree](#)

Enum rgaFilamentXTripMode

Filament external trip configuration

- filamentTrip=0 External trip input trips filaments off
- filamentControl=1 External trip input controls filament on/off state

[Back to Object tree](#)

Enum rgaMeasurementTypes

Measurement types. Allows a generic Measurement to be queried for it's real type, can be useful for verifying the type of measurement before querying for specific measurement interface or in late bound languages to ensure that the measurement is the correct type to run the code coming up.

- measurementAnalog=0 Measurement is an analog type
- measurementBarchart=1 Measurement is a barchart type
- measurementPeakJump=2 Measurement is a peak jump type
- measurementSinglePeak=3 Measurement is a single peak type

[Back to Object tree](#)

Enum rgaFilterModes

Modes used to filter data in Barchart and PeakJump measurements

- PeakCenter=0 1 reading is taken at the nominal peak center
- PeakMaximum=1 1/2 AMU scanned in 8 points over nominal center and max value reported
- PeakAverage=2 1/4 AMU scanned in 8 points over nominal center and averaged

[Back to Object tree](#)

Enum rgaInfoAvailableStates

Available states of a sensor

- Unavailable=0 Sensor is unavailable
- Unconfigured=1 Sensor needs configuring
- InUse=2 Sensor is in use by someone else
- Available=3 Sensor is available for use

[Back to Object tree](#)

Enum rgaPressureUnits

Pressure units that data can be scaled to

- Pascal=0 Pressure reported as Pascal units
- mbar=1 Pressure reported as mbar units
- Torr=2 Pressure reported as Torr units
- milliTorr=3 Pressure reported as millTorr units

[Back to Object tree](#)

Enum rgaLinkDownReasons

Reasons for loss of link to server/sensor

- linkDownSerial=0 Link has been lost due to serial connection between sensor and server computer
- linkDownVSC=1 Link has been lost due to loss of communication between the sensor and it's VSC
- linkDownTCP=2 Link has been lost to sensor due to loss of the TCP/IP connection

[Back to Object tree](#)

Enum rgaCirrusPumpStates

States of the Cirrus pump

- cirrusPumpOff=0 Pump is off
- cirrusPumpAccelerating=1 Pump is accelerating
- cirrusPumpOn=2 Pump is on

[Back to Object tree](#)

Enum rgaCirrusHeaterStates

States of the Cirrus heater

- cirrusHeaterOff=0 Heater is off
- cirrusHeaterWarm=1 Heater is warming up
- cirrusHeaterBake=2 Heater is at bake temperature

[Back to Object tree](#)

Enum rgaRVCPumpStates

States of the RVC pump

- rvcPumpOff=0 Pump is off
- rvcPumpAccelerating=1 Pump is accelerating
- rvcPumpOn=2 Pump is on and at speed

[Back to Object tree](#)

Enum rgaRVCHeaterStates

States of the RVC heater

- rvcHeaterOff=0 Heater is off
- rvcHeaterOn=1 Heater is on
- rvcHeaterCoolingDown=2 Heater is cooling down

[Back to Object tree](#)

Enum rgaRVCValveModes

Operational valve modes of the RVC

- rvcValveModeAutomatic=0 Valves will operate automatically
- rvcValveModeManual=1 Valves operate manually

[Back to Object tree](#)

Chapter 6

The Java Libraries

The Java libraries are distributed in the form of JAR file. This file will be copied to the chosen directories for Microsoft Windows-based users; however for other users this will need to be done manually. A full set of Java Documents is included in the directory with the JAR file and an example Java Console Application.

The Java libraries adopt a very similar approach to that already mentioned in the ActiveX DLL; for a complete description of the libraries please refer to Chapter 5, since the 'object model' has been kept as language-independent as possible. There are slight variations between the two libraries relating to the method in which properties are accessed. In languages such as Visual Basic a property can be set or read as shown below.

```
`Setting
MyRGA.Filaments.filamentOn = true

`Getting
Dim bFilOn
bFilOn = MyRGA.Filaments.filamentOn
```

However, in the case of the Java Libraries a separate get and set method is required. The same code would therefore look something like

```
//Setting
MyRGA.getFilaments().setFilamentOn(true);

//Getting
Bool bFilOn = false;
bFilOn = MyRGA.getFilaments().IsFilamentOn();
```

In both the ActiveX SDK and the Java SDK a number of events are fired. These events are fired at the same time under the same circumstances for both of the SDKs. The Java SDK has one more event than the ActiveX SDK. This added event is fired when an error occurs, for example if an incorrect parameter is passed. At present the Java Libraries lack some of the features offered by the Active X libraries. A list of events, methods and properties available can be seen by looking at the various JavaDoc files. On Windows computers this documentation is accessible direct from the Start Menu.

An example of how to use an MKS RGA to perform a simple Peak Jump measurement is shown in Chapter 11.

Chapter 7

Redistributing the MKS RGA Libraries

When you have used the SDK to write your own RGA applications you will want to redistribute the necessary MKS libraries with your own software.

The license terms that MKS grants for this purpose are explained in the license document on the CD-Rom. In particular, you should note that you must not redistribute the complete SDK CD.

If you are targeting a non-Windows computer with the Java library, you only need to ensure that the JAR file is part of your installation package.

If you are targeting a Windows system, there is a re-distributable installation package contained within the `mksrgasdk_redistnnn.exe` file in the 'Redistributables' folder of the SDK CD. This package offers a choice of installable components.

- | | |
|---------------------|--|
| RGA Device Manager | The full software to manage RGA's on the target computer, optionally including support for serial RS232 instruments connected locally. |
| The ActiveX Library | This is installed to the <code>Windows\System</code> folder. |
| The Java Library | You can choose the location for this item. |

Chapter 8

An Introduction to the Examples

The SDK contains examples for use with the most popular programming languages – Visual Basic, C++, Java and LabView.

Two of the examples, Visual Basic and LabView, provide fully working implementations of a working RGA. The C++ example is a console application that shows a simple means of implementing a working RGA. The Java example provides a template that shows how to connect, control and acquire data from the RGA but leaves you to flesh out the capabilities depending on how you wish to integrate the RGA with your own software.

Whilst MKS Instruments, Spectra Group fully supports the ActiveX and Java libraries, the examples are supplied purely as examples. Although they may seem to build into applications that behave as a fully functioning RGA, they are supplied as-is; MKS cannot modify them to suit your purposes; neither can they be supported as applications in their own right.

MKS supplies its own Process Eye and EasyView software, which is a fully supported implementation of a recipe-driven or interactive RGA package.

Chapter 9

The VB Example

The example for VB6 is a complete implementation of a working RGA. You can include the forms and modules directly in a project of your own and connect to and control an RGA and get access to the data acquired.

The forms and modules are

frmMain This form simulates the world of your complete application into which you would like to embed the RGA example. This is the only form that would never become part of your application.

frmRGA This is the main RGA form that carries the buttons for controlling the RGA and displays a rolling buffer of data points. Barchart, peak jump and leak check data acquisition modes are provided and the accuracy, electronic gain and detector can be changed on-the-fly.

frmSettings This is a dialog form that allows you to change the data acquisition settings like the masses scanned and the units of pressure used.

frmAboutMyRGA A simple About box

modMyRGA Globally available helper functions

clsEventLog The event log class holds all the messages that were raised by the RGA software

clsHugeArrayRow The data stored in one row of the Huge Array

clsHugeArray The data buffer as an array of rows.

When the application runs the three main screens look like this

Your Main Application [X]

This Form Window represents your main application. It shows how you can control the RGA, find out about its state and manipulate the data it returns.

When the Form loads it creates an RGA window. By including this frmRGA in your main project you will have instant access to a wide range of functionality.

Connection

Specify the IP address of the instrument to connect to. If it is a serial instrument, the address will be the address of the PC to which it connects.

Connect to IP Address:

Instrument name (where more than one exists at the specified IP address):

Open the RGA window minimized

Connect automatically when this Form loads

Initial scan mode:

Questions

These buttons demonstrate

Checking the 'state' of the RGA Form:



Querying the RGA object:

Checking data in the Huge Array:

Checking the Event Log:

Answers

Scanning
Scan mode is Barchart



The *frmRGA* captures the following events

Connected
Controlling
MeasurementCreated
FilamentChange

LinkDown

It runs a **State Machine** that has the following enumerated states

```
TCPIPLinkDown
Offline
Connecting
ConnectedToServer
ConnectedToInstrument
ConnectedToUnavailableInstrument
SerialLinkDown
TakingControl
Controlling
FailedMeasurements
ReadyToScan
Scanning
```

The control unit has three measurements downloaded to it – a barchart, a peak jump measurement and a leak check measurement. The accuracy, electronic gain and detector can be changed interactively. Even the scanning mode can be changed, for example from peak jump to barchart. Properties such as first mass, last mass and scanned masses in the peak jump measurement are adjustable through a settings dialog, which is only accessible before connecting to the instrument.

As data becomes available it is stored in a huge array, in which the scan number, scan time and data values form the columns and successive scans make up the rows. In addition the sum of the readings is stored in column 3 of the array. One of the available settings is to set the buffer size for the huge array. When the stored row count exceeds twice the size set as the minimum the newest data is shuffled into the start of the array and the array is truncated back to its minimum size.

The RGA form supports a few registry settings that control the way that it behaves. These registry settings are stored in

HKCU\Software\MKS Instruments\VBRGA

Value name	Value Data
Autoconnect: Connect string	The IP address to connect to. If there is more than one RGA at this address, the name is appended after a separating semi-colon. For example '192.168.0.55' or '192.168.0.55;friendly name'
Autoconnect: Initial mode	One of 'idle', 'barchart', 'peakjump' or 'leakcheck'
Autoconnect: start minimized	One of 'yes' or 'no'

If the RGA form finds these entries in the registry it will automatically connect to the specified location and optionally start scanning – although it will not turn the filaments on.

The main window that represents the world of your application shows examples of how you can access the data and public methods on *frmRGA* from your own code on other forms. Click each of the buttons to discover

- The state of the state machine
- Whether the filaments are on or off
- The most recent data and scan number
- The most recent message held in the event log

○

You can also work directly in the RGA window. Certain controls are disabled while connected, whilst others are disabled until connected.

All the code is fully commented and you are encouraged to study it to find out more about the implementation of an RGA in VB6.


Extending the example code

There are several ways in which you could extend the example code

- You could add methods to access data from the huge array by scan number or by nearest time
- You may need to add functionality to access the analog or digital inputs on the control unit
- You may need to add code to change the Factor property of an individual DetectorSetting item so that the RGA is calibrated to an external reference gauge. You can ratio the external reference to the sum of the peaks in column 3 of the grid and apply this ratio to the Factor property. Note that smaller values for Factor give rise to larger peaks and larger values give rise to smaller peaks.
- If your control unit has an associated Remote Vacuum Controller (RVC) you will have to add code that implements some of the RVC object methods.
- You can easily add code to store the data acquired to a data file in a format of your choice.
- You could enable interactivity of the masses acquired but you need to remember that changing the number of columns of data stored in the huge array will take some careful implementation.
- You can certainly change the source settings interactively – thereby switching from standard electron energy to low electron energy, for example.

Chapter 10

The C++ Example



The C++ example is a small console based application that works by using the mouse or keyboard to control buttons and an edit box, in this respect it is very similar to a true windows GUI application but demonstrates just the key pieces of code to use the control.

Using the sample

To use the sample simply enter the IP address of the sensor you want to communicate with into the box at the top left of the screen. Click the blue button labelled CONNECT to attempt to connect to the sensor. If the address is correct then the sensor will allow the connection and the CONTROL button at the top of the window will turn from grey to blue indicating that it is enabled. At this point the other areas of the display should update to reflect the current state of the sensor, this information includes the active filament and it's on/off state as well as information about any external digital or analog I/O configuration.

Having connected to the sensor and viewed the information about the current state and configuration you must take control in order to do more. Only one application can be in control of a sensor at any time so applications must be prepared for any attempt to take control to fail. Click the CONTROL button to attempt to take control, if successful then the remaining buttons should become active and you can click them to control filaments and run scans. While in control of the sensor any analog/digital I/O changes will be displayed.

The program is a Visual Studio 6 project made up of the following files:

CPPSample.cpp	The main file for the sample
RGACConnection.cpp/RGACConnection.h	A class that handles creation of the RGACConnection COM component and hooking of it's events. It is designed to be a base class from which derived classes only need to provide implementation for the events that are fired by the control. You may want to use this as the basis for your own applications using the component.
Display.cpp/Display.h	These files contain all the code to handle the user interface of the application and is unimportant in terms of understanding how to use the ActiveX component.

The example uses Visual C++'s ability to import an active-x type library and create wrapper classes that simplify working with COM interfaces. For detailed information on how COM works and how to program it from C++ we would recommend Microsofts MSDN website (<http://msdn.microsoft.com>).

The CRGACConnection class

The CRGACConnection class contains the code necessary to hook and unhook events from the Active-X component. To handle the events we must implement the _IRGAEvents interface which is defined inside the objects type library. It is essentially an IDispatch interface but with ID's and parameter lists assigned for the various events. CRGACConnection takes care of instantiating an instance of the RGACConnection object which is accessible through the m_IRGACConnection member variable. Having created the class instance it tries to hook the events which are fired from the object by calling the Invoke() method. The implementation of Invoke() then calls the appropriate handler function which are all defined as pure virtual functions in CRGACConnection.

To use CRGACConnection all that is required is to derive a class from it and provide overrides for the following event handlers:

```
virtual void Connected(bool bConnected)=0;
virtual void Controlling(bool bInControl)=0;
virtual void MeasurementCreated(bool
    bCreated, RGACCommsLib::IMeasurementPtr theMeasurement)=0;
virtual void StartingScan(RGACCommsLib::IScanPtr theScan)=0;
virtual void EndOfScan(RGACCommsLib::IScanPtr theScan)=0;
virtual void FilamentChange(RGACCommsLib::IFilamentsPtr theFilaments)=0;
virtual void FilamentOnTime(RGACCommsLib::IFilamentsPtr theFilaments)=0;
virtual void DigitalInputChange(RGACCommsLib::IDigitalPortPtr
    thePort)=0;
virtual void InletChange(RGACCommsLib::IInletPtr ActiveInlet)=0;
virtual void LinkDown(RGACCommsLib::rgaLinkDownReasons Reason)=0;
virtual void AnalogInputReading(RGACCommsLib::IAnalogInputPtr
    theInput)=0;
virtual void TotalPressureReading(RGACCommsLib::ITotalPressurePtr
    theTotalPressure)=0;
```

It is then straightforward to create instances using the C++ new operator. Due to the reference counting implementation the C++ delete operator should not be called directly to destroy an instance, instead you should call the ShutDown() method which will unhook from the objects events and release our initial reference, if this is the final reference then the class will be destroyed, otherwise when the final reference is released the class will be destroyed. In use this means that if you are trying to destroy the object from within an event handler the final release wont be called on the object until you return from the event handler, however if you call ShutDown() from anywhere else you will be the one with the final reference and the object will be destroyed immediately.

Notice that all the ActiveX components objects, methods, properties and enumerations are defined inside the RGACommsLib namespace. If you would rather not prefix items with RGACommsLib:: in your code then you can always put the using namespace RGACommsLib declaration at the top of your source code.

CPPSample.cpp

This is the main program file and the one that you should focus on to see how to use the active-x component and when the various events get fired.

At the bottom of the file you will find the main() function where everything starts from. You will see that the first thing the code does is to call CoInitialize(NULL) to get the COM libraries set up. After that we enter a loop that simply waits for events to happen and then processes them:

```
while(1)
{
    switch(MsgWaitForMultipleObjects(1, &hConsole, FALSE, INFINITE, QS_ALLINP
UT))
    {
        case WAIT_OBJECT_0:
            // Some console event(s) has occurred, process them
            conin.ProcessEvents();
            break;
        case WAIT_OBJECT_0+1:
            // Window messages pending
            while(PeekMessage(&msg, NULL, 0, 0, PM_REMOVE))
            {
                DispatchMessage(&msg);
            }
            break;
    }
}
```

The active-x dll creates a separate thread for handling all tcp/ip communications for instances of the RGAConnection class. When key events need to be handled by the application that created the object the communication thread posts messages to it. For this mechanism to work your application must process windows messages, for most windows GUI applications this is normal and often conveniently hidden away by frameworks such as MFC. This example is a console application to show how these messages might be processed by a non GUI application that might be common in an embedded type system.

Moving back up to the top of the file we derive the MyRGA class from the CRGAConnection class so that we can implement the event handler functions that the active-x component will call as things happen. This application is only designed to

communicate with one sensor at any time so we have a global variable called pRGA that is a pointer to a MyRGA class instance. Initially this pointer is NULL and will only become valid while we have a valid connection to a sensor.

In the display.h and display.cpp files there is a class called UI that holds class instances that manage the screen display, buttons and address edit box. Each button is given a callback function to call when the button is clicked and each of these callbacks is a static member function of the UI class. These functions are implemented in CPPSample.cpp so that all the code that manipulates the active-x object model is in one location for the sample. The UI functions are as follows:

```
void UI::ClickConnect(Button* btn)
```

This code handles connecting or disconnecting from the sensor. If the pRGA variable is NULL then we create one and try to connect. This will result in the Controlling() event firing to indicate success/failure of the connection operation.

```
void UI::ClickControl(Button* btn)
```

This button is only enabled when we have a valid connection to a sensor. The code here uses the SelectedSensor's InControl property to see if we are currently in control of the sensor or not and if we are then it sets InControl to false, otherwise we call the Control() method passing our application name and version to the sensor. When the InControl state changes the Controlling() event will be fired indicating the new state.

```
void UI::ClickFilamentSelect(Button* btn)
```

When we are in control of the sensor clicking this button toggles the selected filament between 1 and 2. All changes to the filament state are then reported through the FilamentChange() event being called.

```
void UI::ClickFilament(Button* btn)
```

This attempts to toggle filament state between on and off, once again all changes to the filament state are reported through the FilamentChange() event being called.

```
void UI::ClickBarchart(Button* btn)  
void UI::ClickAnalog(Button* btn)  
void UI::ClickSinglePeak(Button* btn)  
void UI::ClickPeakJump(Button* btn)
```

These buttons start or stop the appropriate scan measurement. It starts by stopping any scan that might be currently running. Next it uses the state of the button to determine whether we are starting or stopping the scan mode. If we are starting then the measurement with the name Analog, Barchart, PeakJump or SinglePeak is located in the sensors measurements collection and added to the scan measurements collection, the display is set up for the data that will be acquired and the scan is started.

Once the scan is started the StartingScan() and EndOfScan() events will fire as scans progress.

Note that in this example only one measurement is added to the scan but there is no real restriction on the number or types of measurement that can be added to the scan and acquired.

Further down in this file are the implementations of each of the event handler functions that the active-x control can call:

```
void MyRGA::Connected(bool bConnected)
```

This event code is called when the initial connection has been made or failed. In this example this event is called after clicking the CONNECT/DISCONNECT button. When a successful connection has been made to the sensor the UI is updated to show the currently selected filament and analog/digital I/O values.

```
void MyRGA::Controlling(bool bInControl)
```

After clicking the CONTROL/RELEASE button this event should fire indicating the updated controlling state. Upon a successful attempt to take control of the sensor the code creates measurements and enables any available analog inputs.

Each attempt to create a measurement results in the MeasurementCreated() event being called with the success/failure of the operation.

Note that there is no restriction that forces you to create measurements during this event but for this example it is convenient to create the measurements as soon as we take control of the instrument.

```
void MyRGA::MeasurementCreated(bool bCreated, IMeasurementPtr  
theMeasurement)
```

Following a call to the AddAnalog(), AddBarchart(), AddPeakJump() or AddSinglePeak() methods of the selected sensor's measurements collection a valid measurement instance is returned but it is in an incomplete state until the sensor either confirms or denies creation of the measurement. This event will be called to indicate the completion of each of the measurements that are created in the Controlling() event handler. The code uses this opportunity to enable the appropriate scan mode buttons for each measurement as they are created and in the case of the PeakJump measurement some masses are added for it to acquire.

```
void MyRGA::StartingScan(IScanPtr theScan)
```

This event is called at the beginning of each scan. It allows a program to update any user interface state that might show the current scan number or time. It is also an ideal place to use the Scan's Resume() method to ensure that the sensor continues to acquire and send more data. In this example for the Analog, Barchart, and Peak Jump scans we tell the sensor to do one scan at a time because the time between receiving this event and the end of the scan is easily enough to inform the sensor to scan again. However, for the single peak measurement it is more likely that the sensor will send one reading and then be left waiting before receiving the message to scan again so in this case we start the scan asking for 128 to be done by the sensor and then in this event we check if the number of remaining scans is less than 64, if it is then we call the scans Resume() method asking for a further 128 scans to be done.

Note that whilst this is a convenient place to call the Scan->Resume() method it may be called from elsewhere quite safely. For example, if you had some external event source that indicates when to take data you might choose to call Resume() from there to ensure that data is synchronized in some way.

```
void MyRGA::EndOfScan(IScanPtr theScan)
```

This event is called at the end of each scan and allows the application to access all of the acquired data for each measurement that is part of the scan. In this example

we scroll the screen display data up by one line and then write out the new data in the bottom line. A real application might log data to a file, database and/or check alarm limits for example.

```
void MyRGA::FilamentChange(IFilamentsPtr theFilaments)
```

This event is called whenever filament state changes in some way. The example just updates the user interface text for the selected filament and it's on/off state.

```
void MyRGA::FilamentOnTime(IFilamentsPtr theFilaments)
```

If the sensor is configured to have a maximum time that the filaments stay on for then this will be called at approximately 2 second intervals to let the application know how much time remains. In this example we do nothing with this event.

```
void MyRGA::DigitalInputChange(IDigitalPortPtr thePort)
```

Whenever digital ports change in some way this event will be called. We simply update the display to show the new value.

```
void MyRGA::InletChange(IInletPtr ActiveInlet)
```

If the sensor is configured to have an RVC or VSC then it will have more than one inlet. When valve configuration changes the active inlet will change and this event will be called. All data that is acquired is automatically scaled by the appropriate factors for the active inlet but an application might use this event to log the fact that the inlet changed or scale some display to the new limits expected for this inlet. In this example we do nothing in this event.

```
void MyRGA::LinkDown(rgaLinkDownReasons Reason)
```

If the link to a sensor is lost then this event is called. In this example we indicate that the link was lost on the status bar at the bottom of the display and then release the active-x component.

```
void MyRGA::AnalogInputReading(IAnalogInputPtr theInput)
```

Whenever analog input channels are enabled they will return readings periodically. This event is called when such readings are received. For this example we just update the screen display with the current value.

```
void MyRGA::TotalPressureReading(ITotalPressurePtr theTotalPressure)
```

If the sensor is configured to have a total pressure gauge fitted then this event is called when a new total pressure reading is measured by the sensor. For this example we simply update the total pressure text display.

Chapter 11

The Java Example

The following example illustrates the basic methods required to perform a simple but common use of an MKS RGA. It covers the following topics:

- Connection
- Information Retrieval
- Control
- Filament Control
- Measurement Creation
- Data Acquisition
- Releasing Control

Connection

The first step in the example is to create an instance of the CRGAConnection class. This is done with the line:

```
RGA = new JRGACConnection(this);
```

It can be seen that the constructor for this class requires a parameter being passed. This parameter must be a reference to a class that implements the IEvents interface. In the example, Class CRGAExample does just this, and in doing so creates the methods which will be called in order handle the appropriate events that are generated by the MKS RGA.

After the JRGACConnection Object is created, Connect is called with the IP Address of the MKS RGA Server, in the case of the example this line is:

```
RGA.Connect("10.10.54.100");
```

It is at this point that the true event driven process is entered. When the server is connected, the OnConnected event will be fired. Once this has occurred there are two possible courses of action that can be taken. For old style RS232 type sensors, one RGA server may have several sensors, however with the modern IP range of sensors, there will be only one sensor per server. Therefore, in the case of RS232 sensors it is necessary to use the following line:

```
RGA.Select(0);
```

This will result in the OnConnected event being fired one more time however, this time the SelectedSensor property will no longer be null. This process ensures that the first sensor in the sensors collection has been selected for use during data acquisition. This is not required for the IP based sensors.

Information Retrieval

After the OnConnected event has been received and the selected sensor is no longer null a number of pieces of information are available to the user.

This information is displayed on the console using a standard `System.out.println` call. All of the information available is described in Chapter 5.

Control

Before any data acquisition is possible it is necessary for the user to gain control of the sensor. This is done using the following command:

```
Sensor.Control("Demo Console App","1.0");
```

On completion of this method call the `OnControlling` event will be fired. This event contains a property which indicates whether the control attempt was successful. If it is successful the user can then begin to prepare the sensor for data acquisition.

Filament Control

Before any valid data can be retrieved it is necessary to turn on the filaments of sensor. This is done using the following line:

```
pSensor.getFilaments().setFilamentOn(true);
```

This will cause the `OnFilamentChange` event to be fired. The parameter passed into this event handler contains the relevant information concerning the filaments.

Measurement Creation

In order to gain any useful partial pressure readings it is necessary to create a measurement within the sensor. This example uses a Peak Jump Measurement to gather values for four masses. To create a Peak Jump Measurement the following command is used:

```
pSensor.getMeasurements().AddPeakJump("MyMeasurement1",  
    rgaFilterModes.PeakAverage,5,1,0,0);
```

Once this command has been received and processed by the sensor the `OnMeasurementCreated` event will be fired. If the measurement was successfully created the individual masses can be added to it (Note this is only necessary for the Peak Jump Measurement).

After adding the appropriate masses to the measurement, it is then necessary to add the measurement to the task list or Scan. This is done using:

```
pSensor.getScan().getMeasurements().Add(theMeasurement);
```

Once the measurement has been added to the Scan, the Scan can then be started. This causes the data acquisition process to begin.

Data Acquisition

Upon a Scan starting the `OnStartingScan` event is fired. The event passes a Scan object containing important information such as remaining number of scans left etc. It should be noted that when starting a scan it is possible to specify the number scans to complete before finishing, this enables the client application to control the amount of data being returned at anyone time. Therefore this examples only ever does one scan at a time and then uses:

```
theScan.Resume(1);
```

to continue the scanning process. Once a scan has completed the `OnEndOfScan` event is fired. It is at this point that user can retrieve the partial pressure values for the appropriate mass. This is done using:

```
JReadings pReadings = theScan.getMeasurements().getItem(0).getData();  
  
System.out.println("\tMass 18 = "+pReadings.getReading(0)+" Pascal");  
System.out.println("\tMass 28 = "+pReadings.getReading(1)+" Pascal");  
System.out.println("\tMass 40 = "+pReadings.getReading(2)+" Pascal");
```

Releasing Control

Once the client application has finished, the control of the sensor must be relinquished. This can be done using:


```
pSensor.setInControl(false);
```


All of the above method calls are done using the `RGAComms.jar` file. It necessary to import this library into your own project before compiling the example code, there are various development suites available to do this, each containing the relevant information of how it should be achieved.

Chapter 12

The LabView Example

The LabView version 6 example has much in common with the VB example as the following screenshots will illustrate.





The connection panel is similar to the VB example. You must enter the **IP address** where the server is located (127.0.0.1 is the local computer). If the server has multiple instruments, you specify the specific instrument by name. The **Connect** button is used to connect or disconnect from the server. The VI also has an **event log** that displays a list of the most recent messages thrown up by the VI in latest-at-the-top order. Once connected, the unit can be in one of **three modes** – idle, barchart and peak jump. However this example exploits the advantages of LabView in displaying the acquired data graphically. However, compared with the VB example there is no leak check mode and the five masses scanned in peak jump mode are hard coded in the VI diagrams. The same five masses trended in peak jump mode are also trended automatically in barchart mode. The **accuracy** and **electronic gain** are adjustable while scanning. Like the VB example there is a **settings dialog**, but this dialog only supports the selection of the pressure units for the display.

The coding also follows a similar philosophy to all the other examples, but since LabView is multi-threaded and data driven there are detail differences between this example and all the others.

There are essentially three components to the MyRGA VI

- The event handler to capture the events raised by the ActiveX DLL
- A polled loop that captures the user input on the front panel
- An occurrence driven data display handler that puts new data into the charts as it arrives.

A number of the most important properties of the example are stored in one of three global variable VI's. You can access these from your own top level VI's.

- Event Globals
- RGA Globals. This includes the current state.
- Data Globals. This holds the most recent scan data.

The VI's are fully self-documented and there is a description for all the major features, so you will want to enable the context help window as you browse through the code.

Extending the example code

There are several ways in which you could extend the example code

- You could hold more than one scan of data in the data buffer
- You could enable more interactivity by offering selection of different source and or detector settings
- You could enable interactivity of the masses acquired. However as soon as you permit a free choice of the peak jump masses acquired you need to think about how they will be displayed on the stripchart.
- You could easily add data file storage using the LabView 'SaveToSpreadsheetFile' VI's.

Chapter 13

Using Virtual RGA Systems

Version 1.10 of the SDK introduces the concept of Virtual RGA's. A Virtual RGA provides application developers with a great simulation platform for recipe code or customised software developed using the ASCII SDK.

What is a Virtual RGA System?

A virtual RGA system is a complete PC-memory-based implementation of the RGA hardware along with any peripheral hardware such as an RVC, as well as a model of the vacuum system itself. A virtual RGA system appears to the application software as a fully functioning system in every respect. Here is just a summary of the features of a virtual RGA system.

The control unit and analyser are modelled, including

- Ion source parameters such as emission and electron energy
- Filter parameters including alignment and resolution
- Detector behaviour including a typical multiplier gain/voltage characteristic and optional peak noise
- Filament behaviour including filament trips and RF trip
- Analog and digital I/O
- External hardware such as an RVC or a Cirrus that controls external pumping

The vacuum environment is modelled such that

- The individual peak height of any mass can be set
- One of a selection of standard spectra can be displayed
- A dynamic profile can be 're-played' from a file
- Random variations in selected peaks can be incorporated.

All virtual systems are managed from a central control panel and changes to an individual system can be made interactively. On the other hand an ActiveX control is provided so that any input to a system can be programmatically changed as part of a complete simulation environment.

How are Virtual RGA's Installed?

You must select the 'Virtual RGA Support' feature when selecting the options to install.

How do Virtual RGA's Work?

A virtual system is a memory-mapped image of a control unit and its associated hardware. Version 5.20 of the RGA Server application has been modified to 'connect' to this memory mapping in just the same way as it connects to a 'real' instrument. Having connected to this virtual instrument it presents the same interface to the applications that connect to it, using the standard MKS ASCII protocol.

Thus, provided virtual RGA's have been enabled in the administrative setup, all the applications in the suite that communicate with the server as their source of data, behave exactly as though they were connected to a real RGA. You can write or run recipes for them, see them in RGA Device Manager, configure them with Hardware set-up, degas them and run diagnostics, calibrate them, save and recall data from them and even simulate what happens when the serial or network cable is pulled out.

You can control virtual RGA's using Process Eye, EasyView, the Developer's SDK or, at the most basic level, by using ASCII commands over telnet.

Are there any limitations to Virtual RGA's?

Very few. You can create any number of virtual RGA's of a variety of types, including Microvision Plus, HPQ2, HPQ2S, e-Vision, e-Vision+, Cirrus, Vision 1000P, 1000C and 1000E and 300mm ResistTorr. The only limitations are

- There is no simulation of browsing to a network-enabled instrument
- There is no simulation of a VSC
- There is no simulation of rollover on standard analysers. If the vacuum on a standard analyser is set at $1e^{-3}$ torr, the resultant peak height will be full scale $1e^{-4}$ torr, rather than what would happen in practice where the peaks would be significantly lower than this.
- Likewise there is no internal simulation of the rollover correction procedure in an HPQ2S – that is to say that there is no attempt to simulate the actual (uncorrected) ion current at a given high pressure and then apply rollover correction to that signal.
- There is no core or app download.
- The filaments do not burn out when they are over-pressurised!

e-Vision control units are modelled as though they were serial instruments. They have all the same configuration as a real e-Vision and behave in EasyView (for example) exactly as a real e-Vision would, but they cannot be configured because they cannot be browsed to. However, the factory choices of 40eV or 70eV can be selected when the virtual e-Vision is created.

All the information in a virtual RGA system persists when you log off or close down your computer. This is obviously important for objects like an RVC (which must continue monitoring the state of vacuum all the time) and it also means that, like a real RGA, the source, filter and detector settings also persist through shut-downs and re-starts.

Using a Virtual System

Getting Started

Run the SDK Setup or re-run it using Control Panel and Add / Remove Programs, selecting the 'Modify' option. On a subsequent page you can select whether you wish to have the serial server connect to real instruments in addition to virtual ones. Complete the set-up.


Log off as the administrator. Do this even if you run your software as the administrator.

Log on as the regular RGA user. You will notice a new icon in the system tray.



This is the icon for the Virtual RGA Control Panel. (If you are wondering, it is the symbol for a vacuum ion gauge.)

Double click on the 'Virtual RGA Control Panel' (referred to from here on simply as the control panel). Until you create your first virtual instrument, it will have an empty list. Move the window down to the bottom of the screen to a convenient location.



Now, in the control panel, click the 'Add' button (the top one). From the drop down list choose the system you would like to create. The default is an LM76 Microvision Plus, but you can choose any one of a total of nine virtual systems. The control panel will suggest an unused serial number and name for your instrument, so you generally do not need to change what is proposed, unless you are creating a simulation of a real factory system. Click Ok and your first virtual system will be created. If you want further instruments, repeat the procedure again. Close the Virtual RGA Control Panel with the Close button.

Now run the regular user set-up. If you have just upgraded, run through to the end in the normal way, to register the new version. Then run it again. At the hardware page select the first option to search for and configure the installed hardware. This, as usual, will launch RGA Device Manager. Your instrument(s) will be in the list. Right click on the instrument and select the 'Access...' menu to set 'exclusive use' of this instrument.


Close RGA Device Manager.

You are now ready to run your software and acquire virtual data.

The Virtual RGA Control Panel

Run your software. Double click the control panel icon and check that the position is somewhere convenient on your screen. The position you select will be remembered the next time you use it. Double click an entry in the control panel and a new multi-tabbed window will appear that is your interface to that virtual system. Position this window somewhere convenient. This position will also be preserved. You can click the Close box on any of these virtual RGA windows. The effect will merely be to hide the window. You restore a particular window by right mouse clicking in the system

tray on the  icon. This will display the main control panel menu.



The remaining buttons on the control panel are, in turn



The **'trash'** button. This removes the virtual system from your computer.



The **'disconnect / re-connect'** button. This simulates a link down between the computer and the control unit. Click the button once to invoke a link down. The status in the control panel will show 'link down'. Click the button a second time to re-connect with the server.




The **'properties'** button. This opens the multi-tabbed window that is the user interface to the selected virtual system. Clicking this button is exactly the same as double clicking the entry in the list.

The items on the preferences menu are global settings for the application as a whole.

Running the RGA Control Panel is the equivalent of powering up your control units. If you run the server but do not run the control panel, this is the exact equivalent of de-powering all your control units. Equally, if you run the control panel but close the server, this is the exact equivalent of plugging all your control units in to the 24V supply, but with no server running their COMMS LED's will be blinking slowly.

All the remaining interaction with a specific virtual RGA is carried out by activating its individual System Properties window as will now be described.

System Properties – The Control Unit Tab




The first tab of the system properties window is a simulation of the standard control unit itself. You will notice it has the same LED's as a Microvision Plus/HPQ, which function in exactly the same way as a real control unit. If the control panel is running but the RGAServer is not, the comms light will be slowly blinking. If the server is running the comms light will be solid green.

There are three options on this page.

- You can set the external trip open circuit. The effect of this depends on how the external trip has been configured.
- You can set the RF to fail above mass 30. If you do this you will see peaks below mass 30 but none above.
- You can set the filament emission to be 'bad' when filaments are switched on. Again the exact effect of this depends on how the unit has been configured.

The 'apply changes immediately' checkbox (which is checked by default) means that any change in the check state of the three main checkboxes is applied immediately. If this box is not checked you can make more than one change and no change will take effect until you click the 'Apply' button.


System Properties – The Digital I/O Tab



If the control unit that you have created has available digital I/O, you will have a digital I/O tab which will show the state of the 16 digital I/O lines. Once again the exact content of this page depends on how the control unit is configured. Hardware set-up can of course be used to change the I/O configuration, or it can be set when the virtual instrument is created, using the Options... button.

Inputs are shown green when TTL high. Outputs are shown red when TTL high. You can change the state of any input simply by clicking it. You cannot change outputs. You have the choice of 'applying the changes immediately' or making more than one change and then clicking the 'Apply' button.


System Properties – The Cirrus Tab



If the virtual system is a Cirrus, you will have a Cirrus tab. This shows the same LED's that are found on the front of the real Cirrus hardware. The only interactive 'input' that is available is the 'turbo at speed' setting. The behaviour of the filaments mirrors the real system. Filaments can only be on when the turbo is at speed and the system pressure is less than $5e-5$ mbar.

If the Cirrus is configured to have a multi-valve inlet, the current valve position will be displayed.

System Properties – The RVC Tab




If the system has an RVC, you will have an RVC tab, which is a full virtualisation of a real RVC. You will see exactly the same front panel as a real RVC and you have the ability to mimic the external input functionality using the user interface. Once again you can set the turbo 'at speed' state, the two general purpose digital inputs, the maintenance key switch and the settings for the pressure switch that governs the Status1 and Status2 lines. Like a real RVC, the Pump and Bake buttons can be used to toggle the state of the relevant item.

If you change the pressure switch set point pressures or the general purpose digital inputs you must click the 'Apply' button when you have the settings that you want.

In common with all virtual outputs, valves show red when open.

On a Vision 1000P system, which only has one valve, both the Status1 and Status2 lines are and-ed together and the two valve outputs work as one. This is equivalent to the standard behaviour.

System Properties – The Analog Input Tab




If the control unit has analog inputs, you can set their value on the analog input tab. In all cases of setting a numeric value such as the voltage where you have the option of entering the text or using the spin buttons, you make a change either by

- Editing the text and clicking Apply
- Editing the text and pressing the tab key to move to another control, then clicking Apply
- Using the spin button(s), then clicking Apply

If the control unit has a reserved analog input for a pressure gauge, you can choose whether the simulator automatically sets the value to the appropriate value for the particular gauge. A Cirrus has a gauge on Input 1, whereas an HPQ2S may have a gauge on Input 4. When you link the total pressure to the analog input in this way you cannot change that channel manually.

System Properties – The Vacuum Tab



This tab is in many ways the most interesting since it permits interactive changes to the data displayed by the virtual RGA. There are two choices.

The first choice is to pick one of a range of standard spectra at a pressure of your choice. Be aware that the pressure specified is the actual chamber pressure. If your system has a special inlet, the inlet factor will be applied to this value. If the pressure specified results in an individual partial pressure exceeding the maximum measurable pressure for the gain conditions in force at the time, the displayed peak height will be full scale. On the other hand, if the peak is so small that it cannot be measured, the peak will be 'in the noise'.

You can elect whether the simulation should add appropriate measurement 'noise' to the returned peak height.

Any change to the selected spectrum or pressure will take effect when you click the Apply button.

The information at the bottom of the window shows the actual sum of peaks currently measured – this will often differ just slightly from the requested value at the top.

When you are displaying a standard spectrum you can change the partial pressure of an individual mass using the 'Set mass to ...' and clicking the Set button. The new setting will apply for the remainder of the current 'replay scan' (whose length you can determine); it will then revert to the normal setting for that spectrum. The progress bar shows how much time remains for the scan. Some of the spectra are 'dynamic' which means that some constituents such as helium, water or air change on a random basis. Each time the progress bar starts over, a new composition will be applied.

You can use the Restart button to restart the scan prior to clicking the Set button to change an individual mass.

The duration of the standard replay scan is an option that can be set on the preferences menu.

The second choice is to 'replay' a vacuum profile from a file. This file can have one or many individual 'scans'. The file is typically generated using your favourite spreadsheet. If you want to 'replay' real live data from a Process Eye or EasyView installation, you use Recall to save your data as a VVP file (virtual vacuum profile). This can be read directly by the virtual RGA.

The format of the file is explained in appendix 1.

The Restart button is used to synchronise the data back to the start of the first scan.

A replay file can be 'real' or entirely artificial. A real replay file is one derived from real data where each 'scan' has a similar duration, for example every 4 or 5 seconds. Replaying this file would produce a trend view that looked very similar to the one stored in the original data file. An 'artificial' replay file is typically created by a recipe developer to test the behaviour of the recipe in alarm conditions. The scans in this replay file may have widely varying durations for reasons that will become apparent later.

Once again, you can change an individual mass peak using the Set button, but remember that the change will only last until the next replay file scan.

Please understand that it is not the replay scan that **drives** the scan in your software. There is no mechanism for 'asking for' a particular reading at a particular time. The replay scan is the provider of data values during that replay scan's lifetime, but there is **no synchronisation** between a scan in your software and the start of a replay scan. This means that if the data in the file was acquired at accuracy 5 but it is replayed at accuracy 3, the replay duration will remain the same but you will get more scans than in the original. Likewise, if it is replayed at accuracy 8, you will get fewer scans than the original. Where synchronisation is critical (usually only to developers rather than software evaluators) a synchronisation tool is provided and is described below.

If you are particularly interested in what is happening at low partial pressures, you will probably choose to check the 'Remove 'noise' from the spectra' checkbox because the simulator will be adding more noise to the noise already stored in the file.

More advanced information about the use of replay files is provided in a later section.

General Information about Partial Pressure Simulation

When the simulation generates a partial pressure reading, it takes note of a number of virtual RGA configuration settings. The mass alignment is modelled exactly and the resolution is also modelled, particularly on the over-resolved side.

Electron emission is modelled so that the peak height is proportional to the emission up to a maximum of 10% over the recommended value. Electron energy is modelled to the degree that mass 4 is dependent on it and mass 20 is very small if the 20/40 ratio looks like argon.

The multiplier detector has a recognised voltage / gain characteristic.

For these reasons, please be aware that you will only get the partial pressure reading recorded in the file or standard spectrum when

- The alignment and resolution settings are all set to 32767
- The filament emission current is set to the default value for the source setting
- The electron energy is at its appropriate default for the source setting.
- The multiplier gain/voltage matches a two-part graph of 30/-650 to 300/-750 or 300/-750 to 3000/-900.
- The sensitivity (amps per unit of pressure) for the active source setting is the default value for that setting.

If any of these items are not the case, the peak heights will not match the value stated on the Vacuum tab of the virtual RGA window.

Noise is simulated in three parts – ADC noise, ion source noise and detector noise.

Calibrating a Virtual RGA


It is perfectly possible to calibrate a virtual RGA using MKS' EasyView application or indeed any programmatic method of your own. The Vacuum tab of the virtual RGA window will show what is the current sum of the partial pressures, but you are free to calibrate the height of any peak to any value you choose.

For example, suppose you are displaying the argon spectrum, so that the simulation peak height of mass 40 is $1e^{-5}$ torr (as determined by the vacuum tab).

Now run a calibration on mass 40 but state that you want to calibrate mass 40 to $1e^{-6}$ torr. The calibration will run, and when you run a bar chart afterwards the height of mass 40 will indeed be $1e^{-6}$ torr and the sensitivity (amps per unit of pressure) of the analyser will be 10 times higher than the 'standard' value.

Modifying the Initial Configuration of a Virtual RGA

When you create a new virtual RGA using the 'Add' button, the 'Add New Instrument' dialog has an 'Options...' button which allows you to set a number of general options that will apply to all new systems that you create.



In the example shown here, the following defaults have been selected.

- Where possible, analysers are assumed to be 200 amu, and when an HPQ2S is created it will have a Granville Phillips ion gauge with Torr units attached.
- e-Vision control units will be set to the 40 eV option.
- The default trip actions are as shown.
- Where a control unit allows it, the digital I/O will be set to have the first 4 bits as inputs and the next 2 bits as outputs on port A, with no inputs or outputs on port B. General purpose digital outputs will be initialised TTL high. Where a Cirrus system is created, it will be configured to have 4 valves.
- When an RVC is created, the cool down parameters will be set so that there are only 10 seconds of bake before cool down will be enforced, and the cool down period will last for 30 seconds. The valve mode for V3 behaviour is the standard 1000C/E nitrogen purge. The alternatives are to select the 300mm ResistTorr behaviour of admitting a calibration gas under scripted recipe control, or to have the 1000P behaviour where V3 is not connected.

Things to do with a Virtual RGA

Once you have created your virtual RGA systems you might like to try some of the following using the EasyView recipe or installation.

- Set up your software so that it takes up the full height of your screen but leave space on the right hand side so that you can see the System Properties window for your installed system(s). (If the System Properties window is not visible, right-click on the Virtual RGA Control Panel icon and select the system from the pop-up menu.)
- Switch on the filament and check that the filament LED goes green on the virtual UI.
- Start a barchart scan and check that you see an air spectrum. Using the UI select the vacuum tab and change the system pressure, then click the Apply button – watch the peak heights change.
- Select a different spectrum from the drop-down list and click 'Apply'.
- On the Control Unit tab, check the 'X-trip' check box and watch the filaments turn off and the peaks disappear. Uncheck the box and manually turn the filaments back on.
- Check the 'RF fail above mass 30' check box and watch the spectrum.
- Do the same with the 'bad emission' check box. The behaviour depends on how the unit is configured.
- Use the user interface in the EasyView recipe to change the gain range and the detector. Notice how, as you use a higher gain you see peaks that you did not see at lower gains. Can you see any helium? What about Krypton? Any other unexpected peaks?
- Set up a spectrum of argon and see what you get on both 'standard' and 'low' electron energy.
- Now run an analog RGA spectrum and select the tuning tab. Experiment with the low and high mass resolution and alignment.
- Using the source settings tab, experiment with the filament emission current and see the effect on the peaks.
- See if you can observe any effect on the air spectrum of changes to electron energy.

Now create a Cirrus benchtop system. (See if you can create one with a 4-valve multi-inlet.)

- Run the EasyView recipe (or run the EasyView application) and select the Cirrus button in the recipe UI. Select the Cirrus tab in the Virtual System UI.

- Turn on the Cirrus pump – use the UI to get the turbo at speed. Notice that the gauge pressure readout is activated.
- Turn on the RGA filaments and run a barchart.
- Notice how the virtual UI LED's mimic exactly the ones shown on a real Cirrus unit.
- Turn on the capillary heater and/or the main heater.
- Can you find the effect of checking the 'Chamber hotter than 90C' checkbox? (Clue: run a scan using the multiplier).
- What happens if you turn off the pump?
- What happens if, after turning the pump back on, you increase the system pressure on the Vacuum tab of the virtual UI to about 6 atmospheres? (Actually make the Cirrus internal gauge read greater than $5e-5$ mbar or $3.75e-5$ torr). Don't forget to click the 'Apply' button when you make system pressure changes.
- Check that when you change the inlet valve that the virtual UI recognises the fact.

Now create a virtual system that uses an RVC – maybe a Vision 1000P or a Vision 1000C.

Start EasyView again and use this new virtual system.

- Check that you can start and stop the pumps using the EasyView UI or the virtual system UI. The virtual UI RVC front panel is laid out exactly like the real one.
- Experiment with the RVC key switch toggling between maintenance and manual / automatic mode.
- Do you know the steps required to get the valves open in automatic mode? You need to make sure that the system vacuum is low enough to get the status1 and status2 lines green. Then you need the filaments on and the Quad OK light to be green on the RVC. The behaviour is subtly different between a 1000P and a 1000C.
- Run this system and acquire a barchart. What happens to the spectrum when the valves close? Actually the answer depends on whether you created the virtual system with a purge or a cal gas or a direct connection in the virtual RGA UI options dialog.
- Try baking the chamber. What effect does this have on the scan? Check that if you leave the bake for long enough, the system goes into cooldown.

- Check that the general purpose digital inputs work correctly – you need to click ‘Apply’ for these.
- You can change the Status1 and Status2 set points on the virtual UI.

Back on the standard control unit with standard general purpose digital I/O

- Check that when you change an input on the virtual UI that you see the corresponding change in Process or EasyView.
- Try running diagnostics on the control unit.
- Try and Degas the ion source.

Close Process / EasyView and run RGA Device Manager. Right mouse click on the control unit and choose ‘Configure or Install....’. You will be asked for a password – any non-empty password works for a virtual RGA. You can change any behaviour that you wish to – then run the control unit again in EasyView to see the effect. A new inlet factor? A new trip behaviour?

Try calibrating your virtual control unit. If the UI Vacuum tab says that the pressure is $1e^{-6}$, try calibrating at a pressure of $1e^{-5}$. Then run the barchart again. Is the peak height $1e^{-5}$ or $1e^{-6}$?

Bring up the main Virtual RGA Control Panel window and click the ‘Link-down’ button. (The third one down). You can restore the link by clicking the button again. What happened?

If you think you have made a mess of these virtual control units, just send them to the trash-can (button two) and create some new ones!

If you stored data files using EasyView, you can Recall them using the Recall application.

Synchronising a Virtual System with a Full System Simulation

This is an advanced topic that will mainly be of interest to developers.

You have already seen that the scans of vacuum data in a replay file are not directly linked to your scanning mechanism. Some background concerning data acquisition using a real RGA may be helpful here. The control unit is (within limits which we will ignore here) free to acquire data asynchronously at a rate determined by the accuracy setting. The server (Windows or IP) continuously goes round the loop of asking the control unit if it has new data, processing that data and sending it to Your software and then asking the control unit again. Sometimes the control unit has no data at all (at high accuracies) because new values only come every one hundred milliseconds or so. On other occasions (low accuracy leak check, for example) the control unit may send back many items of data in one response.

The server pushes its calibrated data to your software, which is a passive recipient. Items of data are 'labelled' so you know where the new items belong.

In a virtual RGA the data server effectively talks to itself to get new data. This is implemented using a 100-millisecond timer. Each time this goes off, the server, in conjunction with the virtual UI, determines how many readings could have been acquired since the previous report, allowing for any unused time from the previous period. These readings are, of course partial pressures but also analog inputs from virtual gauges and other external 'hardware'. The partial pressures for the next masses in the scan are calculated from the currently applicable partial pressures in the replay scan. The number of readings reported depends on the accuracy and the zero and detector delays that are part of the hardware setup.

Recipe application developers will need much tighter synchronisation of the replay file to the other drivers of the system, such as SECS messages or even digital inputs to the control unit. You can take control of all the virtual system inputs programmatically, using the supplied ActiveX control. A simple Visual Basic 6 project is available that demonstrates the use of this control.

The Synchronization ActiveX Control

You can use this ActiveX control in a VB project or in a scripted recipe.

The object model contains two read only properties, one write-only property and several (write only) methods. None of the methods either cause or return any errors. If you invoke a method that cannot be executed, for example because the particular virtual instrument is off-line, the method does nothing.

You invoke the ActiveX control either by adding a project reference to 'MKS Virtual RGA Synchronization v1.0', then coding

```
Dim sync As VRGASYNCLib.VrgaSynchronize      `VB declaration  
  
Set sync = CreateObject("vrga.synchronize")
```

Or, for VB script

```
Dim sync  
  
Set sync = CreateObject("vrga.synchronize")
```

This is a brief description of the properties and methods

Properties

IsInstalled

returns 0 if the administrator has disabled virtual RGA's or nonzero if they are enabled. This property is read-only.

CanConnect (Name As String)

Returns non-zero if the virtual instrument named '**Name**' is powered up and connected to the server. Returns 0 if virtualRGA.exe is not running or the control unit is in a link-down state.

Note that if this property returns zero, none of the following methods will succeed, but they will not generate errors.

This property is read-only.

StandardScanDuration (TimeInSeconds as Long)

Sets the duration for standard spectra scans. Must be greater than 1. Does not have to be one of the standard preferences on the main menu.

This property is write-only.

Methods

AnalogInput

AnalogInput(Name As String, Port As Byte, Value As Single)
Sets analog input **port** on virtual instrument **name** to a value of **Value**.

Port is 0-based.

```
AnalogInput "MV Plus", 0, -5.00
```

sets analog input 1 to -5.00 volts

DigitalInput

DigitalInput(Name As String, Port As Byte, FromBit As Byte, ToBit As Byte, SetTo As Byte)

Sets digital port **Port** on virtual instrument **name** to a value specified by **FromBit**, **ToBit** and **SetTo**. **Port**, **FromBit** and **ToBit** are all 0-based. Unconnected inputs are ignored, as are all outputs.

```
DigitalInput "MV Plus", 1, 2, 4, 1
```

sets digital inputs as follows: PB2=1, PB3=0 and PB4=0

PeakHeight

PeakHeight(Name As String, Mass As Integer, Value As Single)

Sets the peak height of mass **Mass** to a value of **Value** Pascal on virtual instrument **Name**

```
PeakHeight "MV Plus", 28, 1e-4
```

sets the peak height of mass 28 on "MV Plus" to 1-e^4 Pa **for the duration of the current scan**.

Restart

Restart(Name as String)

Performs the equivalent of clicking the Restart button on the Vacuum page. Restarts the scan or file. Functionally equivalent to setting ScanNumber to 1.

```
Restart "MV Plus"
```

RVCGeneralPurposeDIIs

RVCGeneralPurposeDIIs(Name As String, SetTo As Byte)

Sets the two general purpose digital inputs on the RVC attached to virtual instrument **Name** to a bit pattern **SetTo**. **SetTo** has a value between 0 and 3 where bit 0 is the state of input 1 and bit 1 is the state of input 2.

```
RVCGeneralPurposeDIIs "MV Plus", 2
```

sets input 1 low and input 2 high

ScanNumber

```
ScanNumber(Name As String, ScanNumber As Long)
```

Sets the current scan number (1-based) being 'played' by virtual instrument **Name**.

```
ScanNumber "MV Plus", 19
```

sets the active scan to 19.

ScanTime

```
ScanTime(Name As String, HH As Integer, MM As Integer, SS As Integer)
```

Sets the current scan to the one relevant to the specified time. HH can be greater than 23.

```
ScanTime "MV Plus" 2, 30, 0
```

sets the elapsed time in the file to be 2 hours 30 minutes. Note that the specified time does not have to explicitly match any particular elapsed time in the replay file.

StandardSpectrum

```
StandardSpectrum(Name As String, SpectrumID As Integer, TP As Single)
```

Sets the active standard spectrum for virtual instrument **Name** to be the one corresponding to the 0-based list index **SpectrumID** at a total pressure **TP** Pascal.

```
StandardSpectrum "MV Plus", 2, 1.333e-3
```

sets the standard spectrum to air on virtual RGA "MV Plus" at a total pressure of 1e-5 torr.

TurboAtSpeed

```
TurboAtSpeed(Name As String, bAtSpeed As Byte)
```

Sets the turbo 'at speed' on virtual instrument **Name** if the **bAtSpeed** is non-zero, otherwise sets it to 'not at speed'. This method works equally well for an RVC system and a Cirrus system.

```
TurboAtSpeed "V1000P", 1
```

sets the turbo 'at speed'.

TurboPumpOn

TurboPumpOn (Name As String, bOn As Byte)

Sets the pump on or off on an RVC system (it has no effect on a Cirrus).

```
TurboPumpOn "V1000P", 0
```

sets the pump off on virtual instrument "V1000P".

You will remember that there is never any error when calling any of these functions, so on a real system without virtual RGA support installed they become do-nothing commands and can therefore safely exist in your shipping recipes. You just need to protect against failing to create the Synchronize object, which would occur if the recipe were used on a system running a version prior to v5.20.

```
Dim sync
Set Sync = Nothing
If App.Version >= 5.20 Then
    Set Sync = CreateObject("vrga.synchronize")
End If

If Not Sync Is Nothing Then Sync.TurboPumpOn "V1000P", 1
```

Appendix 1

Vacuum Profile 'Replay' File Format

Basic Format

The vacuum profile or 'replay' file can be either a csv file or a tab-delimited text file, either of which can be created directly by Excel. Alternatively you can use a Recall data file to create a .vvp (virtual vacuum profile) file for you.

The rules for the file are as follows.

The rules for the file are as follows.

1. You can have any number of lines of text at the top of the file for description and comments.
2. The first critical line in the file is a line that contains [UNITS] in column A. Column B must then contain one of pascal, torr, mbar or millitorr. None of the text in a vacuum profile file is case-sensitive.
3. The next critical line is one that contains [DATA] in column A. Columns B onwards then contain the masses that are going to have specified peak heights. These masses must be in ascending order but there can be gaps, i.e. you can specify 2, 14, 17, 18, 28, 32. You can specify any masses in the range 1 to 300. You can optionally prefix the numeric value with the word 'Mass' (without quotes).
4. Every succeeding row must contain partial pressure data. There can be no blank rows, nor rows with comment text, although you can write comments in columns beyond the highest mass. Each row has an elapsed time in column A using the hh:mm:ss format. Columns B onwards contain the partial peak height for the mass specified above in the units specified. You must specify the partial pressures in the format x.yyE-*nn*. It does not matter what the decimal separator is but there must be exactly 3 significant figures to the mantissa. The exponent does not have to have a + sign. The exponent can have 1, 2 or 3 digits with or without leading zeros. If you follow these rules you can work with any international settings.

When the file is 'played' the simulator executes the first row (scan 1) until its elapsed time is completed. The simulator then moves on to scan 2 until the total elapsed time exceeds the time in column A of scan 2. Finally, when the total elapsed time exceeds the column A entry of the last row, the whole thing is replayed again.

Examples

Here is an example of a tab delimited file

Any text here on multiple lines as a description
Of what the file contains. Then a line with [UNITS]

```
[UNITS]      mbar
[DATA]       14          28          29          32
0:00:15     1.12e-7     2.34e-5     5.67e-8     8.99e-6
0:00:30     1.12e-7     2.34e-5     5.67e-8     9.23e-6
0:00:45     1.12e-7     2.34e-5     5.67e-8     9.57e-6     a comment
0:01:00     1.12e-7     2.34e-5     5.67e-8     9.81e-6
```

or

```
[UNITS]      mbar
[DATA]       Mass 14     Mass 28     Mass 29     Mass 32
0:00:15     1.12e-7     2.34e-5     5.67e-8     8.99e-6
0:00:30     1.12e-7     2.34e-5     5.67e-8     9.23e-6
0:00:45     1.12e-7     2.34e-5     5.67e-8     9.57e-6     a comment
0:01:00     1.12e-7     2.34e-5     5.67e-8     9.81e-6
```

This plays a file for 1 minute showing a gradually increasing 32 level, while masses 14, 28 and 29 stay constant. After 1 minute the simulation repeats itself.

There is no limit to the number of 'scans' that your file contains – until the computer runs out of memory.

That was a simple example. A more complex example would be when you want to construct a replay file that models an entire semiconductor processing fab.

Here is how you might integrate with your fab simulation. Define the first hour of the replay period to correspond to idle. So from 0:00:00 to 1:00:00 define a scan profile that may change peak heights every 20 seconds for the first 20 minutes and then has a scan line that lasts until 1:00:00.

Now define that the period in the file from 1:00:00 to 1:30:00 is wafer processing and 1:30:00 to 2:00:00 is inter-wafer. Set up your desired profile for each of these, but only use scans every 2 or 3 seconds for the first minute, then add the 1:30 and 2:00 lines to use up all the remaining time.

Now in your VB simulation, when the time comes to idle, simply call

```
Sync.ScanTime MyChamberName, 0, 0, 0
```

or

```
Sync.Restart MyChamberName
```

And when you send a wafer start message to Process call

```
Sync.ScanTime MyChamberName, 1, 0, 0
```

And a wafer end message, call

```
Sync.ScanTime MyChamberName, 1, 30, 0
```

Then all you have to do is ensure that your idle simulation duration never normally exceeds 20 minutes and certainly never exceeds one hour (or your RGA data will run into your wafer processing profile). Likewise ensure that your wafer process and inter-wafer time last no more than one minute.

Random Behaviour

You can introduce random behaviour by setting a peak height to be negative. If you set a value for, say mass 18 to be $-1.00e^{-8}$ torr, the simulator will generate a new random number at the time that the scan becomes active and display a partial pressure somewhere between $1.00e^{-8}$ and $1.00e^{-13}$, ie somewhere evenly logarithmically distributed over 5 decades.

Different random numbers are generated for mass 4, mass 17/18, 14/16/28/32 and the remainder of the peaks. Over the course of a number of scans this means that sometimes helium will be critically high, whereas on other occasions air, water, nitrogen, oxygen or general contamination will be critically high. This allows you to evaluate many different alarm conditions.

If the ratio of 17/18 is between 0.2 and 0.3, 17 and 18 are considered to be water and get modified by the random water factor.

If the ratio of 32/28 is between 0.18 and 0.28 then 14/16/28 and 32 are air and get modified by the air random factor

Otherwise, if 16/32 is between 0.03 and 0.04 then 16 and 32 are Oxygen and get modified by the random oxygen factor

Otherwise, if 14/28 is between 0.046 and 0.056, then 14 and 28 are nitrogen and get modified by the random nitrogen factor

Otherwise 14, 16, 28 and 32 get modified by the general random factor.

Mass 4 is modified by the random helium factor and all other masses are modified by the general random factor.

Chapter 13 of this document described how, by using the ActiveX control, you can programmatically position the scan number or current elapsed time within the file mapping.

Using a Recall File

The steps in converting a Recall file to a Replay file are:

1. Open the file in Recall
2. Choose File | Save As... and be sure to select 'Virtual vacuum profile file (*.vvp)' from the File Type drop-down list. Click Save.
3. At the options dialog, the options on the left hand side are all pre-selected for you.
4. If there is more than one measurement listed on the right hand side, you must choose just one – usually the main barchart. You cannot use a user measurement or an analog scan but you can use a peak jump measurement. Recall will automatically sort peak jump measurements into ascending mass order.
5. Save the file.

Noise Behaviour

Finally – and this is only really of interest to those who enjoy 'splitting hairs' – can you see what happens to file noise when it is replayed with no simulator noise added? A typical barchart file has a number of masses where no ion peaks exist so the data for that column is a mixture of positive and negative values around 'zero'. Take for example a traditional ¼ inch rod RGA on Faraday at low electronic gain. The range of values might be, say, from $1.2e^{-9}$ torr to $-1.2e^{-9}$ torr, representing +/- 100 ADC units. Now you already know that positive values are displayed as they are, but negative values are displayed as somewhere between the full positive value and 100,000 times less. So, in this example, and bearing in mind the logarithmic distribution of random values, two out of five occurrences of $-1.2e^{-9}$ torr will actually be converted to a displayable positive reading. Likewise two out of fifty readings of $-1.2e^{-10}$ torr will be displayed as a positive value.

Alternatively, at the dialog that follows the 'File | Save As' menu selection in Recall, you can clip all the data values at a very low positive value. This should have the effect of making the replay spectrum the same as the original because the clip value will usually be less than 1 ADC unit since the clip value is actually $1e^{-16}$ mbar.

Appendix 2

Background Information

The ActiveX and Java libraries provide a wrapper around the communications between the computer and the data server or control unit. This enables you as a programmer to work at a higher level than the basic communications layer.

The communications between computer and server is over TCP/IP and is all based on human readable formatted plain ASCII text. The SDK contains full documentation on this communications protocol, but you can see it in action simply by opening a telnet-type connection to the server as follows.

In Windows open a DOS window and at the DOS prompt type the following, substituting the actual IP address that you wish to connect to.

```
telnet 127.0.0.1 10014
```

This is the syntax for making a connection to 127.0.0.1 using port 10014.

If you connect successfully you will see a response such as

```
MKSRGA Multi
  Protocol_Revision 1.3
  Min_Compatibility 1.2
```

Now you type any one of the ASCII commands (you will not see what you type echoed to the screen) and press ENTER. The response from the unit will be displayed.

So, for example, type

```
Sensors
```

And press the ENTER key. The response will be something like

```
Sensors OK
State SerialNumber Name
Ready LM76-00101001 "Friendly Name"
```

Successful responses have the command echoed back with OK, followed by a row of heading text followed by one or more rows of further data.

If you make a typing error, hit ENTER and start again. You cannot use the backspace key.

Type (substituting your instrument serial number)

```
Select lm76-00101001
```

And press ENTER. The response will be

```
Select OK
SerialNumber LM76-00101001
State Ready
```

Now type

Info

And press ENTER. A long response will be generated that details all the available information about the control unit.

You can obtain all this information even when someone else is connected to and controlling the RGA. It can provide a quick and easy way of confirming that your application is working with the correct information.