

## **Inline TTL Shutter**

**INLINE-TTL-S**



## **Installation and Operation Manual**

**Document Number 000-10000-090-02-0505**

**Offices:** **Ocean Optics, Inc.**  
830 Douglas Ave., Dunedin, FL, USA 34698  
Phone 727.733.2447  
Fax 727.733.3962  
8 a.m.– 8 p.m. (Mon-Thu), 8 a.m.– 6 p.m. (Fri) EST

**Ocean Optics B.V. (Europe)**  
Geograaf 24, 6921 EW DUIVEN, The Netherlands  
Phone 31-(0)26-3190500  
Fax 31-(0)26-3190505

---

**E-mail:** **Info@OceanOptics.com** (General sales inquiries)  
**Info@OceanOpticsBV.com** (European sales inquiries)  
**Orders@OceanOptics.com** (Questions about orders)  
**TechSupport@OceanOptics.com** (Technical support)



Ocean Optics offers the most comprehensive, innovative and high-quality line of modular spectroscopy tools in the world. Mikropack, a leading supplier of spectroscopy and thin film components, is an essential and valuable partner in this enterprise. We have partnered with Mikropack because they are committed to the same goals of innovation and quality that inspire us here at Ocean Optics. As always, Ocean Optics conducts its business in an open, honest and technically available fashion. We invite you to contact us at Ocean Optics, Inc. (see front cover for contact information) or Mikropack GmbH with any technical questions, comments, or applications inquiries. Mikropack GmbH can be contacted at the following location:

MIKROPACK GmbH  
Maybachstraße 11  
D-73760 Ostfildern  
Germany  
Tel.: +49 (0)711 34 16 96-0 • Fax.: +49 (0)711 34 16 96-85  
e-mail: [info@mikropack.de](mailto:info@mikropack.de)  
internet: [www.mikropack.de](http://www.mikropack.de)

**Copyright © 2001-2005 Ocean Optics, Inc.**

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission from Ocean Optics, Inc.

This manual is sold as part of an order and subject to the condition that it shall not, by way of trade or otherwise, be lent, re-sold, hired out or otherwise circulated without the prior consent of Ocean Optics, Inc. in any form of binding or cover other than that in which it is published.

**Trademarks**

Microsoft, Windows, Windows 95, Windows 98, Windows Me, Windows NT, Windows 2000, Windows XP and Excel are either registered trademarks or trademarks of Microsoft Corporation.

**Limit of Liability**

Every effort has been made to make this manual as complete and as accurate as possible, but no warranty or fitness is implied. The information provided is on an “as is” basis. Ocean Optics, Inc. shall have neither liability nor responsibility to any person or entity with respect to any loss or damages arising from the information contained in this manual

---

# Table of Contents

About This Manual .....	iii
Document Purpose and Intended Audience .....	iii
What's New in this Document.....	iii
Document Summary .....	iii
Product-Related Documentation.....	iii
Upgrades .....	iii
<b>Chapter 1: Setup .....</b>	<b>1</b>
Overview.....	1
Unpacking .....	2
Contents .....	2
<b>Chapter 2: INLINE-TTL-S Specifications.....</b>	<b>3</b>
Operating Environment.....	3
Specifications .....	3
Pinout Information .....	4
Pinout Diagram .....	4
<b>Index.....</b>	<b>5</b>



---

# About This Manual

## Document Purpose and Intended Audience

This document provides you with an installation section to get your system up and running.

## What's New in this Document

This version of the *Inline TTL Shutter Inline-TTL-S Installation and Operation Manual* adds our partnership agreement.

## Document Summary

Chapter	Description
Chapter 1: <a href="#">Setup</a>	Contains instructions for setting up your Inline TTL Shutter.
Chapter 2: <a href="#">INLINE-TTL-S Specifications</a>	Contains operating specifications and pinout information.

## Product-Related Documentation

You can access documentation for Ocean Optics products by visiting our website at <http://www.oceanoptics.com>. Select *Technical* → *Operating Instructions*, then choose the appropriate document from the available drop-down lists. Or, use the **Search by Model Number** field at the bottom of the web page.

You can also access operating instructions for Ocean Optics products on the *Software and Technical Resources* CD included with the system.

Engineering-level documentation is located on our website at *Technical* → *Engineering Docs*.

## Upgrades

Occasionally, you may find that you need Ocean Optics to make a change or an upgrade to your system. To facilitate these changes, you must first contact Customer Support and obtain a Return Merchandise Authorization (RMA) number. Please contact an Ocean Optics for specific instructions when returning a product.



---

# Chapter 1

# Setup

## Overview

INLINE-TTL-S Inline TTL Shutter is a small mechanical laser-cut shutter that is positioned between two UV/VIS-optical components. The shutter is controlled by a small circuit board that is powered by a 12 VDC signal. The circuitry in the INLINE-TTL-S monitors Pin 13 of the 15-pin D-Sub connector for a 5-volt TTL signal and triggers the shutter when the 5-volt signal is received.

The maximum frequency at which the INLINE-TTL-S can operate is 5Hz (5 shutter cycles per second).



**Inline TTL Shutter (15-pin input component and control circuitry not shown)**



**Inline TTL Shutter (bottom) with 15-pin input component and control circuitry (top)**

# Unpacking

## ► Procedure

1. Unpack your lamp assembly carefully. Although the deuterium lamp is rigidly mounted, dropping this instrument can cause permanent damage.
2. Inspect the outside of the instrument and make sure that there is no damage. Do not use the instrument if damage is present. Contact your dealer for repair or replacement information, if necessary.
3. Use this instrument in a clean laboratory environment (see [Operating Environment](#)).

# Contents

Your package should contain the following:

- ❑ Inline-TTL-S Inline TTL Shutter
- ❑ One IC-DB15-2 interface cable for shutter operation

Additionally, you will need an approved power supply for the LS-475, such as the Mikropack PS-12V/1.25A power supply.

## ► Procedure

To set up your Inline TTL Shutter,

1. Connect a UV or VIS optic to the SMA 905 Connectors on each end of the INLINE-TTL-S shutter.
2. Connect the 12 VDC power supply to the control box of the INLINE-TTL-S.
3. Connect the trigger signal for your application to Pin 13 of the 15-pin input connector on the INLINE-TTL-S control box (see [Pinout Information](#)).

You have now configured the INLINE-TTL-S for use. Ensure that during your acquisition process, the 5-volt signal reaches Pin 13 of the INLINE-TTL-S and that the shutter operates on command.



# INLINE-TTL-S Specifications

This section provides information on the operating environment and physical specifications of the INLINE-TTL-S. It also provides pinouts for the 15-pin connector.

## Operating Environment

The following table provides information on optimizing the operating environment of your INLINE-TTL-S.

Operating Environment	The INLINE-TTL-S Unit . . .
Moisture	Is designed for operation in dry rooms only.
Ventilation	Should be situated so that its location or position does not interfere with proper ventilation.
Heat	Should be situated away from any device that emits excessive heat.
Object and Liquid Entry	Should be positioned so that objects do not fall on top of the unit. Additionally, ensure that no liquids are spilled into the enclosure through openings.
Power Sources	Should be connected to an approved power supply, such as the Mikropack 12 VDC 1250mA analog regulated power supply (PS-12V/1.25A)

## Specifications


Specifications	Criteria
Spectral Range	UV-VIS
Shutter Input	TTL maximum 5 Hz
Power Requirements	12 VDC
Power Consumption	Maximum 100 mA
Weight	Approximately 600 g
Size	140 x 50 x 50 mm

# Pinout Information

The following table contains pinout information for the INLINE-TTL-S:

Pin	Description
1	na
2	na
3	na
4	na
5	na
6	na
7	na
8	na
9	na
10	Ground
11	na
12	na
13	TTL Signal
14	na
15	na
na = not applicable	

## Pinout Diagram



---

# Index

## D

document  
  audience, iii  
  purpose, iii  
  summary, iii

## O

operating environment, 3

## P

package contents, 2  
pinouts, 4  
product-related documentation, iii

## S

setup, 1  
specifications, 3

## U

unpacking procedure, 2  
upgrades, iii

## W

what's new, iii

