

## Optical Multiplexer

MPM-2000



## Installation and Operation Manual

Document Number 000-10000-120-02-0505

**Offices:**

**Ocean Optics, Inc.**

830 Douglas Ave., Dunedin, FL, USA 34698

Phone 727.733.2447

Fax 727.733.3962

8 a.m.– 8 p.m. (Mon-Thu), 8 a.m.– 6 p.m. (Fri) EST

**Ocean Optics B.V. (Europe)**

Geograaf 24, 6921 EW DUIVEN, The Netherlands

Phone 31-(0)26-3190500

Fax 31-(0)26-3190505

---

**E-mail:**    [Info@OceanOptics.com](mailto:Info@OceanOptics.com)                    (General sales inquiries)  
              [Info@OceanOpticsBV.com](mailto:Info@OceanOpticsBV.com)                (European sales inquiries)  
              [Orders@OceanOptics.com](mailto:Orders@OceanOptics.com)                (Questions about orders)  
              [TechSupport@OceanOptics.com](mailto:TechSupport@OceanOptics.com)        (Technical support)



Ocean Optics offers the most comprehensive, innovative and high-quality line of modular spectroscopy tools in the world. Mikropack, a leading supplier of spectroscopy and thin film components, is an essential and valuable partner in this enterprise. We have partnered with Mikropack because they are committed to the same goals of innovation and quality that inspire us here at Ocean Optics. As always, Ocean Optics conducts its business in an open, honest and technically available fashion. We invite you to contact us at Ocean Optics, Inc. (see front cover for contact information) or Mikropack GmbH with any technical questions, comments, or applications inquiries. Mikropack GmbH can be contacted at the following location:

MIKROPACK GmbH  
Maybachstraße 11  
D-73760 Ostfildern  
Germany  
Tel.: +49 (0)711 34 16 96-0 • Fax.: +49 (0)711 34 16 96-85  
e-mail: [info@mikropack.de](mailto:info@mikropack.de)  
internet: [www.mikropack.de](http://www.mikropack.de)

**Copyright © 2001-2005 Ocean Optics, Inc.**

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission from Ocean Optics, Inc.

This manual is sold as part of an order and subject to the condition that it shall not, by way of trade or otherwise, be lent, re-sold, hired out or otherwise circulated without the prior consent of Ocean Optics, Inc. in any form of binding or cover other than that in which it is published.

**Trademarks**

Microsoft, Windows, Windows 95, Windows 98, Windows Me, Windows NT, Windows 2000, Windows XP and Excel are either registered trademarks or trademarks of Microsoft Corporation.

**Limit of Liability**

Every effort has been made to make this manual as complete and as accurate as possible, but no warranty or fitness is implied. The information provided is on an “as is” basis. Ocean Optics, Inc. shall have neither liability nor responsibility to any person or entity with respect to any loss or damages arising from the information contained in this manual

---

# Important Safety Notices

1. Read all safety notices and operating instructions before operating this unit.
2. Inspect the item for transport damage before using the 24VDC power supply for the first time.
3. Adhere to all warning stickers on the unit and all warnings contained in this manual.

## Warranty

Mikropack GmbH warrants to the original user of this instrument that it shall be free of any defects resulting from faulty manufacture of this instrument for a period of 12 months from the original date of shipment.

**This instrument should not be used for any Clinical or Diagnostic purposes.** Data generated in these areas is not warranted in any way by Mikropack GmbH. Any defects covered by this Warranty shall be corrected either by repair or by replacement, as determined by Mikropack GmbH.

There are no warranties that extend beyond the description herein.

This Warranty is in lieu of, and excludes, any and all other warranties or representations expressed, implied, or statutory, including merchantability and fitness, as well as any and all other obligations or liabilities of Mikropack GmbH including, but not limited to, special or consequential damages. No person, firm, or corporation is authorized to assume for Mikropack GmbH. Any additional obligation or liability not expressed provided for herein except in writing duly executed by an officer of Mikropack GmbH:

MIKROPACK GmbH  
Maybachstraße 11  
D-73760 Ostfildern

Tel.: +49 (0)711 3428088 • Fax.: +49 (0)711 3428085  
e-mail: [info@mikropack.de](mailto:info@mikropack.de)

## Warranty Handling

### ► *Procedure*

Follow the procedure below to process a warranty claim:

1. Determine the problem or fault with your local distributor.
2. If a problem is evident, obtain an RMA number from your local distributor.

### **Important Safety Notices**

---

3. Send the equipment to the local distributor for repair. If the item is under warranty, shipping will be free-of-charge both ways.
4. Contact your distributor for repair and delivery time. If the item is out of warranty, your distributor will provide a repair cost to you. In this situation, the distributor will not proceed with the repair until you order it.

Your system will be shipped back to you free of charge with insurance (if under warranty).

---

# Table of Contents

About This Manual .....	iii
Document Purpose and Intended Audience .....	iii
What's New in this Document.....	iii
Document Summary .....	iii
Product-Related Documentation.....	iii
Upgrades .....	iii
<b>Chapter 1: Setup .....</b>	<b>1</b>
Overview.....	1
Unpacking the MPM-2000 .....	2
Contents .....	2
Typical System Configuration.....	3
<b>Chapter 2: MPM-2000 Specifications.....</b>	<b>5</b>
Operating Environment.....	5
MPM-2000 Components.....	6
Front Panel .....	6
Rear Panel.....	7
<b>Chapter 3: Operating Instructions .....</b>	<b>9</b>
Operating Software .....	9
Channel Mode .....	10
Program Mode.....	11
Small Mode.....	12
Main Toolbar.....	13
Menu Options .....	13
Options Dialog Box.....	14
Calibrating the MPM-2000.....	15
MPM-2000-Calibration Options Menu .....	16
Performing Calibration.....	18
MPM-2000.drv .....	19

**Index.....21**

---

# About This Manual

## Document Purpose and Intended Audience

This document provides you with an installation section to get your system up and running.

## What's New in this Document

This version of the *Optical Multiplexer MPM-2000 Installation and Operation Manual* adds our partnership agreement.

## Document Summary

Chapter	Description
Chapter 1: <a href="#">Setup</a>	Contains a list of package contents and unpacking instructions.
Chapter 2: <a href="#">MPM-2000 Specifications</a>	Contains operating environment specifications, as well as other physical details of the product.
Chapter 3: <a href="#">Operating Instructions</a>	Provides instructions for operating Windows-based software that comes with the MPM-2000.

## Product-Related Documentation

You can access documentation for Ocean Optics products by visiting our website at <http://www.oceanoptics.com>. Select *Technical* → *Operating Instructions*, then choose the appropriate document from the available drop-down lists. Or, use the **Search by Model Number** field at the bottom of the web page.

You can also access operating instructions for Ocean Optics products on the *Software and Technical Resources* CD included with the system.

Engineering-level documentation is located on our website at *Technical* → *Engineering Docs*.

## Upgrades

Occasionally, you may find that you need Ocean Optics to make a change or an upgrade to your system. To facilitate these changes, you must first contact Customer Support and obtain a Return Merchandise

## About This Manual

---

Authorization (RMA) number. Please contact an Ocean Optics Application Scientist for specific instructions when returning a product.



---

# Chapter 1

# Setup

## Overview

The MPM-2000 Fiber Optic Multiplexer is a compact instrument that can optically couple one input channel to 16 different output channels. The MPM-2000 consists of a precisely controlled DC-motor with encoder and end switch on a rotator block. The optical path couples via UV collimating lenses.

The Multiplexer features a CNC controller working in positioning mode, which eliminates step loss. The MPM-2000 is software controlled and interfaces to a PC through the RS-232 protocol. It automatically saves all calibrated positions in the memory of the controlling unit, and easy-to-use software facilitates calibration. The software enables full control of the switching order and delay time.

The MPM-2000 is ideal for applications in process industry, where you need to measure multiple locations with multiple probes but with a single spectrometer channel and/or light source. The MPM-2000 operates as a stand-alone unit and includes a sample program to guide you in setting up your application. The CD contains a driver for Windows, and you can configure the MPM-2000 in Ocean Optics OOIbase32 Platinum software.

The following sections provide instructions on unpacking and setting up your MPM-2000 Optical Multiplexer.



# Unpacking the MPM-2000

## ► *Procedure*

1. Unpack the Optical Multiplexer carefully. Dropping this instrument can cause permanent damage.
2. Inspect the outside of the instrument and make sure that there is no damage. Do not use the instrument if damage is present. Contact your dealer for repair or replacement information, if necessary.
3. Use this instrument in a clean laboratory environment.
4. Submit the Registration Card to for warranty and support purposes.

## Contents

Your MPM-2000 Optical Multiplexer package should contain the following:

- MPM-2000 Main System
- MPM-2000 24VDC Power Supply

---


### **Caution**

Before using the power supply of the MPM-2000 for the first time, inspect the item for transport damage. Be sure to adhere to all warnings on the unit and in this operational manual.

---

- Software CD (contains MPM-2000 software and Calibration software)
- PC-MPM-2000 Serial Cable

# Typical System Configuration





---

## Chapter 2

# MPM-2000 Specifications

This section provides information on the operating environment and physical controls of the MPM-2000.

## Operating Environment


The following table provides information on optimizing the operating environment of your MPM-2000.

Operating Environment	The MPM-2000 Unit . . .
Moisture	Is designed for operation in dry rooms only.
Ventilation	Should be situated so that its location or position does not interfere with proper ventilation.
Heat	Should be situated away from any device that emits excessive heat.
Object and Liquid Entry	Should be positioned so that objects do not fall on top of the unit. Additionally, ensure that no liquids are spilled into the enclosure through openings.
Power Sources	The unit should be connected to a power supply only of the type described in the operating instructions or as marked on the unit.

# MPM-2000 Components

The following sections describe the components located on the front and rear of the MPM-2000 unit.

## Front Panel




---

### Note

The front panel of both the 1 x 16 and 2 x 8 versions of the MPM-2000 are identical. Only the rear panel will differ (two inputs on the 2 x 8 version).

---

## Rear Panel








# Operating Instructions

## Operating Software

The MPM-2000 comes with easy-to-use Windows-based software that allows you to control the MPM-2000 Optical Multiplexer manually or create sequence programs (scripts) with preconfigured channel settings and delay times. The graphical software interface allows you to choose from three modes of operation:

- [\*Channel Mode\*](#)
- [\*Program Mode\*](#)
- [\*Small Mode\*](#)


## Channel Mode



Channel Mode allows you to change the active channel on the MPM-2000. You can click on any of the blue buttons associated with the SMA 905 Connectors on the front of the MPM-2000 to activate that channel.

Once you activate a channel, the button for that channel turns red until you activate another channel.

## Program Mode



Program Mode allows you to create or change sequence programs used by the MPM-2000. The Program Mode interface contains the following items:











Item	Description
Step Column	Identifies the individual step number in the sequence program
Channel	Identifies the channel that the MPM-2000 will switch to in this step
Time to wait	Identifies the delay time before the MPM-2000 proceeds to the next step
Time unit	Defines the unit of time for the Time to wait column
Continuous	When checked, instructs the software to run the sequence program a specified number of times (specified in the <i>no of runs</i> box)
No of runs	Specifies the number of times the software will execute the sequence program (a value of $-1$ instructs the software to run the program until manually stopped)

## Small Mode










Small Mode provides all functionality of the software in a condensed menu bar. When using another software package, the Always on Top feature of Small Mode eliminates the need to consistently switch back and forth between applications. This mode is only accessible via the **View** Menu.

The Small Mode screen contains the following options:

Item	Function
	Load an existing program sequence
	Activate the loaded program sequence
	Stop the currently running program sequence
	Enable continuous mode
	Specifies the number of times the program sequence will execute (if continuous mode is active)
	Increase the active channel on the MPM-2000 by one (+ 1 channel)
	Decrease the active channel on the MPM-2000 by one (- 1 channel)
	Specify a channel on the MPM-2000
	Switch to the channel specified in the text box to the left of this button
	Exit Small Mode and return to Channel Mode

## Main Toolbar

The Main Toolbar is only visible in Channel Mode and Program Mode. It contains the following buttons and features:

Button	Function
	Switch to Channel Mode
	Switch to Program Mode
	Open new program sequence
	Open an existing program sequence
	Save the current program sequence
	Start the currently loaded program sequence
	Stop the currently running program sequence

## Menu Options

The MPM-2000 software menu bar has four options:

- [File Menu](#)
- [View Menu](#)
- [Tools Menu](#)
- ? (help)

The sections that follow detail each of these options.

### File Menu

The File Menu contains the following options:

Menu Option	Description
New	Clear the program grid and open a new program sequence
Open	Open an existing program sequence
Save	Save the program sequence that you are currently editing
Save As	Save the program sequence that you are currently editing with a new filename
Options	Opens the Options dialog box
Exit	Exits the MPM-2000 software

## View Menu

The View Menu allows you to toggle between the three different application modes. You can only access Small Mode using this option.

## Tools Menu

The Tools Menu contains the following options:


Menu Option	Description
Connect	Connects to the MPM-2000
Disconnect	Disconnects from the MPM-2000
Options	Opens the Options dialog box

## Options Dialog Box

The Options dialog box is accessible from the File Menu or Tools Menu. It contains the following two tabs:

- [Comport Tab](#)
- [Settings Tab](#)


### Comport Tab



Select the COM port on the PC that you will use to communicate to the MPM-2000, and then click the **OK** button.

You must ensure that this option is correct before the MPM-2000 will communicate with the PC.

## Settings Tab



Click the **Always On Top** check box to position the software window above any other running applications, and then click the **OK** button to save your changes.


To specify the default mode that the MPM-2000 software will enter upon startup, select one of the three radio button options under the **Startup View** section of the **Settings** tab, and then click the **OK** button.

## Calibrating the MPM-2000

### Caution

Do not execute this software unless MPM-2000 calibration is required. If used when not needed, this software can alter the correct settings of the MPM-2000 and require you to service or re-calibrate the unit.









The MPM-2000 comes with software that allows you to calibrate the unit to ensure proper operation.



### 3: Operating Instructions

To launch the calibration software, you must execute the MPM-2000-Calibration.exe file from a command line or the **Start | Run** option from the Windows toolbar.

The MPM-2000-Calibration software contains the following options:

Option	Description
	Displays the serial number of the multiplexer
	Displays the Version of the MPM-2000-Calibration software
	Specifies the active channel number. Select a different channel from the drop-down menu to switch to that channel
	Displays the intensity of light (when used with an Ocean Optics spectrometer)
	Change the sampled wavelength position for the chosen channel. Each button represents large (<<< >>>), medium (<< >>), and small steps (< >), respectively.
	Saves the actual multiplexer position in the software
	Writes the saved multiplexer positions to the memory of the MPM-2000
	Exits the program (saved values will be lost unless saved)

## MPM-2000-Calibration Options Menu


The Options Menu contains the following items:

Option	Description
Use with OOI Spectrometer	Check this box if you want to use an Ocean Optics spectrometer for the calibration process
Spectrometer Settings	Select the type of spectrometer (if using an Ocean Optics spectrometer)
Acquisition Parameter	Change the acquisition parameters (if using an Ocean Optics spectrometer)
COMport	Select the COM port on which the MPM-2000 will communicate



## Acquisition Parameter Options

Select the Acquisition Parameter option from the Options menu to view the Acquisition Parameters Options dialog box:



Option	Description
Integration time:	Sets the integration time in OOIBase32 (when using an Ocean Optics spectrometer)
Boxcar width:	Sets the boxcar width in OOIBase32 (when using an Ocean Optics spectrometer)
Averages:	Specifies the number of scans that OOIBase32 will average (when using an Ocean Optics spectrometer)
Pixel no.:	Specifies the pixel number to sample (when using an Ocean Optics spectrometer)
Master	Applies the displayed options to the Master spectrometer channel (when selected)
Slave 1	Applies the displayed options to the Slave 1 spectrometer channel (when selected and available)
Slave 2	Applies the displayed options to the Slave 2 spectrometer channel (when selected and available)
Slave 3	Applies the displayed options to the Slave 3 spectrometer channel (when selected and available)

When using an Ocean Optics spectrometer, you can enter the acquisition parameters directly into this dialog box. Click the **OK** button when you have set the parameters you wish to use for your experiment.

## Comport Options

Select the Comport option from the Options menu to view the **Comport Options** dialog box:



Select the COM port you wish to use from the drop-down menu, and then click the **OK** button to enable the COM port change.

## Performing Calibration

### ► Procedure

Follow the steps below to perform a calibration on the MPM-2000:

1. Establish the connection between the MPM-2000 and the PC, ensuring that the MPM-2000 is connected to the power supply.

---

### Note:

Go to Step 5 if you are not using an Ocean Optics spectrometer

---

2. Open the MPM-2000-Calibration software.
3. Open the Options Menu and check the **Use with OOI Spectrometer** option.
4. Select the correct spectrometer and A/D Converter from the **Select Spectrometer** option.
5. Tweak the acquisition parameters to your environment, ensuring that the spectrometer is not saturated.
6. Activate Channel 1 on the MPM-2000.
7. Optimize the transmission by modifying the position of the multiplexer until the best transmission is obtained.
8. Save the position for Channel 1 in the software.
9. Activate the next channel, if needed. Otherwise, proceed to Step 11.
10. Repeat Steps 5-8 until all desired channels are optimized.

11. Click the **Write to MPM** button to save the data to the multiplexer. After writing the data to the MPM-2000, the calibration completes.
  
12. Close the program with **Exit** command from the Options menu.

## MPM-2000.drv

The MPM-2000.DRV is an interface DLL for use with Visual Basic or C/C++. It allows you to integrate the MPM-2000 in your own custom application.

The MPM-2000.DRV package contains a sample program written in Visual Basic® and a sample program written in Visual C++®.

For further details on how to use the MPM-2000.drv, consult the online documentation included with the MPM-2000.DRV.

### 3: Operating Instructions

---

---

# Index

## A

Acquisition Parameter Options, 17

## C

calibration, 15  
  Options menu, 16  
  performing, 18  
Channel Mode, 10  
components, 6  
  front panel, 6  
  rear panel, 7  
Comport Options, 18  
Comport Tab, 14  
configuration  
  typical system, 3

## D

document  
  audience, iii  
  purpose, iii  
  summary, iii

## F

File Menu, 13

## M

Main Toolbar, 13  
menu  
  File, 13  
  Tools, 14  
  View, 14  
Menu Options, 13  
mode  
  Channel, 10  
  Program, 11  
  Small, 12  
MPM-2000.drv, 19

## O

operating environment, 5  
operating instructions, 9  
options  
  Acquisition Parameter, 17  
  Comport, 18  
Options Dialog Box, 14

## P

package contents, 2  
product-related documentation, iii  
Program Mode, 11

## S

Settings Tab, 15  
setup, 1  
Small Mode, 12  
software, 9  
specifications, 5

## T

tab  
  Comport, 14  
  Settings, 15  
Tools Menu, 14

## U

unpacking procedure, 2  
upgrades, iii

## V

View Menu, 14

## W

warranty, A  
what's new, iii

