

# Collimation Tester Instructions

## Description

Use shear-plate collimation testers to examine and adjust the collimation of laser light, or to measure the wavefront curvature and divergence/convergence magnitude of large-radius optical components.


The setting up of a laser beam expander to give a collimated or parallel beam can be a problem if the full capability of any precision laser optics is to be achieved. Methods using autocollimation or measuring beam diameter over some distance are inadequate to obtain diffraction-limited performance unless tedious measurements are made. To accomplish this task, Ocean Optics has an interferometer that is extremely simple to use, requires almost no adjustment and produces results in a few moments.

## Method of Use

The method of use is so simple it may be hard to realize that the device is a precision interferometer.

### ► Procedure

1. Set up the beam expander with the lenses spaced according to their focal lengths.
2. Place the collimation tester in the beam so as to illuminate the front face as fully as possible, and to reflect the beam by about  $90^\circ$  to a screen placed at right angles to the reflected beam (see figure below).



## Collimation Tester Instructions

---

3. The maximum sensitivity of the collimation tester is achieved when the reflected angle is  $108^\circ$ , but sensitivity is reduced by only 1% at  $90^\circ$ , so the angle is not critical. Two overlapping images will be seen on the screen, one from each face of the collimation tester. Fringes (parallel patterns of light and dark) will be seen in the overlap region.
  4. Adjust the lens spacing until the fringes are parallel to the shadow of the cursor wire of the tester. The beam is now collimated.
- 

### Tips

If the fringes are not straight, check to see that the lenses of the beam expander are centered, square to the beam and with the proper face towards the collimated light (if these adjustments are appropriate). Any residual fringe curvature or wiggles indicates aberrations or errors that cannot be focused out.

If no fringes are seen in the initial setup, it may be that the beam is so far from collimation that the fringes are too fine to be seen. In that case, reduce the reflected beam angle to  $45^\circ$  or even  $30^\circ$ , and move the screen accordingly. This reduces the sensitivity of the tester. Make the initial adjustment for collimation and then return to  $90^\circ$  for the final tuning.

---

## USEFUL RANGE, APERTURE, WAVELENGTH

Each size of collimation tester is designed to have 5 – 6 horizontal fringes across the aperture when illuminated by a parallel beam. In theory, to have at least one fringe to set against, the smallest beam that can be tested is about 1/5 full aperture. In practice, it is found that the fringe slope changes from positive to negative when going through collimation and that the fringe spacing decreases away from collimation so that smaller aperture beams can be collimated by finding the center position between two defocused positions. When testing smaller apertures, it may also be necessary to reduce the reflected beam angle to get an appropriate amount of overlap of the two beams. The fringes only appear in the overlap area.

The choice of 5 – 6 fringes across the aperture is a trade-off between sensitivity and utility. For sensitivity, it is desirable to have only one full fringe over the aperture. The actual design allows an individual tester to be usable over a range of beam sizes of at least 5 to 1. If a collimation tester is to be used at a specific wavelength for a single aperture with critical demands on collimation, you can either select a larger aperture such that one fringe will be seen across your aperture or a custom design can be made for optimum sensitivity. Collimation testers are usable in wavelength over the transmission range of the glass, that is, from about 350 nm to about 2500 nm. To observe the fringes in the UV, a fluorescent screen may be necessary while in the IR, an image converter is suitable out to 1200 nm and an IR phosphor card can be used further out.

COLLIMATION Tester	Aperture Size (in millimeters)
CT-10	10
CT-20	20
CT-50	50

COLLIMATION Tester	Aperture Size (in millimeters)
CT-75	75
CT-100	100
CT-125	125
CT-150	150
CT-200	200

## Using the Collimation Tester for Quantative Measurements

Ocean Optics' collimation testers are shearing interferometers. There are many shearing interferometers but what they have in common is that the wavefront divided in amplitude, some shift applied to one or both portions, which are then recombined to form an interference pattern. In the collimation tester, the division arises from reflection at two faces of a slightly wedged plate. The shift is due to the displacement from the plate thickness and the recombination occurs in the overlapped images on the screen. Since the wavefront is being compared to itself rather than a reference flat wavefront, interpretation of the fringes differs from the usual interferometer where the departure from planarity is indicated by a fringe displacement. In the collimation tester, departure from planarity is indicated by fringe slope. A treatment of the complete interpretation and analysis of shearing interferograms is given by Rimmer(1).

## Wavefront Curvature

The collimation tester was designed as a null instrument to indicate when the incident beam is parallel. However, with large radius, spherical wavefronts, the fringes remain linear and can be analyzed to determine the wavefront curvature (2)(3). With two, separated, spherical wavefronts, the interference pattern can be treated in the same manner as the 2-pinhole experiment.

In the collimation tester, there are two separations:


- A vertical angle due to the wedge, and
- A horizontal angle proportional to the shear (beam displacement) and inversely proportional to the wavefront radius.

The figure below illustrates the measurements required,

Where:  $d$  = Fringe spacing perpendicular to the fringes

$A$  = Angle of the fringes with respect to the cursor

$S$  = Shear



The wavefront radius, at the observation screen, is given by the following equation, with the wavelength,  $\lambda$ , in the same units as S and d:

$$R = sd/\lambda \sin A$$

► **Procedure**

To find the wavefront radius at the lens plane,

1. Add the distance of the lens from the screen for a converging wavefront, or subtract the distance for a diverging wavefront.
2. Measure the values of shear, fringe spacing and angle.
3. Divide the distance between several fringes by the number of fringes to get the average spacing.

Where the shear is appreciable, it can be measured by placing a thin wire in the beam and measuring the separation of the two images. Shear can be calculated by the following relation,

Where:

I = Angle of incidence

N = Index of refraction

T = Plate thickness

$$S = \frac{T \sin 2I}{(N^2 - \sin^2 I)^{1/2}}$$


The index of the glass used is:

- 1.530 at 405 nm,
- 1.522 at 489 nm,
- 1.515 at 633 nm, and
- 1.507 at 1060 nm.

To find the divergence of the laser beam, divide the width of the lens or beam by the radius of the curvature.

## Component Testing

The collimation tester is suitable for measuring mirrors with a large radius of curvature. Murty (4) has examined various techniques and has found the shearing interferometer one of the more accurate methods. The setup is shown in the figure below:



The wavefront radius at the mirror is calculated by adding or subtracting the mirror screen distance to the radius measured at the screen. The mirror radius is twice the wavefront radius. The collimation tester can also be used to measure optical material homogeneity as suggested by De Vany (5). It provides a simple method of obtaining single pass effects without the complications of a Mach-Zehnder interferometer.

The collimation tester can be used to do more complex wavefront analysis including lens aberration testing. For these applications, the reader is referred to the classical literature on interferometers as well as the references cited.

### References

- (1) M.P. Rimmer, "Method For Evaluating Lateral Shearing Interferograms:?" Applied Optics, 13,3,1974. P 623
- (2) F.M. Dickey and T.M. Harder, "Shearing Plate Optical Alignment". Optical Engineering, 17,3, 1978,P 295
- (3) M.E. Riley and M.A. Gusinow, "Laser Beam Divergence Using A Lateral Shearing Interferometer". Applied Optics; 16. 10. 1977,P 2753
- (4) M.V.RK Murty and R.P. Shukla, "Measurement of Long Radius of Curvature", Optical Engineering, 22, 2,1983, P 231
- (5) A.S. De Vany, Applied Optics,10, 6. 1971, P 1459

