

---

# Engineering Note

Topic: Reprogramming a USB Device

Products Affected: ADC1000-USB, HR2000, HR4000, USB2000

Date Issued: 6/18/2002

## ► Procedure

Follow the steps below to reprogram the Ocean Optics USB Device:

1. Download and install the USB EEPROM Programmer utility from our web site at <http://www.OceanOptics.com/Technical/SoftwareDownloads.asp>. This utility installs the latest version of the microcode for all Ocean Optics devices.
2. Plug in the USB device to the PC.
3. Run the USB Programmer by selecting **Start | Programs | Ocean Optics | USB Programmer | USB Programmer**. This utility displays all of Ocean Optics USB devices that are plugged in the left-hand pane. Once a device is selected, its detailed information is displayed in the right-hand pane.
4. To reprogram a device, select the device in the left-hand pane. The current firmware version is displayed as the first entry in the right-hand pane.
5. Press the PROGRAM button and select the appropriate EEPROM update file (\*.iic file) for the Ocean Optics USB Device. The EEPROM update files are installed in the same location as the utility. The default location is C:\Program Files \ Ocean Optics \ USB Programmer.
6. The program immediately starts the reprogramming process. The estimated time for the reprogramming process is displayed in the bottom status bar.

---

### Caution

While programming, do NOT remove the cable, turn off the computer, or use the computer in any manner.

---

7. Once the programming process has completed, unplug the USB Device from the PC and then plug it back into the PC for the new program to take effect.

