

# Solo PIV

## Nd:YAG Laser Systems



### Reliable Operation

- Thermally compensated resonator assures stable operation.
- Rugged I-Beam resonator design on Solo 120XT and Solo 200XT
- Requires minimal maintenance, increasing system up-time.
- Field-proven reliability permits users to concentrate on their applications, rather than on system upkeep.

### Exceptional Performance

- Superior, proven design provides stable, high-energy output with excellent beam quality and pulse-to-pulse stability.
- Compact resonator design provides excellent beam pointing and energy stability.
- Predictable, high performance ensures that your work gets done faster.
- Solo 120XT and Solo 200XT design allows change over from 532 nm to 266 nm operation.

**S**olo PIV is a compact, dual laser-head system designed to provide a highly stable green light source for Particle Image Velocimetry (PIV) applications. It is ideally suited for most liquid and many air-based PIV experiments, and its small size provides excellent flexibility in setting-up such experiments.

### Features

- **Small laser head requires minimum space**
- **Single power supply simplifies setup and enhances mobility**
- **High output energy**
  - 15 - 200 mJ at 532 nm
- **Highly flexible design with repetition rates**
  - From 1 to 15, 30, or 50 Hz, depending on model selected
- **Operating convenience provided through multiple triggering capabilities**
  - Continuous internal trigger
  - External TTL trigger
  - Single input pulse activating laser lamp and Q-switch
  - Separate pulses to control lamp & Q-switch independently for precise laser pulse timing control
- **Easy set up:**
  - Single power supply features internal, closed-loop cooling system
  - Operates on 95-240 VAC single phase source
- **Convenient operation made possible with:**
  - Remote positioning of a single power supply - saves valuable lab space
  - Local control panel on power supply with all system controls, including optional optical attenuator
  - Hi/Lo power switch permits energy reduction during optics alignment



**NEW WAVE**<sup>™</sup>  
**R E S E A R C H**

## Solo PIV Products

	Solo I-15	Solo II-15	Solo II-30	Solo III-15	Solo IV-50	Solo 120XT	Solo 200XT
<b>Repetition Rate (Hz)</b>	15	15	30	15	50	15	15
<b>Energy<sup>1,2</sup> (mJ)</b>	532 nm	15	30	30	50	50	120
	355 nm <sup>5</sup>	NA	NA	NA	NA	NA	30
	266 nm	NA	NA	NA	NA	NA	20
<b>Energy Stability<sup>2</sup> (±%)</b>	532 nm	4	4	4	4	6	4
	355 nm <sup>5</sup>	NA	NA	NA	NA	NA	7
	266 nm	NA	NA	NA	NA	NA	9
<b>Beam Diameter (mm)</b>	3	3	3	4	4	5	6
<b>Pulse Width<sup>3</sup> (ns)</b>	3-5	3-5	3-5	3-5	3-5	3-5	3-5
<b>Divergence<sup>4</sup> (mrad)</b>	< 3	<3	<3	<4	<5	<3	<4
<b>Beam Pointing (urad)</b>	<100	<100	<100	<100	<200	<100	<100
<b>Jitter (±ns)<sup>6</sup></b>	1	1	1	1	1	1	1

- Optical losses due to optional attenuator will reduce maximum energy by 10%.
- Energy and pulse-to-pulse stability for 98% of shots after 30 minute warm up.
- Full width half maximum.

- Full angle for 86% of the energy, at 1/e<sup>2</sup> point.
- For single-head operation. Only one laser head may be optimized for 355 nm.
- From Q-switch synch out to light pulse, 98% of 1000 shots.

## Physical Characteristics

	Laser Head*			Power Supply		
	Solo I, II, III	Solo IV Solo 120XT	Solo 200XT	Solo I, II, III	Solo 120XT	Solo IV Solo 200XT
<b>Length</b>	13.775"/350 mm	22.4"/569mm	24.4"/620 mm	18.15" / 461 mm	19.0"/483 mm	21.2" / 538 mm
<b>Width</b>	7.0"/178 mm	9.24"/235 mm	9.24"/235 mm	7.77" / 194 mm	8.6"/218 mm	10.6" / 269 mm
<b>Height</b>	3.187"/81 mm	4.86"/123 mm	4.86"/123 mm	14.32" / 363 mm	15.0"/381 mm	15.16" / 385 mm
<b>Weight</b>	10 lbs./4.5 kg	34 lbs./15.5 kg	40 lbs./18.2 kg	48 lbs. / 22 kg	53 lbs/24 kg	53 lbs. / 24 kg
<b>Length Umbilical</b>	8 ft / 2.4 m	10 ft / 3 m	10ft / 3 m			

\* Width and height include mounting plate


## Operating Requirements

Temperature	50° - 86° F (10° - 30° C)	
Relative Humidity	20—80% non-condensing	
Voltage	95—240 V, 50/60 Hz	
Power	Solo I, II, III	15 Hz-800 watts; 30 Hz-900 watts
	Solo 120XT	1000 watts
	Solo IV, Solo 200XT	1500 watts

This document is for informational purposes only and does not set forth any warranty, expressed or implied, concerning any hardware and software, feature, or service offered or to be offered by New Wave™ Research. Specifications and product offering subject to change without notice.

## Warranty

One year warranty, details provided upon request.



**USA**  
New Wave™ Research, Inc.  
48660 Kato Road  
Fremont CA. 94538-7339  
Tel: 510-249-1550  
Tel: 800-566-1743  
Fax: 510-249-1551  
Email: [Lasers@new-wave.com](mailto:Lasers@new-wave.com)

**Japan**  
New Wave™ Research, KK  
5F Chojiya Building, 1-36-4,  
Shinjuku-ku, Shinjuku  
Tokyo, 160-0022 Japan  
Tel: +81-3-3351-0131  
Fax: +81-3-3351-0121  
Email: [NewWaveKK@new-wave.com](mailto:NewWaveKK@new-wave.com)

**Taiwan**  
New Wave™ Research G. C. Co., Ltd.  
No. 58-3, XingShan Rd.,  
Neihu Dist.,  
Taipei 114 Taiwan  
Tel: 886-2-8792-7585  
Fax: 886-2-8792-7584  
Email: [NewWaveGC@new-wave.com](mailto:NewWaveGC@new-wave.com)

**Europe**  
New Wave™ Research Co., Ltd.  
8 Avro Court, Emine Business Park  
Huntingdon, Cambs  
England PE29 6XS, UK  
Tel: 44-(0)1480 456566  
Fax: 44-(0)1480 456545  
Email: [NewWaveEU@new-wave.com](mailto:NewWaveEU@new-wave.com)

**China**  
New Wave™ Research (China) Co., Ltd.  
Rm. 1102, 11/F, Jiaxing Mansion,  
No. 877, Dongfang Rd.  
Shanghai 200122, China  
Tel: 86-21-5058-7785  
Fax: 86-21-5058-7786  
Email: [NewWaveCN@new-wave.com](mailto:NewWaveCN@new-wave.com)